



Daresbury Proteins

## Product description

**Name:** Mature BMP4 Human

**Synonyms:** Bone morphogenetic protein 4; Bone morphogenetic protein 2B; BMP2B

**Species:** Human

**Source:** HEK293

**Amino Acids:** 293-408

**Tag:** None

**Predicted Molecular Weight:** 13,128.92 Da (with glycosylation giving 25kDa)

**Protein ID:** P12644

**Sequence:**

**Mature BMP4 Amino acid Sequence**

```


      10      20      30      40      50      60
SPKHHSQRAR KKNKNCRRHS LYVDFSDVGW NDWIVAPPGY QAFYCHGDCP FPLADHLNST

      70      80      90     100     110
NHAIVQTLVN SVNSSIPKAC CVPTELSAIS MLYLDEYDKV VLKNYQEMVV EGCGCR
  
```

## Product specifications

**Estimated Molecular Weight, SDS-PAGE:** ≈ 25 kDa

**Grade & Purity:** >95% by SDS-PAGE stained with Instant Blue Stain (Expedeon).



1. Mw (NEB P7719) 5µl
2. BMP4 2µg
3. BMP4 4µg

175V for 45minutes (Biorad Mini-Protean TGX 4-20%, 10well 50ul)  
Gel samples were made 32.5ul sample + 6.5ul 6x reducing Laemmli

Daresbury Proteins Ltd. A company registered in England, UK. Company number 10835544.

Address: Daresbury Labs, Keckwick Lane, Warrington WA4 4AD, United Kingdom.

Web address: [www.daresburyproteins.co.uk](http://www.daresburyproteins.co.uk) Tel: +44 7398 623734 Email: [myprotein@daresburyproteins.co.uk](mailto:myprotein@daresburyproteins.co.uk)

**Amount sent:** 1 vial containing X $\mu$ l (X $\mu$ g) of Xmg/ml, SDS PAGE analysis and Nanodrop estimation at 280 nm

**Batch/Lot number:** #DSB-BMP4-XXXX

**Formulation:** 20mM histidine pH5.5, 0.2M NaCl, 0.6M arginine.

### Shipping

Product is shipped on wet or dry ice or frozen packs. Upon receipt keep refrigerated if for immediate use or frozen at -70<sup>0</sup>/80<sup>0</sup>C for medium to long-term storage.

### Product application and Storage

**Storage:** Store at 4<sup>0</sup>C for short term. -70<sup>0</sup>/80<sup>0</sup>C for medium and long term.


**Stability:** 3 weeks at 2<sup>0</sup>C to 8<sup>0</sup>C. Longer storage might be possible without loss of activity but has not been tested.

**Application Note:** For research purposes only. Not for use in humans.



### Activity data

#### I. Cardiomyocytes differentiation


##### Cardiac differentiation protocol (from Jeziorowska et al, 2015):



##### Embroid Body Differentiation




Still photos of the differentiated beating EBs composed of cardiac muscle cells expressing GFP marker post BMP4 treatment:




BMP4 (R&D Systems #314-BP, 25ng/ml), Daresbury Proteins manufactured BMP4 (11ng/ml)

## II. Trophoblast differentiation

### *Differentiation of Human Embryonic Stem Cells into Trophoblast Lineage – day 5*



**qRT-PCR analysis of trophoblast markers in MAN13 HESCs differentiated into trophoblast lineages with BMP4 –day 5**



Co1 = MAN13 HESCs in mTeSR media

CoTBd5 = MAN13 in mTeSR + BMP4 standard, R&D Systems #314-BP (10ng/ml)

dpTBd5 = MAN13 in mTeSR + BMP4, Daresbury Proteins (10ng/ml)

**Protocol references:**

**Comparison of Four Protocols for In Vitro Differentiation of Human Embryonic Stem Cells into Trophoblast Lineages by BMP4 and Dual Inhibition of Activin/Nodal and FGF2 Signaling**

Anvar Z, Chakchouk I, Sharif M, Mahadevan S, Su L, Anikar S, Naini FA, Utama AB, Van den Veyver IB.

Reprod Sci. 2024 Jan;31(1):173-189. doi: 10.1007/s43032-023-01334-5. Epub 2023 Sep 1. PMID: 37658178; PMCID: PMC10784360.

**BMP4-directed trophoblast differentiation of human embryonic stem cells is mediated through a  $\Delta$ Np63<sup>+</sup> cytotrophoblast stem cell state**

Li Y, Moretto-Zita M, Soncin F, Wakeland A, Wolfe L, Leon-Garcia S, Pandian R, Pizzo D, Cui L, Nazor K, Loring JF, Crum CP, Laurent LC, Parast MM.

Development. 2013 Oct;140(19):3965-76. doi: 10.1242/dev.092155. Epub 2013 Sep 4. PMID: 24004950; PMCID: PMC3775414.