Restek ProFLOW 6000 Electronic Flowmeter

Version 5.6
(Cat.# 22656)

Available Languages

Instruction Manual	English	Page 3
Gebrauchsanweisung	Deutsch	Seite 16
Notice d'utilisation	Français	Page 29
Manuale di istruzioni	Italiano	Pagina 42
Manual de instrucciones	Español	Página 55
使用说明书	Chinese	页 68
取扱説明書	Japanese	ページ 81

Table of Contents

1.0	Introduction	4
2.0	Specifications.....	5
3.0	Installing the batteries	5
4.0	Battery power consumption	6
4.1	Battery lifetime	6
4.2	Battery charge indicator	6
5.0	Operating instructions	6
6.0	Interpreting results	7
6.1	Flow range display	7
7.0	Data collection on the PC	8
8.0	ProFLOW 6000 flowmeter menu structure.....	10
8.1	Unit power up/power down.....	10
8.2	USB activation.....	10
8.3	Adjust LCD image backlight.....	10
8.4	Show battery charge indicator.....	11
8.5	Firmware version information	11
8.6	Adjust auto shutoff duration.....	12
8.7	Display calibration date	12
8.8	Display unit serial number	13
9.0	Troubleshooting	13
10.0	Product back label legend.....	14
11.0	Volumetric vs. mass flow measurements	15
12.0	Bubble flowmeters	15
13.0	Product case specifications	15
14.0	Calibration and service	15
	Declaration of Conformity	94

Restek ProFLOW 6000 Electronic Flowmeter

1.0 Introduction

Restek's ProFLOW 6000 flowmeter is specifically designed for use with gas chromatography (GC) systems. The probe is applied directly to the gas flow stream, and the measured flow rate is presented on the LCD screen. Units of flow are measured in mL/min.


This unit provides continuous real-time measurements of gas streams ranging from 0.50 mL/min to 500 mL/min. Because the technology uses volumetric flow measurement, the unit is compatible with all laboratory gases.

If this instrument is used in any manner other than described in the manual, the CE declaration is void.

CAUTION: Do NOT exceed maximum operating flow rates. Recalibration may be required if the unit has been subjected to extreme flow rates.

Always use appropriate laboratory safety practices when operating this device. Wear laboratory safety goggles when operating this unit.

Figure 1: Product diagrams

NOTE: There are no serviceable parts in this unit. Opening the device—other than to change the batteries—or tampering with the internal parts will void the factory warranty.

NOTE: To ensure accurate measurements and effective clearance of the flow gas from the unit, DO NOT obstruct the vent holes on the back of the unit.

NOTE: Modifying the length of tubing can cause a shift in calibration.

2.0 Specifications

Table I:

Type of measurement	Volumetric flow
Accuracy of measurements	Accuracy of ± 2.00% of flow reading or ± 0.200 mL/min, whichever is greater.
Power requirements	2 AA alkaline batteries 1.5 V DC each/3 V DC 200 ma
Operating flow range	0.50 to 500 mL/min
Operating temperature range	32 °-120 °F (0 °- 48 °C)
Available communication	USB type B data port
Warranty	One year (excludes recalibration)
Calibration	NIST traceable Yearly recalibration is recommended
Certification/Compliance	CE, Ex, Canadian ICES-003, WEEE, RoHS 3, China RoHS 2, UKCA (see section 10.0)

3.0 Installing the batteries

This unit uses two AA alkaline batteries.

To install batteries, extend the wire stand. Open the cover. Insert the batteries with the polarity (⊕ and ⊖) correctly aligned. Close the cover. (Figures 2 and 3.)

Precautions for battery replacement:

- Load the new batteries with their polarity (⊕ and ⊖) aligned correctly.
- **Do not use rechargeable batteries.**


Figure 2: Extend the wire stand before opening battery door.


Figure 3
Insert the batteries as marked.

4.0 Battery power consumption

4.1 Battery lifetime

The battery lifetime is dependent on the number of options the user has enabled.
The unit is shipped with the most energy-demanding options disabled (Table II).

The power saving functions can be changed.

- ➔ See Section 8.2: USB Activation
- ➔ See Section 8.3: Adjust LCD Image Backlight
- ➔ See Section 8.6: Adjust Auto Shutoff Duration

Table II: Default settings for the ProFLOW 6000 flowmeter

Auto shutoff duration	6 minutes
LCD backlight	0 (off)
USB port	disabled

NOTE: Store your ProFLOW 6000 flowmeter in its protective storage case following use. Do not store this manual or any other items on top of the ProFLOW 6000 flowmeter or the unit may turn on when the storage case lid is closed.

4.2 Battery charge indicator

The unit includes a battery charge indicator. Replace batteries as needed.

- ➔ See Section 8.4: Show Battery Charge Indicator

5.0 Operating instructions

 **CAUTION: Do not exceed maximum operating flow rates.**
Recalibration may be required if the unit has been subjected to extreme flow rates.

Connect the white probe end tip to the output of the gas flow line to be measured. Be sure the probe tip connection is completely sealed around the flow source outlet and is free of leaks (Figure 4). (White tubing ID 0.125"; OD 0.250")

Press and hold the  (Power) button until the unit responds with a regular clicking sound. The Pro-FLOW will immediately begin to provide flow measurements (Figure 5). Wait for the measured values to stabilize. It takes a few seconds for the unit to reach a steady state with the gas flow line.

To power down the unit, press and hold the  (Power) button until the unit stops clicking.

The unit is equipped with a timed auto shutoff option. (Default: 6 minutes.)

- ➔ See Section 8.6: Adjust Auto Shutoff Duration

Figure 4: Probe connected to a GC gas outlet.

Figure 5: LED displays the measured flow value.

6.0 Interpreting results

The unit has an operating range of 0.50 mL/min to 500 mL/min (Figure 7). If the flow is less than 0.50 mL/min, the display will read “under range.” If the flow exceeds 515 mL/min, the display will read “over range.” Excessively high flow rates may damage this unit.

NOTE: units of mL/min are equivalent to ccm.

6.1 Flow range display

The unit automatically adjusts the resolution of the display depending on the flow range being measured. Table III shows the resolution of the flow ranges.


Table III: Display resolution vs. flow range.

Flow range	Display resolution (mL/min)
0.50–9.99	0.01
10.0–99.9	0.1
100–500	1

Figure 6: Example flow value.

Figure 7: Description of flow ranges.

7.0 Data collection on the PC

WARNING: ONLY connect USB cable to USB port while unit is OFF.

The ProFLOW 6000 flowmeter provides you with a data stream of real-time flow values via the USB port (Figure 1, p. 4). In order to use this feature, you must first install the appropriate FTDI Virtual Com Port (VCP) Driver available at:

<http://www.ftdichip.com/Drivers/VCP.htm>


The VCP driver will cause the ProFLOW 6000 flowmeter to appear as a standard RS-232 port. This will work on any operating system for which there is an FTDI VCP driver. After installing the driver, connecting the device, and determining which port it creates, you can access the data stream through any programmatic means, or by using any serial terminal software.

For Windows systems:

To determine which port the ProFLOW 6000 flowmeter is using, go to the Control Panel and open System. Go to the Hardware tab and click the Device Manager button. Expand the Ports (COM & LPT) entry. Make sure the VCP driver is installed and then connect a powered ProFLOW 6000 flowmeter to the USB port. You will see the new COM port appear. Open your serial terminal.

➔ See Section 8.2: USB Activation

Figure 8: Screen capture of data collection.

If you do not already have serial terminal software, free, open-source options are available online (i.e., RealTerm, etc.). You can download RealTerm* software from:

<http://realterm.sourceforge.net/>

After installation, click the Port tab and set the following:


- Baud: 115200
- Port: The appropriate VCP for your ProFLOW 6000 flowmeter
- Parity: None
- Data Bits: 8 bits
- Stop Bits: 1 bit
- HardwareFlow Control: None

Go to the USB menu entry in the ProFLOW 6000 flowmeter and turn the transmission on. Finally, click **Open** on the Port tab in RealTerm and you will see the serial data stream begin in the terminal window (Figure 8).

If you would like to log the flow data, this can be done by clicking on the Capture tab. Set File to the name and location of the log file that you would like to save and click either the **Start: Overwrite** or **Start: Append** buttons appropriately (Figure 9).

*This software is not supplied or supported by Restek. User assumes all responsibility for the downloading and use of the program.

Figure 9: Screen capture of data collection.

8.0 ProFLOW 6000 flowmeter menu structure

8.1 Unit power up/power down

Press the  (Power) button:

The LCD screen will display the device intro screen:

Followed by calibrated date:

01/10/2012

Followed by measured flow data:

50.5 ml/min

To power off, press and hold the  (Power) button.

Other messages encountered at power up

under range

If the unit is hooked up to a flow stream with a flow rate less than 0.50 mL/min, the unit will report an “under range” status for the flow. This message will appear until the flow rate exceeds 0.50 mL/min.

8.2 USB activation

To enable the USB, press the  (Menu) button.

Use the   (arrow) keys to select the USB menu:

►USB

Press the  (menu) button again to enter the value select screen.

Use the   (arrow) keys to toggle between USB “on” and “off.”

To return to the main menu screen, press the  (menu) button.

To exit and return to measuring flow, press the  (menu) button again.

➔ See Section 7.0: Data Collection on the PC

8.3 Adjust LCD image backlight

Press the  (Menu) button.

Use the   (arrow) keys to select the Backlight menu:

►Backlight

Press the  (Menu) button again to enter the value select screen.

Use the   (arrow) keys to select the backlight value.

Backlight Values: 0 (off) and 5 (maximum).


To return to the main menu screen, press the  (Menu) button.

To exit and return to measuring flow, press the  (Menu) button again.

8.4 Show battery charge indicator

Press the  (Menu) button.

Use the  (arrow) keys to select the Batteries menu:

Press the  (Menu) button again.

The battery life is displayed.


To return to the main menu screen, press the  (Menu) button.

To exit and return to measuring flow, press the  (Menu) button again.

8.5 Firmware version information


Press the  (Menu) button.

Use the  (arrow) keys to select the Firmware menu:

Press the  (Menu) button again.

The most recent version of firmware is displayed.


To return to the main menu screen, press the  (Menu) button.

To exit and return to measuring flow, press the  (Menu) button again.

8.6 Adjust auto shutoff duration

To conserve battery life, the unit automatically turns off after 6 minutes. To customize the auto shutoff setting, press the  (Menu) button.

Use the   (arrow) keys to select the Power Time menu:

Press the  (Menu) button again to enter the value select screen.

Use the   (arrow) keys to select the auto shutoff setting.

Values: 1–59 minutes or “constant on” (max.)


To return to the main menu screen, press the  (Menu) button.

To exit and return to measuring flow, press the  (Menu) button again.


8.7 Calibration Date

To display the calibration date, press the  (Menu) button.

Use the   (arrow) keys to select the Cal Date menu:

Press the  (Menu) button again to display the calibration date.


To return to the main menu screen, press the  (Menu) button.

To exit and return to measuring flow, press the  (Menu) button again.


8.8 Serial Number

To display the flowmeter serial number, press the  (Menu) button.

Use the  (arrow) keys to select the Serial Num menu:

Press the  (menu) button again to display the serial number.

To return to the main menu screen, press the  (Menu) button.

To exit and return to measuring flow, press the  (Menu) button again.

9.0 Troubleshooting


Problem	Possible Cause(s)	Suggested Solution(s)
Multiple readings are not giving reproducible results	• Unit is out of calibration • Value is being compared to a bubble flowmeter	• Return the unit to Restek for recalibration* • See Section 12.0 for a discussion of the weaknesses of bubble flowmeters
Unit does not power up	• Dead batteries	• Replace with 2 new AA alkaline batteries
Flow value display is erratic/jumpy	• The ProFLOW 6000 flowmeter is very sensitive to small changes in flow	• Allow more time for flow to stabilize
Display showing under range.	• No flow or flow below 0.5 mL/min applied to the flowmeter. • Batteries have low power.	• Check that the correct flow below 500 mL/min and above 0.5 mL/min are established before applying the flowmeter. • Check the battery level following the instructions in section 8.4 and/or change the batteries.
Display constantly cycling through the startup sequence.	• Batteries have low power.	• Check the battery level following the instructions in section 8.4 and/or change the batteries.

*Contact returns@restek.com for return instructions for servicing a damaged unit. Additional charges may apply if the warranty has expired or the unit is damaged due to misuse.

Call Restek Customer Service at 1-800-356-1688 or 1-814-353-1300, ext. 3 (or your Restek representative) if you have any questions about this product or any other Restek product.

10.0 Product back label legend

#	Description
1	Serial number
2	Company name
3	This unit conforms to EU/EMC Directive 2004/108/EC; standards to which conformity is declared include 61326:1997 w/A3 Class A.
4	This Class A digital apparatus complies with Canadian ICES-003.
5	This complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
6	Country of origin
7	Product catalog number
8	Company address
9	This unit is WEEE compliant.
10	This unit is UKCA compliant.
11	This unit is RoHS 3 and China RoHS 2 compliant. 
12	
Ex nL	EN60079-0: 2006; Electrical apparatus for explosive gas atmospheres—Part 0: General Requirements. EN60079-15: 2005; Electrical apparatus for explosive gas atmospheres—Part 15: Construction, test, and marking of type of protection “nL” energy limited apparatus.
IIC	Group II applies to areas above ground environments. Gas Group IIC relates to hydrogen and related gas types.
3G	Category 3 relating to gas analysis; normal safety measure. Sufficient safety during normal operation. Normal operation described as measuring flows of flammable or explosive gases in a nonflammable environment.
T6	During testing, neither internal nor external elements exceed 85 °C.
X	Additional information: Operating range: 32 °F ≤ Tamb ≤ 120 °F 0 °C ≤ Tamb ≤ 48 °C Not intended for outdoor use or wet locations.
2010	Year of product design release.

11.0 Volumetric vs. mass flow measurements


The Restek ProFLOW 6000 flowmeter is a volumetric flow measurement device. Volumetric flow is the measurement of the volume of gas through a conveyance per quantity of time. Standard units of measure for this parameter are given in mL/min. The advantage of measuring volumetric flow is its independence to the composition of the flow gas. It is not necessary to correct the flow values based on the gas composition, as is required for mass flow devices.

Mass flow measures the weight of the gas flowing through the instrument per quantity of time. Mass flow units of measure are commonly g/sec.

12.0 Bubble flowmeter measurements

If you employ bubble flowmeters in your laboratory, you may find that they give slightly different flow rate values than the Restek ProFLOW 6000 flowmeter. This error is due to technology limitations inherent in the bubble flowmeter device, error from variances in air humidity within the bubble chamber, and its direct contribution to the measured flow rate. In the event a bubble flowmeter is used to measure flow gas where the gas is at elevated temperatures, the error due to humidity contributions can be extreme. For the most accurate measurement of laboratory gas flow rates, we recommend using the Restek ProFLOW 6000 flowmeter over bubble flowmeters.

13.0 Product case specifications

14.0 Calibration and service

The Restek ProFLOW 6000 flowmeter comes factory calibrated and carries a one-year warranty (excluding recalibration) from time of purchase. All units are calibrated to NIST traceable standards.

Recommended schedule for recalibration is once every year from time of purchase. Customers will need to return the unit to Restek for recalibration. At that time, preventative maintenance services can also be performed. A fee will be charged for recalibration and servicing of the unit. Prolonged failure to recalibrate the instrument may result in increased error.

Please have the serial number available when calling Restek with any concerns you may have.

Additional charges may apply if the warranty is expired or the damage is due to misuse. This manual is also available in electronic format at www.restek.com.

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MODEL / TYPE: Restek ProFLOW 6000 (22656)

DESCRIPTION: Electronic Flowmeter

SERIAL NUMBER: TBD

REPORTS: F2P30377B-01S, F2P30377B -C1-01E, F2P30377B-02RS

UK LEGISLATION: The Electrical Equipment (Safety) Regulations 2016, Electromagnetic Compatibility Regulations 2016, The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

STANDARDS CONSIDERED, FULL OR IN PART:

- EN 61010-1:2010+A1:2019
- EN IEC 61326-1:2021
- EN IEC 63000:2018

THIRD PARTY TESTING:

www.f2labs.com

26501 Ridge Road
Damascus, MD 20872 USA
Phone: 301.253.4500
Fax: 301.253.5179
Toll Free: 877.405.1580

16740 Peters Road
Middlefield, OH 44062 USA
Phone: 440.632.5541
Fax: 440.632.5542

8583 Zionsville Road
Indianapolis, IN 46268 USA
Phone: 317.610.0611
Fax: 317.610.0525

Restek Corporation declares under its sole responsibility that the **Electronic Flowmeter, model Restek ProFLOW 6000 (22656)**, is in conformity with The Electrical Equipment (Safety) Regulations 2016, Electromagnetic Compatibility Regulations 2016 and The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012.

Authorized by: Michael Straw

(signature)

Date: July 11, 2025

Name: Michael Straw

Title: Director Regulatory Compliance and EHS

Location: 110 Benner Circle, Bellefonte PA, 16823

**FC Supplier's Declaration of Conformity 47 CFR § 2.1077
Compliance Information**

PRODUCT IDENTIFICATION: Electronic Flowmeter, model Restek ProFLOW 6000

RESPONSIBLE PARTY (U.S. Contact Information required):

Restek Corporation
110 Benner Circle
Bellefonte, Pennsylvania 16823 USA

CONTACT NAME: Michael Straw

TELEPHONE NUMBER AND/OR E-MAIL ADDRESS: 814.353.1300 / regulatorycompliance@restek.com

FCC COMPLIANCE STATEMENT:

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

THIRD PARTY TESTING:

16740 Peters Road
Middlefield, OH 44062 USA
Phone: 440.632.5541
Fax: 440.632.5542
Toll Free: 877.405.1580

RESPONSIBLE PARTY SIGNATURE:

Michael Straw

(signature)

Date: July 11, 2025

Name: Michael Straw

Title: Director Regulatory Compliance and EHS

Location: 110 Benner Circle, Bellefonte PA, 16823

EU DECLARATION OF CONFORMITY

MANUFACTURER: Restek Corporation
110 Benner Circle
Bellefonte, Pennsylvania 16823 USA
Phone: 814.353.1300
URL: <https://www.restek.com/>

MODEL / TYPE: Restek ProFLOW 6000 (22656)

DESCRIPTION: Electronic Flowmeter

SERIAL NUMBER: TBD

REPORTS: F2P30377B-01S, F2P30377B-C1-01E, F2P30377B-02RS

DIRECTIVES: Low Voltage Directive 2014/35/EU, EMC Directive 2014/30/EU, RoHS Directive 2011/65/EU+(EU)2015/863

STANDARDS CONSIDERED, FULL OR IN PART:

- EN 61010-1:2010+A1:2019
- EN IEC 61326-1:2021
- EN IEC 63000:2018

THIRD PARTY TESTING:

www.f2labs.com

26501 Ridge Road
Damascus, MD 20872 USA
Phone: 301.253.4500
Fax: 301.253.5179
Toll Free: 877.405.1580

16740 Peters Road
Middlefield, OH 44062 USA
Phone: 440.632.5541
Fax: 440.632.5542

8583 Zionsville Road
Indianapolis, IN 46268 USA
Phone: 317.610.0611
Fax: 317.610.0525

Restek Corporation declares under its sole responsibility that the **Electronic Flowmeter, model Restek ProFLOW 6000 (22656)**, is in conformity with the Low Voltage Directive 2014/35/EU, the EMC Directive 2014/30/EU, and the RoHS Directive 2011/65/EU+(EU)2015/863.

Authorized by: Michael Straw

(signature)

Date: July 11, 2025

Name: Michael Straw

Title: Director Regulatory Compliance and EHS

Location: 110 Benner Circle, Bellefonte PA, 16823

CERTIFICATE OF TEST

Restek Corporation

110 Benner Circle
Bellefonte, Pennsylvania 16823 USA

Product Name: Electronic Flowmeter

Model(s): Restek ProFLOW 6000

Report No(s): F2P30377B-C1-01E, F2P30377B-C1-02E, F2P30377B-C1-03E,
F2P30377B-C1-04E, F2P30377B-C1-05E

The above-listed product was tested and found to comply with the following standard(s):*

Standard: EN IEC 61326-1:2021
Standard: Federal Register CFR 47, Part 15, subpart B, Class A
Standard: ICES-003, Issue 7:2020
Standard: VCCI-CISPR 32:2016
Standard: IEC 61326-1:2020
Standard: KS C IEC 61326-1:2021

A handwritten signature in black ink, appearing to read "Ben Coolbear".

Ben Coolbear
EMC Technical Manager
Issue Date: 2025-07-11

Page 1 of 1

**This Certificate is based on the evaluation of a representative sample(s) of the above-mentioned product(s) on a voluntary basis. This is to confirm that the tested sample(s) were tested to all or part of the standard(s) listed herein, as detailed in the applicable test report(s). It does not imply the assessment of the production of the product or approval of the manufacturing process or facility. The holder is authorized to use the Certificate to provide evidence of conformity as detailed above with the authority of F2 Labs.*

F2 Labs - Maryland
26501 Ridge Rd
Damascus, MD 20872
(301) 253-4500

F2 Labs - Ohio
16740 Peters Rd
Middlefield, OH 44062
(440) 632-5541

F2 Labs - Indiana
8583 Zionsville Rd
Indianapolis, IN 46268
(877) 405-1580

Notes:

Notes:

For information on Restek patents and trademarks, visit www.restek.com/patents-trademarks To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2025 Restek Corporation. All rights reserved.

www.restek.com

#204-04-002 Rev. date: 11/25