

The Applications of Empore™ Membrane SPE for LC mobile Phase Clean-up

May 6, 2020


CDS Analytical

Outline


1. Empore SPE Introduction
2. Application in LC Mobile Phase Clean-up

Empore™ Technology Highlights



- 3M™ patented SPE particle-loaded membrane




- Best Uniformity



- High Density




- Least Elution Volume



- Reduce Particle Fines

Empore™ Portfolio



• Empore™ by Package Style

- Disks
- Cartridges
- Plates
- StageTips

• Empore™ by Sorbent Style

Reverses Phases: C8, C18, UR, SDB-XC, Activated Carbon, Oil & Grease

Mixed Phases: MPC, UR, SDB-RPS

Ion Exchange: Cation, Anion, Chelating, SDB-RPS



Empore™ Unique Features



- Ultra-fast flow rate up to 700ml/min – **2X** faster than other leading brands
- Least elution volume – **1/3** of other leading brands and **1/10** of loose-packed SPE
- Highest consistency & reproducibility – **10-15% higher** than Waters Oasis plates
- Least particle fines in eluates -**1/10** of other leading brands to increase efficiency (reducing tube clogging & system downtime)

Recommended in Dozens of EPA Methods:

- 1664 (Rev. A) - N- Hexane Extrachable Material (HEM; Oil and Grease)**
- 506 - Phthalate and Adipate Esters in Drinking Water**
- 507 - Nitrogen- and Phosphorous-Containing Pesticides in Water**
- 508.1 - Chlorinated Pesticides, Herbicides, and Organohalides in Water**
- 512.2 - Chlorinated Acids in Water**
- 525.3 - Organic Compounds in Drinking Water**
- 549.1 - Diquat and Paraquat in Drinking Water**
- 550.1 - Polycyclic Aromatic Hydrocarbons in Drinking Water**
- 552.1 - Haloacetic Acids and Dalapon in Drinking Water**
- 553 - Benzidines and Nitrogen-Containing Pesticides in Water**
- 1613 (Rev. B) - Tetra-Through Octa- Chlorinated Dioxins and Furans by Isotope Dilution e.g. in Water**
- SW846 method 3535 – Test Methods for TCLP Leachates**
- QTM – Aqueous Phases Quick Turnaround Methods**
 - PAH**
 - Phenols**
 - Pesticides & PCBs**

Typical Customers:



- Environmental



- Food and Agricultural



- Pharmaceutical/Clinical



- Research



CDS has rebuilt Empore production line in a new-construct, GMP-compliant, clean room facility at Oxford, PA.




Empore production line at 3M




Empore new production line at CDS


Empore Disks with Improved Quality at CDS than at 3M



EPA 8141.B Organophosphorous Pesticides Recoveries with CDS Empore SDB-RPS Disk



Part 2. Application in LC Mobile Phase Clean-up



- **Ghost peaks or contaminant peaks sources:**
- airborne plasticizers particles, organics, colloids, ions, impure additives, dissolved gases, organic solvent impurities, accidental contamination from pH meters and laboratory glassware, plasticware, etc.

Comparison of Different Clean-up Methods

Mobile phase treatment	Advantages	Disadvantages
SPE cartridge (off-line)	remove both non-polar and polar contaminants	time-consuming, small volumes, channeling and cavitation-caused inconsistency. Subsequent contamination from airborne plasticizers, dust, and microbes
Empore disks (off-line)	fast, consistent, effectiveness	need to use combined disks to remove both polar and non-polar contaminants.
PDMS SBSE (off-line or in-line)	simple, continuously	small volume, special devices.
TFC Trap (SDB or C18, in-line)	good capacity	system contamination, inconvenient operation
DS micro-column (in-line)	compatible with HPLC and UPLC system	limited capacity, contaminant breakthrough
C18 semi-prep column (recirculating)	remove 95% contaminant. High capacity	extra system, high cost of prep columns

SBSE: Stir-Bar Sorptive Extraction; **TFC:** Turbulent Flow Chromatography

GSK Application Case: SDB-XC Disks for Mobile Phase Clean-up



4 Methods at Different pH Conditions

Table I: Summary of HPLC gradient methods used in investigation

Method	Aqueous Component	Gradient*				Column	Flow Rate (mL/min)	UV Detection Wavelength (nm)
		Mobile Phase A (v/v Aqueous–Acetonitrile)	Mobile Phase B (v/v Aqueous–Acetonitrile)	Time (min)	% B			
Method A	0.044% Trifluoroacetic acid (pH 2.3)	90:10	5:95	0; 3.4; 10	10; 13; 85	100 mm × 4.6 mm, 3.5- μ m d_p Waters Symmetry C18	1.5	226
Method B	0.05 M monobasic ammonium phosphate (pH 2.9)	70:30	22:78	0; 60	0; 100	250 × 4.6 mm, 5- μ m d_p Inertsil ODS2	1.0	228
Method C	0.05 M monobasic potassium phosphate (pH 7.1)	98:2	78:22	0; 10; 15	0; 100; 100	100 mm × 4.6 mm, 3- μ m d_p Hypersil HyPurity Elite C18	1.5	230
Method D	0.020 M ammonium acetate (pH 10.0)	95:5	5:95	0; 2; 32	0; 0; 100	100 mm × 4.6 mm, 3.5- μ m d_p Waters Xterra MS C18	2.0	282


* The gradient program could include high-organic solvent flush and reequilibration steps that are not shown.

Empore SDB-XC Disk has the best performance on all 4 testing conditions than other clean-up methods





Effect on the system suitability:

Empore disks can selectively remove typical mobile phase contaminants and don't affect the target analytes' retention times and resolutions.



- M. Ringo, J. Allen, and D. Mattocks, *LCGC North Am.* 21, 167-170 (2003).

Trifluoroacetic acid removal by Empore disks: passive adsorption. Will affect analytes' retention times and resolutions.



• M. Ringo, J. Allen, and D. Mattocks, *LCGC North Am.* 21, 167–170 (2003).

Method to Compensate TFA Removal by Empore SDB-XC Disk:

2 L of dilute trifluoroacetic acid (0.044%) might be necessary to equilibrate the SDB-XC disk before use for systems that are sensitive to trifluoroacetic acid concentration.

Biopharmaceutical Company Customers:



remdesivir



Empore™ C18, SDB-RPS and Carbon Disks for in-house LC-mobile phase clean-up.

Summary



Empore membrane technique is a powerful tool for LC mobile phase clean-up to improve HPLC method development efficiency.