

# FT-IR Purge Gas Generators

- ▲ Eliminates the need for costly, dangerous, inconvenient nitrogen cylinders in the laboratory
- ▲ Compact design frees up valuable laboratory floor space
- ▲ Improves signal-to-noise ratio even on non-purge systems
- ▲ Increases FT-IR sample thru-put and maximizes up-time
- ▲ Recommended and used by all major FT-IR manufacturers



Models 75-52NA, 75-62NA, and 75-45NA

**The Parker Balston® FT-IR Purge Gas Generator** is specifically designed for use with FT-IR Spectrometers to provide a purified purge and air bearing gas from compressed air. The generator supplies carbon dioxide-free air at less than -100°F (-73°C) dew point with no suspended impurities larger than 0.01 µm. The unit is designed to operate continuously 24 hours/day, 7 days/week. The Parker Balston Purge Gas Generator completely eliminates the inconvenience and the high costs of nitrogen cylinders and dewars, and significantly reduces the costs of operating FT-IR instrumentation. The Parker Balston unit offers cleaner back-

ground spectra in a shorter period of time and more accurate analysis by improving the signal-to-noise ratio. The typical payback period is less than one year. The generator is also ideally suited for use with CO<sub>2</sub> Analyzers and Matrix GC's in addition to supplying gas to other laboratory instruments.

The generators are quiet, reliable, and easy to install - simply attach the inlet and outlet air lines (at least 60 psig and 1/4 inch pipe), plug the power cord into a wall outlet, and enjoy trouble-free operation.

#### Here's what your colleagues say:


"A Parker Balston® FT-IR Purge Gas Generator and Self Contained Lab Gas Generator were used in conjunction with the Society for Applied Spectroscopy Fourier Transform Infrared Spectrometry Workshop at the University of Georgia, June 2000 (organized by Dr. James A de Haseth and Dr. Peter R. Griffiths). The Self-Contained Lab Gas Generator provided excellent purge for six spectrometers. The organizers were so pleased with the performance of the Parker Balston® systems, they have requested that Parker Hannifin Corporation, Inc. participate in future workshops."

- Dr. James A. de Haseth and  
Dr. Peter R. Griffiths


# FT-IR Purge Gas Generators

## Comparative Spectral Analysis in Purging an FT-IR Sample Chamber

### No Purge Gas



### 2 Minutes Parker Balston®



This spectra comparison illustrates that a Parker Balston® FT-IR Purge Gas Generator allows an FT-IR to be purged at a much higher flow rate than is practical with nitrogen.

The sample chamber purged by the Parker Balston unit is free of CO<sub>2</sub> and water faster than the sample chamber purged by nitrogen.

## Principal Specifications

|                                      |         |  |
|--------------------------------------|---------|--|
| Flow Rate for Specified Dew Point    |         |  |
| Inlet Pressure 125 psig              | 75-45NA | 36 scfh (17 lpm)                           |
| Inlet Pressure 60 psig               |         | 18 scfh (9 lpm)                            |
| Inlet Pressure 125 psig              | 75-52NA | 72 scfh (34 lpm)                           |
| Inlet Pressure 60 psig               |         | 36 scfh (17 lpm)                           |
| Inlet Pressure 125 psig              | 75-62NA | 216 scfh (102 lpm)                         |
| Inlet Pressure 60 psig               |         | 120 scfh (57 lpm)                          |
| CO <sub>2</sub> Concentration        |         | < 1 ppm                                    |
| Dew Point                            |         | -100°F (-73°C)                             |
| Min/Max Inlet Air Pressure           |         | 60 psig/125 psig                           |
| Max Inlet Air Temperature (1)        |         | 78°F (25°C)                                |
| Air Consumption for regeneration (2) | 75-45NA | 30 scfh (14 lpm)                           |
|                                      | 75-52NA | 60 scfh (28 lpm)                           |
|                                      | 75-62NA | 120 scfh (57 lpm)                          |
| Inlet/Outlet Port Size               |         | 1/4" NPT (female)                          |
| Electrical Requirements (3)          |         | 120 VAC/60 Hz/10 watts                     |
| Dimensions                           | 75-45NA | 7" w x 13" h x 6" d (18cm x 33cm x 15cm)   |
|                                      | 75-52NA | 13" w x 28" h x 9" d (32cm x 71cm x 23cm)  |
|                                      | 75-62NA | 13" w x 42" h x 9" d (32cm x 102cm x 23cm) |
| Shipping Weight                      | 75-45NA | 25 lbs (11 kg)                             |
|                                      | 75-52NA | 40 lbs (20 kg)                             |
|                                      | 75-62NA | 80 lbs (36 kg)                             |

## Ordering Information for assistance, call 800-343-4048, 8 to 5 Eastern Time

| Description                             |                                 | Model Number              |
|---|---------------------------------|---------------------------|
| FT-IR Purge Gas Generator               |                                 | 75-45NA, 75-52NA, 75-62NA |
| Annual Maintenance Kit                  | 75-45NA                         | MK7505                    |
|   | 75-52NA                         | MK7552                    |
|   | 75-62NA                         | MK7520                    |
| Installation Kit for All Models         |                                 | IK7572                    |
| Preventative Maintenance Plan           | 75-45NA                         | 75-45-PM                  |
|   | 75-52NA                         | 75-52-PM                  |
|   | 75-62NA                         | 75-62-PM                  |
| Extended Support with 24 Month Warranty | 75-45-DN2, 75-52-DN2, 75-62-DN2 |                           |

### Notes

- 1 Outlet dew point will increase at higher inlet compressed air temperatures.
- 2 Total air consumption = regeneration flow + flow demand.
- 3 Refer to voltage appendix for electrical and plug configurations for outside North America.

# Self-Contained FT-IR Purge Gas Generator

- ▲ Less expensive and more convenient than nitrogen cylinders and dewars
- ▲ Includes state-of-the-art, oil-less compressor
- ▲ Compact, portable design is ideal for mobile labs
- ▲ Improves signal-to-noise ratio even on non-purge systems
- ▲ Increases FT-IR sample thru-put and maximizes up-time
- ▲ Special sound insulation design ensures quiet operation

**The Parker Balston® Model 74-5041NA FT-IR Purge Gas Generator** is specifically designed for use with FT-IR spectrometers to provide a purified purge and air bearing gas supply from compressed air. The Parker Balston model 74-5041NA provides instruments with CO<sub>2</sub>-free compressed air at less than -100°F (-73°C) dew point with no suspended impurities larger than 0.01 micron 24 hours/day, 7 days/week. The Parker Balston Self-Contained FT-IR Purge Gas Generator completely eliminates the inconvenience and the high costs of nitrogen cylinders and Dewars, and significantly reduces the costs of operating FT-IR instruments.

The Parker Balston unit generates cleaner background spectra in a shorter period of time and more accurate analysis by improving the signal-to-noise ratio. The typical payback period is less than one year.

The generator is quiet, very reliable, and easy to install - simply attach the outlet air line, plug the electrical cord into a wall outlet, and the unit is ready for trouble-free operation.




Model 74-5041NA

#### Here's what your colleagues say:

"A Parker Balston® FT-IR Purge Gas Generator and Self Contained Lab Gas Generator were used in conjunction with the Society for Applied Spectroscopy Fourier Transform Infrared Spectrometry Workshop at the University of Georgia, June 2000 (organized by Dr. James A de Haseth and Dr. Peter R. Griffiths). The Self-Contained Lab Gas Generator provided excellent purge for six spectrometers. The organizers were so pleased with the performance of the Parker Balston® systems, they have requested that Parker Hannifin Corporation, Inc. participate in future workshops."

- Dr. James A. de Haseth and  
Dr. Peter R. Griffiths

# Self-Contained FT-IR Purge Gas Generator



This spectra comparison illustrates that a Parker Balston FT-IR Purge Gas Generator allows an FT-IR to be purged at a much higher flow rate than is practical with nitrogen. The sample chamber purged by the Parker Balston unit is free of CO<sub>2</sub> and water faster than the sample chamber purged by nitrogen.

## Principal Specifications

|   |   |
|---|---|
| <b>Self-Contained FT-IR Purge Gas Generator</b> | <b>74-5041NA</b>                              |
| Maximum Flow Rate (at 80 psig)                  | 60 SCFH (28 lpm)                              |
| Maximum Output Pressure                         | 80 psig                                       |
| CO <sub>2</sub> Concentration                   | < 1 ppm                                       |
| Dew Point                                       | -100°F (-73°C)                                |
| Outlet Port Size                                | 1/4" NPT (female)                             |
| Min/Max Ambient Temperature                     | 30°F/90°F (-1°C/32°C)                         |
| Electrical Requirements (single phase)          | 120 VAC/60 Hz, 20 amps (1)                    |
| Compressor                                      | 3/4 hp  |
| Dimensions                                      | 18"w x 31"h x 32"d<br>(46 cm x 76 cm x 81 cm) |
| Shipping Weight                                 | 250 lbs. (114 kg)                             |

(1) Refer to voltage appendix for electrical and plug configurations for outside North America.

## Ordering Information for assistance, call 800-343-4048, 8 to 5 Eastern Time

| Description                             | Model Number |
|---|--------------|
| FT-IR Purge Gas Generator               | 74-5041NA    |
| Annual Maintenance Kit                  | 74065        |
| Replacement Compressor                  | 74156        |
| Preventative Maintenance Plan           | 74-5041-PM   |
| Extended Support with 24 Month Warranty | 74-5041-DN2  |

# Ultra Dry Gas Generator

- ▲ Supplies ultra-dry, purified compressed air to NMR Spectrometers and other analytical instruments
- ▲ Ideal gas supply for spindle and automatic sample changer
- ▲ Completely eliminates costly, inconvenient nitrogen dewars - never pay for or change out another dewar
- ▲ Compact design frees up valuable laboratory floor space
- ▲ Completely automatic - plug it in and forget about it



Model UDA-300NA

**The Parker Balston® Model UDA-300NA Compressed Air Dryer** provides ultra-dry, purified compressed air to analytical instruments. The model UDA-300 reduces the dewpoint to -100°F (-73°C) without operator attention.

Each system is delivered complete, and ready for easy installation. A high efficiency prefiltration system, automatic drains, a 0.01µm final filter, a moisture indicator, and pretested controls are integral to the design of each dryer.

To install, simply connect your house compressed air supply (at least 60 psig and 1/4 inch pipe) to the dryer inlet port, and connect the dryer outlet port to your instruments. Plug the electrical cord into a wall outlet - no electrician required - and the unit is ready for trouble-free operation.

Designed specifically for NMR instrumentation, the generator is completely automatic, and virtually maintenance free. It is ideal for injecting, spinning, and lifting opera-

tions. It is recommended by major NMR instrument manufacturers and is currently installed in several thousand locations.

## Principal Specifications

### Model UDA-300NA Compressed Air Dryer

|                               |   |
|-------------------------------|---|
| Dew Point                     | -100°F (-73°C)                          |
| Flow Rate at 60 psig          | 390 scfh (184 lpm)                      |
| Flow Rate at 125 psig         | 720 scfh (340 lpm)                      |
| Min/Max Inlet Air Pressure    | 60 psig/125 psig                        |
| Max Inlet Air Temperature (1) | 78°F (25°C)                             |
| Inlet/Outlet Port Size        | 1/4" NPT (female)                       |
| Electrical Requirements (2)   | 120 VAC/60 Hz, 10 Watts                 |
| Dimensions                    | 41"h x 15"w x 8"d (104cm x 38cm x 20cm) |
| Shipping Weight               | 50 lbs (23 kg)                          |

## Ordering Information for assistance, call 800-343-4048, 8 to 5 Eastern Time

| Description                             | Model Number |
|---|--------------|
| Compressed Air Dryer                    | UDA-300NA    |
| Inlet Pressure Regulator                | 72-130-V883  |
| Annual Maintenance Kit                  | MK7525       |
| Preventative Maintenance Plan           | UDA-300-PM   |
| Extended Support with 24 Month Warranty | UDA-300-DN2  |

Notes:

- 1 Outlet dew point will increase at higher inlet compressed air temperatures
- 2 (1) Refer to voltage appendix for electrical and plug configurations for outside North America.

# ICP Spectrometer Nitrogen Generator

- ▲ Produces a continuous supply of ultra high purity nitrogen gas from existing compressed air
- ▲ Eliminates the need for costly, dangerous, inconvenient nitrogen cylinders or dewars in the laboratory
- ▲ Extends ICP Analysis into far-UV range below 170 (nm)
- ▲ Compact design frees up valuable laboratory floor space
- ▲ Offers long term cost stability - uncontrollable vendor price increases, contract negotiations, long term commitments and tank rentals are no longer a concern



Model 76-98NA  
Nitrogen Generator

**The Parker Balston® 76-97NA and 76-98NA UHP Nitrogen Generators** can produce 5-12 lpm of ultra high purity nitrogen gas. These systems are completely engineered to transform standard compressed air into 99.9999% of 99.995% pure nitrogen, exceeding the specification of UHP cylinder gas and dewars. Nitrogen is produced by utilizing a combination of state-of-the-art purification technologies and high efficiency filtration. Pressure swing absorption is utilized for the removal of O<sub>2</sub>, CO<sub>2</sub>, and water vapor. A catalyst module is incorporated in the 76-98NA to oxidize hydrocarbons from the inlet air supply. The generators also have a combination of high efficiency prefilters and a

0.01 micron (absolute) membrane filter incorporated into their design. The Parker Balston UHP Nitrogen Generators are engineered and packaged in a laboratory cabinet to fit nearly any laboratory. The systems eliminate the needs for costly, inconvenient high pressure nitrogen cylinders or dewars. The 76-97NA and 76-98NA are ideal for ICP Purge gas applications.

## Applications

Other applications include high flow GC carrier gas needs, DNA Synthesis and Sequencing Equipment, Mocon Moisture Analyzers, Circular Dichroism and Gel Permeation needs.

# ICP Spectrometer Nitrogen Generator

## Flow Table@ 99.9999% Purity

| Inlet Air Pressure (psig)<br>Models 76-97NA and 76-98NA | Max Outlet Flow (lpm) | Max Outlet Pressure (psig) |
|---|-----------------------|----------------------------|
| 120   | 5                     | 83                         |
| 110   | 5                     | 73                         |
| 100   | 5                     | 63                         |
| 90  | 4                     | 62                         |
| 80  | 4                     | 51                         |
| 70  | 2                     | 50                         |
| 60  | 2                     | 42                         |

## Flow Table@ 99.995% Purity

| Inlet Air Pressure (psig)<br>Models 76-97NA and 76-98NA | Max Outlet Flow (lpm) | Max Outlet Pressure (psig) |
|---|-----------------------|----------------------------|
| 120   | 12                    | 60                         |
| 110   | 12                    | 55                         |
| 100   | 12                    | 45                         |
| 90  | 10                    | 45                         |
| 80  | 8                     | 40                         |
| 70  | 8                     | 35                         |
| 60  | 6                     | 33                         |

## Principal Specifications

| Model                          | 76-97NA/76-98NA                             |
|--------------------------------|---|
| Nitrogen Purity                | 99.995% and 99.9999%                        |
| Max Nitrogen Output Pressure   | See Table                                   |
| CO Concentration               | < 1 ppm                                     |
| CO <sub>2</sub> Concentration  | < 1 ppm                                     |
| O <sub>2</sub> Concentration   | < 1 ppm                                     |
| H <sub>2</sub> O Concentration | < 2 ppm                                     |
| Hydrocarbon Concentration (1)  | < 0.1 ppm                                   |
| Argon Concentration (2)        | 0.9%  |
| Min/Max Inlet Pressure         | 60 psig/120 psig                            |
| Recommended Inlet Temperature  | 78°F (25°C)                                 |
| Ambient Operating Temperature  | 60°F-100°F (16°C-38°C)                      |
| Average Air Consumption        | 3.0 scfm                                    |
| Inlet Connection               | 1/4" NPT                                    |
| Outlet Connection              | 1/8" NPT, convertible to 1/4" NPT           |
| Electrical Requirements (3, 4) | 120 VAC/60 Hz                               |
| Dimensions                     | 41"h x 25"w x 25"d<br>(104cm x 64cm x 64cm) |
| Shipping Weight                | 500 lbs (227 kg)                            |

## Ordering Information

| Model Numbers        | Description                             |
|----------------------|---|
| 76-97NA and 76-98NA  | Ultra High Purity Nitrogen Generator    |
| 76-97-PM, 76-98-PM   | Preventive Maintenance Plan             |
| 76-97-DN2, 76-98-DN2 | Extended Support with 24 Month Warranty |

Notes:

1 Model 76-97NA does not remove hydrocarbons.

2 Purity specification for Nitrogen does not include Argon concentration.

3 Power Consumption is as follows:  
Model 76-97NA = 10 Watts, Model 76-98NA = 1 KW

4 Refer to voltage appendix for electrical and plug configurations for outside North America.