

Application Note

Automated purification of His-tagged β-galactosidase

Category	Bioanalytics
Matrix	Enzymes
Method	FPLC
Keywords	automation, IMAC, desalting, β -galactosidase
Analytes	β-galactosidase
ID	VBS0026N



Summary

Introduction

This application note describes an automated method using BioFox 40 IDAlow loaded with Ni²⁺ and a desalting column for a two-step automated isolation of a recombinant β -galactosidase. The high metal-ion capacity of the media and the automation of the method make it well-suited for a wide range of research laboratories and for the industry.

The hydrolase enzyme β -galactosidase catalyzes the hydrolysis of β -galactosides into monosaccharides. It is an essential enzyme in the human body and commonly used in molecular biology as a reporter marker in various applications.

Several IMAC applications have been developed for the purification of the recombinant β -galactosidase¹. Ni-NTA agarose is the most common used resin for IMAC preparations. In comparison to BioFox 40 IDAlow loaded with Ni²⁺, Ni-NTA has a lower pressure resistance followed by a slower filling process and a poorer resolution.

BioFox 40 IDAlow agarose is produced from highly purified agarose using a proprietary cross-linking method that results in a highly porous and physically stable agarose matrix. The ligand IDA is comparable to NTA, but the capacity of binding metal-ions is much lower for IDAlow because of the final spherical interaction between the proteins of interest.

In this study, the resin BioFox 40 IDAlow loaded with Ni²⁺ was pressure-filled in a Bioline HR glass column and connected to a Bioline protein purification system as well as a desalting column. In an automated run, the isolation of His-tagged β -galactosidase will be shown.

Kluyveromyces lactis cells were cultured and the genomic DNA isolated from the harvested cells. The *LAC4*-coding DNA fragment was amplified from *K. lactis* genomic DNA by PCR. *K. lactis* β -galactosidase gene was expressed in *E.coli* as a His-tagged recombinant enzyme.


Experimental: Sample preparation

HPLC · SMB · Osmometry



Method parameters

Column	BioFox 40 IDAlow					
Eluent A	50 mM Na/K-phosphate buffer (pH 6.5) + 500 mM NaCl					
Eluent B	50 mM Na/K-phosphate buffer (pH 6.5) + 500 mM NaCl					
	+ 500 mM imidazole					
Eluent C	20 mM Tris-buffer (pH 7.2)					
Gradient	t [min]	% A	% B	% C		
	0.0	96	4	0		
	20.0	96	4	0		
	20.1	96	4	0		
	35.0	96	4	0		
	35.1	0	100	0		
	42.5	0	100	0		
	42.6	0	0	100		
	120.0	0	0	100		
Flow rate	0.8 ml/min					
m [injection]	200 mg total protein					
T [column]	RT					
Detection	UV at 280 nm					
Run time	100 min					



Using a KNAUER Bioline biochromatography system and an HR glass column filled with BioFox 40 IDAlow loaded with Ni²⁺, recombinant His-tagged β -galactosidase was successfully separated from *E.coli* in an automated run (Fig. 1). After the desalting step, the analysis of the SDS-PAGE with Comassie Blue staining allowed for the identification of β -galactosidase using a molecular weight marker (lane 1). The yields of the separated product before desalting (lane 3) and after desalting (lane 4) were compared to the extract and showed a increase of β -galactosidase after desalting (lane 4).

Fig. 1

Purification of recombinant, His_6 -tagged β -galactosidase from Kluyveromyces lactis.

www.knauer.net

HPLC · SMB · Osmometry



Conclusion	The recombinant β -galactosidase expressed in <i>E.coli</i> is a His-tagged fusion protein and was purified by Ni- BioFox 40 IDAlow agarose affinity chromatography. The separation on an IMAC medium involving a final desalting step, was run on a Bioline protein purification system with an HR glass column. This fully automated method produced in a highly-purified protein fraction.		
References	1 Chang Sup Kim et al. Biotechnology Letters, Vol. 25, 1769-1774 (2003)		
Authors	Dr. Melanie Naether , Bioanalytics Department, KNAUER Timo Stressler, Thomas Eisele, Lutz Fischer, University of Hohenheim, Institute for Food		
	Science and Biotechnology, Department for Biotechnology		
Physical properties of recommended column			
	Stationary phase	BioFox 40 IDAlow [loaded with Ni ²⁺]	
	Specific surface area	320 m ² /g	
	Particle size	40 μm (32-60 μm]	
	Metal-ion capacity	10-20 μegu Ni ²⁺ /ml]	
	Order number resin	Y4025	
	HR Bioline glass column	10 x 300 mm	
	Order number	A71110	
Recommended instrumentation	The biochromatography run was performed on a KNAUER Bioline system, equipped with a gradient pump \$1050, degasser unit \$5050, UV detector, conductivity detector, and a fraction collector.		
1 t	Description	Order No.	
	Smartline Pump 1050	A50351-1	
	Smartline Manager 5050	A5331-1	
RR III	SmartMix 350	A5351	
	Bioline Assistant 6000	A5004V316	
	Smartline UV Detector 2520	A5150	
	3 mm flow cell	A4132	
	Smartline Conductivity Monitor 2900	A70090	
	Smartline Fraction Collector	A53712	
	Bioline affinity module	A7004V322	
	Bioline MPLC Rack	A70190	

Contact information

Wissenschaftliche Gerätebau Dr. Ing. Herbert Knauer GmbH Hegauer Weg 38 14163 Berlin, Germany

Bioline Thermostate

ChromGate software

 Tel:
 +49 (0)30 / 809727-0

 Fax:
 +49 (0)30 / 8015010

 E-Mail:
 info@knauer.net

 Internet:
 www.knauer.net

A1456-8/A1470-8/A1512-8

A70050

