

CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

Basic HPLC Troubleshooting

Technical Note CS06

The most common problems and solutions observed in HPLC separations are outlined in this technical note.

Basic troubleshooting involves systematically isolating the problem and seeing if it is reproducible. Once the problem is isolated, the next step is to change only one variable at a time to see if it resolves the problem. Always document the problem and solution for future reference, to reduce the amount of time spent troubleshooting.

Problem: Extra or ghost peaks

Probable Causes

Extend analysis run time or flush column with stronger eluent at end of run time. Peaks from previous injection

Column contamination Improve sample clean-up by adding either a filtration or solid phase extraction step to sample

> preparation prior to analysis. High back pressure can be an indication of contaminant buildup at the head of the column. Disconnect the column from the detector reverse the direction of flow through the column and backflush. It maybe necessary to remove strongly retained contaminants with a series of miscible solvents. An example of a solvent rinse for heavily contaminated reverse phase columns is water, acetonitrile, tetrahydrofuran, methylene chloride, hexane and reverse the order until back to water.

Refer to the column manufacturer's instructions for rinsing heavily contaminated columns.

Contaminated solvents or water Use HPLC grades solvents and water for mobile phase preparation.

Filter all mobile phases through a $0.45\mu m$ or $0.20\mu m$ prior to use.

Use solvent reservoir and in-line filters to trap particulate.

Store mobile phase in glass containers to avoid contamination from plasticizers.

Problem: Tailing Peaks

Probable Causes

Basic analytes interacting

with surface silanols

Solution

For reverse phase columns, use an endcapped column that has been treated to reduce free surface silanol

Triethylamine (TEA) can be added to MP to irreversibly react with free surface silanols.

Metals in silica Use high purity silica based columns with low metal content.

Problem: Broad Peaks

Probable Causes

Sample overload

Solution Reduce injection volume or sample concentration.

Extra column volume Small area peaks that are quite broad often indicate excessive extra column volume between the injector

and the detector. Use zero dead volume connections and the minimum tubing diameter/length required.

Too long a run time If possible use gradient elution program or stronger mobile phase.

Poor column efficiency Use a column with a smaller particle size to increase column efficiency.

Blocked column frit Disconnect column from detector, reverse flow and rinse to remove blockage.

Problem: Split Peaks

Probable Causes

Solution Sample volume too large Reduce injection volume or sample concentration.

Injection solvent too strong The sample should be diluted in mobile phase when possible to avoid peak shape distortion

when the sample diluent is stronger than the mobile phase.

Column void or channeling Replace column.

Blocked column frit Reverse column flow and rinse column. Installation of $0.2-0.5\mu$ m in-line filter after the injector can

reduce the amount of particulate contaminating the analytical column.

Change in mobile phase If only one or two peaks are splitting, resolution may have been reduced due to changes in mobile

composition

phase composition. Prepare fresh mobile phase.

Contamination Additional sample clean up with filtration or solid phase extraction may be required to remove

interfering peak.

Channels present in column Replace analytical column.

> P.O. Bag 1150, 300 Laurier Blvd., Brockville, ON K6V 5W1 Toll Free: 1-800-267-8103 • Fax: 613-342-1144 • sales@chromspec.com

Problem: Changing Retention Times

<u>Probable Causes</u> <u>Solution</u>

Fluctuations in column Use a column heater or chiller to maintain temperature stability

temperature

Insufficient column equilibration

during gradient run

Increase equilibration time between injections. Allow at least 10-15 column volumes of mobile phase to

pass through the column prior to sample injection.

Contamination build-up Improve sample clean-up by adding either a filtration or solid phase extraction step to sample

preparation prior to analysis. High back pressure can be an indication of contaminant buildup at the head of the column. Disconnect the column from the detector reverse the direction of flow through the column

and backflush.

Changes in mobile phase

composition

Keep solvent reservoirs covered to prevent evaporation.

Buffered mobile phases should be prepared fresh.

Helium sparging could increase evaporation of volatile solvents. Use of in-line solvent degasser can

reduce evaporation caused by helium sparging.

Column aging Column life can be extended with the use of a guard column.

Replace analytical column.

Problem: Baseline Noise

Probable Causes

<u>Solution</u>

Contamination Random base line noise or baseline drift often indicates contamination. Flush column especially after

using buffered mobile phases, backflush column if necessary.

Improve sample clean-up prior to injection by filtering sample through $0.45\mu m$ or $0.20\mu m$ syringe filter

or adding a solid phase extraction step.

Use HPLC grade solvents, reservoir filters and inline filters to trap contaminants.

Detector lamp deterioration Continuous baseline noise often indicates lamp deterioration, replace lamp. Most UV detector lamps

last approximately 1000 hours.

Air bubbles If the baseline noise stops when the flow is turned off, there is likely an air bubble in the detector.

A back-pressure regulator on the outlet of the flow cell can help remove air bubbles.

Oscillating back pressure is frequently caused by inadequate degassing of the mobile phase(s).

Electrical interference Occasional baseline noise could be due to electrical interference by water baths and other

instrumentation. Use of a voltage stabilizer can reduce noise due to electrical interference

Problem: High Column Back Pressure

<u>Probable Causes</u> <u>Solution</u>

Column blockage or blocked frit Reverse solvent flow and backflush column while disconnected from detector.

Particle size too small Use a column with a larger particle size.

Mobile phase too viscous

Use lower viscosity solvents or increase column temperature.

Microbial growth in mobile phase Prepare aqueous buffers fresh daily, or add >10% organic solvent to prevent microbial growth.

Problem: Drop in system pressure or leaks
Probable Causes Solution

Loose fitting If there is precipitate around the fittings, re-tighten fitting or replace ferrule and tubing if necessary.

Pump seal failure Replace pump seals.

Worn injection valve rotor Replace rotor.

Clogged solvent reservoir filters Clean or replace solvent reservoir filter.