

# EZ No-Vent<sup>™</sup> GC-MS Connector

**Change Columns in Minutes, Without Venting** 







### EZ No-Vent™ GC Column-Mass Spectrometer Connector


For Agilent GCs with 5971/5972 or 5973 GC/MS or Varian Saturn 2000 Series GC/MSs

- Change columns in minutes without venting—100 $\mu$ m transfer line throttles vacuum and prevents MS venting.
- Easy to install and maintain—no special tools or plumbing required.
- Gold-plated body for inertness.
- Deactivated transfer line keeps analytes focused; high-temperature polyimide ferrules eliminate leaks at the problematic transfer line fitting.
- · Lower cost than other "no-vent" fittings.

We designed the EZ No-Vent™ GC column-mass spectrometer connector to be simple and easy to use. After studying user feedback concerning our EZ-Vent™ 2000 connector, we re-engineered the connector fitting for even better performance.

A critical orifice in the EZ No-Vent™ connector minimizes the amount of oxygen allowed into the MS source, eliminating the need for purge gas and enabling you to skip the lengthy vent and pump-down cycle otherwise required when you make a column change. Save nearly a day of downtime with each column change by using an EZ No-Vent™ connector. The EZ No-Vent™ connector easily attaches to the MS source without special tools or extra plumbing. The diagram below shows the internal structure of the connector.


The EZ No-Vent<sup>™</sup> connector incorporates a laser machined orifice for optimal flow through the fitting. Coupled with the 100µm ID transfer line, it minimizes the amount of oxygen entering the MS, which eliminates the need for purge gases.



restek innovation!


### The EZ No-Vent™ Connector Does Not Add Dead Volume

We tested the EZ No-Vent<sup>™</sup> connector, using highly volatile gaseous sample components that are very susceptible to tailing in the presence of dead volume. We operated the system as a purge-and-trap GC/MS system, initially with a split at the injection port and the column inserted directly into the MS interface. Subsequently, we included the EZ No-Vent<sup>™</sup> connector at the MS interface. We anticipated that any dead volume in this fitting would produce significant tailing and broader peaks, relative to the direct connection. Peak shape was excellent using the new fitting, as shown below.



# The EZ No-Vent™ Connector Allows Reproducibility

Similarly, we used an application-specific test to evaluate the ability of the MS (equipped with a turbomolecular pump) to stabilize after we changed columns without venting. Again we used a purge-and-trap system, halogenated volatiles as the sample, and an Rtx®-624 column to separate the analytes. We acquired the first chromatogram at 08:12 AM, then changed the column. We acquired the second chromatogram 76 minutes later, at 09:28 AM. Note the excellent peak shapes and responses. In the interval between the two analyses we verified MSD tuning, and the system passed bromofluorobenzene (BFB) criteria. Subsequent detailed investigations have established that the EZ No-Vent<sup>TM</sup> connector will allow several column changes in a single day, with no harm to the MS or loss of data quality.



# EZ No-Vent™ GC Column-Mass Spectrometer Connector

For Agilent GCs with 5971/5972 or 5973 GC/MS or Varian Saturn 2000 Series GC/MSs

- · Change columns in minutes without venting.
- Easy to install and maintain.
- Gold-plated body for inertness.
- Keeps analytes focused; no leaks at the transfer line fitting.
- Lower cost than other "no-vent" fittings.

### **EZ No-Vent™ Connector Kit for Agilent GCs**


Description	qty.	cat.#
EZ No-Vent™ Connector Kit for Agilent 5971/5972 and 5973 GC/MS		
Kit includes: EZ No-Vent™ Connector, two 0.4mm ID ferrules for capillary		
column, two 0.4mm ID ferrules for transfer line, 100 $\mu$ m deactivated		
transfer line (3 ft.), column plug, column nut.	kit	21323

### **EZ No-Vent™ Connector Kit for Varian GCs**

Description	qty.	cat.#
EZ No-Vent™ Connector Kit for Varian Saturn 2000 Series MSs		
Kit includes: EZ No-Vent™ Connector, two 0.4mm ID ferrules for capillary		
column, two 0.4mm ID ferrules for transfer line, $100\mu$ m deactivated		
transfer line (3 ft.), column plug, column nut.	kit	22423

# Replacement Parts for EZ No-Vent™ Connector for Agilent or **Varian GCs**

qty.	cat.#
2-pk.	21015
2-pk.	21016
2-pk.	21043
3 ft.	21018
5-pk.	21900
2-pk.	21915
2-pk.	20110
	2-pk. 2-pk. 2-pk. 3 ft. 5-pk. 2-pk.





The EZ No-Vent<sup>™</sup> connector kit installs easily, without special tools or plumbing.

Restek Trademarks: EZ No-Vent, "Turning Visions into Reality," Restek logo



Lit. Cat.# 59498A-INT © 2005 Restek Corporation.

