GC Troubleshooting Tips

BASIC STEPS

Follow these basic troubleshooting steps to isolate problems related to the sample, injector, detector, and column. Check the obvious explanations first and change only one thing at a time until you identify and resolve the problem.

Check the Obvious:

- Power supply
- Gas purity

Gas flows

Electrical connections

Signal connections

- Temperature settings
- Syringe condition
- Sample preparation
- Analytical conditions

Identify the Cause:


- Define the problem clearly; for example, "Over the last 4 days, only the phenols in my sample have been tailing."
- Review sample and maintenance records to identify trends in the data or problem indicators, such as area counts decreasing over time or injector maintenance not being performed as scheduled.
- Use a logical sequence of steps to isolate possible causes.

Document Work and Verify System Performance:

- Document all troubleshooting steps and results; this may help you identify and solve the next problem faster.
- Always inject a test mix and compare to previous data to ensure restored performance.

EXAMPLE TROUBLESHOOTING SEQUENCE

An analyst observed that no peaks appeared during a GC-FID analysis. The flowchart below shows a logical progression of steps that can be used to identify the cause and correct the problem.





For an in-depth discussion on how to choose the right column and improve your chromatographic results, visit **www.restek.com**

- Download our *Guide to GC Column Selection and Optimizing Separations* (lit. cat.# GNBR1724- UNV)
- Request our Simplifying Column Selection poster (GNWC1612-UNV)

SYMPTOMS AND SOLUTIONS

Good chromatography is critical to obtaining accurate, reproducible results. Coelutions, asymmetric peaks, baseline noise, and other issues are common challenges in the GC laboratory. These analytical problems and others can be overcome by troubleshooting your separations using the tips below.

Solutions
Choose appropriate stationary phase and column dimensions.
 Optimize carrier gas linear velocity and GC oven temperature program.
 Adjust sample concentration or amount on column.
Verify temperature program, flow rates,

Poor Retention Time Reproducibility

Causes	Solutions
	Leak check injector and press-fit connections
Leaks	 Replace critical seals (i.e., septa, O-rings, inlet disc, etc.)
	• Maintain inlet liner and GC column.
Analyte adsorption	 Use properly deactivated liners, seals, and columns.
Resolution/integration issues	Avoid sample overload.
Incorrect column/oven temperature program	 Verify column temperature and oven temperature program.
Incorrect or variable carrier gas flow rate/linear velocity	 Verify the carrier gas flow and linear velocity. Repair or replace parts if neccessary.
Poor control of oven temperature programming	Confirm GC oven program falls within instrument manufacturer's recommendation.
Incorrect oven equilibration time	• Extend GC oven equilibration time.
If manual injection, delay between pushing start and actual injection	 Use autosampler or standardize manual injection procedure.

Fronting Peaks



............

.....

CausesSolutionsIncompatible stationary• Choose appropriate stationary phase.

phase

Reduce amount injected, dilute sample.
 Increase column inner diameter and/or

Tailing Peaks

lin
ace
(sqi

film thickness.

Verify that the column is cut properly (square).Verify correct installation distances.

Sp	lit	Pea	ks

	Causes	Solutions
M	Mismatched solvent/stationary phase polarity	 Adjust solvent or stationary phase to allow wetting.
	Incomplete vaporization	 Add surface area, such as wool, to the inlet liner to enhance vaporization. Use proper injector temperature.
	Sample loading capacity exceeded	 Inject less sample (dilute, use split injection, reduce injection volume).
	Fast autosampler injection into open liner	• Use wool or slow injection speed.

ттт


Carryover/Ghost Peaks			
Т		Causes	Solutions
		Contaminated syringe or rinse solvent	Replace rinse solvent.
	Injection 1		Rinse or replace syringe.
			Inject a smaller amount.
			• Use a liner with a large internal diameter.
······	Backflash (sample volume exceeds liner volume)	 Increase head pressure (i.e., flowrate) to contain the vapor cloud. 	
		Use slower injection rate.	
		exceeds liner volume,	Lower inlet temperature.
			Increase split flow.
	Injection 2		Use liner with packing.
			Use pressure-pulse injection.
		Last analysis ended too soon	• Extend analysis time to allow all components and/or matrix interferences to elute.

High Bleed			
	Causes	Solutions	
	Improper column conditioning	 Increase conditioning time and/or temperature. 	
+++++++++++++++++++++++++++++++++++++++	Contamination	 Trim column and/or heat to maximum temperature to remove contaminants. 	
		 Replace carrier gas and/or detector gas filters. Clean injector and detector. 	
	Leak in system and oxidation of stationary phase	 Check for oxygen leaks across the entire system and replace seals and/or filters. Replace column. 	

Unstable Baseline (Spiking, Noise, Drift)

		Causes	Solutions
Spiking	Spiking	Carrier gas leak or contamination	 Leak check connections and replace seals if needed.
			Replace carrier gas and/or detector gas filters.
		Injector or detector contamination	Clean system and perform regular maintenance.
TT		Column contamination or stationary phase bleed	Condition, trim, and rinse column.
			Replace septum.
	Noise	Septum coring/bleed	 Inspect inlet liner for septa particles and replace liner if needed.
		Loose cable or circuit board connections	Clean and repair electrical connections.
		Variable carrier gas or detector gas flows	 Verify flow rates are steady and reproducible; may need to replace or repair flow controller. Leak check system.
Drift			· Leak check system.
		Detector not ready	• Allow enough time for detector temperatures and flows to equilibrate.

Response Variation




Causes	Solutions
	Check sample concentration.
Sample issues	Check sample preparation procedure.
	Check sample decomposition/shelf life.
Syringe problems	Replace syringe.
Synnge problems	Check autosampler operation.
Flectronics	• Verify signal settings and adjust if needed.
Licetonies	Repair or replace cables or boards.
Dirty or damaged	Perform detector maintenance or
detector	replace parts.
Flow/temperature	Verify steady flow rates and temperatures, then
settings wrong or variable	adjust settings and/or replace parts if needed.
Adsorption/reactivity	 Remove contamination and use properly deactivated liner, seal, and column.
Leaks	Check for leaks at all connections and repair connections as needed.
Change in sample	 Verify injection technique and change back to original technique.
introduction/injection method	Check that split ratio is correct.
incurou	• Verify that the splitless hold time is correct.

No Peaks

Causes	Solutions
Injection problems	Plugged syringe; clean or replace syringe.
	 Verify there is sample in the syringe.
	 Injecting into wrong inlet; reset autosampler.
	• Verify carrier gas is flowing.
Broken column	• Replace column.
Column installed into wrong injector or detector	• Re-install column.
Detector problems	 Signal not recorded; check detector cables and verify that detector is turned on.
	 Detector gas turned off or wrong flow rates used; turn detector on and/or adjust flow rates



Broad Peaks



Causes	Solutions
High dead volume	 Minimize dead volume in the GC system; verify proper column installation, proper connectors, proper liners, etc.
Low flow rates	 Verify injector and detector flow rates and adjust if needed.
	• Verify make-up gas flow and adjust if needed.
Slow GC oven program	Increase GC oven programming rate.
Poor analyte/solvent focusing	Lower GC oven start temperature.
Column film is too thick	 Reduce retention of compounds by decreasing film thickness and length.
Sample carryover	See Carryover/Ghost Peaks solutions.

Lit. Cat.# GNWC1723-UNV



