

Solutions for Your Petro Analyses


| | | | |
|------------------------------------|---|-------------------------------------|----|
| Simulated Distillation..... | 2 | Options for Biodiesel | 10 |
| Detailed Hydrocarbon Analysis..... | 4 | Reference Standards..... | 11 |
| Bonded PLOT Columns | 5 | Sample Cylinders & Accessories..... | 14 |
| D3606 Analysis | 8 | Petro Essentials..... | 16 |
| Permanent Gases ... | 9 | | |

Visit us at www.restek.com/petro

RESTEK

Innovative Chromatography Solutions

www.restek.com



**CHROMATOGRAPHIC
SPECIALTIES INC.**

www.chromspec.com

1-800-267-8103 • sales@chromspec.com • tech@chromspec.com




Improve Productivity—Get More Runs from Your SimDist Setup Using Stabilized MXT®-1HT Columns

- Long lasting, low bleed columns are stable up to 450 °C.
- Individually tested for guaranteed performance.
- Meet all ASTM D6352 and D7500 specifications.
- Unbreakable columns offer security when using hydrogen.

Accurate determination of the boiling range distribution of medium and heavy fractions using GC simulated distillation requires columns and phase polymers that are robust enough to withstand the high method temperatures without significant degradation. Metal columns are a much better alternative than fused silica, and the MXT®-1HT SimDist columns, with stabilized dimethyl polysiloxane polymer, are rugged enough to meet critical method parameters, including bleed and C50/C52 resolution (Figures 1 and 2). Field testing of the MXT®-1HT SimDist column shows excellent performance, even under faster analytical conditions than those in the published method. Bleed is minimal, even at 430 °C, which is essential for precise time slices and accurate final boiling point determination. MXT®-1HT SimDist columns are the lowest bleed column on the market, which translates directly into more analyses per calibration and longer column lifetimes. In addition, all Restek® SimDist columns are individually tested up to C90 and thus require only minimal conditioning. This gives you guaranteed method performance while saving you time and money.

Figure 1: Low-bleed, high-efficiency MXT®-1HT SimDist columns are individually tested for guaranteed performance (ASTM D6352 conditions).



Column: MXT®-1HT SimDist, 5 m, 0.53 mm ID, 0.10 µm (cat.# 70112); **Sample:** C8-C100 hydrocarbons (Separation Systems SD-ss3e-05), Diluent: CS₂, Conc.: 1.0%; **Injection:** Inj. Vol.: 0.5 µL cold on-column, Temp. Program: 53 °C to 430 °C at 10 °C/min (hold 5 min); **Oven:** Oven Temp.: 50 °C to 430 °C at 10 °C/min (hold 5 min); **Carrier Gas:** He, constant flow, Flow Rate: 18 mL/min; **Detector:** FID @ 430 °C, Make-up Gas Flow Rate: 24 mL/min, Make-up Gas Type: N₂; **Instrument:** Shimadzu 2010 GC; **Note:** Inlet temperature program tracks oven temperature program.



 Find everything you need for each method in our
ASTM Product Guide at www.restek.com/petro



Figure 2: Superior resolution and peak shape on MXT®-1HT SimDist columns result in more accurate final boiling point determinations.



Column: MXT®-1HT SimDist, 5 m, 0.53 mm ID, 0.20 µm (cat.# 70115); **Sample:** C5-C100, Diluent: carbon disulfide, Conc.: 1 %; **Injection:** Inj. Vol.: 1 µL cold on-column, Inlet Temp. Program: 53 °C to 430 °C at 10 °C/min (hold 5 min); **Oven:** Oven Temp.: 50 °C to 430 °C at 10 °C/min (hold 5 min); **Carrier Gas:** He, constant flow; Flow Rate: 18 mL/min; **Detector:** FID @ 430 °C; **Instrument:** Shimadzu 2010 GC; **Note:** Inlet temperature program tracks oven temperature program.

SimDist Products

Method Recommended Columns

| ASTM Method | Hydrocarbon Range | Dimensions | cat. # |
|----------------------|-------------------|-------------------------|--------|
| D2887 | C5–C44 | 5 m x 0.53 mm, 0.88 µm | 70131 |
| | | 10 m x 0.53 mm, 2.65 µm | 70132 |
| | | 5m x 0.53 mm, 0.88 µm | 70131 |
| D7213 (D2887-ext) | C5–C60 | 5 m x 0.53 mm, 0.20 µm | 70115 |
| | | 5 m x 0.53 mm, 0.10 µm | 70112 |
| | | 5 m x 0.53 mm, 0.20 µm | 70115 |
| D5307 | crude up to C42 | 5 m x 0.53 mm, 0.20 µm | 70115 |
| D6352 | C10–C90 | 5 m x 0.53 mm, 0.10 µm | 70112 |
| | | 5 m x 0.53 mm, 0.20 µm | 70115 |
| D7500 | C7–C110 | 5 m x 0.53 mm, 0.10 µm | 70112 |
| | | 5 m x 0.53 mm, 0.20 µm | 70115 |
| D7169 | C5–C100 | 5 m x 0.53 mm, 0.10 µm | 70112 |
| | | 5 m x 0.53 mm, 0.20 µm | 70115 |

MXT®-1HT SimDist Column

(Siltek®-treated stainless steel) (nonpolar phases)

| Description | temp. limits | qty. | cat.# |
|---------------------------|-------------------|------|-------|
| 5 m, 0.53 mm ID, 0.10 µm | -60 to 430/450 °C | ea. | 70112 |
| 5 m, 0.53 mm ID, 0.20 µm | -60 to 400/430 °C | ea. | 70115 |
| 10 m, 0.53 mm ID, 0.21 µm | -60 to 400/430 °C | ea. | 70118 |
| 5 m, 0.53 mm ID, 0.88 µm | -60 to 380/430 °C | ea. | 70131 |
| 10 m, 0.53 mm ID, 0.88 µm | -60 to 400/430 °C | ea. | 70134 |
| 10 m, 0.53 mm ID, 1.00 µm | -60 to 380/400 °C | ea. | 70130 |
| 10 m, 0.53 mm ID, 1.20 µm | -60 to 380/380 °C | ea. | 70119 |
| 10 m, 0.53 mm ID, 2.65 µm | -60 to 360/400 °C | ea. | 70132 |
| 10 m, 0.53 mm ID, 5.00 µm | -60 to 360/400 °C | ea. | 70133 |




Get More Accurate DHA Analysis Including Alcohols Using Rtx®-DHA Columns

- Columns individually QC tested to meet or exceed all ASTM and CAN/CSGB method guidelines.
- Excellent responses and peak symmetry for polar oxygenates.
- Guaranteed column-to-column reproducibility for resolution, asymmetry, k, and low bleed.
- Columns designed for ASTM Methods D6729, D6730, D6733, D5134, and D5501.

Gasolines are complex mixtures of hundreds of compounds. Information about concentrations of the individual components is important for evaluating raw materials and for controlling refinery processes. ASTM D6730 outlines a high-resolution GC method for detailed hydrocarbon analysis (DHA) of gasolines. Rtx®-DHA columns are ideal for DHA methods and easily meet or exceed both ASTM D6730 and Canadian General Standards Board CAN/CSGB 14.3-99 requirements (Figure 3). Every Rtx®-DHA column is tested for resolution, asymmetry, k, and bleed—guaranteeing reproducible column-to-column performance.

Figure 3: Critical pairs of gasoline components resolved per ASTM specifications on an Rtx®-DHA column.



DHA Products

did you know?

Using hydrogen instead of helium can cut analysis time in half! Visit www.restek.com/petro for complete analytical details.

Rtx®-DHA Columns (fused silica)
(Crossbond® 100% dimethyl polysiloxane—optimized for hydrocarbon analysis)

| Description | temp. limits | qty. | cat.# |
|---------------------------------|-------------------|------|-------|
| 50 m, 0.20 mm ID, 0.50 μ m | -60 to 300/340 °C | ea. | 10147 |
| 100 m, 0.25 mm ID, 0.50 μ m | -60 to 300/340 °C | ea. | 10148 |
| 150 m, 0.25 mm ID, 1.00 μ m | -60 to 280/340 °C | ea. | 10149 |

Rtx®-5 DHA Tuning Column (fused silica)
(Crossbond® 5% diphenyl/95% dimethyl polysiloxane—optimized for hydrocarbon analysis)

| Description | temp. limits | qty. | cat.# |
|-------------------------------|-------------------|------|-------|
| 5 m, 0.25 mm ID, 1.00 μ m | -60 to 325/350 °C | ea. | 10165 |

Method Recommended Columns

| ASTM Method | Column | Dimensions | cat. # |
|-------------|---------------------------------------|--|---------------|
| D6729 | Rtx-DHA-100 | 100 m x 0.25 mm, 0.50 μ m | 10148 |
| D6730 | Rtx-DHA-100 & Rtx-5 DHA Tuning Column | 100 m x 0.25 mm, 0.50 μ m w/ precolumn | 10148 & 10165 |
| D6733 | Rtx-DHA-50 | 50 m x 0.20 mm, 0.50 μ m | 10147 |
| D5501 | Rtx-DHA-150 | 150 m x 0.25 mm, 1.0 μ m | 10149 |

similar phases

Petrocol DH, DB-Petro, HP-PONA

Note: Rtx®-1 PONA columns have been renamed as Rtx®-DHA columns. There are no changes in manufacturing process or column performance.



Protect Equipment and Assure Predictable Retention Times Using Bonded PLOT Column Technology


- Bonding process minimizes particle release, reducing column blockage and protecting instrument parts.
- More consistent flow means stable retention times; ideal for Deans and related flow switching techniques.
- Outstanding peak symmetry improves accuracy of impurity analysis for gases, solvents, and hydrocarbons.

PLOT columns are widely used in the petroleum industry because their retention and selectivity characteristics allow gases and volatiles to be separated with high resolution at above ambient temperatures. However, the utility of PLOT columns is significantly limited by the mechanical instability of the particle layers. Routine vibrations or changes in gas pressure can cause a release of particles that clog the column, resulting in highly variable flow behaviors over time and between columns.

Restek has developed a bonding process that established a new benchmark for PLOT column stabilization. This innovative process significantly reduces particle release and column blockage, resulting in highly stable flows and retention times, both run-to-run and column-to-column (Figure 4). By replacing conventional PLOT columns with bonded PLOT columns manufactured using this Restek® technology, labs can increase the accuracy of impurities analysis and make better process decisions.

This stabilization technology is especially beneficial to flow switching applications and is available for porous polymer (4 selectivities), alumina (MAPD, CFC, KCl, and Na₂SO₄ deactivations), and molecular sieve PLOT columns. Our newest PLOT column, the Rt®-Silica BOND PLOT column, is so stable it is virtually particle-free and can be used without pressure conditioning.

Figure 4: Restek® PLOT column technology assures stable retention times, even after 500 pressure cycles.



PLOT Columns

Fused Silica PLOT Columns



Rt®-Silica BOND Columns (fused silica PLOT)


- Versatile column ideal for analysis of light hydrocarbons, sulfur gases, halocarbons, and carbon dioxide.
- Individually QC tested with sensitive C4 probes to ensure consistent selectivity.
- Proprietary manufacturing process practically eliminates particle release, reducing downtime due to clogged FID jets.
- Bonded silica stationary phase minimizes impact of water, resulting in reproducible retention times for water-containing samples.

similar phases

GS-GasPro, CP-SilicaPLOT

| ID | temp. limits | 15-Meter cat.# | 30-Meter cat.# | 60-Meter cat.# |
|---------|---------------|----------------|----------------|----------------|
| 0.32 mm | -80 to 260 °C | 19784 | 19785 | 19786 |

Saturated and Unsaturated Hydrocarbons on Rt®-Silica BOND PLOT Column



Rt®-Q-BOND Columns (fused silica PLOT)

100% divinylbenzene

- Nonpolar PLOT column incorporating 100% divinylbenzene.
- Excellent for analysis of C1 to C3 isomers and alkanes up to C12.
- High retention for CO₂ simplifies gas analysis; CO₂ and methane separated from O₂/N₂/CO (Note: O₂/N₂/CO not separated at room temperature).
- Use for analysis of oxygenated compounds and solvents.
- Maximum temperature of 300 °C.

| ID | df | temp. limits | 15-Meter cat.# | 30-Meter cat.# |
|---------|-------|---------------|----------------|----------------|
| 0.25 mm | 8 µm | to 280/300 °C | 19764 | 19765 |
| 0.32 mm | 10 µm | to 280/300 °C | 19743 | 19744 |
| 0.53 mm | 20 µm | to 280/300 °C | 19741 | 19742 |

did you know?

Restek® chemists developed a new process for manufacturing PLOT columns. This process bonds the particles to the walls of the tubing, so there is virtually no particle generation. Reduced particle generation assures reproducible selectivity and flow on every run and every column.

Rt®-QS-BOND Columns (fused silica PLOT)

porous divinylbenzene homopolymer

- Intermediate polarity PLOT column incorporating low 4-vinylpyridine.
- Separates ethane, ethylene, and acetylene to baseline.
- Designed for the best possible separation between all C2 isomers.

| ID | df | temp. limits | 15-Meter cat.# | 30-Meter cat.# |
|---------|-------|--------------|----------------|----------------|
| 0.25 mm | 8 µm | to 250 °C | 19767 | 19768 |
| 0.32 mm | 10 µm | to 250 °C | 19739 | 19740 |
| 0.53 mm | 20 µm | to 250 °C | 19737 | 19738 |

Rt®-S-BOND Columns (fused silica PLOT)

porous divinylbenzene homopolymer

- Midpolarity PLOT column, incorporating high 4-vinylpyridine.
- Use for the analysis of nonpolar and polar compounds.

| ID | df | temp. limits | 15-Meter cat.# | 30-Meter cat.# |
|---------|-------|--------------|----------------|----------------|
| 0.25 mm | 8 µm | to 250 °C | 19769 | 19770 |
| 0.32 mm | 10 µm | to 250 °C | 19747 | 19748 |
| 0.53 mm | 20 µm | to 250 °C | 19745 | 19746 |



Rt®-U-BOND Columns (fused silica PLOT)

divinylbenzene ethylene glycol/dimethylacrylate

- Restek's highest polarity porous polymer column.
- Polar PLOT column, incorporating divinylbenzene ethylene glycol/dimethylacrylate.
- Highly inert for the analysis of polar and nonpolar compounds.

| ID | df | temp. limits | 15-Meter cat.# | 30-Meter cat.# |
|---------|-------|--------------|-------------------|-------------------|
| 0.25 mm | 8 µm | to 190 °C | 19771 | 19772 |
| 0.25 mm | 12 µm | to 190 °C | 19782 | — |
| 0.32 mm | 10 µm | to 190 °C | 19751 | 19752 |
| 0.53 mm | 20 µm | to 190 °C | 19749 | 19750 |

Rt®-Alumina BOND/KCI Columns (fused silica PLOT)

(KCI deactivation)

- Restek's lowest polarity alumina column.
- Low moisture sensitivity reduces the need for frequent regeneration.
- Acetylene elutes before *n*-butane.
- Methyl acetylene (impurity in 1,3-butadiene) elutes before 1,3-butadiene.

| ID | df | temp. limits | 30-Meter cat.# | 50-Meter cat.# |
|---------|-------|--------------|-------------------|-------------------|
| 0.25 mm | 4 µm | to 200 °C | 19776 | — |
| 0.32 mm | 5 µm | to 200 °C | 19761 | 19762 |
| 0.53 mm | 10 µm | to 200 °C | 19759 | 19760 |

Rt®-Alumina BOND/Na₂SO₄ Columns (fused silica PLOT)

(Na₂SO₄ deactivation)

- Acetylene and propadiene elute after butanes.
- Best separation for butene isomers (impurities in butene streams).
- Methyl acetylene elutes after 1,3-butadiene.
- Cyclopropane (impurity in propylene) elutes well before propylene.

| ID | ID | temp. limits | 30-Meter cat.# | 50-Meter cat.# |
|---------|-------|--------------|-------------------|-------------------|
| 0.25 mm | 4 µm | to 200 °C | 19775 | — |
| 0.32 mm | 5 µm | to 200 °C | 19757 | 19758 |
| 0.53 mm | 10 µm | to 200 °C | 19755 | 19756 |

Rt®-Alumina BOND/CFC Columns (fused silica PLOT)

- Improved inertness for chlorofluorocarbon (CFC) compounds.
- Highly selective alumina-based column, separates most CFCs.
- High retention and capacity for CFCs.

| ID | df | temp. limits | 30-Meter cat.# |
|---------|-------|--------------|-------------------|
| 0.53 mm | 10 µm | to 200 °C | 19763 |

Rt®-Alumina BOND/MAPD Columns (fused silica PLOT)

- Optimized deactivation produces maximum response when analyzing trace levels of acetylene, methyl acetylene, and propadiene.
- Stable response factors make this column ideal for process-type applications where recalibration must be minimized.
- High loadability reduces peak tailing and improves separations.
- Extended temperature range up to 250 °C for fast elution of high molecular weight (HMW) hydrocarbons and accelerated column regeneration following exposure to water.

| ID | df | temp. limits | 30-Meter cat.# | 50-Meter cat.# |
|---------|-------|--------------|-------------------|-------------------|
| 0.25 mm | 4 µm | to 250 °C | 19781 | — |
| 0.32 mm | 5 µm | to 250 °C | 19779 | 19780 |
| 0.53 mm | 10 µm | to 250 °C | 19777 | 19778 |

Rt®-Msieve 5A Columns (fused silica PLOT)

- Fast separation of all inert/permanent gases, including argon and oxygen at temperatures above ambient.
- Special deactivation ensures carbon monoxide elutes as a sharp peak allowing ppm level quantification.
- Unique immobilization process ensures minimal particle generation—even after continuous valve cycling.
- Stable up to 300 °C for fast regeneration following exposure to water or carbon dioxide.
- Available in flexible 0.53, 0.32, and 0.25 mm ID fused silica tubing for easy installation, cutting, and coupling.

| ID | df | temp. limits | 15-Meter cat.# | 30-Meter cat.# |
|---------|-------|--------------|-------------------|-------------------|
| 0.25 mm | 20 µm | to 300 °C | 19773 | — |
| 0.32 mm | 30 µm | to 300 °C | 19720 | 19722 |
| 0.53 mm | 50 µm | to 300 °C | 19721 | 19723 |

Metal MXT® PLOT Columns

Advantages of metal MXT® PLOT columns include:

- Can be made in small coil diameters—perfect for tight spaces.
- Rugged material withstands rough handling and shock.
- Designed for robust performance in process GCs and field instruments.
- Available in 3.5"-coil diameter or 7"-diameter, 11-pin cage.

| ID | df | temp. limits | 3.5" coil | 7" diameter | 3.5" coil | 7" diameter |
|--|-------|---------------|-------------------|-------------------|-------------------|-------------------|
| | | | 15-Meter cat.# | 15-Meter cat.# | 30-Meter cat.# | 30-Meter cat.# |
| MXT-Msieve 5A | | | | | | |
| 0.25 mm | 20 µm | to 300 °C | 79717-273 | 79717 | — | — |
| 0.53 mm | 50 µm | to 300 °C | — | — | 79723-273 | 79723 |
| MXT-Alumina BOND/Na ₂ SO ₄ | | | | | | |
| 0.53 mm | 10 µm | to 200 °C | — | — | 79714-273 | 79714 |
| MXT-Alumina BOND/MAPD | | | | | | |
| 0.53 mm | 10 µm | to 250 °C | — | — | 79728-273 | 79728 |
| MXT-Q-BOND | | | | | | |
| 0.25 mm | 8 µm | to 300 °C | 79718-273 | 79718 | — | — |
| 0.53 mm | 20 µm | to 280/300 °C | — | — | 79716-273 | 79716 |
| MXT-S-BOND | | | | | | |
| 0.53 mm | 20 µm | to 250 °C | — | — | 79712-273 | 79712 |

Visit www.restek.com/petro for more PLOT column applications.




Ensure Accurate Results for Benzene and Toluene in Spark Ignition Fuels Using the D3606 Column Set

- Complete baseline resolution makes accurate quantification easy and reliable.
- Guaranteed performance: every set is tested for method application and includes chromatogram.
- Fully conditioned, ready-to-use column set.

Laboratories analyzing aviation fuel or reformulated spark ignition fuels that contain ethanol for the determination of benzene and toluene must use a modified ASTM D3606 method to prevent the coelution of ethanol and benzene. The primary challenge in this analysis is twofold: the tailing of the ethanol peak into benzene and the retention time shift of the aromatics, particularly benzene, toward the ethanol peak. Restek solved these issues by developing the D3606 column set, specifically for the modified ASTM D3606 application. Using this column set, which is specified in the Method D3606 Addendum, the aromatic compounds are fully resolved and can easily be quantified using the internal standard 2-butanol (Figure 5). This column set is fully conditioned and ready to use out of the box with only a 10-minute carrier gas purge at ambient temperature, followed by a 30-minute hold at 165 °C. It can be used for the analysis of gasoline with or without ethanol, using either helium or nitrogen carrier gas.

Figure 5: Complete resolution of benzene from ethanol using a D3606 column set for modified ASTM D3606.



Column: D3606 application column (2 column set). Column 1: 6' (1.8 m), 1/8" OD, 2.0 mm ID, nonpolar Rtx®-1; Column 2: 16' (4.9 m), 1/8" OD, 2.0 mm ID, proprietary packing material (cat.# 83606-800); **Sample:** Ethanol-containing gasoline with internal standard (IS). Diluent: ; **Injection:** Direct, Sample Loop Vol.: 1.5 µL, Valve Temp.: 135 °C; **Oven:** Oven Temp: 135 °C (hold 12 min); **Carrier Gas:** He, constant flow, Flow Rate: 20.0 mL/min; **Detector:** TCD @ 200 °C; **Notes:** 2.38 minute backflush (must be determined for each GC system).

D3606 Products

D3606 Application Column (2 column set)

| Description | cat.#* |
|---|--------|
| D3606 Application Column (2 column set)** | |
| Column 1: 6' (1.8 m), 1/8" OD, 2.0 mm ID, nonpolar Rtx-1 | 83606- |
| Column 2: 16' (4.9 m), 1/8" OD, 2.0 mm ID, proprietary packing material | |

*Please add column instrument configuration suffix number to cat.# when ordering. See chart on next page.

**The column set is designed to accommodate both valve injection and/or syringe injection. Column 1 is configured with a 2" inlet void to facilitate on-column injection. The inlet is identified on both column 1 and column 2. Note: The inlet of column 2 is identified for proper orientation for connection to the valve.

ordering note

www.restek.com/petro for a complete list of D3606 standards.



Single Column Separation of CO and CO₂ in Presence of O₂/N₂ Using a ShinCarbon ST Packed/Micropacked Column

- Separate permanent gases, as well as CO/CO₂, without cryogenic cooling.
- Excellent compatibility with most GC detectors—low bleed and minimal baseline rise.
- Preconditioned; takes less than 30 minutes to stabilize.

Analyze oxygen, nitrogen, methane, carbon monoxide, and carbon dioxide with one column and at room temperature. Restek's ShinCarbon ST material is a high surface area carbon molecular sieve (~1,500 m²/g) that is ideal for separating gases and highly volatile compounds by gas-solid chromatography (GSC). ShinCarbon ST is an exceptionally stable material with good loadability. Its 300 °C upper temperature limit minimizes bleed and baseline rise during temperature programming, making it compatible with most detection systems, including TCD or HID (Figure 6). All ShinCarbon ST columns are fully conditioned in an oxygen/moisture-free environment to prevent contamination. This minimizes stabilization time (less than 30 minutes) when installing a new column, which, in turn, minimizes downtime. Custom-made ShinCarbon ST columns are available on request.

Restek® innovation!

ShinCarbon ST is an ideal packing material for permanent gases, low molecular weight hydrocarbons, sulfur dioxide, and Freon® gases.

a plus 1 story

"Being one of the first labs to utilize the ShinCarbon column in a real working environment, I was pleased to find that I was able to do all my permanent gas analysis on one column instead of the customary two. The peaks were sharper than I had experienced in the past and run time was significantly reduced. We are extremely pleased with the performance of the ShinCarbon column and will continue to find even more applications for it."

Bruce Nasser, Quality Control Chemist, Oxygen Service Spec Lab

ordering note

Installation kits available at www.restek.com

Packed Column Instrument Configurations



General Configuration
Suffix -800



Agilent 5880, 5890,
5987, 6890, 7890:
Suffix -810*



Varian 3700,
Vista Series, FID:
Suffix -820



PE 900-3920,
Sigma 1,2,3:
Suffix -830




PE Auto System 8300,
8400, 8700:
Suffix -840

Visit www.restek.com
for additional configurations.

Note: Initial 2" of column will be empty to accommodate a needle. For a completely filled column (not on-column) add suffix -901.

*-810 suffix also includes 1 1/2" void on detector side.

Figure 6: Separate permanent gases in 10 min, without cryogenics.



ShinCarbon Products

ShinCarbon ST Columns (packed) (SilcoSmooth® Stainless Steel)*

| OD | ID | Mesh | 2-Meter cat.#* |
|------------------|--------|--------|----------------|
| 1/8" Silcosmooth | 2.0 mm | 80/100 | 80486- |

ShinCarbon ST Columns (micropacked) (SilcoSmooth® Stainless Steel)**

| OD | ID | Mesh | 1-Meter cat.# | 2-Meter cat.# |
|---------|---------|---------|---------------|---------------|
| 1/16" | 1.0 mm | 100/120 | 19809 | 19808 |
| 0.95 mm | 0.75 mm | 100/120 | 19810 | — |
| 0.74 mm | 0.53 mm | 80/100 | 19045 | 19043 |

*Please add column instrument configuration suffix number to cat.# when ordering.

**Does not include column nuts and ferrules. Optional installation kits can be ordered separately.

Visit www.restek.com/packed to order a custom packed column.




Accurate, Reliable Analysis of Glycerides in Biodiesel Using Either Metal or High Temperature Fused Silica Column Solutions

- Complete resolution of all compounds from interference peaks.
- Columns available in metal and fused silica.
 - Use Rtx®-Biodiesel TG fused silica columns up to 380 °C—low bleed for accurate, reliable results.
 - Use MXT®-Biodiesel TG metal columns up to 430 °C—same performance as fused silica with longer column lifetimes.

One of the biggest challenges in biodiesel analysis is the accurate determination of residual glyceride content. The high temperatures required to elute triglycerides cause most fused silica columns to deteriorate rapidly. Restek offers high-temperature fused silica columns that are stable up to 380 °C, and the metal MXT®-Biodiesel TG column line is an even better solution. MXT®-Biodiesel TG columns offer the same chromatographic performance as fused silica columns, but are designed for maximum heat tolerance, even up to 430 °C. These columns provide good resolution and peak shape for all glycerides, as well as highly reproducible retention times (Figure 7). Additionally, metal MXT® columns are available with Integra-Gap® technology, a built-in retention gap that eliminates connector-related leaks, reducing tailing and minimizing downtime for maintenance.

Figure 7: Robust metal MXT®-Biodiesel TG columns meet resolution criteria and show excellent response for determining glycerides in biodiesel.



Biodiesel Products

MXT®-Biodiesel TG Columns (Siltek®-treated stainless steel)

| Description | temp. limits | 3.5" Coil cat.# | 7" diameter 11-pin cage cat.# |
|--|-------------------|--------------------|-------------------------------------|
| 14 m, 0.53 mm ID, 0.16 µm with 2 m Integra-Gap* | -60 to 380/430 °C | 70289-273 | 70289 |
| 10 m, 0.32 mm ID, 0.10 µm | -60 to 380/430 °C | — | 70292 |
| 10 m, 0.32 mm ID, 0.10 µm with 2 m x 0.53 mm Retention Gap** | -60 to 380/430 °C | 70290-273 | 70290 |
| 15 m, 0.32 mm ID, 0.10 µm | -60 to 380/430 °C | — | 70293 |
| 15 m, 0.32 mm ID, 0.10 µm with 2 m x 0.53 mm Retention Gap** | -60 to 380/430 °C | 70291-273 | 70291 |
| 2 m x 0.53 mm MXT Biodiesel TG Retention Gap | | — | 70294 |

*Total column length = 16 meters.

**Connected with low dead volume MXT connector.

Rtx®-Biodiesel TG Columns (fused silica)

| Description | temp. limits | cat.# |
|---|---------------|-------|
| 10 m, 0.32 mm ID, 0.10 µm | to 330/380 °C | 10292 |
| 10 m, 0.32 mm ID, 0.10 µm with 2 m x 0.53 mm ID Retention Gap** | to 330/380 °C | 10291 |
| 15 m, 0.32 mm ID, 0.10 µm | to 330/380 °C | 10294 |
| 15 m, 0.32 mm ID, 0.10 µm with 2 m x 0.53 mm ID Retention Gap** | to 330/380 °C | 10293 |

**Connected with Alumaseal® connector.

Visit www.restek.com/biodiesel for a complete list of biodiesel standards, products, and applications.



Reference Standards

Petroleum

Gravimetrically prepared, NIST-traceable by weight, and verified by one or more analytical methods.

Sulfur Simulated Distillation Standard (SSDS)

30 ppm total sulfur by weight from ethanethiol
60 ppm total sulfur by weight from 1-propanethiol
30 ppm total sulfur by weight from 1-butanethiol
60 ppm total sulfur by weight from 1-pentanethiol
30 ppm total sulfur by weight from 1-hexanethiol
60 ppm total sulfur by weight from 1-heptanethiol
30 ppm total sulfur by weight from 3,5-dimethylbenzenethiol
60 ppm total sulfur by weight from 1-octanethiol
30 ppm total sulfur by weight from 1-nonanethiol
60 ppm total sulfur by weight from 1-decanethiol
30 ppm total sulfur by weight from 1-pentadecanethiol
60 ppm total sulfur by weight from 1-hexadecanethiol
30 ppm total sulfur by weight from 1-octadecanethiol
Balance: toluene:isooctane (1:15)

1 mL pre-scored amber ampul.

cat.# 33049 (ea.)

SimDist

D2887 Calibration Mix (17 components)

| Compound | Conc. (% w/w) | Compound | Conc. (% w/w) |
|---|------------------|--|------------------|
| <i>n</i> -Hexane (C6) (110-54-3) | 6 | <i>n</i> -Eicosane (C20) (112-95-8) | 2 |
| <i>n</i> -Heptane (C7) (142-82-5) | 6 | <i>n</i> -Tetracosane (C24) (646-31-1) | 2 |
| <i>n</i> -Octane (C8) (111-65-9) | 8 | <i>n</i> -Octacosane (C28) (630-02-4) | 1 |
| <i>n</i> -Nonane (C9) (111-84-2) | 8 | <i>n</i> -Dotriacontane (C32) (544-85-4) | 1 |
| <i>n</i> -Decane (C10) (124-18-5) | 12 | <i>n</i> -Hexatriacontane (C36) (630-06-8) | 1 |
| <i>n</i> -Undecane (C11) (1120-21-4) | 12 | <i>n</i> -Tetracontane (C40) (4181-95-7) | 1 |
| <i>n</i> -Dodecane (C12) (112-40-3) | 12 | <i>n</i> -Tetratetracontane (C44) (7098-22-8) | 1 |
| <i>n</i> -Tetradecane (C14) (629-59-4) | 12 | | |
| <i>n</i> -Hexadecane (C16) (544-76-3) | 10 | | |
| <i>n</i> -Octadecane (C18) (593-45-3) | 5 | | |

Packaged 1 mL/ampul

cat.# 31222 (ea.)

No data pack available.

ASTM D2887-12 Calibration Standard

(20 components)

| | |
|--|--|
| <i>n</i> -Pentane (C5) (109-66-0) | <i>n</i> -Octadecane (C18) (593-45-3) |
| <i>n</i> -Hexane (C6) (110-54-3) | <i>n</i> -Eicosane (C20) (112-95-8) |
| <i>n</i> -Heptane (C7) (142-82-5) | <i>n</i> -Tetracosane (C24) (646-31-1) |
| <i>n</i> -Octane (C8) (111-65-9) | <i>n</i> -Octacosane (C28) (630-02-4) |
| <i>n</i> -Nonane (C9) (111-84-2) | <i>n</i> -Dotriacontane (C32) (544-85-4) |
| <i>n</i> -Decane (C10) (124-18-5) | <i>n</i> -Hexatriacontane (C36) (630-06-8) |
| <i>n</i> -Undecane (C11) (1120-21-4) | <i>n</i> -Tetracontane (C40) (4181-95-7) |
| <i>n</i> -Dodecane (C12) (112-40-3) | <i>n</i> -Tetratetracontane (C44) (7098-22-8) |
| <i>n</i> -Tetradecane (C14) (629-59-4) | |
| <i>n</i> -Pentadecane (C15) (629-62-9) | |
| <i>n</i> -Hexadecane (C16) (544-76-3) | |
| <i>n</i> -Heptadecane (C17) (629-78-7) | |

1% w/w in carbon disulfide, 1 g solution/ampul

cat.# 31674 (ea.)

5% w/w, Neat, 1 g /ampul

cat.# 31675 (ea.)

No data pack available.

Polywax Standards

These high molecular weight hydrocarbon waxes are useful for simulated distillation and other high-temperature GC work.

Volume is 1 mL/ampul.

| Compound | CAS # | Solvent | Conc. | cat.# |
|---------------|-----------|---------|-------|-------|
| Polywax 500 | 9002-88-4 | Neat | 1 g | 36224 |
| Polywax 655 | 9002-88-4 | Neat | 1 g | 36225 |
| Polywax 850 | 9002-88-4 | Neat | 1 g | 36226 |
| Polywax 1,000 | 9002-88-4 | Neat | 1 g | 36227 |

No data pack available.

RESTEK[®] REFINED

Solutions for Your Petro Analyses

- Improved, best-in-class columns
- Standards to support ASTM methods
- Industry experts at your service

www.restek.com/petro



Reference Standards

DHA

DHA PONA Standard (188 components)

This standard is a qualitative mixture of various gasoline and refinery materials prepared to provide nearly every component that may be encountered in feedstock and finished gasolines. Some oxygenates have been added to allow this blend to be used by refineries for Detailed Hydrocarbon Analysis (DHA) method setup and to fulfill requirements such as ASTM Methods D6729, D6730, and D6733.

Neat, 0.15 mL in Autosampler Vial
cat.# 30731 (ea.)

No data pack available.

For a full component list, visit www.restek.com and search for "30731".

Oxy Set-Up Blend (30 components)

Gravimetrically prepared and NIST-traceable.

Visit www.restek.com/petromixes for a component list.

2 mL prescored ampul
cat.# 33034 (ea.)

PiONA Blends

These standards are produced for refineries performing detailed hydrocarbon analysis (DHA) of crude oil feedstocks and fuels. They help with calibrating complex hydrocarbon analyses and provide the greatest number of gravimetrically determined values for quantitative calibration to help customers optimize production and maximize profitability—while fulfilling requirements such as ASTM Methods D6729, D6730, and D6733.

DHA PiONA Standard (133 components)

Neat, 0.15 mL in Autosampler Vial
cat.# 30730 (ea.)

No data pack available.

DHA Paraffins Standard (11 components)

Neat, 0.15 mL in Autosampler Vial
cat.# 30725 (ea.)

No data pack available.

DHA Isoparaffins Standard (33 components)

Neat, 0.15 mL in Autosampler Vial
cat.# 30726 (ea.)

No data pack available.

DHA Olefins Standard (26 components)

Neat, 0.15 mL in Autosampler Vial
cat.# 30727 (ea.)

No data pack available.

DHA Naphthenes Standard (26 components)

Neat, 0.15 mL in Autosampler Vial
cat.# 30728 (ea.)

No data pack available.

DHA Aromatics Standard (37 components)

Neat, 0.15 mL in Autosampler Vial
cat.# 30729 (ea.)

No data pack available. Quantity discounts not available.

D3606

ASTM D3606 Calibration Kit Without Internal Standard

Contains 25 mL each of these mixtures.

30647: ASTM D3606 Calibration Standard #1 Without Internal Standard
30648: ASTM D3606 Calibration Standard #2 Without Internal Standard
30649: ASTM D3606 Calibration Standard #3 Without Internal Standard
30650: ASTM D3606 Calibration Standard #4 Without Internal Standard
30651: ASTM D3606 Calibration Standard #5 Without Internal Standard
30652: ASTM D3606 Calibration Standard #6 Without Internal Standard
30653: ASTM D3606 Calibration Standard #7 Without Internal Standard

Contains 25 mL each of these mixtures.

cat.# 30672 (kit)

kit

ASTM D3606 Calibration Kit With MEK Internal Standard

Contains 1 mL each of these mixtures.

30654: ASTM D3606 Calibration Standard #1 With MEK Internal Standard
30655: ASTM D3606 Calibration Standard #2 With MEK Internal Standard
30656: ASTM D3606 Calibration Standard #3 With MEK Internal Standard
30657: ASTM D3606 Calibration Standard #4 With MEK Internal Standard
30658: ASTM D3606 Calibration Standard #5 With MEK Internal Standard
30659: ASTM D3606 Calibration Standard #6 With MEK Internal Standard
30660: ASTM D3606 Calibration Standard #7 With MEK Internal Standard

Contains 1 mL each of these mixtures.

cat.# 30673 (kit)

kit

ASTM D3606 Calibration Kit With sec-Butanol Internal Standard

Contains 1 mL each of these mixtures.

30661: ASTM D3606 Calibration Standard #1 With sec-Butanol Internal Standard
30662: ASTM D3606 Calibration Standard #2 With sec-Butanol Internal Standard
30663: ASTM D3606 Calibration Standard #3 With sec-Butanol Internal Standard
30664: ASTM D3606 Calibration Standard #4 With sec-Butanol Internal Standard
30665: ASTM D3606 Calibration Standard #5 With sec-Butanol Internal Standard
30666: ASTM D3606 Calibration Standard #6 With sec-Butanol Internal Standard
30667: ASTM D3606 Calibration Standard #7 With sec-Butanol Internal Standard

Contains 1 mL each of these mixtures.

cat.# 30674 (kit)

kit

ASTM D3606 Quality Control Standard with MEK Internal Standard (3 components)

| | |
|----------------------------|-----------------|
| Benzene (71-43-2) | 0.6432% vol/vol |
| 2-Butanone (MEK) (78-93-3) | 4 |
| Toluene (108-88-3) | 4.8 |

In isooctane, 25 mL/ampul

cat.# 30669 (ea.)

ASTM D3606 Quality Control Standard with sec-Butanol Internal Standard (3 components)

| | |
|---|-----------------|
| Benzene (71-43-2) | 0.6432% vol/vol |
| 2-Butanol (sec-butyl alcohol) (78-92-2) | 4 |
| Toluene (108-88-3) | 4.8 |

In isooctane, 25 mL/ampul

cat.# 30670 (ea.)

Visit www.restek.com/petromixes for our complete line of petrochemical reference standards.

For custom standards, visit www.restek.com/solutions



Reference Standards

Natural Gas and Refinery Gas Standards

- Each available in three varying concentrations.
- Mini-regulator designed specifically for these standards.

Natural Gas Standards

Available in three mixes, from lean to rich. Each has an extended list of C6+ components.

| | Natural Gas Standard #1 cat.# 34438, ea. % each compound* | Natural Gas Standard #2 cat.# 34439, ea. % each compound* | Natural Gas Standard #3 cat.# 34440, ea. % each compound* |
|--|---|---|---|
| nitrogen | 1.000 | 2.500 | 5.000 |
| carbon dioxide | 0.500 | 1.000 | 1.500 |
| methane UHP | 94.750 | 85.250 | 70.000 |
| ethane UHP | 2.000 | 5.000 | 9.000 |
| propane | 0.750 | 3.000 | 6.000 |
| isobutane | 0.300 | 1.000 | 3.000 |
| <i>n</i> -butane | 0.300 | 1.000 | 3.000 |
| isopentane | 0.150 | 0.500 | 1.000 |
| <i>n</i> -pentane | 0.150 | 0.500 | 1.000 |
| hexanes plus | 0.100 | 0.250 | 0.500 |
| Concentration | mole | mole | mole |
| Volume | 13.16 L @ 200 psig (1,379 kPa) | 13.16 L @ 200 psig (1,379 kPa) | 5.5 L @ 75 psig (517 kPa) |
| Ideal Heating Value (Dry BTU/SCF) | 1,048 gross | 1,142 gross | 1,317 gross |

Refinery Gas Standards

Available in three mixes with varying C5 unsaturates or extended C6+ components.

| | Refinery Gas Standard #1 cat.# 34441, ea. % each compound* | Refinery Gas Standard #2 cat.# 34442, ea. % each compound* | Refinery Gas Standard #3 cat.# 34443, ea. % each compound* |
|-------------------------|--|--|--|
| hydrogen | 40.750 | 12.500 | 12.500 |
| argon | 0.500 | 1.000 | 1.000 |
| nitrogen | 4.000 | 37.200 | 37.200 |
| carbon monoxide | 1.000 | 1.000 | 1.000 |
| carbon dioxide | 3.000 | 3.000 | 3.000 |
| methane | 8.500 | 5.000 | 5.000 |
| ethane | 6.000 | 4.000 | 4.000 |
| ethylene | 2.000 | 2.000 | 2.000 |
| acetylene | — | 1.000 | 1.000 |
| propane | 7.000 | 6.000 | 6.000 |
| propylene | 3.000 | 3.000 | 3.000 |
| propadiene | 0.850 | 1.000 | 1.000 |
| cyclopropane | — | 0.040 | — |
| isobutane | 6.000 | 5.000 | 5.000 |
| <i>n</i> -butane | 4.000 | 4.000 | 4.000 |
| isobutylene | 2.000 | 1.000 | 1.000 |
| 1,3 butadiene | 3.000 | 3.000 | 3.000 |
| <i>cis</i> -2-butene | 2.000 | 2.000 | 2.000 |
| <i>trans</i> -2-butene | 2.000 | 3.000 | 3.000 |
| 1-butene | 2.000 | 2.000 | 2.000 |
| 2-methyl-2-butene | — | 0.200 | 0.200 |
| isopentane | 1.000 | 1.000 | 1.000 |
| <i>n</i> -pentane | 1.000 | 1.000 | 1.000 |
| <i>cis</i> -2-pentene | — | 0.400 | 0.400 |
| <i>trans</i> -2-pentene | — | 0.160 | 0.200 |
| pentene-1 | — | 0.400 | 0.400 |
| <i>n</i> -hexane | 0.500 | 0.100 | — |
| hexanes plus | — | — | 0.100 |
| Concentration | mole | mole | mole |
| Volume | 5.2 L @ 70 psig (483 kPa) | 4.9 L @ 60 psig (414 kPa) | 4.6 L @ 60 psig (414 kPa) |

*Precise concentrations are provided on the data sheet included with each cylinder and may vary slightly from those listed here.



cylinder design

DCG Partnership Cylinders:

Size: 7.6 x 24 cm
CGA-170/110 connection
U.S. DOT Specs:
DOT-4B-240ET

Please note: This cylinder is not approved for use in Canada.

Need a regulator?

Request
cat.# 22032.



Please note: Gas standards on this page are not available in Pi-marked (TPED-rated) cylinders required for EU countries.




Sample Cylinders, Valves & Accessories

Sulfinert®-Treated High-Pressure Sample Cylinders

- Sulfinert® coating provides stable storage of sulfur and mercury at ppb levels in petroleum samples.
- Inert coating doesn't flake; more durable than PTFE.
- TPED compliant cylinders now available for shipping into EU countries.
- All cylinders have 1/4" female NPT threads on both ends.

Sulfur compounds are stable in Sulfinert®-treated stainless steel systems.



also available

Certificates are available upon request.

Sample Cylinders, High Pressure

(Stainless Steel & Sulfinert® Treated)

| Size | 1,800 psig (12,411 kPa), 304L SS | | TPED, 1,450 psig (9,997 kPa), 304L SS | |
|----------|----------------------------------|-------------------|---------------------------------------|-------------------|
| | Stainless Steel | Sulfinert Treated | Stainless Steel | Sulfinert Treated |
| 75 cc | cat.# 22921 | cat.# 24130 | cat.# 22921-PI | cat.# 24130-PI |
| 150 cc | cat.# 22922 | cat.# 24131 | cat.# 22922-PI | cat.# 24131-PI |
| 300 cc | cat.# 22923 | cat.# 24132 | cat.# 22923-PI | cat.# 24132-PI |
| 500 cc | cat.# 22924 | cat.# 24133 | cat.# 22924-PI | cat.# 24133-PI |
| 1,000 cc | cat.# 22925 | cat.# 24134 | cat.# 22925-PI | cat.# 24134-PI |
| 2,250 cc | cat.# 22926 | cat.# 21394 | cat.# 22926-PI | cat.# 21394-PI |

Sample Cylinders, Ultra-High Pressure

(Stainless Steel & Sulfinert® Treated)

| Size | 5,000 psig (34,474 kPa), 316L SS | | TPED, 4,350 psig (29,992 kPa), 316L SS | |
|--------|----------------------------------|-------------------|--|-------------------|
| | Stainless Steel | Sulfinert Treated | Stainless Steel | Sulfinert Treated |
| 150 cc | cat.# 22927 | cat.# 22111 | cat.# 22927-PI | cat.# 22111-PI |
| 300 cc | cat.# 22928 | cat.# 22112 | cat.# 22928-PI | cat.# 22112-PI |
| 500 cc | cat.# 22929 | cat.# 22113 | cat.# 22929-PI | cat.# 22113-PI |

Metering Control Valves (Stainless Steel & Sulfinert® Treated)

- Reduces pressure between sample cylinder and GC injector.
- Maintains fine metering control.
- Contains Kel-F® seat.

| Description | Stainless Steel cat.# | Sulfinert Treated cat.# |
|--|-----------------------|-------------------------|
| 3,500 psig (24,132 kPa) DOT Pressure Rating Metering Control Valve, 1/4" Male NPT x 1/4" Male NPT | 26326 | 26327 |



Sample Cylinder Valves (Stainless Steel & Sulfinert® Treated)

- Multiple valve configurations, including dip tube and rupture disks.
- Large, durable, Kel-F® seat ensures leak-free operation.
- Temperature range: -40 °C to 120 °C

| Description | Stainless Steel cat.# | Sulfinert Treated cat.# |
|---|-----------------------|-------------------------|
| 3,500 psig (24,132 kPa) DOT Pressure Rating | | |
| 1/4" Male NPT x 1/4" Male NPT | 26297 | 21400 |
| 1/4" Male NPT x 1/4" Female NPT | 26298 | 26299 |
| 1/4" Male NPT x 1/4" Male Compression | 26300 | 21401 |
| 1/4" Male NPT x 1/4" Male NPT w/5.25" Dip Tube* | 26301 | 21402* |
| 1/4" Male NPT x 1/4" Male NPT w/1,800 psi (12,411 kPa) Rupture Disc | 26302 | 26303 |
| 1/4" Male NPT x 1/4" Female NPT w/1,800 psi (12,411 kPa) Rupture Disc | 26304 | 26305 |
| Replacement Rupture Disc, 1,800 psig (12,411 kPa) | 26320 | — |
| 5,000 psig (34,474 kPa) DOT Pressure Rating | | |
| 1/4" Male NPT x 1/4" Male NPT | 26306 | 26307 |
| 1/4" Male NPT x 1/4" Female NPT | 26308 | 26309 |
| 1/4" Male NPT x 1/4" Male Compression | 26310 | 26311 |
| 1/4" Male NPT x 1/4" Male NPT w/5.25" Dip Tube* | 26312 | 26313 |
| 1/4" Male NPT x 1/4" Male NPT w/2,850 psi (19,650 kPa) Rupture Disc | 26314 | 26315 |
| 1/4" Male NPT x 1/4" Female NPT w/2,850 psi (19,650 kPa) Rupture Disc | 26316 | 26317 |
| Replacement Rupture Disc, 2,850 psig (19,650 kPa) | 26324 | — |

*To order a sample cylinder valve with dip tube, please call Customer Service at 800-356-1688, ext. 3, or contact your Restek representative. Specify dip tube length or % outage when ordering (maximum length = 5.25" / 13.3 cm). Note: End of part will not be treated after cutting tube to length.



Sample Cylinders, Valves & Accessories

Gas Sampling Valves and Sample Loops (Sulfinert® Treated)

- Ideal for samples containing low concentrations of sulfur or other active compounds.
- Sample loop sizes from 5 µL to 5 mL.

Sulfinert® treatment eliminates active sites in the valve or loop for better recovery of active compounds.

Gas Sampling Valves & Replacement Rotors (Sulfinert® Treated)

(1/16" Fittings, 0.40 mm Port Diameter; "W Type" Valve)

| Description | qty. | cat.# |
|---------------------------------------|------|-------|
| Sulfinert Gas Sampling Valve; 4-Port | ea. | 20584 |
| Sulfinert Gas Sampling Valve; 6-Port | ea. | 20585 |
| Sulfinert Gas Sampling Valve; 10-Port | ea. | 20586 |

Replacement Rotors (Not Coated)

| Description | qty. | cat.# |
|---|------|-------|
| Replacement Rotor (not coated) for 4-Port Sulfinert Gas Sampling Valve | ea. | 20587 |
| Replacement Rotor (not coated) for 6-Port Sulfinert Gas Sampling Valve | ea. | 20588 |
| Replacement Rotor (not coated) for 10-Port Sulfinert Gas Sampling Valve | ea. | 20589 |

Gas Sample Loops (Sulfinert® Treated) (1/16" fittings, for "W Type" valves)

| Description | Size | qty. | cat.# |
|---------------------------------|--------|------|-------|
| Sample Loops, Sulfinert Treated | 5 µL | ea. | 22840 |
| Sample Loops, Sulfinert Treated | 10 µL | ea. | 22841 |
| Sample Loops, Sulfinert Treated | 20 µL | ea. | 22842 |
| Sample Loops, Sulfinert Treated | 25 µL | ea. | 22843 |
| Sample Loops, Sulfinert Treated | 50 µL | ea. | 22844 |
| Sample Loops, Sulfinert Treated | 100 µL | ea. | 22845 |
| Sample Loops, Sulfinert Treated | 250 µL | ea. | 22846 |
| Sample Loops, Sulfinert Treated | 500 µL | ea. | 22847 |
| Sample Loops, Sulfinert Treated | 1 mL | ea. | 22848 |
| Sample Loops, Sulfinert Treated | 2 mL | ea. | 22849 |
| Sample Loops, Sulfinert Treated | 5 mL | ea. | 22850 |

Jumbo Syringe

Clear acrylic syringes, ideal for holding and dispensing large volumes of gas. An adjustable plunger on the O-ring ensures that the syringe is gas-tight over a long period of time. The central port is supplied with a luer lock fitting; the secondary port is supplied with a septum nut. This enables access to the gas sample for adding standards or removing a subsample. The plunger stem is detachable, making sample storage easy.

| Volume | SGE | | Restek | |
|----------|---------------|--------|--------|-------|
| | Model | cat.# | qty. | cat.# |
| 500 mL | 500MAR-LL-GT | 009910 | ea. | 21275 |
| 1,000 mL | 1000MAR-LL-GT | 009920 | ea. | 21276 |
| 2,000 mL | 2000MAR-LL-GT | 009930 | ea. | 21277 |

Syringe O-Rings

| Syringe Volume | SGE | | Restek | |
|----------------|--------|--|--------|-------|
| | cat.# | | qty. | cat.# |
| 500 mL | 032527 | | ea. | 21278 |
| 1,000 mL | 032532 | | ea. | 21279 |



Petro Essentials

Restek® Filter Base Plates

- End fittings available in brass or stainless steel.
- Base plates fit all stand alone Super Clean® gas filters offered.



22025



22026



22027

| Description | Brass | | Stainless Steel | |
|------------------------------------|-------|-------|-----------------|-------|
| | qty. | cat.# | qty. | cat.# |
| Filter Base Plate, Single-Position | ea. | 22025 | ea. | 22344 |
| Filter Base Plate, 2-Position | ea. | 22026 | ea. | 22345 |
| Filter Base Plate, 3-Position | ea. | 22027 | ea. | 22346 |

Each single-position base plate unit measures: 4" x 4" x 1 7/8" (10.2 x 10.2 x 4.8 cm)

Each 2-position base plate unit measures: 8" x 4" x 1 7/8" (10.2 x 10.2 x 4.8 cm)

Each 3-position base plate unit measures: 12" x 4" x 1 7/8" (10.2 x 10.2 x 4.8 cm)

Standard base plate inlet/outlet fittings accept 1/8" tubing. To adapt to 1/4" tubing, order 1/8" to 1/4" tube end reducers.



22020



22022



base plate

Restek® Super Clean® Gas Filter Kits and Replacements

- High-purity output ensures 99.9999% pure gas (at max. flow of 2 L/min).
- “Quick connect” fittings for easy, leak-tight filter cartridge changes.
- Glass inside to prevent diffusion; polycarbonate housing outside for safety.
- All traps measure 10 5/8" x 1 3/4" (27 x 4.4 cm).
- Each base plate unit measures 4" x 4" x 1 7/8" (10.2 x 10.2 x 4.8 cm).

| Description | qty. | cat.# |
|---|------|-------|
| Carrier Gas Cleaning Kit | kit | 22019 |
| Includes: mounting base plate, 1/8" inlet/outlet fittings, and oxygen/moisture/hydrocarbon triple gas filter | | |
| Fuel Gas Purification Kit | kit | 22021 |
| Includes: mounting base plate, 1/8" inlet/outlet fittings, and hydrocarbon/moisture fuel gas filter | | |
| Ultra-High Capacity Hydrocarbon Filter | ea. | 22030 |
| Ultra-High Capacity Moisture Filter | ea. | 22028 |
| Ultra-High Capacity Oxygen Filter | ea. | 22029 |
| Replacement Triple Gas Filter (removes oxygen, moisture, and hydrocarbons) | ea. | 22020 |
| Replacement Fuel Gas Filter (removes moisture and hydrocarbons) | ea. | 22022 |
| Helium-Specific Carrier Gas Cleaning Kit | kit | 21983 |
| Includes: mounting base plate, 1/8" inlet/outlet fittings, and oxygen/moisture/hydrocarbon helium-specific filter | | |
| Replacement Helium-Specific Gas Filter (removes oxygen, moisture, and hydrocarbons) | ea. | 21982 |
| Gas Filter Bundle Kit | kit | 22031 |
| Includes: one triple gas filter (cat.# 22020) and two fuel gas filters (cat.# 22022) | | |



Dynamic Duo (Restek® Leak Detector and ProFLOW 6000 Flowmeter)

Protect your instrument and improve data quality with this powerful pair from Restek. Checking for leaks and verifying flows before you start helps you avoid costly problems later.

| Description | qty. | cat.# |
|--|------|-------|
| Dynamic Duo Combo Pack (Restek Leak Detector and ProFLOW 6000 Flowmeter) | kit | 22654 |
| Related Products and Accessories | | |
| Leak Detector With Hard-Sided Carrying Case and Universal Charger Set (U.S., UK, European, Australian) | ea. | 22655 |
| Small Probe Adaptor for Leak Detector | ea. | 22658 |
| Restek ProFLOW 6000 Electronic Flowmeter With Hard-Sided Carrying Case | ea. | 22656 |
| Soft-Sided Storage Case for Leak Detector or ProFLOW 6000 Flowmeter | ea. | 22657 |

Restek's New Leak Detector

Redesigned and better than ever, our new leak detector is an essential tool for troubleshooting and routine maintenance of your gas chromatograph. Don't risk damaging your system or losing sensitivity; check for leaks often and protect your GC column and instrument with a Restek® leak detector!



Leak Detector Specifications

| | |
|-------------------|--|
| Detectable Gases: | Helium, nitrogen, argon, carbon dioxide, hydrogen |
| Battery: | Rechargeable lithium ion internal battery pack (12 hours normal operation) |

| | |
|------------------------------|---------------------|
| Operating Temperature Range: | 32–120 °F (0–48 °C) |
| Humidity Range: | 0–97% |
| Warranty: | One year |
| Certifications: | CE, Ex, Japan |
| Compliance: | WEEE, RoHS |

Limits of Detection

These gases can be detected with the Restek® electronic leak detector at the following leak rates:

Minimum Detectable Gas Limits and Indicating LED Color:

| |
|---|
| Helium, 1.0×10^{-5} , red LED |
| Hydrogen*, 1.0×10^{-5} , red LED |
| Nitrogen, 1.4×10^{-3} , yellow LED |
| Argon, 1.0×10^{-4} , yellow LED |
| Carbon dioxide, 1.0×10^{-4} , yellow LED |

Gas detection limits measured in atm cc/sec.

ProFLOW 6000 Flowmeter

With its wide range of capabilities, the ProFLOW 6000 flowmeter simplifies gas flow measurement in the lab. Real-time measurements can be made for various types of flow paths, including continually changing gas types.



Flowmeter Specifications:

| | |
|------------------------|---------------------|
| Type of Flowmeter: | Volumetric |
| Battery: | 2-AA |
| Operating Temp. Range: | 32–120 °F (0–48 °C) |
| Warranty: | One year |
| Certifications: | CE, Ex |
| Compliance: | WEEE, RoHS |
| Patented. | |

*Caution: The Restek® electronic leak detector is designed to detect trace amounts of hydrogen in a noncombustible environment. It is NOT designed for determining leaks in a combustible environment. A combustible gas detector should be used for determining combustible gas leaks under any condition. When using it to detect hydrogen, the Restek® electronic leak detector may only be used for determining trace amounts in a GC environment.

Dynamic Duo

Restek's Leak Detector and
ProFLOW 6000 Flowmeter

An Unbeatable Combination!



Sulfinert® and Silcosteel®-CR Treated Tubing and Fittings

Restek sets the standard in tubing for analytical and process applications. Complete your system with precleaned or treated tubing and treated fittings and valves for an inert, corrosion-resistant pathway. Restek offers the following coatings:

Silcosteel®-CR

A corrosion-resistant layer that increases the lifetime of system components in acidic environments containing hydrochloric acid, nitric acid, or seawater.

Sulfinert®

The ultimate passivation of treated surfaces; ideal for automotive exhaust testing, stack gas sampling, or analyzing for parts-per-billion levels of organosulfur compounds.

Minimum Bend Radius for Coated Tubing

| OD | Min. Bend Radius |
|-------|------------------|
| 1/16" | 1" (2.5 cm) |
| 1/8" | 2" (5.1 cm) |
| 1/4" | 4" (10.2 cm) |
| 3/8" | 6" (15.2 cm) |

did you know?

Other lengths and diameters of treated tubing are available on a custom basis.

Call for availability of lengths greater than 1,000 ft

ordering note

Required length in meters x 3.2808 = length in feet.

Treated Welded/Drawn 304 Grade Stainless Steel Tubing

Our most popular grade of tubing. Recommended for:

- Chromatography applications.
- Gas delivery systems.
- Lower pressures.
- Inert applications.

Maximum temperature of 450 °C in an inert atmosphere.

Siltek®/Sulfinert® Treated (Coiled)

| OD | ID | Wall Thickness | cat.# |
|------------------|------------------|----------------|-------|
| 0.022" (0.56 mm) | 0.011" (0.28 mm) | | 22500 |
| 0.029" (0.74 mm) | 0.021" (0.53 mm) | | 22501 |
| 1/16" (1.59 mm) | 0.010" (0.25 mm) | | 22502 |
| 1/16" (1.59 mm) | 0.020" (0.51 mm) | | 22503 |
| 1/16" (1.59 mm) | 0.030" (0.76 mm) | | 22504 |
| 1/16" (1.59 mm) | 0.040" (1.02 mm) | | 22505 |
| 1/8" (3.18 mm) | 0.085" (2.16 mm) | 0.020" | 22506 |
| 1/4" (6.35 mm) | 0.210" (5.33 mm) | 0.020" | 22507 |

Minimum order is 5 ft.

*The availability of long lengths is subject to inventory constraints. Lead times may vary depending on the continuous length needed. Please inquire before ordering. Maximum continuous lengths are: 200** ft (cat.#22500-22501), 2,000 ft (cat.# 22502-22505), 1,500 ft (cat.# 22506), and 750 ft (cat.# 22507).

** Contact us if longer length is needed for cat.# 22500-22501

Treated Seamless 316L Grade Stainless Steel Tubing

High durability tubing. Recommended for:

- Inert applications.
- High temperatures.
- High pressures.
- Corrosive environments.
- Zero bleed.

Siltek®/Sulfinert® Treated (Coiled)

| OD | ID | Wall Thickness | cat.# |
|----------------|------------------|----------------|-------|
| 1/8" (3.18 mm) | 0.055" (1.40 mm) | 0.035" | 22508 |
| 1/4" (6.35 mm) | 0.180" (4.57 mm) | 0.035" | 22509 |
| 3/8" (9.52 mm) | 0.277" (7.04 mm) | 0.049" | 22914 |

Silcosteel®-CR Treated (Coiled)

| OD | ID | Wall Thickness | cat.# |
|----------------|------------------|----------------|-------|
| 1/8" (3.18 mm) | 0.055" (1.40 mm) | 0.035" | 22896 |
| 1/4" (6.35 mm) | 0.180" (4.57 mm) | 0.035" | 22897 |
| 3/8" (9.52 mm) | 0.277" (7.04 mm) | 0.049" | 22915 |

*The availability of long lengths is subject to inventory constraints. Lead times may vary depending on the continuous length needed. Please inquire before ordering.

1/8" OD: 5 ft to 1,500 ft in one continuous coil; 1/4" OD: 5 ft to 750 ft in one continuous coil; 3/8" OD: 5 ft to 250 ft in one continuous coil. Longer lengths will be more than one coil.

Treated Seamless Electropolished 316L Grade Stainless Steel Tubing

Our highest performing tubing. Recommended for:

- Demanding/corrosive environments.
- High temperatures.
- Ultimate inertness.

Siltek®/Sulfinert Treated (Coiled)

| OD | ID | Wall Thickness | cat.# |
|----------------|------------------|----------------|-------|
| 1/8" (3.18 mm) | 0.085" (2.16 mm) | 0.020" | 22538 |
| 1/4" (6.35 mm) | 0.180" (4.57 mm) | 0.035" | 22539 |

Silcosteel®-CR Treated (Coiled)

| OD | ID | Wall Thickness | cat.# |
|----------------|------------------|----------------|-------|
| 1/8" (3.18 mm) | 0.085" (2.16 mm) | 0.020" | 22536 |
| 1/4" (6.35 mm) | 0.180" (4.57 mm) | 0.035" | 22537 |

1/8" OD: 5 ft to 95 ft in one continuous coil; 1/4" OD: 5 ft to 300 ft in one continuous coil. Longer lengths will be more than one coil.

Note: required length in meters x 3.2808 = length in feet.



Swagelok® Fittings (Siltek®/Sulfinert® & Silcosteel®-CR Treated)

- Full line of treated 1/16", 1/8", and 1/4" fittings.
- Silcosteel®-CR treatment enhances corrosion resistance by 10x or more.
- For custom treatment on any Swagelok® fitting or other system parts not listed here, call us or contact your Restek representative.

| | | Siltek/Sulfinert Treated | | Silcosteel-CR Treated | | |
|------------------|---|--------------------------|------|-----------------------|------|-------|
| Fitting Type | Size | Swagelok # | qty. | cat.# | qty. | cat.# |
| Union | ¹ / ₁₆ " | SS-100-6 | ea. | 22540 | ea. | 22575 |
| | ¹ / ₈ " | SS-200-6 | ea. | 22541 | ea. | 22576 |
| | ¹ / ₄ " | SS-400-6 | ea. | 22542 | ea. | 22577 |
| | ³ / ₈ " | SS-600-6 | ea. | 22909 | ea. | 22904 |
| Tee | ¹ / ₁₆ " | SS-100-3 | ea. | 22543 | ea. | 22578 |
| | ¹ / ₈ " | SS-200-3 | ea. | 22544 | ea. | 22579 |
| | ¹ / ₄ " | SS-400-3 | ea. | 22545 | ea. | 22580 |
| | ³ / ₈ " | SS-600-3 | ea. | 22910 | ea. | 22905 |
| Reducing Union | ¹ / ₈ " to ¹ / ₁₆ " | SS-200-6-1 | ea. | 22546 | ea. | 22581 |
| | ¹ / ₄ " to ¹ / ₁₆ " | SS-400-6-1 | ea. | 22547 | ea. | 22582 |
| | ¹ / ₄ " to ¹ / ₈ " | SS-400-6-2 | ea. | 22548 | ea. | 22583 |
| | ³ / ₈ " to ¹ / ₄ " | SS-600-6-4 | ea. | 22911 | ea. | 22906 |
| Elbow | ¹ / ₈ " | SS-200-9 | ea. | 22549 | ea. | 22584 |
| | ¹ / ₄ " | SS-400-9 | ea. | 22550 | ea. | 22585 |
| Plug | ¹ / ₈ " | SS-200-P | ea. | 22573 | ea. | 22620 |
| | ¹ / ₄ " | SS-400-P | ea. | 22574 | ea. | 22597 |
| Cross | ¹ / ₈ " | SS-200-4 | ea. | 22551 | ea. | 22586 |
| | ¹ / ₄ " | SS-400-4 | ea. | 22552 | ea. | 22587 |
| Tube End Reducer | ¹ / ₈ " to ¹ / ₁₆ " | SS-100-R-2 | ea. | 22553 | ea. | 22588 |
| | ¹ / ₄ " to ¹ / ₁₆ " | SS-100-R-4 | ea. | 22554 | ea. | 22589 |
| | ¹ / ₈ " to ¹ / ₄ " | SS-400-R-2 | ea. | 22555 | ea. | 22590 |
| | ¹ / ₄ " to ¹ / ₈ " | SS-200-R-4 | ea. | 22556 | ea. | 22591 |
| Port Connector | ¹ / ₈ " | SS-201-PC | ea. | 22557 | ea. | 22592 |
| | ¹ / ₄ " | SS-401-PC | ea. | 22558 | ea. | 22593 |
| | ¹ / ₈ " to ¹ / ₄ " | SS-401-PC-2 | ea. | 22559 | ea. | 22594 |
| Male Connector | ¹ / ₈ " to ¹ / ₈ " NPT | SS-200-1-2 | ea. | 22561 | ea. | 22595 |
| | ¹ / ₄ " to ¹ / ₄ " NPT | SS-400-1-4 | ea. | 22562 | ea. | 22596 |
| | ¹ / ₁₆ " to ¹ / ₈ " NPT | SS-100-1-2 | ea. | 22563 | ea. | 22610 |
| | ¹ / ₈ " to ¹ / ₄ " NPT | SS-200-1-4 | ea. | 22564 | ea. | 22611 |
| | ¹ / ₄ " to ¹ / ₈ " NPT | SS-400-1-2 | ea. | 22565 | ea. | 22612 |
| | ³ / ₈ " to ³ / ₈ " NPT | SS-600-1-6 | ea. | 22912 | ea. | 22907 |
| Female Connector | ¹ / ₈ " to ¹ / ₈ " NPT | SS-200-7-2 | ea. | 22566 | ea. | 22613 |
| | ¹ / ₄ " to ¹ / ₄ " NPT | SS-400-7-4 | ea. | 22567 | ea. | 22614 |
| | ¹ / ₄ " to ¹ / ₈ " NPT | SS-400-7-2 | ea. | 22568 | ea. | 22615 |
| | ¹ / ₈ " to ¹ / ₄ " NPT | SS-200-7-4 | ea. | 22569 | ea. | 22616 |
| Bulkhead Union | ¹ / ₈ " | SS-200-61 | ea. | 22570 | ea. | 22617 |
| | ¹ / ₄ " | SS-400-61 | ea. | 22571 | ea. | 22618 |



union



tee



reducing union



elbow



plug



cross



tube end reducer



port connector



male connector



female connector



bulkhead union



100%
Satisfaction
Guaranteed

Sky®

Inlet Liners

True Blue Performance

Exceptionally inert Sky® inlet liners with state-of-the-art deactivation improve trace-level analysis—and now come with a **100% satisfaction guarantee!**

* 100% SATISFACTION GUARANTEE: If your Sky® inlet liner does not perform to your expectations for any reason, simply contact Restek® Technical Service or your local Restek® representative and provide a sample chromatogram showing the problem. If our GC experts are not able to quickly and completely resolve the issue to your satisfaction, you will be given an account credit or replacement product (same cat.#) along with instructions for returning any unopened product. (Do not return product prior to receiving authorization.) For additional details about Restek's return policy, visit www.restek.com/warranty



Ask Our Experts How to Improve Your Analyses!

At Restek, we have invested in highly focused product development and a world-class team of industry professionals. We are dedicated to bringing you the innovative chromatography solutions you need to make better process-control decisions.



Jan Pijpelink

Market Research Manager

Jan is a 30-plus-year veteran of the petrochemical industry with extensive international experience in petro laboratories and with process applications throughout North America, Europe, and Asia. Jan leads our scientific collaborations and key industry partnerships in the petro market. Jan also serves as chairman of ASTM Section D02.D0.03

E-mail: jan.pijpelink@restek.com



Jaap de Zeeuw

International GC Specialist

Jaap is a renowned chromatographer with over 30 years of petrochemical experience, including 28 years with Varian/Chrompack where he focused primarily on industrial analysis issues. Jaap is widely published and travels extensively for seminars and international workshops.

E-mail: jaap.dezeeuw@restek.com



Barry Burger

Sr. Petrochemical Innovations Chemist

Barry has more than 30 years of chromatography experience and has been a voting member of the ASTM D2 committee for over 12 years. He specializes in petrochemical applications and the development of capillary, packed, and micropacked GC columns.

E-mail: barry.burger@restek.com



Katarina Oden

Petrochemical Innovations Chemist

Katarina came to Restek with a B.S. in Chemistry from the University of Ljubljana, Slovenia, and research experience in both commercial and university labs. After three years as a QA Analyst, she joined Innovations, where she supports the petroleum and chemical markets.

E-mail: katarina.odan@restek.com

RESTEK REFINED
Solutions for Your Petro Analyses

- Improved, best-in-class columns
- Standards to support ASTM methods
- Industry experts at your service

www.restek.com/petro


PATENTS & TRADEMARKS

Restek® patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks appearing in Restek® literature or on its website are the property of their respective owners. The Restek® registered trademarks used here are registered in the United States and may also be registered in other countries.



Lit. Cat.# PCBR1195D-UNV

© 2014 Restek Corporation. All rights reserved.
Printed in the U.S.A.



CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

1-800-267-8103 • sales@chromspec.com • tech@chromspec.com