

# "New Mechanisms of Solid Phase Extraction to Improve Your Analytical Results"

UCT • 2731 Bartram Road • Bristol • Pennsylvania 19007 • USA • 800.385.3153 • 215.781.9255 • Fax: 215.785.1226 • www.unitedchem.com

## **Types of Sorbent-Analyte Interactions**

- Polar
- Non-polar
- Ion-exchange
- Covalent
- Copolymeric

**Polar Extractions** 

- Also called hydrophilic or normal phase
- Unequal distribution of electrons
- Involves hydrogen bonding, pi-pi and dipole/dipole interactions
- Sorbents silica, diol, diethylamino, cyanopropyl
- Applications lipids, oil additives, carbohydrates, phenols. oil soluble vitamins
- Analytes amines, hydroxyls, carbonyls, aromatic rings. heteroatoms (O. S. N. P)
- Matrix non-polar organic
- . Elution solvents medium to high polarity

#### **Non-Polar Extractions**

- Also called hydrophobic or reverse phase
- Interactions between sorbent C-H bonds and analyte C-H bonds Involves van der Waals / dispersion forces
- Sorbents C2, C3,C4, iC4, tC4, C5, C6, C7, C8, C10, C12,
- C18, C20, C30 phenyl and cyclohexyl
- Applications drugs of abuse, TDM, pesticides
- Analytes protonated / neutral state, aromatics & alkyl chains
- · Matrix biologicals, water, aqueous buffers
- Elution solvents typically non-polar to moderately polar

#### Ion Exchange Mechanisms

- Ionic interactions occur between charged sorbent & analyte of
- · pH is manipulated to ionize analytes functional group
- · Ionic bonds are strong & retain analyte
- · Hydrophobic interferences washed away with organic solvents
- Polar interferences removed with aqueous or weak aqueous /
- Elute solvents containing stronger counterions or by changing pH
- For ionic/hydrophobic analytes, elute by simultaneously disrupting both

### **Cation Exchange Extractions**

- Cation exchange sorbents negatively charged
   Basic analytes manipulated to carry positive charge
   Opposites attract forming strong bonds
- Benzenesulfonic acid (strong)
- Propvlsulfonic acid (strong) Carboxylic acid (weak) Applications include basic drugs, catecholamines,
- pharmaceuticals, herbicides
  - Pyrimidines (cations)
- · Basic elution solvents to neutralize analyte

## **Anion Exchange Extractions**

- Acidic analytes manipulated to carry negative charge
- Opposites attract forming strong bonds
   Sorbents
- 1°. 2° amine Quaternary amine (strong) Aminopropyl (weak)
- · Applications include phosphates, acidic drugs, organic acids
- Analytes
- Carboxylic acids
- Matrix aqueous
- Acidic elution solvents to neutralize analyte

### **Copolymeric Extractions**

Sulfonic acids (cations)

- · Reverse phase sorbent with cation OR anion exchange
- Acidic, basic & neutral analyte applications
- Matrix aqueous
- Selective washes
- · Elution solvents mixture of organics with acid or base
- · Superior sample clean up


# pKa, pH & Ionization

#### % of Compound in Ionic State

Functionality	Ionization State			away fr at pKa		
Acid	Anion (-)	1	9	50	91	99
Base	Cation (+)	99	91	50	9	1

## **Relative Counter ion Selectivity**

Larger numbers reflect greater ability of the ion to displace other






<u>Cations</u>		<u>Anions</u>	<u>Anions</u>		
Ba <sup>2+</sup> Ag <sup>2+</sup> Pb <sup>2+</sup> Cu <sup>+</sup> Kg <sup>2+</sup> Hg <sup>2+</sup> Cd <sup>2+</sup> Cd <sup>2+</sup> Cd <sup>2+</sup> CO <sup>2+</sup> CS <sup>2+</sup> Rb <sup>+</sup> Kt Fe <sup>2+</sup> Mn <sup>2+</sup> NH <sup>4</sup> Ht Li	8.7 7.6 7.5 7.2 5.3 4.9 3.0 2.9 2.8 2.7 2.7 2.6 2.5 2.5 2.5 2.5 2.5 2.19 1.0 0.8	Benzene Sulfonate Citrate I Phenate HSO <sub>4</sub> CIO <sub>3</sub> NO <sub>3</sub> Br CN HSO BrO NO <sub>2</sub> CI HCO <sub>3</sub> IO3, Formate Acetate Propionate F	500 220 175 110 85 74 65 50 28 27 27 27 24 22 6.0 5.5 4.6 3.2 2.6 1.0		
	Standard cati	ion exchange counter ion			

### **SPE Amine Scavenger**

**Purification of Small Molecule Libraries** by Pharmasil® Ion Exchange SPE



UCT Column Part Number: CUBCX156 Sorbent Amount: 500mg Column Volume: 6 ml

#### **SPE TFAA Removal**

**Purification of Small Molecule Libraries** TFAA Removal by Pharmasil® Ion Exchange SPE

$$H_2 H_2 H_2 + C - C - C - N - (CH_3)_3OH^-$$

UCT Column Part Number: CHQAX156 Sorbent Amount: 500mg Column Volume: 6 mL


#### **SPE Metal Removal**

Purification of Small Molecule Libraries Palladium (Pd) Removal by Pharmasil® Ion Exchange SPE

UCT Column Part Number: CUTAX156 Sorbent Amount: 500mg Column Volume: 6 ml

# Benzondiazepines By Polar Reverse Phase





#### **Chiral Solid Phase Extraction**

Joint venture between UCT and ENS


ENS nanostructured sorbent does not require thousands of theoretical plates to get an effective separation

Chiral Enrichment of DL-Norvaline as a Function of pH Using ENS 1 mL Polar Cationic SPE Cartridges


Extraction Protocol Analyte samples were prepared at a concentration of 10 mg/mL (m/v) in 50:50 Ethanol/Water pH adjusted to values 2 through 6 extract a total volume of 1 mL analyte solution through an ENS 1mL Polar Cationic SPE cartridge. Each cartridge was then aspirated with Nitrogen to flush out any residual analyte solution. A centrifugal evaporator was used residual alialyte solution. A cellinitugal evaporator was to concentrate the samples, which were then reconstitut in 70:30 Water/Methanol mobile phase to be analyzed using the Chirobiotic™ TAG column.


Chirobiotic™ TAG 4.6 x 250 mm Mobile Phase: 70:30 Water/Methanol Flow Rate: 1.0 mL/min

HPI C




The chiral enrichment of norvaline is enhanced by lowering the pH of the extraction solvent At a low pH acids on the surface of the sorbent are not charged, and the hydrophobic amin acid side chain drives adsorption and extraction.





#### Chiral Chromatograph - Norvaline before and after enantioenrichment



	Retention Time		Peak Area	EE
Ī	Enantiomer 1	Enantiomer 2	Ratio	
Before SPE	2.2067	4.0192	50/50	0
After SPE	2.2508	4.3708	62/38	24

Note: The difference in peak height for the first enantiomer results from