



ENVIRO




FOOD
SAFETY



UCT

ENVIRO-CLEAN[®]

SOLID PHASE EXTRACTION
APPLICATIONS MANUAL



www.unitedchem.com

Table of Contents

Extraction Hints.....	6
-----------------------	---

ENVIRO-CLEAN® ENVIRONMENTAL APPLICATIONS

Drinking Water Analyses

EPA Method 506 – Phthalate and Adipate Esters in Drinking Water.....	8-10
EPA Method 508.1 - Rev 2.0 Chlorinated Pesticides, Herbicides, and Organohalides.....	11-12
EPA Method 509 - Determination of Ethylene Thiourea (ETU) in Water Using Gas Chromatography with a Nitrogen-Phosphorus Detector	13-15
EPA Method 515.2 – Determination of Chlorinated Acids in Water Using Liquid-Solid Extraction and Gas Chromatography with an Electron Capture Detector.....	16-21
EPA Method 521 - The Determination of Nitrosamines in Drinking Water	22-25
EPA Method 522 - Determination of 1,4-dioxane in drinking water by Solid-Phase Extraction (SPE) and Gas Chromatography/Mass Spectrometry (GC/MS) with Selected Ion Monitoring (SIM)	26-30
EPA Method 523: Determination of Triazine Pesticides and their Degradates in Drinking Water by Gas Chromatography/Mass Spectrometry (GC/MS)	31-33
EPA Method 525.2 - Extraction of Organic Compounds in Drinking Water	34-38
EPA Method 525.3 Version 1.0 Semivolatile Organic Compounds In Drinking Water By Solid-Phase Extraction And Capillary Column (GC/MS)	39-48
EPA Method 528 - Determination of Phenols in Drinking Water by Solid-Phase Extraction (SPE) and Capillary Column Gas Chromatography (GC/MS).	49-50
EPA Method 529 - Determination of Explosives And Related Compounds in Drinking Water by Solid-Phase Extraction and Capillary Column Gas Chromatography/Mass Spectrometry (GC/MS)	51-52
EPA Method 532 - Determination of Phenylurea Compounds in Drinking Water by Solid-Phase Extraction and High Performance Liquid Chromatography with UV Detection	53-60
EPA Method 535 - Measurement of Chloroacetanilide and Other Acetamide Herbicide Degradates in Drinking Water by SPE	61-63
EPA Method 537 – Determination of Selected Perfluorinated Alkyl Acids in Drinking Water by Solid-Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS).....	64-69
EPA Method 548.1 - Endothall.....	70-73
EPA Method 549.2 - Diquat and Paraquat.....	74-77
EPA Method 550.1 - Determination of Polycyclic Aromatic Hydrocarbons in Drinking Water by Liquid-Solid Extraction and HPLC with Coupled Ultraviolet and Fluorescence Detection.....	78-79
EPA Method 552 - Determination of Haloacetic Acids and Dalapon in Drinking Water by SPE and GC/ECD	80-86
EPA Method 553 - Determination of Benzidines and Nitrogen-Containing Pesticides in Water Liquid-Solid Extraction and Reverse Phase High Performance Liquid Chromatography/Particle Beam/Mass Spectrometry	87

Water Analyses

Bisphenol A Analysis in Water by GC/MS	89-91
Bisphenol A Analysis in Water by LC/MS/MS	92-94
Determination of Hormones in Water by Solid-Phase Extraction (SPE) and Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS).....	95-105
Determination of Diesel Range Organics Using GC/FID by Method 8015D.....	106-110
Polycyclic Aromatic Hydrocarbons from a Water Matrix	111-112
Organochlorine Pesticides and Polychlorinated Biphenyls by Solid-Phase Extraction.....	113-118
Solid Phase Extraction of Pesticides in Waters Using Graphitized Carbon Black	119-122
Determination of Carbonyl Compounds in Water by DNPH Derivatization and HPLC/UV	123-127
EPA Method 1664 - Oil and Grease - Hexane Extractable Material (HEM)	128-130
EPA Method 8323 - SPE and Determination of Organotin by Micro-Liquid Chromatography Electrospray Ion Trap MS.....	131-132
EPA Method 8330 - Analysis of Nitroaromatics, Nitramines and Nitrate Ester Compounds-Explosive and Explosive Residue Compounds	133-137
EPA Method 608 ATP	138-141
Extraction of Metals	142
Ion Exchange Sorbents for Metals Extraction Analysis & Sorbent Use Selection Guide.....	143-146

Agricultural & Food Safety

Pesticides in Fatty Matrices Extraction Procedure	148-150
Analysis of Tobacco Alkaloids.....	151-152
Sample Preparation of Glyphosate (N-phosphonomethyl glycine) and Glufosinate (RS)-2-amino-4-(hydroxyl-methyl-phosphoryl) butanoic acid by Solid-Phase Anion Exchange Extraction	153-154
Extraction of Phenoxyacetic Acid Herbicides from Soil by GC-MS or LC/MS/MS	155-157
QuEChERS Extraction and Clean-Up of Pesticides from Olive Oil	158
QuEChERS Multiresidue Pesticide Method for the Determination of Multiple Pesticides in Wines	159-165
Determination of Pesticides in Red Wine by QuEChERS Extraction, Quick QuEChERS Clean-up, and LC/MS/MS Detection.....	166-171
Flukicides / Anthelmintics by QuEChERS	172-173
Streamlined Method for the Determination of More Than 100 Veterinary Drugs in Animal Tissue Using Dispersive-SPE Clean-up and LC-MS/MS Detection.....	174-178
Acrylamide by QuEChERS Extraction with LC/MS/MS Detection.....	179-180
Trichothecene Type A & B Analysis in Wheat and Corn Using the QuEChERS Approach	181-182
Multiresidue Analysis in Cereal Grains Using Modified QuEChERS Method with UPLC-MS	183-188
Extraction of Pesticides from Tomato Using the QuEChERS Approach	189-190
Extraction of Polycyclic Aromatic Hydrocarbons from Fish Using the QuEChERS Approach	191-192
Determination of Polycyclic Aromatic Hydrocarbons (PAHs) in Seafood Using GC/MS.....	193-197
QuEChERS Analysis of Miticides and other Agrochemicals in Honey Bees, Wax or Pollen	198-200
QuEChERS Pesticide Analysis for Fresh Produce by GC/MS/MS.....	201-205
Analysis of Cyromazine in Poultry Feed Using a QuEChERS Approach	206-207
Determination of Anthelmintic Drug Residues in Milk Using Ultra High Performance Liquid Chromatography	208-211
EURL-FV Multiresidue Method Using QuEChERS.....	212-217
Extraction of Pyrethrin and Pyrethroid Pesticides from Fish Using the QuEChERS Approach.....	218-219
QuEChERS-Based LC/MS/MS Method for Multiresidue Pesticide Analysis in Fruits and Vegetables	220-229
A Summary of U.S. F.D.A. LIB4465: Collaboration of the QuEChERS Procedure for the Multiresidue Determination of Pesticides in Raw Agricultural Commodities by LC-MS/MS	230-237
Modified QuEChERS Procedure for the Analysis of Bisphenol A in Canned Food Products	238-250

QuEChERS Sample Preparation for the Analysis of Pesticide Residue in Olives.....	251-255
Determination of Carbendazim in Orange Juice Using QuEChERS with LC/MS/MS Detection	256-260
Determination of Pesticide Residues in Marijuana and Tea by QuEChERS and LC-MS/MS	261-265
Determination of Pesticides in Strawberries by QuEChERS Extraction, Quick QuEChERS Clean-up, and GC/MS Detection	266-268
Multiresidue Pesticide Analysis of Botanical Dietary Supplements Using SPE Clean Up and GC-Triple Quadrupole MS/MS.....	269-272
ChloroFiltr [®] : A Novel Sorbent for Chlorophyll Removal Using QuEChERS.....	273-275
Analysis of 136 Pesticides in Avocado Using a Modified QuEChERS Method with LC-MS/MS and GC-MS/MS.....	276-288
QuEChERS Determination of 2,4- Dichlorophenoxyacetic Acid (2,4-D) and Other Acidic Herbicides by LC-MS/MS.....	289-292

Clean-Up and Fractionation

Clean-Up of Organochlorine Pesticides and PCB Extracts Using Florisil [®]	294
Fractionation of Aliphatic and Aromatic Hydrocarbons Using ENVIRO-CLEAN [®] EPH Silica	295-296
Removal of Sulfur from Environmental Samples Using Copper Beads.....	297
Determination of Phthalates Leached from Toys into Artificial Saliva.....	298-300

Extraction Hints



Extraction Hints

- Verify sample application pH. Analytes that are not in their proper form (i.e., neutral or charged), will not effectively bind to the sorbent and may result in low or erratic recoveries. The pH of deionized water cannot be correctly determined using pH paper. Use of a calibrated pH probe is necessary.
- Always pre-rinse the cartridge with the strongest solution the cartridge will see to ensure the cleanest extraction of your eluate.
- Do not allow the sorbent to completely dry out between conditioning steps or before sample application. To ensure properly solvated cartridges, apply each solvent immediately after the previous solvent. Improperly conditioned cartridges may lead to erratic recoveries.
- Prior to elution, fully dried cartridges will ensure optimal analyte recovery. To confirm cartridge dryness, touch the sides of the cartridge at the sorbent level at full vacuum. Cartridges should feel about room temperature but not cool. If the cartridge feels cool, water is probably present. Continue drying the cartridge unless otherwise specified in the application note.
- Elution rates and soak times specified in the applications are critical for acceptable and consistent recoveries. Hint: When in doubt, slower is always better.
- Always use fresh ammonium hydroxide (NH_4OH). NH_4OH rapidly loses its strength when exposed to air. Weak NH_4OH may lead to erratic recoveries.
- Some analytes are slightly volatile. Closely monitor eluate concentration to prevent loss of analyte. Hint: Higher water bath temperatures and lower nitrogen flow rates usually provide optimal results.

Drinking Water Analyses





Determination of Phthalate and Adipate Esters in Drinking Water by Liquid-Solid Extraction and Gas Chromatography with Photoionization Detection

UCT Part Numbers:

ECUNIC18 (1100 mg C18, 83 mL cartridge) or **EEC181M6G** (1000 mg C18, 6 mL glass cartridge)

ECSS156 (Drying Cartridge 5 g, Na₂SO₄, 6 mL cartridge)

FLORISIL PR® Clean-up Cartridge EUFLS12M15 (2000 mg, 15 mL cartridge)

Alumina Clean-up Cartridge EUALN1M6 (1000 mg, 6 mL cartridge)

ENVIRO-CLEAN® Zero-Blank™ filter cartridge (optional)

EPA Method 506 Revision 1.1*

March 2012

Procedure

1. Cartridge Activation

- a) Place **ECUNIC18** cartridge(s) on the manifold
- b) Add a 10 mL aliquot of methylene chloride
- c) Slowly draw solvent completely through the cartridge
- d) Repeat with a second 10 mL aliquot of methylene chloride
- e) Hold for 1 minute then draw through
- f) Dry cartridge by drawing air at full vacuum for 2-3 minutes

Optional: To reduce background contamination from laboratory air, use of the **ENVIRO-CLEAN® Zero-Blank™** filter cartridge is highly recommended during drying steps 1) f) and 2) c)

- g) Add 10 mL of methanol to the cartridge
- h) Draw through leaving a thin layer on the frit

Note: Do not allow the cartridge to go dry otherwise repeat step i

- i) Add a second 10 mL portion of methanol
- j) Wait 1-2 minutes to activate sorbent then draw through to level of frit
- k) Add 10 mL reagent water to the cartridge and draw through until meniscus reaches the top of the frit
- l) Cartridge is now ready for sample extraction

2. Sample Extraction

- a) Add the 1 liter water sample to the cartridge and draw through over a period of about 20–30 minutes (fast drip)
- b) Rinse sample bottle and cartridge with a small volume of reagent water then add to cartridge
- c) Dry cartridge under full vacuum for 10 minutes (**Optional**: Use **Zero-Blank™** filter cartridge to prevent airborne contamination)

3. Cartridge Elution

- a) Place a collection vial in the vacuum manifold
- b) Add 5 mL of methylene chloride to the sample bottle and swirl
- c) Using a disposable glass pipette transfer the methylene chloride to the cartridge
- d) Collect dropwise
- e) Repeat this procedure with another 5 mL of methylene chloride

4. Eluate Drying and Concentration

- a) Pour eluant through a 3 gram bed of anhydrous sodium sulfate (or use Drying Cartridge **ECSS156**) and collect
- b) Rinse vial and sodium sulfate with a 3 mL aliquot of methylene chloride
- c) Repeat rinse using an additional 3 mL aliquot of methylene chloride
- d) Evaporate with a gentle stream of N₂ to 1 mL
- e) If sample is clean proceed to GC analysis
- f) If extract requires clean-up for phthalates esters proceed to Florisil or Alumina clean-up stages

Extract Clean-up---Florisil or Alumina

Clean-up procedures are not required for clean drinking water. Under certain circumstances for dirty water, a Florisil or Alumina clean-up may be needed.

Florisil Column Clean-up for Phthalate Esters

1. Add a 1 cm layer of anhydrous sodium sulfate to the top frit of Florisil cartridge **EUFLS12M15**
2. Flush cartridge with 20 mL of hexane leaving enough to cover the frit
3. Add sample extract to the cartridge then rinse vial with 2 mL of hexane
4. Add 20 mL of hexane to the cartridge and elute. Discard the hexane solution

5. Elute using 20 mL of 20% diethyl ether in hexane into a 500 mL K-D flask* equipped with a 10 mL concentrator tube. Elute at a rate of about 2 mL/minute
6. No solvent exchange is required
7. Concentrate eluate in hot water bath at 85° C to 1 mL
8. Sample is ready for GC analysis

Alumina Column Clean-up (neutral) for Phthalate Esters

1. Add a 1 cm layer of anhydrous sodium sulfate to the top frit of a Alumina cartridge **EUALN1M6**
2. Flush cartridge with 10 mL of hexane leaving enough to cover sodium sulfate
3. Add sample extract to the cartridge. Rinse vial with 2 mL of hexane and add to the cartridge
4. Add 15 mL of hexane to the cartridge and elute. Discard hexane solution
5. Elute at a rate of about 2 mL/minute using 15 mL of 20% diethyl ether in hexane
6. Collect eluate in a 500 mL K-D flask equipped with a 10 mL concentrator tube
7. No solvent exchange is required
8. Concentrate eluate in hot water bath at 85° C to 1 mL
9. Sample is ready for GC analysis

*The analyst should refer to EPA Method 506 "Determination of Phthalate and Adipate Esters in Drinking Water by Liquid-Solid Extraction and Gas Chromatography with Photoionization Detection", Revision 1.1 Issued 1995, F.K. Kawahara and J.W. Hodgeson, Ed. By D. J. Munch US EPA, National Exposure Research Laboratory, Office of Research and Development, US Environmental Protection Agency, Cincinnati, Ohio 45268

**K-D Flask = Kuderna-Danish Flask



Determination of Chlorinated Pesticides, Herbicides, and Organohalides by Liquid-Solid Extraction and Electron Capture Gas Chromatography*

UCT Product Numbers:

ECUNIC18 (1100mg C18, 83 mL cartridge)

ECSS15M6 (5000 mg Na₂SO₄, 6 mL cartridge)

EPA Method 508.1 Rev 2.0

August 2009

Procedure

1. Condition Cartridge

- a) Wash the cartridge with 5 mL of ethyl acetate (EtOAc) and 5 mL of methylene chloride (MeCl₂). Let the cartridge drain dry after each wash
- b) Add 10 mL of methanol (MeOH) to the cartridge then draw through to top of frit
- c) Rinse the cartridge with 10 mL of reagent water and draw through to top of frit

Note: Do not let the cartridge go dry after addition of methanol otherwise repeat at step b)

2. Sample Extraction

- a) Samples must be dechlorinated and pH adjusted to <2
- b) Add 5 mL of methanol to the 1-liter sample and mix well
- c) Add 50 µL of surrogate compound to the water sample and mix well
- d) Draw the water sample through the cartridge at a fast drip
- e) Dry the cartridge by drawing air or N₂ through for about 10 minutes

3. Extract Elution

- a) Insert an eluate collection tube in the vacuum manifold
- b) Rinse the inside walls of the sample bottle using 10 mL EtOAc then transfer solvent to the cartridge
- c) Collect dropwise
- d) Rinse the inside walls of the sample bottle using 10 mL of MeCl₂ then transfer to the cartridge
- e) Collect dropwise combining eluents

4. Sample Drying

- a) Pour the combined eluents through drying tube **ECSS15M6**
- b) Rinse the drying tube with 2 x 3 mL portions of MeCl_2
- c) Concentrate the extract to about 0.8 mL under a gentle stream of N_2 while in a water bath
- d) Rinse the inside walls of the concentrator tube three times with EtOAc during evaporation
- e) Add internal standard (IS)
- f) Adjust the final volume of the extract to exactly 1.0 mL
- g) Sample is ready for GC analysis

For complete details on Method 508.1 , rev 2.0 "Determination of Chlorinated Pesticides, Herbicides, and Organohalides by Liquid-Solid Extraction and Electron Capture Gas Chromatograph", the analyst is referred to: J. W. Eichelberger rev 1.0, 1994 and J. Munch, rev 2.0, 1995, Environmental Monitoring Systems Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268



Determination of Ethylene Thiourea (ETU) in Water Using Gas Chromatography with a Nitrogen-Phosphorus Detector

UCT Product Numbers:

CLEAN-ELUTE™ (25,000 mg diatomaceous earth, 200 mL cartridge)

EPA Method 509

October 2009

Method 509 is a gas chromatographic method for the determination of ethylene thiourea (ETU) (CAS 96-45-7) a metabolic byproduct of the ethylene bisdithiocarbonate (EBDC) fungicides in water. Toxicological studies indicate that ETU may produce goitrogenic, tumorigenic, and teratogenic effects in laboratory animals, raising the concern that residues may be found in agricultural commodities. The method uses a packed column of diatomaceous earth to capture the analyte before elution with methylene chloride. Confirmation is made using a nitrogen-phosphorous detector or a mass spectrometer.

Method Summary

A 50 mL water sample is adjusted to ionic strength and pH by the addition of ammonium chloride (NH_4Cl) and potassium fluoride (KF). The sample is poured into a UCT CLEAN-ELUTE™ column and the ETU is eluted from the column using 400 mL of methylene chloride. An excess of a free radical scavenger is added to the eluate. The methylene chloride eluant is concentrated to 5 mL after exchange into ethyl acetate. GC analysis with a nitrogen-phosphorous detector or mass spectrometer is used for quantitation.

Safety

- ETU is a suspected carcinogen. Prepare all standards in a fume hood

Sample Collection and Preservation

- Grab samples must be collected in 60 mL glass containers fitted with Teflon-lined crew caps
- Do not pre-rinse with sample before collection
- After collection shake the sample bottle for 1 minute
- ETU may degrade in water even during refrigeration. Mercuric chloride has been used as a preservative but due to its toxicity and harm to the environment is not recommended
- Store sample on ice or in refrigerator at 4°C and protected from light. Extract as soon as possible after collection

Interferences

Method interferences arise from contaminated glassware, solvents, reagents and other laboratory apparatus in which the sample may come in contact. All reagents and glassware must be shown to be free from interferences under analysis conditions.

- Glassware must be scrupulously clean
- Clean glassware by rinsing with the last solvent used followed by hot water and detergent. Rinse with reagent water, dry and heat in an oven at 400°C for one hour. Do not heat volumetric flasks

- Always use high purity reagents and solvents
- Interfering contamination may occur when a low concentration sample is analyzed after a high concentration sample. Complete rinsing of the syringe using ethyl acetate may be required

Procedure

1. Sample Preparation

- a) Pipette 50 mL of the water sample into a clean bottle
- b) Add 1.5 grams ammonium chloride (NH_4Cl)
- c) Add 25 grams potassium fluoride (KF)
- d) Seal bottle and shake until salts are completely dissolved

2. Sample Extraction

- a) Add 5 mL of 1000 g/mL of dithiothreitol (DTT, Cleland's Reagent) in ethyl acetate as a free radical scavenger to a 500 mL Kuderna-Danish K-D concentrator tube
- b) Support a CLEAN-ELUTE™ 200 mL cartridge using a clamp over a (K-D) tube
- c) Add the entire contents of the bottle from step 1) d) above
- d) Do not use vacuum but allow the cartridge to stand for 15 minutes

3. Sample Collection

- a) Add 400 mL of methylene chloride in 50 mL aliquots to the CLEAN-ELUTE™ column
- b) Collect the eluant in the K-D apparatus

4. Extract Concentration

The following steps must be conducted in a fume hood

- a) Add two boiling chips to the K-D apparatus and attach a macro Snyder column
- b) Attach a condenser to the Snyder column to collect solvent
- c) Place the K-D apparatus in a 65-70°C water bath so that the K-D tube is partially submerged in the water
- d) Once liquid volume had been reduced to 5 mL remove from the water bath
- e) Continue to reduce the liquid volume to < 1 mL in an analytical evaporator at 35-40°C under a stream of nitrogen
- f) Dilute sample to 5 mL with ethyl acetate rinsing the walls of the K-D apparatus
- g) Add 50 µL of internal standard and agitate
- h) Transfer to a GC vial
- i) Sample is ready for analysis

5. GC Analysis Conditions

Primary Conditions:

Column: 10 m long x 0.25 mm I.D. DB-Wax bonded fused

Carrier Gas: He @ 30 cm/sec linear velocity

Makeup Gas: He @ 30 mL/min flow

Detector Gases: Air @ 100 mL/min flow; H₂ @ 3 mL/min flow

Injector Temperature: 220°C

Detector Temperature: 230°C

Oven Temperature: 220°C isothermal

Sample: 2 µL splitless; nine second split delay

Detector: Nitrogen-phosphorus

Confirmation Conditions:

Column: 5 m long x 0.25 mm I.D. DB-1701 bonded fused

Carrier Gas: He @ 30 cm/sec linear velocity

Makeup Gas: He @ 30 mL/min flow

Detector Gases: Air @ 100 mL/min flow; H₂ @ 3 mL/min flow

Injector Temperature: 150°C

Detector Temperature: 270°C

Oven Temperature: 150°C isothermal

Sample: 2 µL splitless; nine second split delay

Detector: Nitrogen-phosphorus

Analyte	Primary Column RT (min)	Confirmation Column RT
ETU	3.5	4.5
THP internal standard	5.1	5.0
PTU surrogate standard	2.7	2.2

*The analyst should refer to EPA Method 509 "Determination of Ehtylene Thiourea (ETU) in Water Using Gas Chromatography with a Nitrogen-Phoshorus Detector", Revision 1.0 Issued 1992, By DJ Munch and RL Graves, US EPA, National Exposure Research Laboratory, Office of Research and Development, US Environmental Protection Agency, Cincinnati, Ohio 45268 and TM Engel and ST Champagne, Battelle, Columbus Division



Determination of Chlorinated Acids in Water Using Liquid-Solid Extraction and Gas Chromatography with an Electron Capture Detector

UCT Product Numbers:

ECDVB156 (500 mg, DVB, 6 mL cartridge)

ECUNIDVB500 (500 mg, DVB, 83 mL cartridge)

EPA Method 515.2 Revision 1.1

January 2011

Procedure

1. Hydrolysis and Separation of Interferences

Hydrolysis is required to convert esters to the parent acid prior to extraction

- a) Allow water samples from cold storage to come to room temperature
- b) Add 80 mg of sodium thiosulfate (or 50 mg sodium sulfite) if residual chlorine is present. Shake
- c) Perform the same steps for all LRB's, LFB's LFM's and QCS samples
- d) Measure 250 mL of each water sample with a graduate cylinder and pour into a 500 mL separatory funnel
- e) Add 250 µL of the surrogate **PDS** to each sample
- f) Dissolve 50 g sodium sulfate in each sample
- g) Add 4 mL of 6 N NaOH to each sample, seal, and shake
- h) If the sample pH is less than 12, adjust by adding more 6N NaOH
- i) Allow sample to sit at RT for one hour, occasionally shaking the contents
- j) Use 15 mL methylene chloride to rinse the sample bottle and the graduated cylinder
- k) Add the methylene chloride to the separatory funnel
- l) Extract the sample by vigorously shaking the funnel for 2 minutes
- m) Allow the organic layer to separate for 10 minutes
- n) Discard the methylene chloride phase
- o) Repeat the extraction procedure by addition of a second 15 mL volume of methylene chloride to the separatory funnel
- p) Discard the methylene chloride
- q) Perform a third extraction in the same manner as above discarding the methylene chloride

- r) Drain the aqueous contents of the separatory funnel into a 500 mL beaker
- s) Adjust the pH to 1.0 ± 0.1 with constant stirring by the dropwise addition of concentrated sulfuric acid

2. Sample Extraction

- a) Place **ECDVB156** or **ECUNIDVB500** cartridge(s) on the manifold
- b) With vacuum off add 20 mL of 10% by volume of methanol in methyl-*t*-butyl ether (MTBE) to the top of each cartridge
- c) Allow the solvent to remain for two minutes
- d) Turn on full vacuum and draw the solvent through
- e) Dry in room air for five minutes under full vacuum
- f) Adjust the vacuum to ≤ 5 in. Hg and add the following in series:
 - i. 20 mL methanol
 - ii. 20 mL reagent water
 - iii. sample

Note: Do not allow the cartridge to dry after steps f) a) or b) otherwise repeat

- g) After the sample has drawn through, apply full vacuum for 20 minutes to dry the cartridge(s)
- h) Place culture tubes in the manifold to collect eluates
- i) Add 2 mL aliquot of 10% methanol in MTBE and allow to sit for 1 minute before applying vacuum
- j) Add a second 2 mL aliquot of 10% methanol in MTBE and allow to sit for 1 minute before applying vacuum
- k) Rinse the beaker(s) with 4 mL of pure MTBE and add to the cartridge(s)
- l) Remove the culture tubes from the manifold and cover

3. Extract Preparation

- a) Prepare a acidified sodium sulfate drying tube(s)
- b) Rinse with 2 mL of MTBE
- c) Remove the entire extract with a 5 mL pipette
- d) Drain the lower aqueous layer back into the culture tube
- e) Add the organic layer to the drying tube

- f) Maintain liquid in the drying tube between this step and the following steps
- g) Collect the dried extract in a 15 mL graduated centrifuge tube or a 10 mL Kuderna-Danish tube
- h) Rinse the culture tube with an additional 1 mL of MTBE and repeat at step 3) c)
- i) Repeat step 3) e) and then add a 1 mL aliquot of MTBE to the drying tube before it empties. The final volume should be 6-9 mL
- j) The extract can now be esterified as in Step 4

4. Extract Esterification with Diazomethane

- a) Assemble a diazomethane generator in a hood as shown in diagram below
- b) Add 5 mL of ethyl ether to Tube 1
- c) Add 4 mL of Diazald* solution and 3 mL of 37% KOH solution to the reaction Tube Each charge of the generator should be sufficient to esterify four samples
 *Add 10 g Diazald (Sigma-Aldrich D28000 or other supplier) in 100 mL of a 50:50 by volume mixture of ethyl ether and Carbitol (2-(2-ethoxyethoxy) ethanol). Solution is stable for one month or longer when refrigerated in a Teflon-lined screw cap amber bottle
- d) Immediately place the exit tube into the collection tube containing the sample extract
- e) Turn on N₂ at 10 mL/min to bubble diazomethane through the extract. The appearance of a persistent yellow color indicates that esterification is complete. The first sample should require 30 seconds to one minute and each subsequent sample somewhat longer. The final sample may require two to three minutes
- f) Cover each collection tube and store at room temperature in a hood for 30 minutes
No significant fading of the yellow color should occur during this time
- g) Fortify each sample with 100 µL of the internal PDS solution and reduce the volume to 5.0 mL with an analytical concentrator or equivalent concentration technique using a gentle stream of dry N₂
- h) Cover tubes and store in a refrigerator if further analysis will not be performed immediately

Note: The excess diazomethane is volatilized from the extract during the concentration procedure

- i) Samples are ready to analyze by GC-ECD

Extract Esterification with Trimethylsilyldiazomethane (TMSD)

Alternative Procedure

Note: GC background is significantly increased when TMSD is used as the derivatizing reagent instead of diazomethane. While no method analyte is affected by this increased background, the recommended surrogate, 2, 4-dichloro-phenylacetic acid, is masked by an interfering peak. Any compound found suitable when TMSD is used is acceptable as a surrogate.

- a) Conduct the hydrolysis, clean-up, and extraction of the method analytes as described up to steps 2) g) above
- b) Elute the herbicides from the cartridge by passing 2 x 2 mL aliquots of MTBE through the cartridge into the collection tube
- c) Rinse the sample container with 4 mL of MTBE and pass it through the cartridge into the tube
- d) Pre-wet an **ECSS156** anhydrous sodium sulfate drying tube with 1 mL MTBE
- e) Transfer the MTBE extract from the collection tube discarding any water layer
- f) Before the extract passes completely through the drying tube, rinse with an additional 2 mL of MTBE
- g) Concentrate the dried extract to approximately 4 mL
- h) Add 1 mL of methanol to the extract to yield a 20% (v/v) methanol in MTBE solution
- i) Adjust the volume to 5 mL with MTBE. (TMSD produces the most efficient methylation of the herbicides in a 20% methanol, 80% MTBE solution)
- j) Add 50 μ L of the 2M TMSD solution to each 5 mL sample extract
- k) Place the tube containing the extract into a heating block at 50°C
- l) Heat the extract for one hour then allow the extract to cool to room temperature
- m) Add 100 μ L of 2M acetic acid in methanol to react with any excess TMSD
- n) Fortify the extract with 100 μ L of **ISS** to yield a concentration of 0.020 μ g/mL
- o) Sample(s) are ready for GC/ECD analysis

Stock Standard Solutions (SSS) 1.00-2.00 µg/µL

Stock standard solutions may be purchased as certified solutions or prepared in the lab

Prepare from pure standard materials using the following procedure:

- Accurately weighing approximately 0.0100-0.0200 g of pure material
- Dissolve the material in methanol and dilute to volume in a 10 mL volumetric flask
- If compound purity is certified at 96% or greater, the weight may be used without correction to calculate the concentration of the stock standard
- Transfer the stock standard solutions into 15 mL fluorocarbon-sealed screw cap amber vials. Store at 4°C or less
- Stock standard solutions should be replaced after two months or sooner

Primary Dilution Standards

- Prepare two sets of standards according to the sets labeled A and B in Table 1
- For each set, add approximately 25 mL of methanol to a 50 mL volumetric flask
- Add aliquots of each stock standard in the range of approximately 20-400 µL
- Dilute to volume with methanol
- Individual analyte concentrations will then be in the range of 0.4-8 µg/mL (for a 1.0 mg/mL stock). The minimum concentration would be appropriate for an analyte with strong ECD response, i.e., pentachlorophenol. The maximum concentration is for an analyte with weak response, i.e., 2, 4-DB

Internal Standard Solution (ISS)


- Prepare an **ISS** solution by accurately weighing approximately 0.050 g of pure 4,4'-dibromooctafluorobiphenyl (DBOB). Dissolve the DBOB in methanol and dilute to volume in a 10 mL volumetric flask. Store at room temperature in a fluorocarbon sealed screw cap bottle
- Prepare a **primary dilution standard (PDS)** approximately 1.00 µg/mL by the addition of 20 µL of the stock standard to 100 mL of methanol. Addition of 100 µL of the PDS solution to the final 5 mL of sample extract results in a final internal standard concentration of 0.020 µg/mL

Note: DBOB has been shown to be an effective internal standard for the method analytes

Surrogate Analyte Solution

- Prepare a surrogate analyte stock standard solution by accurately weighing approximately 0.050 g of pure 2, 4-dichlorophenylacetic acid (DCAA). Dissolve in methanol and bring to volume in a 10 mL volumetric flask. Store in a fluorocarbon-sealed screw cap bottle at room temperature
- Prepare a primary dilution standard of approximately 2.0 µg/mL by addition of 40 µL of the stock standard to 100 mL of methanol
- Add 250 µL of the surrogate analyte solution to 250 mL sample prior to extraction to provide a surrogate concentration of 2 µg/L

Diazomethane Generator



Summarized from, Method 515.2, "Determination Of Chlorinated Acids In Water Using Liquid-Solid Extraction And Gas Chromatography With An Electron Capture Detector", Revision 1.1, 1995, National Exposure Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH



Determination of Nitrosamines in Drinking Water by Solid Phase Extraction and Capillary Column Gas Chromatography with Large Volume Injection and Chemical Ionization Tandem Mass Spectrometry (MS/MS)*

UCT Part Numbers:

EU52112M6 (2000 mg activated coconut carbon, 6 mL)

ECSS156 (6 mL Drying Cartridge with 5 grams anhydrous sodium sulfate)

EPA Method 521

November 2012

Activated carbon is used for the determination of various nitrosamines in finished drinking water and untreated source waters using GC/MS/MS.

Nitroaromatics, Nitramines and Nitrate Ester Analytes

Analyte	Abbreviation	% Recovery n=3
N-Nitrosodimethylamine	NDMA	95
N-Nitrosomethyldiethylamine	NMEA	98
N-Nitrosodiethylamine	NDEA	95
N-Nitrosodi-n-propylamine	NDPA	90
N-Nitrosodi-n-butylamine	NDBA	94
N-Nitrosopyrrolidine	NPYR	76
N-Nitrosopiperidine	NPIP	81

*For complete details on Method 521, September 2004, the analyst is referred to: J.W.Munch & M.V.Bassett, "Determination of Nitrosamines in Drinking Water by Solid Phase Extraction and Capillary Column Gas Chromatography with Large Volume Injection and Chemical Ionization Tandem Mass Spectrometry (MS/MS), National Exposure Research Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268

Procedure

1. Cartridge Conditioning

- Add 3 mL of methylene chloride to the cartridge, then slowly draw all solvent through the cartridge
- Add 3 mL of methanol to the cartridge, turn on vacuum and draw through
- Add 3 mL of methanol again and draw through so that the methanol just covers the top of the cartridge frit

Do not let the cartridge go dry after this step otherwise repeat starting at step 1 b)

- Add 3 mL of reagent water and draw through
- Repeat water rinse, step d) **5 additional times**

Proper conditioning of the cartridge is essential for good precision and accuracy

2. Sample Extraction

- a) Adjust the vacuum setting so that the flow rate is 10 mL/minute
- b) After sample extraction, draw air through the cartridge for **10 minutes** at full vacuum
- c) After drying, proceed immediately to cartridge elution

3. Cartridge Elution

- a) Insert a clean collection tube in the manifold
- b) Fill the cartridge with methylene chloride
- c) Partially draw the methylene chloride through at low vacuum and then turn vacuum off
- d) Allow cartridge to soak for 1 minute
- e) Draw the remaining methylene chloride through in dropwise fashion
- f) Continue to add methylene chloride to the cartridge as it is being drawn through until a total of 12-13 mL have been added
- g) Concentrate the methylene chloride to about 0.9 mL in a water bath near room temperature. Do not concentrate less than 0.5 mL as loss of analyte may occur

Note: Small amounts of residual water from the sample container and SPE cartridge may form an immiscible layer with the extract. To eliminate the water a drying column packed with 5 grams of anhydrous sodium sulfate or use **ECSS15M6** for drying. Wet the cartridge with a small volume of methylene chloride before adding extract. Rinse the drying column with 3 mL of methylene chloride.

4. Sample Analysis

- a) Calibrate the MS in EI mode using FC-43
- b) Inject into a GC/MS/MS
- c) Identify the product ion spectrum to a reference spectrum in a user created data base

Mass Spectral Data

Analyte	Retention Time (min)	Precursor Ion (m/z)	Product/Quantitation Ion (m/z)
NDMA	8.43	75	43(56)
NMEA	11.76	89	61(61)
NDEA	14.80	103	75(75)
NPYR	22.34	101	55(55)
NDPA	22.40	131	89(89)
NPIP	24.25	115	69(69)
NDBA	30.09	159	57(103)
NDMA-d6 surrogate	8.34	81	46(59)
NMEA-d10 IS	14.63	113	81(81)
NDPA-d6 IS	22.07	145	97(97)

Injector Program

Temp (°C)	Rate (°C/min)	Time (min)
37	0	0.72
250	100	2.13
250	0	40

Injector Split Vent Program

Time (min)	Split Status
0	Open
0.70	Closed
2.00	Open

GC Oven Temperature Program

Temperature (°C)	Rate (°C/min)	Hold Time (min)
40	0	3.0
170	4.0	0
250	20.0	3.0

Limits and Lowest Concentration Minimum Reporting Levels

Analyte	DL (ng/L)	LCMRL (ng/L)
NDMA	0.28	1.6
NMEA	0.28	1.5
NDEA	0.26	2.1
NPYR	0.35	1.4
NDPA	0.32	1.2
NPIP	0.66	1.4
NDBA	0.36	1.4



Determination Of 1,4-Dioxane In Drinking Water By Solid Phase Extraction (SPE) AND Gas Chromatography/Mass Spectrometry (GC/MS) With Selected Ion Monitoring (SIM)*

UCT Part Numbers:

EU52112M6 (2 grams coconut carbon in a 6 mL cartridge)

EPA Method 522

November 2012

1,4-dioxane (CASRN 123-91-1) referred to as dioxane, is a highly water soluble non-biodegradable ether. This is a gas chromatography/mass spectrometry (GC/MS) method for the determination of this analyte in drinking water. Scan mode may be used if appropriate sensitivity is attained. Option 1 is designated for sample volumes of 500 mL. This product can be used with most automated systems.

Sample Collection & Preservation

- Prepare bottles before sample collection with sodium sulfite (CASRN 7757-83-7)
- Open the tap and allow the system to flush until the water temperature has stabilized (approximately three to five min)
- Fill sample bottles. Do not flush out the dechlorination reagent
- Cover the bottle and agitate by hand until the **sodium sulfite** has dissolved
- Add **sodium bisulfate** (CASRN 7681-38-1). Mix until dissolved

Sample Shipment and Storage

- Samples must be chilled during shipment and must not exceed 10° C during the first 48 hours after collection
- Sample temperature must be confirmed to be at or below 10° C when they are received at the laboratory
- Verification of sample dechlorination – Upon the receipt of samples at the laboratory, verify that field samples were dechlorinated at the time of collection. The absence of total chlorine can be verified with a N, N-DIETHYL-P-HENYLENEDIAMINE (DPD) CHLORINE TEST KIT – (Hach model CN-66; cat. # 2231-01 or equivalent)
- Holding time may be up to 28 days

Interferences:

- All glassware must be meticulously cleaned with detergent and tap water, rinse with tap water, followed by reagent water
- Non-volumetric glassware should be heated in a muffle furnace at 400°C for 2 hours
- Volumetric glassware should be solvent rinsed with DCM or purge and trap grade methanol after washing, and dried in a low temperature oven (<120°C) or air dried
- Interferences may be caused by contaminants in solvents, reagents (including reagent water), sample bottles and caps, and other sample processing hardware. All items must be routinely demonstrated to be free from interferences (less than 1/3 the MRL for the method analyte) under the conditions of the analysis.

Subtracting blank values from sample results is not permitted

- **Purge and trap grade methanol must be used** for all steps where methanol is required. Other grades of methanol may contain low molecular weight compounds which may prohibit accurate identification and quantitation of the analyte, SUR and IS
- Matrix interferences may be caused by contaminants that are co-extracted from the sample. The matrix interferences will vary considerably from source to source
- Analyte carry-over may occur when a relatively “clean” sample is analyzed immediately after a sample containing relatively high concentrations of compounds
- Syringes and splitless injection port liners must be cleaned carefully and replaced as needed. After analysis of a sample containing high concentrations of compounds, a LRB should be analyzed to ensure that accurate values are obtained for the next sample
- During automated GC/MS analyses, extracts with positive results that were analyzed immediately following a sample with high concentrations of the analyte, should be reanalyzed after analyzing an acceptable LRB. If the analyte is not detected in extracts analyzed immediately after a high concentration extract, no reanalysis is necessary
- Many silicone compounds may be leached from punctured autosampler vial septa, particularly when particles of the septa sit in the vial. This can occur after repeated injections from the same autosampler vial. If this method is performed in full scan mode, silicone compounds appear as regular chromatographic peaks with similar fragmentation patterns. They can unnecessarily complicate the total ion chromatogram
- High laboratory background levels of 1,4-dioxane have been reported to be associated with air contamination. Contact UCT technical support for possible solutions to air contamination issues.

Safety:

- The toxicity or carcinogenicity of each reagent used in this method has not been precisely defined
- 1,4-Dioxane is classified as a class B2 or probable human carcinogen
- Sodium bisulfate is used as a sample preservative to inhibit microbial growth and decay of 1,4-dioxane

Option 1—500 mL samples

1. Prepare Cartridge (vacuum manifold use)

- a) Place a **EU52112M6** cartridge on a single or multi-station vacuum manifold
- b) Add 3 mL of dichloromethane (methylene chloride) to the cartridge and draw through
- c) Add 3 mL of methanol and draw through completely
- d) Add 3 mL of DI water and draw through. **Do not let cartridge go dry**
- e) Repeat step d) 5 times without letting the cartridge go dry

2. Extract Sample

- a) Add surrogate 1,4-dioxane-d8 to the 500 mL sample water
- b) Add sample water to the cartridge. Adjust vacuum for a 10 mL/minute flow rate
- c) After 500 mL of sample has drawn through, dry cartridge by drawing air at full vacuum. Cartridge will feel room temperature when dry

3. Cartridge Elution

- a) Place a clean collection vial in the vacuum manifold
- b) Fill each cartridge with dichloromethane (DCM). Soak sorbent for 1 minute then slowly draw through dropwise
- c) Continue to add DCM to the cartridge until the total volume is 9 mL
- d) Remove collection tube and adjust volume to 10 mL with DCM
- e) Add IS. Adding 5 µl of 1000 µg/mL of IS solution to 10 mL DCM extract will yield 500 ng/mL
- f) Mix by vortexing
- g) Dry extract by adding 2 grams anhydrous sodium sulfate (Na₂SO₄)
- h) Mix well
- i) Transfer aliquot to autosampler vials for analysis

NOTE: This extract cannot be reliably concentrated by nitrogen evaporation because of the volatility of dioxane

Instrument & Conditions

Column: Varian CP-Select 624 CB (6% cyanopropyl phenyl, 94% PDMS, 30 m 0.25mm x 1.4 μ m column (or equivalent)

Injector: 200°C (Splitless mode)

Injector Volume: 1 μ L

Helium: Carrier gas 99.999% or better

Flow: 1 mL/minute

Oven: 30° C for 1 minute 90°C at 8° C/minute
200° C at 20° C/minute for 4 minutes

MS: Any type of mass spectrometer may be used (i.e., quadrupole, ion trap, time of flight, etc.), although the SIM option may not provide enhanced sensitivity, or be an available option on all instruments

SIM MODE

Segment 1: m/z 46**, 78, 80

Segment 2: m/z 58, 62, 64, 88**, 96**

Dwell time 100 μ s


Emission current 100 μ A

** quantitation ions

Retention Times and Quantitation Ions (QIs)

Compound	Retention Time (min)	SIM Ions (m/z)
1,4-dioxane	8.85	58, 88
1,4-dioxane- <i>d</i> ₈	8.77	62, 64, 96
THF- <i>d</i> ₈ (IS)	6.68	46 , 78, 80

Reconstructed total ion current chromatogram and mass chromatograms for THF-*d*8 (IS), 1,4-dioxane-*d*8 (SUR), and 1,4-dioxane at 0.5 mg/mL each (the standard is equivalent to an extract of a 10 mg/L aqueous sample). * Peak at 6.845 min is chloroform, a chemical present in DCM



Summarized from Jean W. Munch and Paul E. Grimmer, "Method 522 Determination Of 1,4-Dioxane In Drinking Water By Solid Phase Extraction (SPE) And Gas Chromatography/Mass Spectrometry (GC/MS) With Selected Ion Monitoring (SIM)", Office of Research and Development, National Exposure Laboratory, Cincinnati, OH Version 1.0, September 2008



Method 523: Determination of Triazine Pesticides and their Degradates in Drinking Water by Gas Chromatography/Mass Spectrometry (GC/MS) Version 1.0

UCT Part Numbers:

EC5232506 (250 mg GCB, 6 mL cartridge)

EPA Method 523

August 2012

Method Summary

This is a gas chromatography/mass spectrometry (GC/MS) method for the determination of triazine pesticides and their degradation products in finished drinking waters. Samples are pH adjusted, dechlorinated with ammonium acetate and protected from microbial degradation with 2-chloroacetamide during collection. Analytes are extracted from a **250 mL sample** using 250 mg carbon cartridges.

The following compounds can be determined using this method:

Analyte	CASRN
Atrazine	1912-24-9
Atrazine-desethyl	6190-65-4
Atrazine-desethyl-desisopropyl	3397-62-4
Atrazine-desisopropyl	1007-28-9
Cyanazine	21725-46-2
Propazine	139-40-2
Simazine	122-34-9
Terbuthylazine-desethyl	30125-63-4
Terbuthylazine	5915-41-3
Prometon	1610-18-0
Prometryn	7287-19-6
Ametryn	834-12-8

Procedure

1. Sample Preparation

- a) Allow samples to reach room temperature prior to extraction
- b) Add an aliquot of the Surrogate Primary Dilution Standards (PDS) to each sample
- c) Fortify Laboratory Fortified Blanks, Laboratory Fortified Sample Matrices, or LFSM Duplicates, with an appropriate volume of analyte PDS and the atrazine-desethyl-desisopropyl stock standard
- d) Cap and invert each sample several times to mix
- e) Proceed with sample extraction using SPE carbon cartridges

2. Cartridge Cleaning & Conditioning

- a) Set up extraction cartridges on the SPE vacuum manifold
- b) Using low vacuum (approximately 1 to 2 inches Hg), rinse each cartridge with two 6 mL aliquots of DCM drawing completely through
- c) Rinse each cartridge with a 6 mL aliquot of MeOH
- d) Draw MeOH to the top of the cartridge frit

Note: Do not let the cartridge dry after addition of MeOH

- e) Add a 6 mL aliquot of reagent water (RW) to the cartridge
- f) Draw RW to the top of the cartridge frit

3. Sample Extraction

- a) Add an additional 4 mL of RW to each cartridge
- b) Attach sample transfer lines to the cartridges. The additional volume prevents the SPE cartridge bed from going dry while the dead volume in the transfer lines is being filled
- c) Extract 250 mL of sample at a cartridge flow rate of 10 mL/minute
- d) Dry the cartridges under high vacuum for 10 seconds
- e) Release vacuum, then add a 0.25 mL aliquot of MeOH to each cartridge
- f) Draw the MeOH to waste, then dry cartridge under full vacuum for 10 minutes

4. Sample Elution

- a) Place 15 mL conical tubes into the manifold for collection
- b) Add 2 mL of EtOAc to the cartridge and elute dropwise
- c) Add 2 x 6 mL aliquots of 9:1 DCM/MeOH to cartridge
- d) Allow the cartridge beds to briefly soak in solvent, then draw the solvent through the cartridges
- e) Dry the eluate by passing it through approximately 3 grams of anhydrous Na_2SO_4 collecting it in a 40 mL centrifuge tube. Pre-rinse the Na_2SO_4 with a 1 mL aliquot of 3:1 DCM/EtOAc
- f) Rinse with 1 mL aliquot of 3:1 DCM/EtOAc collecting it in the centrifuge tube
- g) The dried extracts may be stored overnight in the 40 mL tubes at $-10\text{ }^{\circ}\text{C}$
- h) Warm the 40 mL tubes to $35\text{ }^{\circ}\text{C}$ in a water bath under a stream of N_2 and evaporate solvent to less than 1 mL but no less than 0.5 mL
- i) Transfer the concentrated eluate to 1 mL volumetric tubes
- j) Rinse the conical tube with a small volume of EtOAc, and transfer the rinse to the volumetric
- k) Add IS solution and adjust to volume
- l) Transfer the extracts to autosampler vials for analysis or store in a freezer $\leq -10\text{ }^{\circ}\text{C}$

Complete details at Office of Water (MLK 140) EPA Document No. 815-R-11-002 February 2011 <http://www.epa.gov/safewater/>



EPA Method 525.2 Revision 2.0

Determination of Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column Gas Chromatography/Mass Spectrometry

UCTUCT Part Numbers:

ECUNI525 (1500 mg uncapped C18, 83 mL cartridge)

METHOD 525.2

May 2009

Method 525.2 Analytes

The analyte list for this method is comprised of over 120 compounds representative of several classes of pesticides, polynuclear aromatic hydrocarbons, PCBs, phthalates and adipates and other drinking water pollutants. Recovery ranges from 70-130%. Refer to the published method for compound specific MDL's.

The validation data presented herein were determined on independent lots of UCT ENVIRO-CLEAN[®] Universal Cartridges. MDLs were not determined on all analytes as part of this validation. In addition to the listed method analytes, recovery data for an extended list of analytes is also included.


Table of Compounds
Tested using the UCT ENVIRO-CLEAN[®] Universal Cartridge 525

Analyte	Average 3 Replicates % Recovery	Std Dev
Acenaphthene	100	0.0
2,4-dinitrotoluene	83	NA
2,6-dinitrotoluene	78	NA
4,4''-DDE	91	4.2
4,4''-DDT	94	3.5
4,4'-DDD	94	4.5
Acenaphthylene	96.0	0.012
Acenaphthene	99.1	0.013
Acetochlor	115	0.01
Alachlor	99	0.007
Aldrin	77	4.4
Ametryn	95	4.6
Anthracene	80	0.0
Atraton	84	17.3
Atrazine	111	0.011
Benzo(a) anthracene	75.4	0.049
Benzo(a)pyrene	105	9.9
Benzo(b)fluoranthene	184.9	0.022
Benzo(k)fluoranthene	95.7	0.029
Benzo[g,h,i]perylene	83.1	0.05
BHC, alpha	108	6.9
BHC, beta	97	3.1
BHC,delta	109	7.9
BHC,gamma	102	11.9
bis- (2-ethylhexyl) adipate	95.1	0.033
bis 2 ethylhexyl phthalate	104	0.029
Bromacil	126	0.012
Butachlor	113	0.005
Butylate	103	4.6

Butylbenzylphthalate	97.1	0.02
Caffeine	90.0	0
Captan	86.9	0.273
Carboxin	103	12.9
Chlordane, alpha	97	4.6
Chlordane, gamma	94	2.5
Chlordane, trans nonachlor	115	11.0
Chlorneb	113	11.0
Chlorobenzilate	118	10.0
Chlorpropham	130	0
Chlorpyrifos (Dursban)	107	5.0
Chlorothalonil	117	12.1
Chrysene	100	0.012
Cyanazine (Bladex)	126	0.008
Cycloate	111	15.0
Dacthal (DCPA) methyl ester	118	13.1
Diazinon	135	0.031
Dibenzo[a,h]anthracene	77.4	0.051
Dichlorvos (DDVP)	127	9.5
Dieldrin	96	6.8
Diethylphthalate	99.1	0.071
Dimethoate	106	0.008
Dimethylphthlate	78.6	0.022
Di-n-butylphthalate	113	0.12
Diphenamid	119	0.008
Disulfoton	92.1	0.01
Disulfoton Sulfone	108	12.5
Endosulfan I	116	11.1
Endosulfan sulfate	114	6.8
Endrin	88	0.0
Endrin Aldehyde	97	3.6
Endrin Ketone	90	3.8
EPTC	102	0.005
Ethion	112	0.005
Ethoprophos	109	5.8
Etridiazole (terrazole)	97	1.2
Fenarimol	70	0.0
Fluoranthene	100	0.018
Fluorene	99.7	0.012
Heptachlor	79	8.2
Heptachlor Epoxide Iso A	116	16.3
Hexachlorobenzene	94	17.4
Hexachlorocyclopentadiene	82	8.4
Hexazinone (Velpar)	105	8.1
Indeno[1,2,3-cd]pyrene	77.4	0.16
Isophorone	91	NA
Lindane	127	4.8
methoxychlor	123	7.6
Methyl Paraxon (Parathion)	115	5.0
Metolachlor	111	0.004
Metribuzin	109	0.005
Mevinphos (phosdrin)	117	
	12.1	
MGK 264	121	5.8
Molinate	114	0.013
Naphthalene	90.3	0.013
Napropamide (Devrinol)	115	2.3
Nonachlor, trans	116	11.1
Norflurazon	133	6.1
PCNB (carbaryl)	91.4	0.021
Pebulate	101	1.7
Pentachlorophenol	80	0.017
Permethrin, cis	124	2.1
Permethrin, trans	123	3.1
Perylene-d12	119	0.0
Phenanthrene	96.9	0.014
Phenanthrene-d10	99	6.6
Prometon	78.6	0.008
Prometryn	110	0.012
Pronamide (propyzamide)	101	1.2
Propachlor	113	15.0
Propazine	105	4.6

Pyrene	94.6	0.022
Simazine	91.4	0.005
Simetryn	93	4.6
Stirofos (tetrachlorvinphos)	126	6.9
Thiobencarb	112	0.008
Tebuthiuron	85	33.5
Terbacil	120	3.5
Terbutryn	103	2.3
Triademefon	98	6.9
Tricyclazole	107	5.0
Trifluralin	82	9.7
Trifluran	83	9.2
Trithion (carbofenothion)	101	0.004
Terbufos	95	7.0
Vernolate	107	1.2

PCB Congeners	Average	Std Dev
2-chlorobiphenyl	93	2.3
2,3-Dichlorobiphenyl	113	15.0
2,4,5-trichlorobiphenyl	97	3.1
2,2,4,4-tetrachlorobiphenyl	98	5.3
2,2,3,4,6-pentachlorobiphenyl	104	2.0
2,2,4,4,5,6-hexachlorobiphenyl	103	3.1
2,2,3,3,4,4,6-heptachlorobiphenyl	85	1.2
Octachlorobiphenyl (BZ#200)	79	1.2



Procedure

1. Wash the extraction apparatus and cartridge

- Add 10 mL of a 1: 1 mixture of ethyl acetate: methylene chloride (EtOAc: MeCl₂) to the reservoir.
- Draw a small amount through the cartridge with vacuum
- Turn off the vacuum and allow the cartridge to soak for about one minute

- d) Draw the remaining solvent through the cartridge to waste
- e) Allow the cartridge to dry for 3 minutes under full vacuum

2. Condition Cartridge

- a) Add 10 mL of methanol
- b) Draw a small amount through the cartridge
- c) Let soak for about one minute
- d) Draw most of the remaining methanol through the cartridge, leaving 3 to 5 mm of methanol on the surface of the cartridge frit
- e) Immediately add 20 mL of reagent water to the cartridge and draw most of the water through leaving 3 to 5 mm on the top of the cartridge frit

Note: Do not let the cartridge dry out after the addition of water

- f) Add 5 ml of methanol to the water sample and mix well
- g) Add the water sample to the cartridge and under vacuum, filter at a rate of approximately 50 mL per minute
- h) Drain as much water from sample bottle as possible
- i) Dry the cartridge under vacuum for 10 minutes

Note: Exceeding a 10-minute dry time could result in low recoveries. For faster drying, remove the cartridge and tapping the excess moisture from the bottom of the cartridge before continuing vacuum drying

3. Elution

- a) Insert a suitable sample tube for eluate collection
- b) Add 10 mL of EtOAc to the sample bottle
- c) Rinse the sample bottle thoroughly
- d) Transfer the solvent to the cartridge with a disposable pipette, rinsing sides of filtration reservoir
- e) Draw half of solvent through cartridge then release the vacuum. Allow the remaining solvent to soak the cartridge for about one minute
- f) Draw remainder through under vacuum
- g) Repeat the solvent rinse of the sample bottle and apparatus using 10 mL of 1:1 EtOAc:MeCl₂
- h) Using a disposable pipette, rinse down the sides of the cartridge and bottle holder with another 10 mL aliquot of 1: 1 EtOAc:MeCl₂

- i) Add the rinse to the cartridge, then draw through

4. Dry the combined eluant

- a) Use granular anhydrous sodium sulfate
- b) Rinse the collection tube and sodium sulfate with two x 3 mL portions of MeCl₂ and place combined solvent in a concentrator tube
- c) Draw through using vacuum
- d) Concentrate the extract to 1 mL under gentle stream of nitrogen (may be warmed gently) being careful not to spatter the contents

Note: Do not concentrate to <0.5 mL or loss of analytes could occur. Rapid extract concentration could result in loss of low molecular weight analytes

5. Analyze by GC/MS

Revision 2.0, 1995. Method authors: Eichelberger, J. W., Behymer, T. D., Budde, W. L., Munch, J., National Exposure Research Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268

This summary highlights major steps in the 525.2 method. Complete details about the preparation and composition of reagent solutions can be found in method and should be referenced by anyone needing complete details. It is available as a part of Supplement 11 from National Technical Information Service (NTIS), Springfield, VA 22161; publication PB 92 207703. (800) 553-6847 or at www.epa.gov/safewater/methods/methods.html



Semivolatile Organic Compounds In Drinking Water By Solid-Phase Extraction And Capillary Column (GC/MS) EPA Method 525.3 Version 1.0

UCT Part Number:

ECUNI525 (1500 mg unendcapped C18, 83 mL)

October 2012

Method Summary

A 1-liter water sample is fortified with surrogate analytes then extracted using a solid phase extraction (SPE) cartridge (ECUNI525). Analytes are eluted from the solid phase with a small amount of organic solvents. The extract is dried using anhydrous sodium sulfate and concentrated to approximately 0.7 mL using N₂. IS are added and the volume adjusted to 1 mL with ethyl acetate. A splitless injection is made into a GC equipped with a capillary column, interfaced to an MS with either scan, SIM or SIS detection used for analysis. **The GC/MS may be calibrated using standards prepared in solvent or using matrix-matched standards.**

Internal standards are added after the extract concentration step. If the analyte pentachlorophenol is being measured, use IS ¹³C-pentachlorophenol at 1000 µg/mL.

Sample Preservation

Preservation reagents, listed in the table below, are added to each sample bottle as dry solids prior to shipment to the field (or prior to sample collection).

Compound	Amount	Purpose
L-ascorbic acid	0.10 g/L	Dechlorination
Ethylenediaminetetraacetic acid, trisodium salt (EDTA)	0.35 g/L	Inhibit metal-catalyzed hydrolysis of targets
Potassium dihydrogen citrate	9.4 g/L	pH 3.8 buffer, microbial inhibitor

Procedure

1. Cartridge Clean-up

- Assemble an extraction system
- Rinse bottle holders and cartridges with 5 mL 1:1 EtOAc:DCM (ethyl acetate:dichloromethane)
- Draw half the volume through the cartridge and then soak for 1 min
- Draw remaining solvent through the cartridge

- e) Maintain full vacuum for 2 min to dry cartridge

2. Cartridge Conditioning

- a) Add 10 ml of methanol to each cartridge
- b) Soak for 1 minute
- c) Draw through leaving a thin layer of methanol on the cartridge frit

Note: Do not let the cartridge go dry from this point until elution, otherwise recondition

- d) Add 10 mL of reagent water to each cartridge
- e) Draw through leaving a thin layer of water on the cartridge frit

3. Sample Extraction

- a) All field and QC samples, including LRBs and LFBs, must contain preservatives
- b) Ensure that sample pH is ≤ 4 (use a pH meter for reagent water)
- c) Place sample bottle(s) in holder
- d) Adjust vacuum to fast drip flow rate. A flow of 10 mL/min is optimum
- e) After entire sample has extracted, rinse bottle with 10 mL reagent water
- f) Add rinse to cartridge
- g) Rinse cartridge using 10 mL reagent water to remove sample preservatives
- h) Dry cartridge for 10 min under full vacuum or nitrogen positive pressure

4. Cartridge Elution

- a) Insert 40 mL glass vial in manifold
- b) Rinse bottle, holder, and cartridge with 5 mL EtOAc
- c) Pour rinsate into cartridge
- d) Draw $\frac{1}{2}$ volume through cartridge, soak 1 min then draw through completely
- e) Repeat using 5 mL DCM
- f) Repeat using 5 mL EtOAc
- g) Repeat using 5 mL DCM rinse

5. Extract Drying

- a) Pre-rinse a drying tube containing 10 g of anhydrous sodium sulfate with DCM

- b) Quantitatively transfer the eluant through the sodium sulfate tube and collect
- c) Rinse the collection tube 5 mL of DCM
- d) Pass the DCM through the sodium sulfate and collect

6. Extract Concentration

- a) Concentrate extract to about 0.7 mL (not < 0.5 mL) under a gentle stream of N₂ in a water bath at 40 °C
- b) Transfer to a 1 mL volumetric flask, add IS and bring to volume using EtOAc

7. Analyze by GC/MS

Internal Standards	CASRN	Solvent	PDS conc.
acenaphthene- <i>d</i> 10 (IS 1)	15067-26-2	acetone	500 µg/mL
phenanthrene- <i>d</i> 10 (IS 2)	1517-22-2	acetone	500 µg/mL
chrysene- <i>d</i> 12 (IS 3)	1719-03-5	acetone	500 µg/mL
¹³ C-pentachlorophenol (IS 4)	85380-74-1	methanol	1000 µg/mL

Surrogates	CASRN	Solvent	PDS conc.
1,3-dimethyl-2-nitrobenzene (SUR 1)	81-20-9	acetone	500 µg/mL
triphenyl phosphate (SUR 2)	115-86-6	acetone	500 µg/mL
benzo[<i>a</i>]pyrene- <i>d</i> 12 (SUR 3)	63466-71-7	acetone	500 µg/mL

Instrument Conditions for GC Analysis	
Agilent 5975C MSD with 6890N GC	Restek RXI-5sil-MS 30m x 0.25 mm x 0.25 µm column
4-mm i.d. splitless gooseneck injection port liner	UCT#GCLGN4MM
Injection Port	250 °C
Injection Vol	1 µL with 1 min split delay
GC Oven Temp	Initial: 55 °C, hold 1 min Ramp 10 °C/min to 200 °C Ramp 7 °C/min to final T 320 °C Hold 0.36 min

Precision and Accuracy Data Obtained for Method 525.3 Analytes Fortified in Reagent Water at Three Concentrations and Extracted Using UCT 525 Universal Cartridges; N=4; Full Scan GC/MS Analyses^a

Analytes	Fortified Conc. 0.25 µg/L ^b		Fortified Conc. 2.0 µg/L ^c		Fortified Conc. 5.0 µg/L ^d	
	Mean % Recovery	RSD	Mean % Recovery	RSD	Mean % Recovery	RSD
acenaphthylene	101	2.0	93.6	0.51	99.8	0.88
acetochlor	99.0	3.9	93.6	2.1	104	2.1
alachlor	100	7.3	89.8	0.72	92.8	1.0
aldrin	77.0	5.0	78.4	2.9	85.0	3.4
ametryn	105	4.8	93.1	1.3	95.8	1.1
anthracene	106	3.8	92.3	1.0	104	0.71
atraton	112	2.9	90.3	4.1	96.8	2.2
atrazine	111	3.5	96.1	3.2	97.3	1.5
benzo[a]anthracene	112	5.1	99.1	3.8	112	3.6
benzo[a]pyrene	109	5.5	103	1.7	111	1.2
benzo[b]fluoranthene	119	5.0	102	1.2	114	2.4
benzo[g,h,i]perylene	112	2.9	102	4.3	113	2.8
benzo[k]fluoranthene	105	1.9	103	2.4	113	3.1
BHT	ND ^e		ND		ND	
bromacil	102	9.3	98.9	0.86	103	2.1
butachlor	107	3.6	86.3	1.1	99.4	1.4
butylate	85.0	7.1	83.0	2.2	84.0	3.0
butylbenzylphthalate	122	1.9	95.9	3.7	114	3.6
chlordane, cis	98.0	5.3	102	2.5	101	1.3
chlordane, trans	103	1.9	103	2.3	96.6	0.71
chlorfenvinphos	113	1.8	110	2.5	111	3.9
chlorobenzilate	82.0	9.3	99.8	5.1	94.1	1.3
chloroneb	93.0	2.2	92.0	2.9	100	1.4
chlorothalonil	116	2.8	106	3.8	105	1.5
chlorpropham	109	3.5	93.1	2.5	98.6	1.1
chlorpyrifos	102	5.1	93.4	3.3	97.2	2.5
chrysene	117	1.7	97.1	1.7	114	2.1
cyanazine	99.0	3.9	88.1	4.9	106	2.4
cycloate	102	3.9	87.4	1.0	88.8	1.2
dacthal (DCPA)	105	3.6	102	3.5	101	1.7
DDD, 4,4'-	107	3.6	85.8	0.75	105	1.4
DDE, 4,4'-	99.0	3.9	82.3	1.3	101	1.0
DDT, 4,4'-	116	2.8	87.6	2.8	112	0.83
DEET	103	1.9	98.3	3.2	104	2.0
di(2-ethylhexyl)adipate	112	4.1	96.6	3.1	111	1.8
di(2-ethylhexyl)phthalate	137	3.7	97.6	1.3	110	2.4
dibenzo[a,h]anthracene	110	3.6	95.4	2.5	109	1.5
dibutyl phthalate	115	3.3	101	1.5	114	2.6
dichlorvos	104	3.1	91.6	1.7	88.8	2.7
dieldrin	103	1.9	87.4	0.55	98.1	0.39

diethylphthalate	111	1.8	111	2.3	114	1.4
dimethipin	24.0	14	29.5	6.5	24.9	2.5
dimethylphthalate	110	3.6	113	0.25	113	0.78
DIMP	112	5.8	90.0	1.9	93.7	3.0
dinitrotoluene, 2,4-	126	1.8	105	2.6	113	2.5
dinitrotoluene, 2,6-	121	1.7	106	0.71	111	0.67
diphenamid	106	2.2	95.1	0.90	97.8	0.69
disulfoton	79.0	2.5	91.5	8.9	85.3	1.5
endosulfan I	95.0	5.3	88.4	1.3	101	1.0
endosulfan II	103	1.9	89.6	3.5	103	1.0
endosulfan sulfate	112	7.1	96.5	2.4	106	0.75
endrin	89.0	5.7	82.9	3.4	91.3	4.0
EPTC	89.0	2.2	88.0	0.80	85.8	0.60
ethion	106	2.2	100	2.7	108	3.1
ethoprop	110	2.1	91.3	1.6	96.2	1.4
ethyl parathion	117	4.3	97.6	2.3	105	3.5
etridiazole	118	3.4	90.6	2.2	101	1.5
fenarimol	110	4.7	87.1	2.0	91.7	3.4
fluorene	106	3.8	97.5	2.1	101	1.1
fluridone	92.0	5.0	103	4.9	98.6	2.9
HCCPD	92.0	3.5	65.6	1.7	68.0	5.8
HCH, α	101	3.8	92.5	1.3	95.2	0.47
HCH, β	101	3.8	94.0	4.0	102	1.7
HCH, δ	97.0	6.2	96.4	1.3	101	0.49
HCH, γ (lindane)	90.0	4.4	95.6	2.1	97.9	1.8
heptachlor	96.0	3.4	83.1	2.0	86.2	1.2
heptachlor epoxide	104	3.1	86.9	2.0	95.9	1.6
hexachlorobenzene	94.0	5.5	78.4	3.8	93.0	1.9
hexazinone	107	1.9	84.6	1.7	94.6	2.6
indeno[1,2,3-c,d]pyrene	113	4.5	95.0	2.1	112	2.4
isophorone	108	3.0	108	3.2	102	1.1
methoxychlor	122	1.9	89.9	1.5	109	0.72
methyl parathion	129	3.0	103	2.1	112	2.5
metolachlor	109	1.8	93.1	1.1	97.8	0.42
metribuzin	116	2.8	97.3	0.30	106	2.5
mevinphos	115	3.3	96.1	3.1	97.0	1.5
MGK 264(a)	94.0	2.5	75.5	1.9	88.3	4.1
MGK 264(b)	94.0	2.5	82.8	0.35	92.3	0.80
molinate	88.0	3.7	89.4	1.2	88.9	2.1
napropamide	105	3.6	89.9	2.4	99.2	1.5
nitrofen	129	3.0	106	2.7	113	3.9
nonachlor, trans	119	3.2	103	2.3	96.2	1.0
norflurazon	106	2.2	91.9	1.6	102	1.6
oxyfluorfen	129	1.6	93.9	2.9	111	3.5
pebulate	85.0	8.0	84.5	1.7	84.7	2.4
pentachlorophenol	104	4.1	100	1.3	96.0	3.6
permethrin, cis	110	3.6	107	1.3	107	2.2
permethrin, trans	115	3.3	96.1	3.1	97.0	1.5

phenanthrene	94.0	2.5	75.5	1.9	88.3	4.1
phorate	94.0	2.5	82.8	0.35	92.3	0.80
phosphamidon	88.0	3.7	89.4	1.2	88.9	2.1
profenofos	105	3.6	89.9	2.4	99.2	1.5
prometon	129	3.0	106	2.7	113	3.9
prometryn	119	3.2	103	2.3	96.2	1.0
pronamide	106	2.2	91.9	1.6	102	1.6
propachlor	129	1.6	93.9	2.9	111	3.5
propazine	85.0	8.0	84.5	1.7	84.7	2.4
pyrene	104	4.1	100	1.3	96.0	3.6
simazine	110	3.6	107	1.3	107	2.2
simetryn	115	3.3	96.1	3.1	97.0	1.5
tebuconazole	94.0	2.5	75.5	1.9	88.3	4.1
tebuthiuron	94.0	2.5	82.8	0.35	92.3	0.80
terbacil	88.0	3.7	89.4	1.2	88.9	2.1
terbutryn	105	3.6	89.9	2.4	99.2	1.5
tetrachlorvinphos	129	3.0	106	2.7	113	3.9
triadimefon	119	3.2	103	2.3	96.2	1.0
tribufos+merphos	106	2.2	91.9	1.6	102	1.6
trifluralin	129	1.6	93.9	2.9	111	3.5
vernolate	85.0	8.0	84.5	1.7	84.7	2.4
vinclozolin	104	4.1	100	1.3	96.0	3.6
PCB Congeners by UPAC#						
2-chlorobiphenyl (1)	75.0	2.7	81.0	2.6	85.0	1.1
4-chlorobiphenyl (3)	85.0	2.4	84.4	2.8	88.9	1.5
2,4'-dichlorobiphenyl (8)	85.0	2.4	82.5	2.2	87.1	0.51
2,2',5-trichlorobiphenyl (18)	104	3.1	89.3	2.9	95.3	4.0
2,4,4'-trichlorobiphenyl (28)	81.0	2.5	88.6	3.3	92.4	0.48
2,2',3,5'- tetrachlorobiphenyl (44)	85.0	9.7	91.1	3.1	93.4	1.6
2,2',5,5'- tetrachlorobiphenyl (52)	84.0	3.9	90.3	4.4	96.4	1.1
2,3',4',5'- tetrachlorobiphenyl (70)	84.0	3.9	92.4	3.4	97.0	0.94
2,3,3',4',6- pentachlorobiphenyl (110)	81.0	2.5	94.3	3.2	97.1	0.90
2,3',4,4',5- pentachlorobiphenyl (118)	90.0	2.6	94.5	3.7	98.8	1.1
2,2',3,4,4',5'- hexachlorobiphenyl (138)	86.0	4.7	99.0	3.2	105	1.3
2,2',3,4',5',6- hexachlorobiphenyl (149)	87.0	5.8	96.5	4.3	101	1.2
2,2',4,4',5,5'- hexachlorobiphenyl (153)	79.0	2.5	96.9	2.8	101	1.1
2,2',3,4,4',5,5'- heptachlorobiphenyl (180)	101	2.0	93.1	2.9	91.6	2.2

Surrogate Analytes						
1,3-dimethyl-2-nitrobenzene	92.9	4.0	98.9	3.1	88.8	4.0
benzo[a]pyrene- <i>d</i> 12	112	2.1	101	2.6	101	4.9
triphenyl phosphate	107	2.5	97.9	4.1	104	3.2

a. Data obtained on the instrumentation described in Sect. 13.1.1.4

b. Exceptions to the stated concentration are as follows: Surrogate concentrations are 5.0 µg/L, pentachlorophenol is 1.0 µg/L, c-permethrin is 0.13 µg/L, t-permethrin is 0.38 µg/L, MGK 264 (a) is 0.085 µg/L and MGK 264 (b) is 0.17 µg/L.

c. Exceptions to the stated concentration are as follows: Surrogate concentrations are 5.0 µg/L, pentachlorophenol is 8.0 µg/L, c-permethrin is 1.0 µg/L, t-permethrin is 3.0 µg/L, MGK 264 (a) is 0.67 µg/L and MGK 264 (b) is 1.3 µg/L.

d. Exceptions to the stated concentration are as follows: Surrogate concentrations are 5.0 µg/L, pentachlorophenol is 20.0 µg/L, c-permethrin is 2.5 µg/L, and t-permethrin is 7.5 µg/L, MGK 264 (a) is 1.7 µg/L and MGK 264 (b) is 3.3 µg/L.

e. ND = Not determined. Analyte could not be determined because of the low fortified concentration relative to the LRB.

Precision and Accuracy Data Obtained for Method Analytes Fortified into Finished Drinking Waters from Ground and Surface Water Sources, and Extracted Using UCT 525 Universal Cartridges; N=4; Full Scan GC/MS Analyses^a

Analytes	Fortified Conc. (µg/L)	Ground Water ^b		Surface Water ^c	
		Mean % Recovery ^d	RSD	Mean % Recovery ^d	RSD
acenaphthylene	2.0	99.3	5.7	95.1	1.8
acetochlor	2.0	97.4	5.5	106	4.9
alachlor	2.0	93.5	3.5	95.0	3.3
aldrin	2.0	94.9	4.0	81.1	0.59
ametryn	2.0	98.9	5.7	93.4	6.2
anthracene	2.0	102	4.7	101	1.6
atraton	2.0	92.6	5.9	87.3	5.3
atrazine	2.0	100	4.6	95.6	1.5
benzo[a]anthracene	2.0	104	1.5	102	2.4
benzo[a]pyrene	2.0	106	3.2	100	2.7
benzo[b]fluoranthene	2.0	104	3.3	100	3.9
benzo[g,h,i]perylene	2.0	101	3.9	101	5.3
benzo[k]fluoranthene	2.0	103	3.0	101	2.6
BHT	2.0	95.5	1.2	114	1.9
bromacil	2.0	98.9	5.5	104	4.6
butachlor	2.0	96.4	2.9	95.6	3.3
butylate	2.0	87.6	3.2	83.6	1.8
butylbenzylphthalate	2.0	107	4.1	107	2.9
chlordane, cis-	2.0	96.1	6.2	98.0	5.9
chlordane, trans	2.0	94.5	5.9	98.3	7.3
chlorfenvinphos	2.0	93.1	3.9	111	3.4
chlorobenzilate	2.0	101	5.6	97.1	4.3
chloroneb	2.0	104	2.0	108	4.0
chlorothalonil	2.0	108	1.5	110	2.2
chlorpropham	2.0	95.9	5.6	98.1	2.4

chlorpyrifos	2.0	97.8	5.3	103	4.1
chrysene	2.0	106	2.2	100	3.2
cyanazine	2.0	97.4	5.2	91.0	11
cycloate	2.0	92.3	5.6	95.1	1.6
dacthal (DCPA)	2.0	92.1	7.3	107	4.3
DDD, 4,4'-	2.0	93.1	3.7	91.0	2.0
DDE, 4,4'-	2.0	90.0	3.6	85.6	2.7
DDT, 4,4'-	2.0	91.4	3.7	90.4	3.5
DEET	2.0	101	1.8	109	1.1
di(2-ethylhexyl)adipate	2.0	102	5.8	106	3.4
di(2-ethylhexyl)phthalate	2.0	107	2.7	104	2.2
dibenzo[a,h]anthracene	2.0	106	6.2	102	2.2
dibutyl phthalate	2.0	110	3.7	107	1.5
dichlorvos	2.0	90.4	7.3	88.8	2.3
dieldrin	2.0	96.0	5.4	96.3	2.6
diethylphthalate	2.0	110	2.5	107	1.1
dimethipin	2.0	29.4	8.7	38.1	4.2
dimethylphthalate	2.0	111	2.4	111	1.0
DIMP	2.0	89.6	5.9	99.6	5.9
dinitrotoluene, 2,4-	2.0	95.6	7.0	105	4.1
dinitrotoluene, 2,6-	2.0	95.6	6.2	101	5.0
diphenamid	2.0	98.5	5.6	101	2.3
disulfoton	2.0	79.6	5.9	103	6.3
endosulfan I	2.0	96.9	7.0	93.0	4.5
endosulfan II	2.0	98.5	5.7	94.3	4.1
endosulfan sulfate	2.0	102	2.8	100	2.1
endrin	2.0	97.9	2.6	88.1	3.4
EPTC	2.0	89.0	7.2	85.8	3.0
ethion	2.0	95.8	6.9	98.9	4.3
ethoprop	2.0	94.9	5.7	103	3.9
ethyl parathion	2.0	97.6	5.9	103	3.0
etridiazole	2.0	104	2.4	104	3.1
fenarimol	2.0	93.5	4.4	86.6	3.3
fluorene	2.0	101	5.1	104	1.8
fluridone	2.0	111	3.9	97.4	6.0
HCCPD	2.0	64.6	6.7	69.1	2.8
HCH, α	2.0	93.5	5.7	91.8	2.4
HCH, β	2.0	97.9	4.2	97.8	1.8
HCH, δ	2.0	102	6.4	94.3	3.7
HCH, γ (lindane)	2.0	100	6.8	90.0	2.1
heptachlor	2.0	91.0	3.0	86.1	1.5
heptachlor epoxide	2.0	98.3	6.2	91.6	1.6
hexachlorobenzene	2.0	90.4	4.0	89.0	2.7
hexazinone	2.0	93.4	3.1	97.0	8.1
indeno[1,2,3-c,d]pyrene	2.0	107	6.6	103	4.6
isophorone	2.0	104	3.9	101	2.6
methoxychlor	2.0	94.0	5.6	94.0	2.6
methyl parathion	2.0	100	3.7	110	2.5
metolachlor	2.0	97.5	5.6	97.8	2.1
metribuzin	2.0	102	2.6	120	0.77
mevinphos	2.0	89.9	5.3	96.3	4.5
MGK 264(a)	1.3	92.9	2.3	86.4	4.1

MGK 264(b)	0.67	95.9	4.1	92.6	3.9
molinate	2.0	94.9	5.0	89.1	2.3
napropamide	2.0	95.9	4.6	104	3.1
nitrofen	2.0	111	6.2	96.4	5.3
nonachlor, trans	2.0	105	4.0	98.4	5.3
norflurazon	2.0	98.1	3.8	101	6.4
oxyfluorfen	2.0	91.8	7.7	94.9	5.0
pebulate	2.0	89.3	4.4	84.3	1.0
pentachlorophenol	8.0	98.3	2.6	97.3	2.9
permethrin, cis	1.0	92.6	6.0	108	4.3
permethrin, trans	3.0	91.4	7.1	100	2.3
phenanthrene	2.0	106	3.7	104	2.5
phorate	2.0	95.5	2.0	98.9	4.8
phosphamidon	2.0	100	3.3	114	2.8
profenofos	2.0	98.4	4.4	108	1.7
prometon	2.0	96.9	2.1	88.4	5.7
prometryn	2.0	98.8	4.9	94.6	4.9
pronamide	2.0	95.1	5.5	96.6	3.4
propachlor	2.0	106	1.4	110	3.4
propazine	2.0	101	2.8	96.4	2.5
pyrene	2.0	106	1.2	107	2.2
simazine	2.0	101	3.7	96.1	3.0
simetryn	2.0	98.6	4.2	87.9	3.8
tebuconazole	2.0	94.4	5.4	96.0	5.2
tebuthiuron	2.0	88.5	7.1	101	1.8
terbacil	2.0	101	5.6	95.8	10
terbutryn	2.0	97.3	5.1	87.4	5.1
tetrachlorvinphos	2.0	97.5	6.6	104	5.4
triadimefon	2.0	99.3	2.1	101	3.0
tribufos+merphos	4.0	96.6	7.3	107	2.0
trifluralin	2.0	90.8	2.2	90.0	4.5
vernolate	2.0	88.0	6.0	88.3	1.8
vinclozolin	2.0	98.6	7.0	111	1.7
PCB Congeners by IUPAC#					
2-chlorobiphenyl (1)	2.0	90.8	1.7	96.3	5.6
4-chlorobiphenyl (3)	2.0	96.1	0.65	99.9	5.4
2,4'-dichlorobiphenyl (8)	2.0	97.9	0.87	86.6	6.3
2,2',5-trichlorobiphenyl (18)	2.0	101	3.8	90.9	7.5
2,4,4'-trichlorobiphenyl (28)	2.0	101	2.2	86.1	4.9
2,2',3,5'-tetrachlorobiphenyl (44)	2.0	93.0	3.0	88.1	6.9
2,2',5,5'-tetrachlorobiphenyl (52)	2.0	97.1	3.1	87.8	6.6
2,3',4',5-tetrachlorobiphenyl (70)	2.0	107	1.5	88.1	5.2
2,3,3',4',6-pentachlorobiphenyl (110)	2.0	107	1.3	92.0	6.3
2,3',4,4',5-pentachlorobiphenyl (118)	2.0	108	1.2	91.3	6.4
2,2',3,4,4',5'-	2.0	110	2.0	93.8	5.3

hexachlorobiphenyl (138)					
2,2',3,4',5',6-hexachlorobiphenyl (149)	2.0	106	1.4	91.5	6.0
2,2',4,4',5,5'-hexachlorobiphenyl (153)	2.0	108	0.89	91.3	6.3
2,2',3,4,4',5,5'-heptachlorobiphenyl (180)	2.0	99.6	1.0	87.6	7.2
Surrogate Analytes					
1,3-dimethyl-2-nitrobenzene	5.0	91.7	7.5	89.0	6.5
benzo[a]pyrene-<i>d</i>₁₂	5.0	103	2.7	104	3.0
triphenyl phosphate	5.0	104	0.58	112	3.5

- Data obtained on the instrumentation described in Sect. 13.1.1.4.
- Tap water from a ground water source with high mineral content. Tap water hardness was 300mg/L as calcium carbonate
- Tap water from a surface water source. TOC of 2.4 mg/L
- Recoveries have been corrected to reflect the native amount in the unfortified matrix water.

Toxaphene from Fortified Reagent Water Precision and Accuracy: Extracts Analyzed by SIM

Fortified Concentration 10 µg/L N=4	Mean % Recovery	RSD
	111	1.8

Method 525.3, "Determination Of Semivolatile Organic Chemicals In Drinking Water By Solid Phase Extraction And Capillary Column Gas Chromatography/ Mass Spectrometry (GC/MS)," Ver 1.0, February 2012, Jean W. Munch and Paul E. Grimmett (U.S. EPA, Office of Research and Development, National Exposure Research Laboratory), David J. Munch and Steven C. Wendelken (U.S. EPA, Office of Water, Office of Ground Water and Drinking Water, Technical Support Center) Mark M. Domino (Industrial and Environmental Services, LLC) Alan D. Zaffiro and Michael L. Zimmerman (Shaw Environmental and Infrastructure, Inc.), National Exposure Research Laboratory Office Of Research And Development, U. S. Environmental Protection Agency, Cincinnati, Ohio 45268

Complete details at www.epa.gov/safewater/methods/methods.html



Determination of Phenols in Drinking Water by Solid-Phase Extraction (SPE) and Capillary Column Gas Chromatography (GC/MS)

UCT Product Numbers:

ECDVB156 (500 mg DVB, 6 mL cartridge) **or**
ECUNIDVB500 (500 mg DVB, 83 mL cartridge)
ECSS15M6 (5 g Na₂SO₄, 6 mL cartridge)

Method 528 Revision 1.0

February 2009

Procedure

1. Condition Cartridge(s)

- Add 6N HCl solution dropwise to water sample to adjust pH to less than 2
- Rinse the cartridge(s) with 3 mL of methylene chloride then draw to waste
- Repeat** methylene chloride rinse 2 more times
- Rinse the cartridge(s) with 3 mL of methanol then draw to waste
- Repeat** methanol rinse 2 more times
- After the third rinse, leave enough methanol in the cartridge(s) to cover the frit

Do not let the cartridge(s) dry out at this point otherwise start over at step d)

- Rinse the cartridge(s) 3 times with 3 mL of 0.05N HCl and draw to top of frit
- With vacuum off, add 3 mL of 0.05N HCl to the cartridge(s)

2. Sample Addition

- Begin adding the water sample to the cartridge(s)
- Adjust the vacuum so the flow rate is about 20 mL/minute (extraction will require about 50 minutes)
- After extraction allow cartridge(s) to dry under vacuum for 10-15 minutes before proceeding to the next step

A dry cartridge is important for good recoveries

3. Extract Elution

- Rinse the inside of the sample bottle a 10 mL portion of methylene chloride


- b) Add this to the cartridge(s) then draw this through to the collection tube in a dropwise fashion
- c) Add another 3 mL of methylene chloride to the cartridge(s) then slowly draw through to collection tube

4. Eluate Drying

- a) Dry the eluate by passing through a pre-rinsed anhydrous sodium sulfate column **ECSS15M6** and collect eluate in a clean tube
- b) Rinse the sodium sulfate column with a 2 x 3 mL aliquots of methylene chloride and collect in the tube
- c) Concentrate the extract to about 0.9 mL in a warm water bath (40° C) under a gentle stream of N₂
- d) Adjust final volume to 1.0 mL with methylene chloride

5. Analysis

- a) Analyze extract by GC/MS



For complete details on Method 528 the analyst is referred to: Method 528 Revision 1.0, Determination of Phenols in Drinking Water by Solid-Phase Extraction (SPE) and Capillary Column Gas Chromatography (GC/MS), J. W. Munch, April 2000, National Exposure Research Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268



Determination of Explosives and Related Compounds in Drinking Water by Solid-Phase Extraction and Capillary Column Gas Chromatography/Mass Spectrometry (GC/MS)

UCT Product Numbers:

ECDVB156 (500 mg DVB, 6 mL cartridge) or

ECUNIDVB500 (500 mg DVB, 83 mL cartridge)

ECSS15M6 (5 g Na₂SO₄, 6 mL cartridge)

EPA Method 529

March 2012

Procedure

1. Cartridge Conditioning

- Mark the level of liquid on the sample bottle for accurate volume determination
- Verify that the sample pH is between 5 and 7 and free of chlorine
- Add an aliquot of the surrogate fortification solution for a final concentration of 1-5 µg/L and mix (See Notes)
- Place cartridge(s) on vacuum manifold
- Rinse the cartridge with 3 x 5 mL aliquots of ethyl acetate
- Let the cartridge drain completely between each solvent addition
- Rinse the cartridge with 3 x 5 mL aliquots of methanol

Do not let the cartridge go dry after this step otherwise repeat starting at step h)

- Rinse the cartridge with 2 x 10 mL aliquots of reagent water. Leave a visible layer of water covering the cartridge frit

Proper conditioning of the cartridge is essential for good precision and accuracy

2. Sample Extraction

- Add the water sample to the cartridge
- Adjust the vacuum setting so that the flow rate is 10-15 mL/minute and draw through all sample
- After extraction draw air through the cartridge for **10 minutes** at full vacuum

3. Cartridge Elution

- Insert a clean collection tube in the manifold
- Rinse the inside of the sample bottle with 5 mL ethyl acetate and add to cartridge

- c) Add another 5 mL aliquot of ethyl acetate to the cartridge and draw through dropwise to the collection tube
- d) Pass the ethyl acetate eluate through a **ECSS15M6** drying tube or use a bed of anhydrous sodium sulfate
- e) Wash the drying tube using 4 mL of ethyl acetate
- f) Concentrate the eluate under a gentle stream of N₂ in a warm water bath to 0.9 mL
- g) Add internal standards and bring final volume to 1.0 mL for analysis
- h) Sample is ready for analysis

4. Sample Analysis by GC/MS

- a) Inject 1-2 µL

Notes

Surrogate standards (nitrobenzene-*d*₅; 1, 3, 5-trimethyl-2-nitro-benzene; and 1, 2, 4-trimethyl-5-nitro-benzene) are added to all calibration standards, samples, LFBs, LFM, FDs, and LRBs

Internal Standard Stock - Prepare a stock standard of the internal standard, 3, 4-dinitrotoluene (CAS# 610-39-9), at approximately 2.0 mg/mL in ethyl acetate. Store at 0° C or less.

Surrogate Standard Stock Solutions - Prepare or purchase individual stock standards of the following surrogate analytes: 1, 3, 5-trimethyl-2-nitrobenzene commonly known as 2-nitromesitylene [CAS603-71-4], 1, 2, 4-trimethyl-5-nitrobenzene commonly known as 2-nitropseudocumene [CAS# 610-91-3], and nitrobenzene-*d*₅ [CAS# 4165-60-0]. Standards of 2-nitromesitylene and 2-nitropseudocumene should be prepared in methanol at concentrations of 1-2 mg/mL. Nitrobenzene-*d*₅ should be prepared in methylene chloride at a concentration of 1-2 mg/mL. Store at 0° C or less. During method development, these solutions were stable for at least 1 year.

*For complete details on Method 529, September 2002, the analyst is referred to: J.W.Munch, "Determination Of Explosives And Related Compounds In Drinking Water By Solid-Phase Extraction And Capillary Column Gas Chromatography/Mass Spectrometry (GC/MS)", National Exposure Research Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268



Determination of Phenylurea Compounds In Drinking Water By Solid-Phase Extraction And High Performance Liquid Chromatography With UV Detection

UCT Part Numbers:

EUC18156 (6 mL cartridge with 500 mg C18)

EPA Method 532 Revision 1.0

September 2010

This is an HPLC method for the determination of phenylurea pesticides in drinking water. It is applicable to phenylurea compounds that are efficiently extracted from the water using a C18 Solid-Phase cartridge. Accuracy, precision, and method detection limit (MDL) data have been generated for the following compounds* in reagent water and finished ground and surface waters.

Phenylurea Compounds

Analyte	CAS #
Diflubenzuron	35367-38-5
Diuron	330-54-1
Fluometuron	2164-17-2
Linuron	330-55-2
Propanil	709-98-8
Siduron	1982-49-6
Tebuthiuron	34014-18-1
Thidiazuron	51707-55-2

Method Summary

A 500 mL water sample is extracted using an SPE C18 cartridge. The phenylurea pesticides and surrogate compounds are captured by the solid-phase then analytes and surrogates are eluted with methanol. The extract is concentrated to a final volume of 1.0 mL. Analytes are chromatographically separated using an HPLC and detected using UV/Vis. Identification of target and surrogate analytes and quantitation is accomplished by retention times and analyte responses using external standard procedures. Sample extracts with positive results are solvent exchanged and confirmed using a second, dissimilar HPLC column that is also calibrated using external standard procedures.

Safety

- Each chemical compound should be treated as a potential health hazard, and exposure to these chemicals should be minimized using proper handling technique

Sample Collection, Preservation, and Storage

- Grab samples must be collected using 500 mL amber or clear glass bottles fitted with PTFE lined screw caps
- Prior to field shipment, sample bottles must contain, 0.25 g of cupric sulfate and 2.5 g of Trizma crystals to each bottle for each 500 mL of sample collected
- Alternately, the Tris buffer may be prepared by adding 2.35 g of Tris HCl and 0.15 g Tris to the sample bottle in addition to the 0.25 g of cupric sulfate
- When sampling from a cold water tap, remove the aerator so that no air bubbles will be trapped in the sample. Allow the system to flush until the water temperature has stabilized (usually about 3-5 minutes). Collect samples from the flowing system.
- When sampling from an open body of water, fill a 1 quart wide-mouth bottle or 1L beaker with sample. Carefully fill sample bottles from the container
- Sampling equipment, including automatic samplers, must be free of plastic tubing, gaskets, and other parts that may leach interfering analytes into the water sample
- Fill sample bottles, taking care not to over fill bottles so as not to flush out the sample preservation reagents. Samples do not need to be collected headspace free
- After collecting the sample, replace cover then agitate by hand for 1 minute. Keep samples sealed from collection time until extraction

Preservation Reagents

- Cupric Sulfate, $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ (ACS Grade or equivalent) – Added as a biocide to guard against potential degradation of method analytes by microorganisms
- Trizma Preset Crystals, pH 7.0 (Sigma # T 3503 or equivalent) – Reagent grade. A premixed blend of Tris [Tris(hydroxymethyl)aminomethane] and Tris HCL [Tris(hydroxymethyl)aminomethane hydrochloride]
- Alternatively, a mix of the two components with a weight ratio of 15.5/1; Tris HCL/Tris may be used. These blends produce a pH of 7.0 at 25°C in reagent water. Tris functions as a buffer, binding free chlorine in chlorinated finished waters and prevents the formation of a copper precipitate

Sample Shipment and Storage

- All samples should be iced during shipment and must not exceed 10° C during the first 48 hours after collection. Samples should be confirmed to be at or below 10°C when they are received at the laboratory
- Samples stored in the lab must be held at or below 6° C until extraction, but should not be frozen

Sample and Extract Holding Times

- Method analytes are stable for 14 days in water samples that are collected, dechlorinated, preserved, shipped and stored as described in above
- Samples must be extracted within 14 days
- Sample extracts may be stored in methanol at 0°C or less for up to 21 days after extraction. Samples that are exchanged into reagent water/acetonitrile (60/40) for confirmational analysis may be stored 7 days at 0° C or less. The combined extract holding time should not exceed 21 days

Interferences

- All glassware must be meticulously cleaned. Wash glassware with detergent and tap water, rinse with tap water. Follow by a reagent water rinse
- Non-volumetric glassware can be heated in a muffle furnace at 400° C for 2 hours
- Method interferences may also be caused by contaminants in solvents, reagents (including reagent water), sample bottles and caps, and other sample processing hardware. This may lead to discrete artifacts or elevated baselines in the chromatograms
- All solvents and equipment must be routinely demonstrated to be free from interferences under the conditions of the analysis. **Subtracting blank values from sample results is not permitted in this method**
- Matrix interferences may be caused by contaminants that are co-extracted from the sample

Surrogate Analytes (SUR),

Monuron (CAS #150-68-5) & Carbazole (CAS# 86-74-8)

SUR STOCK SOLUTION (5 to 7 mg/mL)

- Accurately weigh 25 to 35 mg of the neat SUR to the nearest 0.1 mg into a tared, 5 mL volumetric flask. Dilute to the mark with the appropriate solvent, methanol for monuron and acetonitrile for cabazole
- Prepare each compound individually. They will be combined in the SUR primary dilution standard

SUR Primary Dilution Standard (500 µg/mL)

- Prepare the SUR Primary Dilution Standard (PDS) by dilution of the SUR stock standards. Add enough of each of the SUR stock standards to a volumetric flask partially filled with methanol to make a 500 µg/mL solution when filled to the mark with methanol

Analyte Stock Standard Solution

- Prepare analyte stock standard solutions for all compounds in methanol except thidiazuron and Diflubenzuron which are prepared in acetone
- Acetone elutes early in the chromatogram and will not interfere with compound quantitation
- Accurately weigh 10 to 12 mg of neat material to the nearest 0.1 mg in volumetric flasks. Thidiazuron is difficult to dissolve but 10 mg of pure material should dissolve in a 10 mL of acetone. Sonication may be used to help dissolve these compounds.

Analyte Primary Dilution Standard (PDS, 200 µg/mL and 10 µg/mL)

- Prepare the Analyte PDS by dilution of the stock standards. Add enough of each stock standard to a volumetric flask partially filled with methanol to make a 200 µg/mL solution when filled to the mark with methanol
- Once prepared, a dilution of the 200 µg/mL solution may be used to prepare a 10 µg/mL solution used for low concentration spiking. The PDS's can be used for preparation of calibration and fortification solutions
- Use at least 5 calibration concentrations to prepare the initial calibration curve
- The lowest concentration of calibration standard must be at or below the MRL
- In this method, 500 mL of an aqueous sample is concentrated to a 1 mL final extract volume

The calibration standards for the development of this method are prepared as below

Preparation of Calibration Curve Standards					
CAL Level	PDS Conc $\mu\text{g/mL}$	Volume PDS Standard μL	Final Vol CAL Standard mL	Final Conc CAL Standard $\mu\text{g/mL}$	Equivalent Conc in 500 mL sample $\mu\text{g/L}$
1	10	25	1	0.25	0.50
2	10	50	1	0.50	1.00
3	200	5.0	1	1.00	2.00
4	200	25	1	5.00	10.0
5	200	50	1	10.0	20.0
6	200	75	1	15.0	30.0

Procedure

1. Cartridge Conditioning

- Rinse the cartridge(s) using 5 mL of methanol
- Slowly draw through leaving a layer covering the frit
- Repeat using a second 5 mL aliquot and slowly draw through. Leave enough methanol to cover the frit

Note: Do not let the cartridge dry out after addition of methanol otherwise repeat

- Add 5 mL of reagent water to the cartridge and draw through to the level of the frit
- Repeat with a second 5 mL aliquot of reagent water

2. Sample Extraction

- Add the water sample to the cartridge. Adjust vacuum setting for a flow rate of about 20 mL/minute
- Rinse the sample container with reagent water and add to the cartridge
- Once the sample has been extracted, turn vacuum setting to high and draw air through for 15 minutes to dry cartridge

3. Sample Elution

- Insert a clean collection tube in the vacuum manifold
- Add about 3-5 mL of methanol to each cartridge
- Draw through to the top of the frit then turn off vacuum
- Allow to soak for 30 seconds

- e) Add a second 2-3 mL of methanol and draw through cartridge as a drip
- f) Repeat methanol addition with a third 1-2 mL portion

4. Extract Concentration

- a) Concentrate the extract to about 0.5 mL in a warm water bath 4° C under a gentle stream of N₂
- b) Transfer to a 1 mL volumetric flask rinsing the collection tube with methanol
- c) Bring to volume with methanol
- d) Sample is ready for analysis

5. Solvent Exchange for Confirmation Analysis

- a) Samples to be confirmed must be exchanged into reagent water/acetonitrile (60/40)
- b) Transfer the remaining 980 µL of the extract to a 1 mL volumetric (or other appropriate collection tube)
- c) Mark the sample volume then take the extract to dryness in a warm water bath (at ~ 40° C) under a gentle stream of N₂
- d) Reconstitute the residue with a mixture of reagent water/acetonitrile (60/40) to the mark made before the extract was taken to dryness. Redissolve the film as thoroughly as possible. Use of a vortex mixer is recommended. Transfer to an appropriate autosampler vial
- e) The best recoveries are obtained when extracts are immediately reconstituted

6. HPLC Analysis—Primary Analysis

Primary HPLC Column: Symmetry 4.6 x 150 mm with 3.5 μ m C18 phase or equivalent

HPLC Conditions	
Solvent A	25 mM phosphate buffer
Solvent B	Acetonitrile
40% B	0-9.5 minutes
Linear gradient 40-50% B	9.5-10.0 minutes
Linear gradient 50-60% B	10.0-14.0 minutes
Linear gradient 60-40% B	14.0-15.0 minutes
Flow Rate	1.5 mL/minute
Wavelength	245 nm
Equilibration time prior to next injection = 15 minutes	

Chromatographic Retention Time Data for the Primary Column

Peak Number	Analyte	Retention Time (min)
1	Tebuthiuron	2.03
2	Thidiazuron	2.48
3	Monuron (SUR)	2.80
4	Fluometuron	4.45
5	Diuron	5.17
6	Propanil	8.53
7	Siduron A	8.91
8	Siduron B	9.76
9	Linuron	11.0
10	Carbazole (SUR)	12.8
11	Diflubenzuron	13.9

7. HPLC Analysis--Confirmation Analysis

Confirmation Column: Supelcosil 4.6 x 150 mm packed with 5 μ m cyanopropyl stationary phase or equivalent

Chromatographic Retention Time Data for the Confirmation Column

Peak Number	Analyte	Retention Time (min)
1	Tebuthiuron	2.56
2	Thidiazuron	3.98
3	Monuron (SUR)	4.93
4	Fluometuron	5.94
5	Diuron	7.67
6	Propanil	9.53
7	Siduron A	10.1
8	Siduron B	10.8
9	Linuron	12.2
10	Carbazole (SUR)	14.3
11	Diflubenzuron	15.2

*The analyst should refer to EPA Method 532 "Determination of Phenylurea Compounds in Drinking Water by Solid-Phase Extraction and High Performance Liquid Chromatography with UV Detection", Revision 1.0 Issued June 2000, by M. V. Bassett, S.C. Wendelken, T.A. Dattilio, and B.V. Pepich (IT Corporation) D.J. Munch (US EPA, Office of Ground Water and Drinking Water) EPA, National Exposure Research Laboratory, Office of Research and Development, US Environmental Protection Agency, Cincinnati, Ohio 45268

Measurement of Chloroacetanilide and Other Acetamide Herbicide Degradates in Drinking Water by Solid-Phase Extraction and Liquid Chromatography/ Tandem Mass Spectrometry (LC/MS/MS)

UCT Part Number:

EC535156 (500 mg graphitized carbon black (GCB) 90 m²/g, 6 mL)

EPA METHOD 535

February 2009

Method Summary

A 250 mL water sample is drawn through and captured on a **EC535156** cartridge containing 0.5 grams of nonporous graphitized carbon. Acetanilide and acetamide compounds are eluted from the cartridge using a small quantity of methanol containing 10 mM ammonium acetate. The methanol extract is concentrated to dryness by blow down with N₂ in a water bath at 65°C then reconstituted with 1 mL of water containing 5 mM ammonium acetate. A 100 µL portion of the aqueous reconstitution is injected into an HPLC fitted with a C18 reverse phase analytical column. Detection occurs by tandem mass spectrometry and is compared to internal standards. A surrogate analyte of known concentration is measured with the same internal standard calibration procedure.

Interferences

Humic and/ or fulvic acid material, if present in the water source, is co-extracted with this method. High concentrations of these compounds can cause enhancement or suppression of the in the electrospray ionization source or low recoveries on the carbon SPE. Total organic carbon (TOC) is a good indicator of these interferences if present in the water sample.

1. Condition Cartridge

- a) Rinse the cartridge with 20 mL of 10 mM ammonium acetate/methanol solution
- b) Rinse cartridge with 30 mL of reagent water. Do not let water drop below level of cartridge packing
- c) Add about 3 mL of reagent water to the top of the cartridge

Note: Do not let the cartridge go dry during any step otherwise start over

2. Sample Addition

- a) Add sample water to the cartridge and adjust vacuum so the flow is about 10-15 mL/minute
- b) Rinse cartridge with 5 mL of reagent water
- c) Draw air or N₂ through the cartridge at high vacuum (10-15 in/Hg) for 3 minutes

3. Extract Elution

All glassware must be meticulously washed to avoid contamination

- a) Insert a clean collection tube into the extraction manifold
- b) Use 15 mL of 10 mM ammonium acetate/methanol and adjust vacuum to draw through at 5 mL/minute. Solvent will exit the cartridge in a dropwise fashion at this vacuum setting

4. Eluate Drying

- a) Concentrate the extract to dryness under a gentle stream of N₂ in a heated water bath at 60°-70° C to remove all of the ammonium acetate/methanol
- b) Reconstitute the dried eluate by adding 1 mL of 5 mM ammonium acetate/methanol solution

5. Extract Analysis

- a) Establish operating conditions for the liquid chromatograph and mass spectrometer according to Tables 1-4 in Section 17. See Table A below for RT and precursor ions
- b) If the analyte peak area exceed the range of the initial calibration curve, the extract may be diluted with 5 mM ammonium acetate/reagent water and adjusting internal standards to compensate for this dilution

Table A
Triple Quadrupole MS/MS Method Conditions

Analyte	Retention Time	Precursor Ion	Product Energy	Collision Energy
Propachlor OA	7.33	206	134	8
Flufenacet OA	8.67	224	152	10
Propachlor ESA	10.01	256	80	25
Flufenacet ESA	10.81	274	80	25
Dimethenamid OA	13.25	270	198	10
Dimethenamid ESA	14.87&15.11	320	80	25
Alachlor OA	15.86	264	160	10
Acetochlor OA	16.34	264	146	10
Alachlor ESA	18.46	314	80	25
Metolachlor OA	18.60	278	206	8
Acetochlor ESA	19.12	314	80	30
Metolachlor ESA	20.95	328	80	25
Dimethachlor ESA (sur)	12.18	300	80	25
Butachlor ESA (IS)	36.95	356	80	25

*For complete details on Method 535 Version 1.1 the analyst is referred to: J. A. Shoemaker and M. V. Bassett, April 2005, National Exposure Research Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268



Determination of Selected Perfluorinated Alkyl Acids in Drinking Water by Solid-Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

UCT Product:

ECDVB156P (500 mg DVB, PE frits, 6 mL cartridge)

EPA Method 537 Version 1.1

January 2011

Method 537 Analytes

Analyte	Acronym	CASRN
N-ethyl perfluorooctanesulfonamidoacetic acid	NEtFOSAA	--
N-methyl perfluorooctanesulfonamidoacetic acid	NMeFOSAA	--
Perfluorobutanesulfonic acid	PFBS	375-73-5
Perfluorodecanoic acid	PFDA	335-76-2
Perfluorododecanoic acid	PFDoA	307-55-1
Perfluoroheptanoic acid	PFHpA	375-85-9
Perfluorohexanesulfonic acid	PFHxS	355-46-4
Perfluorohexanoic acid	PFHxA	307-24-4
Perfluorononanoic acid	PFNA	375-95-1
Perfluorooctanesulfonic acid	PFOS	1763-23-1
Perfluorooctanoic acid	PFOA	335-67-1
Perfluorotetradecanoic acid	PFTA	376-06-7
Perfluorotridecanoic acid	PFTTrDA	72629-94-8
Perfluoroundecanoic acid	PFUnA	2058-94-8

Interferences

- Do not cover glassware with aluminum foil because PFAAs can be potentially transferred from the aluminum foil
- PFAA standards, extracts and samples should not come in contact with any glass containers or pipettes as these analytes can potentially adsorb to glass
- PFAA analyte, IS and SUR standards commercially purchased in glass ampoules are acceptable however all subsequent dilutions must be prepared and stored in polypropylene containers

- The method analytes in this method can also be found in many common laboratory supplies and equipment, such as PTFE (polytetrafluoroethylene) products, LC solvent lines, methanol, aluminum foil and SPE sample transfer lines. These items must be routinely demonstrated to be free from interferences (less than 1/3 the MRL for each method analyte) under the conditions of the analysis by analyzing laboratory reagent blanks
- Humic and/or fulvic material can be co-extracted during SPE and high levels can cause enhancement and/or suppression in the electrospray ionization source or low recoveries.
Total organic carbon (TOC) is a good indicator of humic content of the sample

Preservation Reagent

Compound	Amount	Purpose
Trizma®*	5.0 g/L	Buffer & free Cl removal

* Synonym: TRIS HCl (Tris(hydroxymethyl)aminomethane hydrochloride) CASRN 1185-53-1

Internal Standards and Surrogates

Internal Standards Table

Internal Standards	Acronym
Perfluoro-[1,2- ¹³ C ₂]octanoic acid	¹³ C-PFOA
Sodium perfluoro-1-[1,2,3,4- ¹³ C ₄]octanesulfonate	¹³ C-PFOS
N-deuteriomethylperfluoro-1-octanesulfonamidoacetic acid	d ₃ -NMeFOSAA

Internal Standard (IS) Primary Dilution Table

IS	Conc. of IS Stock (µg/mL)	Vol. Of IS Stock (µL)	Final Vol. of IS PDS (µL)	Final Conc. of IS PDS (ng/µL)
¹³ C-PFOA	1000	5.0	5000	1.0
¹³ C-PFOS	50	300.0	5000	3.0
d ₃ -NMeFOSAA	50	400.0	5000	4.0

Surrogates (SUR) Table

Surrogates	Acronym
Perfluoro-n-[1,2- ¹³ C ₂]hexanoic acid	¹³ C-PFHxA
Perfluoro-n-[1,2- ¹³ C ₂]decanoic acid	¹³ C-PFDA
N-deuterioethylperfluoro-1-octanesulfonamidoacetic acid	d ₅ -NEtFOSAA

Surrogate Primary Dilution Table

SUR	Conc. Of SUR Stock (µg/mL)	Vol. of SUR Stock (µL)	Final Vol. of SUR PDS (µL)	Final Conc. Of SUR PDS (ng/µL)
¹³ C-PFHxA	50	100.0	5000	1.0
¹³ C-PFDA	50	100.0	5000	1.0
d ₅ -NEtFOSAA	50	400.0	5000	4.0

Analyte Solvent Dilution Table

Analyte	Analyte Stock Solvent
PFHxA	96:4% (vol/vol) methanol:water
PFHpA	96:4% (vol/vol) methanol:water
PFOA	96:4% (vol/vol) methanol:water
PFNA	96:4% (vol/vol) methanol:water
PFDA	96:4% (vol/vol) methanol:water
PFUnA	96:4% (vol/vol) methanol:water
PFDaA	96:4% (vol/vol) methanol:water
PFTTrDA	100% ethyl acetate
PFTA	100% ethyl acetate
PFBS	100% methanol
PFHxS	100% methanol
PFOS	100% methanol
NEtFOSAA	100% methanol
NMeFOSAA	100% methanol

Internal Standard Primary Dilution (IS PDS) Standard

Prepare (or purchase) IS PDS at a suggested concentration of 1-4 ng/µL in 96:4% (vol/vol) methanol:water. Use 10 µL of this 1-4 ng/µL solution to fortify the final 1 mL extracts. This will yield a concentration of 10-40 pg/µL of each IS in the 1 mL extracts.

Procedure

1. Cartridge Preparation

- a) Insert a **ECDVB156P** cartridge in a vacuum manifold or automated extraction system
- b) Add 15 mL of methanol to the cartridge and slowly draw through under vacuum

Note: Do not let the cartridge dry out after addition of methanol otherwise start over

- c) Add 18 mL of reagent water
- d) Draw through under vacuum but do not let water level drop below cartridge frit

2. Sample Extraction

- a) Adjust vacuum setting to achieve a flow rate of about 10-15 mL/min
- b) Draw water through sample cartridge
- c) After sample extraction, rinse sample bottles and reservoir with 2 x 7.5 mL aliquots of reagent water and add to cartridge
- d) Dry cartridge by drawing air through it for 5 minutes at high vacuum

3. Sample Elution

- a) Insert a polypropylene (PP) collection tube in the vacuum manifold
- b) Rinse sample bottle and reservoir with 4 mL of methanol and add to cartridge
- c) Draw through cartridge in dropwise manner
- d) Rinse sample bottle and reservoir with another 4 mL of methanol and add to the cartridge
- e) Draw through cartridge in dropwise manner

4. Extract Concentration

- a) Concentrate the extract to dryness using a gentle stream of N₂ in a heated water bath 60-65° C
- b) Add the appropriate amount of 96:4% (vol/vol) methanol:water solution and the IS PDS to the collection vial
- c) Bring the volume to 1 mL and vortex
- d) Transfer a small aliquot with a plastic pipette to a polypropylene autosampler vial

NOTE: Do not transfer the entire 1 mL aliquot to the autosampler vial because the polypropylene autosampler caps do not reseal after injection. Do not store the extracts in the autosampler vials as evaporation losses can occur. Extracts can be stored in 15 mL centrifuge tubes

5. HPLC Analysis Conditions

LC Method Conditions

Time (min)	% 20 mM	% Methanol
Initial	60.0	40.0
1.0	60.0	40.0
25.0	10.0	90.0
32.0	10.0	90.0
32.1	60.0	40.0
37.0	60.0	40.0
Flow rate of 0.3 mL/min		10 µL injection
Waters Atlantis® dC ₁₈ 2.1 x 150 mm packed with 5.0 µm C ₁₈ stationary phase or equivalent		

ESI-MS Conditions

Polarity	Negative Ions
Capillary Needle Voltage	-3 kV
Cone Gas Flow	98 L/hr
Nitrogen Desolvation Gas	1100 L/hr
Desolvation Gas Temperature	350° C

MS/MS Method Conditions

Segment	Analyte	Precursor Ion (m/z)	Product Ion (m/z)	Cone Voltage	Collision Energy (v)
1	PFBS	299	80	40	25
2	PFHxA	313	269	15	10
3	PFHpA	363	319	12	10

3	PFHxS	399	80	40	40
4	PFOA	413	369	15	10
4	PFNA	463	419	12	10
4	PFOS	499	80	40	40
5	PFDA	513	469	15	10
5	NMeFOSAA	570	419	25	20
5	NEtFOSAA	584	419	25	20
5	PFUnA	563	519	15	10
5	PFDoA	613	569	15	10
6	PFTA	663	619	15	10
6	PFTA	713	669	15	10
2	¹³ C-PFHxA	315	270	15	10
5	¹³ C-PFDA	515	470	12	12
5	d ₅ -NEtFOSAA	589	419	25	20
4	¹³ C-PFOA	415	370	15	10
4	¹³ C-PFOS	503	80	40	40
5	d ₃ -NMeFOSAA	573	419	25	20

Summarized 15 from: Shoemaker, J.A., Grimmer P.E., Boutin, B.K., Method 537, Determination Of Selected Perfluorinated Alkyl Acids In Drinking Water By Solid-Phase Extraction And Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS), Version 1.1, September 2009, National Exposure Research Laboratory, Office Of Research And Development, U. S. Environmental Protection Agency, Cincinnati, Ohio 45268

Listing of instrument manufacturers does not constitute endorsement by UCT



Determination of Endothall in Drinking Water by Ion-Exchange Extraction, Acid Methanol Methylation and Gas Chromatography/Mass Spectrometry

UCT Product:
EC548006 (Bio-Rex® 5, 6 mL)

EPA Method 548.1 Revision 1.0

September 2009

Procedure

1. Condition Cartridge

- a) Remove the seal caps on each end of the **EC548006** cartridge(s) and place on vacuum manifold
- b) Draw each of the following reagents through the cartridge at a rate of 10 mL/minute, slow drip
- c) Add 10 mL of methanol and draw through to top of frit
- d) Add 10 mL of reagent water and draw through to top of frit
- e) Add 10 mL of 10% H₂SO₄ in methanol and draw through to top of frit
- f) Add 10 mL of reagent water and draw through to top of frit
- g) Add 20 mL of 1 N NaOH and draw through to top of frit
- h) Add 20 mL of reagent water and draw through to top of frit

Note: Do not allow the cartridge to become dry between steps otherwise repeat steps starting with c)

It is critical that the extraction steps be followed exactly in order for the cartridge to effectively function in sample clean up and extraction

2. Sample Addition

- a) Fill the reservoir with 60 mL of sample and adjust vacuum for a sample flow rate of 3 mL/minute
- b) Add the remaining sample to keep the reservoir filled

- c) After the sample has been drawn through the cartridge add 10 mL of methanol and draw through
- d) Dry cartridge for 5 minutes under 10-20 in Hg vacuum
- e) Place a culture tube inside the manifold to collect the eluant

3. Extract Elution

- a) Elute the cartridge with 8 mL of 10% H₂SO₄ in methanol
- b) Follow with 6 mL of methylene chloride (CH₂Cl₂).
- c) Elute over a 1 minute period

4. Sample Derivatization and Partition

- a) Place a cap on the culture tube and hold for 1 hour at 50° C
- b) Pour the contents of the culture tube into a 125 mL separatory funnel rinsing the tube with 2 x 0.5 mL aliquots of methylene chloride
- c) Add the rinse to the separatory funnel
- d) Add 20 mL of 10% Na₂SO₄ in reagent water to the separatory funnel. Vigorously shake the separatory funnel several times venting the funnel each time
- e) Allow the phases to separate then drain the organic layer into a 15 mL graduated centrifuge tube
- f) Repeat the above extraction with 2 additional 2 mL aliquots of methylene chloride. Add this to the methylene chloride in the centrifuge tube
- g) Sample is ready for GC/MS analysis

5. Analysis

- a) Analyze the extract by injecting 2 µL of the concentrated extract into a GC/MS
- b) Identify endothall by comparison of its mass spectrum to a reference sample


Retention Times and Method Detection Limits

Compound	Retention Time (min)			Method Detection Limits		
	Column A	Column B	Column C	GC/MS	2 µg/L spike	FID
Endothall	16.02	19.85	18.32	1.79		0.7
D10-Acenaphthene	14.69					

Column A: DB-5 fused silica capillary for GC/MS, 30 m x 0.25 mm, 0.25 micron film
 MS inlet temperature = 200°C
 Injector temperature = 200°C
 Temperature Program: Hold five minutes at 80°C, increase to 260°C at 10°/min, hold 10 minutes

Column B: FID primary column, RTX Volatiles, 30 m x 0.53 mm I.D., 2 micron film
 Detector temperature = 280°C
 Injector Temperature = 200°C
 Carrier gas velocity = 50 cm/sec
 Temperature program: Same as Column A.

Column C: FID confirmation column, DB-5, 30 m x 0.32 mm ID, 0.25 micron film.
 Carrier Gas velocity = 27 cm/sec
 Same injector, detector, and temperature program as Column A



Endothall Recovery Using EC548006 Cartridges

Interferences

- Major potential interferences in this ion-exchange procedure are other naturally occurring ions such as calcium, magnesium and sulfate. Calcium and magnesium (>100 mg/L) can complex with the endothall anion and make it unavailable for capture as an anion
- Sulfate anions (>250 mg/L) can act as a counter ion displacing anionic endothall on the ion exchange column. Elevated levels of these ions may contribute to reduced recovery of the primary analyte

One or both of the following remedies may be used reduce these interferences:

- Sample dilution to reduce the concentration of these ions (10:1)

- Ethylenediamine tetraacetic acid (EDTA) addition to complex the cations (186 mg/100 mL sample)

For samples containing moderately high levels of these ions, add 186 mg of EDTA per 100 mL sample (0.005 M). For samples containing higher levels of sulfate, sample dilution may be required in addition to the EDTA.

Using western surface water as an example (2000 mg/L sulfate) it was successfully analyzed after dilution by a factor of 10 and the addition of 75 mg EDTA per 100 mL of the diluted sample (0.002 M). Samples containing intermediate levels of sulfate can be analyzed with smaller dilution factors. Guidelines on dilution factors and EDTA addition are shown below.

Sulfate mg/L	Dilution Factor	Added EDTA mg/100 mL
<250	1:1	186
250-500	1:2	125
500-1250	1:5	75
>1250	1:10	75

Note: Dilution should not be employed if adequate recovery is attained by the addition of EDTA alone

*For complete details on Method 548.1 "Determination of Endothall in Drinking Water by Ion-Exchange Extraction, Acid Methanol Methylation and Gas Chromatography/Mass Spectrometry", the analyst is referred to: J. W. Hodgeson, August 1992, Environmental Monitoring Systems Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268



Determination of Diquat and Paraquat in Drinking Water by Liquid-Solid Extraction and High Performance Liquid Chromatography with Ultraviolet Detection*

UCT Product Number:

EEC08156 (500 mg C8, 6 mL cartridge) or

ECUNI549 (500 mg C8, 83 mL Universal cartridge)

EPA Method 549.2 Revision 1.0

July 2011

Analyte	CASRN
Diquat 1,1'-ethylene-2,2'-bipyridium dibromide salt	85-00-7
Paraquat 1,1'-dimethyl-4,4'-bipyridium dichloride salt	1910-42-5

Initial Preparation

- Since diquat and paraquat are ionic analytes there is the potential for adsorption on glass surfaces
- **Use only plastic labware.** Labware must be thoroughly washed and dried before use
- Adjust a 250 mL of sample to pH 7 - 9 with 10 % aqueous sodium hydroxide or 10% aqueous hydrochloric acid solution depending upon initial pH
- Assemble a C8 extraction cartridge in an appropriate manifold apparatus
- If the sample contains particulates, filter through 0.45 µm Nylon membrane filter
- **Ammonium hydroxide is volatile. Make fresh solutions daily from relatively new ammonium hydroxide stock**

Sample Clean-up

- Clean-up procedures may not be necessary for a relatively clean sample matrix
- If the sample contains particulates the entire sample should be passed through a 0.45 mm Nylon or PTFE membrane filter into a plastic container before starting extraction

Stock Standard Solutions

Diquat dibromide and Paraquat dichloride Stock Solutions (1000 mg/L)

1. Dry diquat (diquat dibromide monohydrate) and paraquat (paraquat dichloride tetrahydrate) salts in an oven at 110°C for three hours. Cool in a desiccator
2. Repeat process to a constant weight.
3. Weigh 0.1968 g of dried diquat salt and 0.1770 g of dried paraquat salt
4. Transfer to a silanized glass or polypropylene 100 mL volumetric flask. Add approximately 50 mL of deionized water then dilute to the mark with deionized water

Calibration

In order to closely match calibration standards to samples, process standards by the following method:

- Condition a cartridge according to section 1 below.
- Pass 250 mL of reagent water through the cartridge and discard the water.
- Dry the cartridge by passing 5 mL of methanol through it. Discard the methanol.
- Pass 4.0 mL of the eluting solution through the cartridge and catch in a 5 mL silanized volumetric flask.
- Fortify the eluted solution with 100 µL of the ion-pair concentrate and with 500 µL of the stock standard and dilute to the mark with eluting solution. This provides a 10:1 dilution of the stock.
- Use serial dilution of the calibration standard by the same method to achieve lower concentration standards.

Procedure

The cartridge must be conditioned properly before extraction

1. Condition Cartridge

- a) Place C8 cartridge(s) on a vacuum manifold system
- b) Draw the following solutions through the cartridge in the stated order. The flow rate through the cartridge should be approximately 10 mL/min

Note: Do not let the cartridge go dry once starting the addition of solutions

- c) Add 5 mL of reagent water to the cartridge and draw through to waste
- d) Add 5 mL of methanol to the cartridge and soak for about one minute

- e) Apply vacuum to draw most of the methanol through the cartridge. Leave a thin layer on top of the frit
- f) Add 5 ml reagent water to the cartridge
- g) Apply vacuum and draw most of the water through the cartridge. Leave a thin layer of water on the frit
- h) Apply 5 mL of conditioning **Solution A** to the cartridge

Solution A: Dissolve 0.500 grams cetyl trimethyl ammonium bromide and 5 mL of ammonium hydroxide in 500 mL of reagent water. Dilute to 1000 mL

- i) Draw a small amount through the cartridge leaving a thin layer on the frit
- j) Soak for one minute
- k) Use 5 mL of reagent grade water to rinse the **Solution A** from the cartridge. Allow a thin layer of water to remain on the cartridge frit
- l) Rinse the cartridge with 10 mL of methanol
- m) Rinse the cartridge with 5 mL of reagent grade water
- n) Condition the cartridge with 20 mL of **Solution B**
- o) **Solution B:** Dissolve 10 g 1-hexanesulfonic acid sodium salt and 10 mL of ammonium hydroxide in 250 mL of DI water then dilute to 500 mL
- p) Retain **Solution B** in the cartridge to keep it activated. **Do Not Rinse**

2. Sample Extraction

- a) Determine the pH of the sample. Adjust to 7.0 – 9.0 with 10% NaOH or 10% v/v HCl before extracting
- b) Using a volumetric flask add 250 mL of the water sample to the reservoir and start the vacuum at a rate of 3 to 6 mL per minute
- c) Draw the sample through the cartridge draining as much water from the sample bottle as possible
- d) Rinse the cartridge with 5 ml of HPLC grade methanol
- e) Draw vacuum through the cartridge for 1 minute to dry
- f) Remove the filtration assembly and insert a silanized 5 mL volumetric (plastic vessel is preferred) flask for collection of the eluate

3. Cartridge Elution

- a) Add 4.5 ml of **Cartridge Eluting Solution** to the cartridge
- b) Allow to soak for one minute

Cartridge Elution Solution: Dissolve 13.5 mL of orthophosphoric acid and 10.3 mL of diethylamine in 500 mL of DI water, then dilute to 1 liter

- c) Elute at 1-2 mL (drop by drop) per minute drawing all solution through the cartridge
- d) Using cartridge **ion-pair solution**, add 100 µL to the flask

Ion-pair Concentrate: Dissolve 3.75 grams of 1-hexanesulfonic acid in 15 mL of the **Cartridge Elution Solution** and dilute to 25 mL in a volumetric flask with additional **Cartridge Elution Solution**

- e) Bring the eluate to a known volume of 5 mL using **Cartridge Elution Solution**
- f) The extract is now ready for HPLC analysis as shown below

4. HPLC Analysis

Mobile Phase – Prepare mobile phase by adding reagents 1-4 to 500 mL DI water:

- a) 13.5 mL of orthophosphoric acid
- b) 10.3 mL of diethylamine
- c) 3.0 g of 1-hexanesulfonic acid, sodium salt
- d) Mix and bring to a final volume of 1 L with DI water

HPLC Conditions

Column: Phenomenex Spherisorb, 3F, 4.6 mm x 100 mm or equivalent

Column Temperature: 35° C

Flow Rate: 2.0 mL/min., Ion-Pair Mobile Phase

Injection Volume: 200 µL

Photodiode Array Detector Settings:

Wavelength Range: 210 - 370 nm

Sample Rate: 1 scan/sec.

Wavelength Step: 1 nm

Integration Time: 1 sec.

Run Time: 5.0 min.

Quantitation Wavelengths: Diquat 308 nm, Paraquat 257 nm

*EPA Method 549.2 Revision 1.0, Determination of Diquat and Paraquat in Drinking Water by Liquid-Solid Extraction and High Performance Liquid Chromatography with Ultraviolet Detection, J.W. Munch (USEPA) and W.J. Bashe (DynCorp/TAI) - Method 549.2, Revision 1.0 (1997), National Exposure Research Laboratory, Office Of Research And Development, U.S. Environmental Protection Agency, Cincinnati, Ohio 45268



Determination of Polycyclic Aromatic Hydrocarbons in Drinking Water by Liquid-Solid Extraction and HPLC with Coupled Ultraviolet and Fluorescence Detection*

UCT Product Number:

ECUNIPAH (2000 mg unendcapped C18, 83 mL cartridge)

ECSS15M6 (5000 mg Na₂SO₄, 6 mL cartridge)

EPA Method 550.1

September 2010

Procedure

1. Cartridge Preparation

- Wash with four (4) x 10 mL aliquots of methylene chloride (MeCl₂)
- Wash with four (4) x 10 mL aliquots of methanol (MeOH)
- Wash with two (2) x 10 mL aliquots of reagent water

Do not let the cartridge dry out after step 1) c. otherwise repeat starting at 1) b.

2. Sample Extraction

- Adjust the vacuum setting for a flow rate of 10-15 mL per minute
- Add the 1 liter sample to the cartridge
- Rinse sample bottle with reagent water, add to cartridge and draw through
- Dry cartridge by drawing full vacuum for 10 minutes

3. Sample Elution and Drying

- Elute the cartridge dropwise by using two (2) x 5 mL aliquots of MeCl₂ and collect
- Rinse sample container with 5 mL of MeCl₂, add to cartridge and draw through
- Prepare a sodium sulfate drying column **ECSS15M6** by rinsing with 10 mL of MeCl₂ and discard
- Add the eluate to the drying column, draw through and collect
- Rinse the drying column with a 2 mL aliquot of MeCl₂ and collect

4. Sample Evaporation

- a) Evaporate the extract using a gentle stream of N₂ with a water bath or heating block temperature of 40°C. Evaporate to 1.0 mL
- b) Add 3.0 mL of acetonitrile (ACN)
- c) Concentrate to a final volume of 0.5 mL

5. Sample Analysis

- a) Inject 5 - 100 µL into an HPLC

*See "Determination of Polycyclic Aromatic Hydrocarbons in Drinking Water by Liquid-Solid Extraction and HPLC with Coupled Ultraviolet and Fluorescence Detection", W. J. Bashe & T.V. Baker (Technology Applications, Inc, Environmental Monitoring Systems Laboratory, US Environmental Protection Agency, Cincinnati, OH



Determination of Haloacetic Acids and Dalapon in Drinking Water by SPE and GC/ECD*

UCT Product Numbers:

EUQAX156 (quaternary amine with Cl⁻ counter ion, 6 mL cartridge)**

CLTTP050 (CLEAN-THRU[®] Tips)

AD0000AS (cartridge adapter)

RFV0075P (75 mL reservoir)

ECSS25K (sodium sulfate anhydrous, ACS, 60 mesh)

Method 552.1 Revision 1.0

April 2011

Summary

Method 552.1 is a gas chromatographic method for determination of haloacetic acid compounds, including Dalapon in drinking water, ground water, raw water and other waters at any intermediate stage of treatment. This method is used over the concentration ranges typically found in drinking water ranging from less than 1 to over 50 µg/L. The observed MDL's may vary according to the particular matrix analyzed and the specific instrumentation employed.

In this procedure, a 100 mL sample of water is adjusted to pH 5 and then extracted using an anion exchange cartridge. Method 552.1 suggests the use of polymeric-based anion exchange resins. This application uses silica-based anion exchange, therefore the conditioning steps have been modified to accommodate this silica-based sorbent. Exceptional recoveries are reported as shown in the data below.

Haloacetic Acids

The following RCRA compounds are determined by method 552.1 in water

Analyte	Abbreviation	CASRN
Monochloroacetic acid	MCAA	79-11-8
Dichloroacetic acid	DCAA	79-03-9
Trichloroacetic acid	TCAA	66-03-9
Monobromoacetic acid	MBAA	79-08-3
Bromochloroacetic acid	BCAA	5589-96-8
Dibromoacetic acid	DBAA	631-64-1
Dalapon		75-99-0

Safety

- The toxicity or carcinogenicity of each reagent used in this method has not been defined

Sample Collection, Preservation, and Storage

- Collect samples in all glass amber bottles with PTFE lined screw caps
- Prior to sample collection add ammonium chloride (NH_4Cl) to the bottles to produce a concentration of 100 mg/L
- Alternatively, add 1.0 mL of a 10 mg/mL aqueous solution of NH_4Cl during collection then shake well
- Refrigerate at 4°C or keep on ice and away from light until extraction
- Analyze within 48 hours after collection

Interferences

Care and attention to potential interferences will result in better analysis

- Method interferences may be caused by contamination of solvents, reagents, glassware or other lab apparatus
- Glassware must be scrupulously clean
- Heat all glassware except volumetric flasks in a muffle oven at 400°C. Thermally stable compounds such as PCB's may not be eliminated by this treatment. Thorough rinsing with reagent grade acetone may be substituted for the heating step
- Glassware must be rinsed with 1:9 HCl prior to use
- Reduced analyte recoveries may be observed in high ionic strength matrices particularly waters containing sulfate concentrations above 200 mg/L
- Improved recoveries may be obtained by sample dilution but at the expense of higher MDL's
- Tribromoacetic acid has not been included because of problems associated with stability and chromatography
- Organic acids and phenols, especially chlorinated compounds, are the most direct potential interferences with this method
- For each set of samples analyzed, the reagent blank concentration values exceeding 0.1 µg/L can be subtracted from the sample concentration
- Interfering contamination may occur when a sample with a low concentration of analytes is analyzed following a sample of relatively high concentration. Routine between-sample rinsing of the sample syringe and associated equipment with methyl-*t*-butyl ether (MTBE) can minimize sample cross-contamination

Standard Stock Solution

Prepare a stock solution of each analyte of interest at a concentration of 1-5 mg/mL in MTBE. Method analytes may be obtained as neat materials or ampulized solutions (>99% purity) from commercial suppliers

Prepare Primary Dilution Standard (PPDS)

Prepare PPDS solutions by combining and diluting stock standard solutions with methanol. Prepare stock standard solutions in the 1-2 mg/mL range for all analytes including the surrogate. Aliquots of each stock standard solution (approximately 50-250 µL) are added to 100 mL methanol to yield a primary dilution standard containing the following approximate concentrations of analytes

Analyte	Concentration, µg/mL
Monochloroacetic Acid (MCAA)	3
Monobromoacetic acid (MBAA)	2
Dalapon	2
Dichloroacetic acid (DCAA)	3
2-bromopropionic acid ^b	1
Trichloroacetic acid (TCAA)	1
Bromochloroacetic acid	2
Dibromoacetic acid	1

The primary dilution standards are used to prepare calibration standards of at least three concentration levels (optimally five) of each analyte with the lowest standard being at or near the MDL of each analyte. The concentrations of the other standards should define a range containing the expected sample concentrations or the working range of the detector


Calibration Standards

- Perform by extracting calibration standards in fortified reagent water using the SPE cartridge
- Prepare a five-point calibration curve 100- mL reagent water samples at pH 5 with 20, 50, 100, 250, and 500 µL of the primary dilution standard
- Analyze each calibration standard in triplicate including a reagent water blank
- Include the surrogate 2-bromopropionic acid with the calibration standards in the table

Internal Standard Fortifying Solution

- Prepare a solution of 1,2,3-trichloropropane at 1 mg/mL by adding 36 µL of the neat material to 50 mL of MTBE. From this stock standard solution, prepare a primary dilution standard at 10 mg/L by the addition of 1-100 mL MTBE

Analyte Recovery at Various Spike Concentrations (N=7)



Procedure

1. Cartridge Conditioning

- Assemble a glass block manifold system
- Place **EUQAX156** cartridge(s) in the bulkhead fittings of the vacuum manifold. If nylon fittings are used, Clean-Thru[®] tips are recommended to prevent damage to the fittings by sulfuric acid
- Attach adapters and reservoirs to the cartridges
- Add a 10 mL aliquot of methanol (MeOH) to the cartridge
- Draw slowly through at 2 mL/min flow
- Add 10 mL of reagent water then slowly draw through

Note: After the conditioning steps do not let the cartridge go dry before addition of sample

2. Sample Extraction

- Check the sample pH (preferably with a pH meter) to assure pH 5 ± 0.5
- Add 250 µL of the surrogate primary dilution standard to water sample
- Add the 100 mL sample to the cartridge while drawing through under vacuum at 2 mL/minute

- d) After the sample has passed through the cartridge, draw 10 mL of MeOH through the cartridge to dry the sorbent

3. Cartridge Elution

- a) Disassemble the vacuum manifold and insert a screw cap culture tube
- b) Add 2 mL of 10% H₂SO₄/MeOH to the cartridge
- c) Elute at a rate of 1.5 mL/minute (slow drip)
- d) Repeat with a second aliquot of 2 mL of 10% H₂SO₄/MeOH

4. Derivatization and Solvent Partition

- a) To each of the elutant samples from 3) d) add 2.5 mL of MTBE
- b) Agitate in a vortex mixer for 5 seconds
- c) Place the culture tubes containing the MTBE in a heating block at 50° C for 1 hour
- d) Remove from heating block and add to each tube 10 mL of a 10% solution of sodium sulfate (NaSO₄)
- e) Agitate for 10 seconds on a vortex mixer
- f) Allow the phases to separate for 5 minutes
- g) Transfer the upper MTBE layer to a 15 mL graduated conical centrifuge tube
- h) Add another 1 mL portion of MTBE
- i) Agitate using a vortex mixer, allow to settle then transfer with a Pasteur pipette to the other MTBE collected in step g)
- j) Add 200 µL IS to each extract then additional MTBE to bring the final volume to 5 mL
- k) Transfer a portion of the extract to capped GC vials
- l) Sample is ready for analysis by GC-ECD

5. Sample Analysis—Gas Chromatography

- a) Inject 2 µL using the recommended conditions listed below

GC/ECD Conditions

GC: Agilent 6890N coupled with ECD, equipped with 7683 auto sampler and Chemstation software for data acquisition and analysis. Other instruments may be used

GC capillary column: Restek Rtx®-1701, 30m x 0.25mm x 0.25µm

Injector: 2 µL splitless injection at 200 °C, with a split delay of 0.5 min.

Liner: 4 mm splitless gooseneck, 4mm ID x 6.5mm OD x 78.5mm (UCT **GCLGN4MM**)

Oven temperature program:

- Initial oven temperature of 55 °C
- hold for 5 minutes
- ramp at 7 °C/min to 115 °C
- ramp at 40 °C/min to a final temperature of 280 °C
- hold for 2.3 minutes
- Total run time 20 minutes


Carrier gas: He flow of 1.5 mL/min

ECD temperature: 280 °C

Make up: N₂ at 30 mL/min

Date rate: 20 Hz, save data from 6 to 14 minutes

Chromatogram Showing Elution of Haloacetic Acids



Analyte Elution Order	RSD @ 0.5-1.5 µg/L %	MDL µg/L @ 0.5-1.5 µg/L	RSD @ 2.5-7.5 µg/L %	RSD @ 5.0-15.0 µg/L %
MCAA	4.6	0.22	1.9	2.2
MBAA	3.2	0.10	1.6	2.0
Dalapon	1.3	0.04	1.0	1.0
DCAA	1.5	0.06	1.1	1.5
2-bromopropionic acid (surrogate)	2.6	0.04	0.6	1.7
TCAA	1.9	0.03	2.2	1.0
1,2,3-trichloropropane (IS)				
BCAA	10	0.25	2.4	3.4
DBAA	10	0.12	2.1	7.6

*For complete details on Method 552.1, August 1992, the analyst is referred to: Hodgeson, Jimmie W, Becker, David, (Technology Applications, Inc., "Determination of Haloacetic acids And Dalapon In Drinking Water by Ion-Exchange Liquid-Solid Extraction and Gas Chromatography with an Electron Capture Detector", Environmental Monitoring Systems Laboratory,, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268

**Other cartridge and sorbent masses are available from UCT



Determination of Benzidines and Nitrogen-Containing Pesticides in Water Liquid-Solid Extraction And Reverse Phase High Performance Liquid Chromatography/Particle Beam/Mass Spectrometry*

UCT Product Number:

ECDVB156 (500 mg, DVB, 6 mL cartridge)

ECUNIDVB500 (500 mg, DVB, 83 mL cartridge)

EPA Method 553 Revision 1.1

September 2010

Procedure

1. Cartridge Preparation

- a) Rinse the cartridge(s) with a 10 mL aliquot of methanol
- b) Slowly draw methanol to the top of the frit
- c) Add a second 10 mL aliquot of methanol and draw through to top of frit
- d) Add 10 mL of reagent water to the cartridge and draw through to top of frit

Note: Do not let the cartridge go dry after addition of methanol otherwise repeat at 1) c.

2. Sample Extraction

- a) Adjust sample pH to 7 using either 1N NaOH or HCl
- b) Add water sample to the cartridge and draw through at 20 mL/minute
- c) Rinse the sample container with 10 mL of reagent water and add to the cartridge
- d) Dry sorbent by drawing full vacuum for 10 minutes
- e) Place a clean collection vial in the vacuum manifold

3. Cartridge Elution


- a) Add 15 mL of methanol to the sample container, swirl then add to the cartridge
- b) Elute the cartridge with 2 x 7.5 mL aliquots of methanol dropwise

4. Extract Concentration

- a) Concentrate the extract under a gentle stream of N₂
- b) Bring to a known volume after concentration step
- c) Sample is now ready for analysis

*The analyst should refer to EPA Method 553 "Determination Of Benzidines and Nitrogen-Containing Pesticides in Water By Liquid-Liquid Extraction or Liquid-Solid Extraction and Reverse Phase High Performance Liquid Chromatography/Particle Beam/Mass Spectrometry", Revision 1.1 Issued August 1992, by Thomas D. Behymer, Thomas A. Bellar, James S. Ho, William L. Budde, US EPA, Office of Ground Water and Drinking Water, EPA, National Exposure Research Laboratory, Office of Research and Development, US Environmental Protection Agency, Cincinnati, Ohio 45268

Water Analyses





Bisphenol A Analysis in Water by GC/MS Using an ENVIRO-CLEAN[®] 200 mg C18 Extraction Cartridge

UCTUCT Part Numbers:
EEC1812Z (Endcapped C18 - 200 mg/10 mL)

February 2009

1. Prepare Sample

- a) Using 100 mL of sample water, adjust the pH to 7 or less using 100mM acetic acid
- b) Add internal standard to water sample

Bisphenol A has a pKa value of approximately 9.5

2. Condition Cartridge

- a) Place a cartridge(s) on a multistation vacuum manifold or automated extraction system
- b) Condition the cartridge by adding 3 mL of methanol
- c) Partially draw the methanol through until the surface of the liquid reaches the top of the cartridge frit
- d) Wait 1 minute then add 3 mL of DI water to the cartridge.
- e) Add 1 mL of 100 mM acetic acid
- f) Draw liquid through until it touches the top of the frit
- g) Cartridge is now ready for sample extraction

Note: Do not allow the sorbent to completely dry out after the addition of methanol, otherwise repeat procedure

3. Apply Sample

- a) Add sample to cartridge at a rate of approximately 5 mL/minute by adjusting vacuum

4. Wash Cartridge

- a) Wash by drawing through 5 mL of deionized water
- b) Dry sorbent (5 minutes at > 10 inches Hg)

5. Elute

- Insert a collection vial in the vacuum manifold
- Rinse sample bottle with 3 mL of methanol
- Add the methanol to the cartridge
- Elute at 5 mL per minute
- Add 3 mL of methanol to the cartridge
- Elute at 5 mL per minute

6. Evaporate

- Evaporate methanol eluate using gentle N₂ (< 40°C) to dryness
- Add 50 µL of ethyl acetate to dissolve
- Add 50 µL of reagent MTBSTFA* or BSTFA** to derivatize. Vortex
- Heat mixture for 20-30 minutes @ 70°C
- Cool. Sample is now ready for GC injection

*MTBSTFA-- N-(t-butyldimethylsilyl)-N-methyltrifluoroacetamide

**BSTFA-- N,O -Bis(trimethylsilyl)trifluoroacetamide

7. Instrument Conditions

- Column:** Rtx-5MS, 30m, 0.25 mm ID, 0.50 µm df (5% diphenyl/95% dimethyl polysiloxane)
- Injector Temperature:** 250°C
- Detector Temperature:** 250° C
- Oven Program:** Initial 70°C, ramp @ 20°C/minute to 320°C, hold 3.0 min
- Purge Flow:** Initially Off. On at 0.75 minutes
- Split Flow:** 30 mL/minute
- Inject:** 2 µL


8. Quantitate

MS in EI (+) mode:

BSTFA	Primary ion	Secondary ion	Tertiary ion
BPA	357	373	207
BPA-D16	368	386	217

MTBSTFA	Primary ion	Secondary ion	Tertiary ion
BPA	441	457	207
BPA-D16	452	470	217

Chromatogram Showing Retention Time for Bisphenol A in Water



Spike Concentration: 0.15 µg/L Calc Concentration: 0.157 µg/L Recovery: 105%



Bisphenol A Analysis in Water by LC/MS/MS Using an ENVIRO-CLEAN® 200 mg C18 Extraction Cartridge

UCT Part number:

EEC1812Z (Endcapped C18 - 200 mg/10 mL)

February 2009

Procedure

1. Prepare Sample

- a) Using 100 mL of sample water, adjust the pH to 7 or less using 100mM acetic acid
- b) Add internal standard to water sample

Note: Bisphenol A has a pKa value of approximately 9.5

2. Condition Cartridge

- a) Place a cartridge(s) on a multistation vacuum manifold or automated extraction system
- b) Condition the cartridge by adding 3 mL of methanol
- c) Partially draw the methanol through until the surface of the liquid reaches the top of the cartridge frit
- d) Wait 1 minute then add 3 mL of DI water to the cartridge
- e) Add 1 mL of 100 mM acetic acid
- f) Draw liquid through until it touches the top of the frit
- g) Cartridge is now ready for sample extraction

Note: Do not allow the sorbent to completely dry out after the addition of methanol, otherwise repeat procedure.

3. Apply Sample

- a) Add sample to cartridge at a rate of approximately 5 mL/minute by adjusting vacuum

4. Wash Cartridge

- a) Wash by drawing through 5 mL of deionized water
- b) Dry sorbent (5 minutes at > 10 inches Hg)

5. Elute

- a) Insert a collection vial in the vacuum manifold
- b) Rinse sample bottle with 3 mL of methanol
- c) Add the methanol to the cartridge
- d) Elute at 5 mL per minute
- e) Add 3 mL of methanol to the cartridge
- f) Elute at 5 mL per minute

6. Evaporate

- a) Evaporate methanol eluate using gentle N₂ (< 40°C) to less than 500 µL
- b) Bring sample volume to 500 µL
- c) Sample is ready for injection

7. Quantitate

- a) LC-MS/MS MRM transition (negative ion mode)

Instrumental & Conditions:

Column: 100 x 2.1 (3 µm) Selectra® Phenyl, UCT, INC

Instrument: Applied Biosystems Triplequad LC/MS/MS (other systems may be used)

- Detector: API3200 QTrap
- bisphenol A Precursor ion mass 226.9, product ion mass 109.1
- *bisphenol A standard deuterated A-D16. Precursor ion mass 241.1, Product ion mass 223.1

Inject: 5-10 µL

Mobile Phase: acetonitrile/0.1% formic acid


Flowrate: 0.5 mL/ minute

Flow Program

Time in minutes	% Acetonitrile	% 0.1% Formic Acid
0	30	70
3.0	90	10
3.5	30	70
5.0	30	70

Chromatogram Showing Retention Time for Bisphenol A in Water

Spike concentration: 15 µg/L, Calc Recovery: 15.70 µg/L, Recovery: 105%





Determination of Hormones in Water by Solid-Phase Extraction (SPE) and Liquid Chromatography – Tandem Mass Spectrometry (LC-MS/MS)*

UCT Part Numbers:

EUC18156 (500 mg unendcapped C18, 6 mL)

November 2012

Analytes Determined Using This Method

Analyte	CASRN
Estriol 16 α -Hydroxyestradiol	50-27-1
17β-Estradiol	50-28-2
17α-Ethinylestradiol	57-63-6
Testosterone	58-22-0
Estrone 3-hydroxy-13-methyl-6,7,8,9,11,12,13,14,15,16-decahydrocyclopenta[a]phenanthren-17-one	53-16-7
4-Androstene-3,17-dione	63-05-8
Equilin 3-hydroxyestra-1,3,5,7-tetraen-17-one	474-86-2

Method Summary

Water samples are dechlorinated with sodium thiosulfate and protected from microbial degradation with 2-mercaptopyridine-1-oxide sodium salt during collection. Samples are fortified with surrogates and extracted using C18 cartridges before elution with methanol. The extract is concentrated to dryness with N₂ before adjusting to a 1 mL volume with 50:50 methanol:water. An aliquot is injected into an LC equipped with a C18 column interfaced to a MS/MS. Detection limits of 0.02-0.37 ng/L can be obtained using this application from samples fortified at 0.25-0.875 ng/mL.

Safety

- The toxicity and carcinogenicity of each reagent has not been defined
- Each chemical should be treated as a potential health hazard and exposure minimized by the use of PPE (personal protective equipment) such as gloves, respirators and other personal safety equipment
- Pure standard materials and stock standard solutions of the method compounds should be handled with suitable protection
- Ammonium hydroxide, used during method development as a pH modifier for the HPLC mobile phase, should be fresh* and handled in a fume hood

*NH₄OH loses its strength after opening due to evaporation. A fresh bottle should be used

Sample Collection, Preservation, and Storage

- Use one-liter amber glass bottles with PTFE-lined screw caps
- If smaller sample sizes are used, adjust the amount of preservatives and surrogate/analyte fortification levels according to the sample size
- Fill sample bottles taking care not to flush out the preservatives. Analytes are not volatile so bottles need not be headspace-free
- When sampling from a cold water tap, remove the aerator and allow the system to flush until the water temperature has stabilized. Invert the bottles several times to mix the sample with the preservation reagents
- Samples should be chilled during shipment and should not exceed 10° C during the first 48 hours after collection
- In the laboratory, store samples at or below 6° C and protect from light until analysis. **Do not freeze unprocessed samples**
- All compounds listed in the method have adequate stability for 28 days when collected, preserved, shipped and stored as described
- Samples should be extracted as soon as possible for best results
- After sample preparation, extracted samples should be stored at 0° C or less and analyzed within 28 days after extraction

Sample Preservatives

Compound	Amount	Purpose
Sodium thiosulfate	80 mg/L	Removes free chlorine
2-mercaptopyridine-1-oxide sodium salt	65 mg/L	Microbial inhibitor

Note: Preservation reagents listed above should be added to each sample bottle prior to shipment to the field or prior to sample collection

Interferences

- Matrix interferences may be caused by contaminants that are co-extracted from the sample and will vary from source to source
- Humic and/or fulvic material in environmental samples may be co-extracted during SPE and can cause enhancement and/or suppression in the electrospray ionization source. Total organic carbon (TOC) is an indicator of the humic content of a sample
- Use only high purity analogs. Depending on the source and purity, labeled analogs used as internal standards may contain a small percentage of the corresponding native analyte which may be significant when attempting to determine LCMRLs and DLs
- Nitrile gloves should be worn at all times. Handling clean glassware may be a potential source of interference
- It may be appropriate to include a Field Blank with the sampling bottles depending on the sampling site. The Field Blank is analyzed along with the samples to ensure that no human hormones were introduced into the samples during the collection and handling process

Internal Standards

Internal Standard	CASRN	Neat Material Cat #	Solution Standard Cat #
16α-Hydroxyestradiol-<i>d</i>2 (Estriol-<i>d</i>2)	53866-32-3	C/D/N Isotopes Cat. No. D-5279	N/A
¹³C6-Estradiol	None	None	Cambridge Isotope Labs 100 μ g/mL in Methanol Cat. No. CLM-7936-1.2
¹³C2-Ethynylestradiol	None	None	Cambridge Isotope Labs 100 μ g/mL in Acetonitrile Cat. No. CLM-3375-1.2
Testosterone-<i>d</i>3	77546-39-5	None	Sigma Drug Std., 100 μ g/mL in dimethoxyethane Cat. No. T5536

Internal Standard Stock Standards

ISSSS 500 μ g/mL - Weigh 5 mg of α -hydroxyestradiol-*d*2 (estriol-*d*2) into a tared 10 mL volumetric flask and dilute to volume with methanol. The remaining internal standards can be purchased as 100 μ g/mL solutions

Internal Standard Primary Dilution Standard (IS-PDS)

IS PDS 1.0 – 4.0 µg/mL: The table below can be used as a guide for preparing the IS PDS. The IS PDS is prepared in acetonitrile and is stable for about six months if stored at a temperature < 6° C. Use 5 µL of the 1.0 – 4.0 µg/mL IS PDS to fortify the final 1 mL extracts. This will yield a final concentration of 5.0 – 20 ng/mL of each IS.

Internal Standard Concentrations

Internal Standard	Conc IS Stock µg/mL	Volume of IS Stock, µL	Final Volume of IS PDS,	Final Conc. of IS PDS (µg/mL)
16α-Hydroxyestradiol-<i>d</i>2 (Estriol-<i>d</i>2)	500	40	10	2.0
¹³C₆-Estradiol	100	400	10	4.0
¹³C₂- Ethinylestradiol	100	400	10	4.0
Testosterone-<i>d</i>3	100	100	10	1.0

Internal Standard Fortification Standard

5 µL of the IS PDS is difficult to accurately and reproducibly transfer into HPLC vials. Therefore, a 1-in-10 dilution of the IS PDS is prepared by transferring 1 mL of the 1.0 – 4.0 µg/mL IS PDS into a 10 mL volumetric flasks and diluting to volume with methanol. This results in a 0.1-0.4 µg/mL IS solution

50 µL of the prepared ISPDS is added to the samples prior to analysis

Surrogate Analytes

Two isotopically labeled surrogates are listed below. Select the surrogate that performs best under the LC-MS/MS conditions employed for the analysis

For this application, Bisphenol A-*d*16 was chosen as surrogate analyte

Surrogate Analytes

Surrogate Analyte	CASRN	Neat Materials Catalog No.
Ethynylestradiol- <i>d</i> ₄	350820-06-3	C/D/N Isotopes, Cat. No. D-4319
Bisphenol A- <i>d</i> ₁₆	96210-87-6	Sigma, Cat. No. 451835

Surrogate Stock Standards

1000 µg/mL - Prepare individual solutions of the surrogate standards by weighing 10 mg of the solid material into tared 10 mL volumetric flasks and diluting to volume with methanol

Surrogate Analyte Primary Dilution Standard

SUR PDS 2.5 – 7.0 µg/mL

Use the table below as a guide for preparation of the **SUR PDS** in methanol.

Use 10 µL of **SUR PDS** to fortify 1-L samples yielding a final concentration of 70 ng/mL ethynylestradiol-*d*₄ or 25 ng/mL bisphenol A-*d*₁₆ in the 1 mL extracts

Surrogate Analyte Solutions

Surrogate Analyte	Conc of SUR Stock (µg/mL)	Volume of SUR Stock (µg/mL)	Final Volume of SUR PDS (mL)	Final Conc. Of SUR PDS (µ g/mL)
Ethynylestradiol- <i>d</i> ₄	1000	70	10	7.0
Bisphenol A- <i>d</i> ₁₆	1000	25	10	2.5

Surrogate Fortification Standard (SUR PDS)

10 µL of the **SUR PDS** is difficult to accurately and reproducibly transfer into samples. Therefore, a 1-in-10 dilution of the **SUR PDS** is prepared by transferring 1 mL of the 2.5 µg/mL IS PDS into a 10 mL volumetric flask and diluting to volume with methanol, resulting in a 0.25 µg/mL IS solution. Use 100 µL of this fortification standard to fortify 1-L samples yielding a final concentration of 25 ng/mL bisphenol A-*d*₁₆ in the 1 mL extracts

Analyte Stock Standard Solution

Each concentration equals 1000 µg/mL - Obtain the analytes listed in Table 1 above as standard solutions or as neat materials. Prepare stock standards individually by weighing 10 mg of the solid standards into tared 10 mL volumetric flasks and diluting to volume with methanol

Analyte Primary Dilution Standard

Prepare concentrations of the stock standard solutions between 1.0 – 3.5 µg/mL – Dilute the Analyte Stock Standard solutions into 50% methanol in reagent water. The concentrations vary based on the instrumental sensitivity. The Analyte PDS is used to prepare calibration standards, and to fortify LFBs, LFSMs, and LFSMDs with the method analyte

Analyte Stock	Stock Concentration µg/mL	Stock Volume µL	Final Volume (mL, 50% MeOH)	Analyte PDS Concentration µg/L
Estriol, 16α- Hydroxyestradiol	100 0	20	10	2.0
Estrone	100 0	20		2.0
17β-Estradiol	100 0	25		2.5
17α-Ethynylestradiol	100 0	35		3.5
Equilin	100 0	20		2.0
Androstenedione	100 0	10		1.0
Testosterone	100 0	10		1.0

Analyte Fortification Standard

10 µL of the Analyte PDS is difficult to accurately and reproducibly pipette into samples. Therefore, a 1-in-10 dilution of the Analyte PDS is prepared by transferring 1 mL of the 1.0 – 3.5 µg/mL Analyte PDS into a 10 mL volumetric flask and diluting to volume with 50% methanol, resulting in a 0.1 – 0.35 µg/mL solution. Use 100 µL of this fortification standard to fortify 1-L samples yielding a final concentration of 10 – 35 ng/mL in the 1 mL extracts.

An additional 1-in-10 dilution of the 0.1 – 0.35 µg/mL fortification solution is prepared in a 10 mL volumetric flask using 50% methanol. This yields a 0.01 – 0.035 µg/mL solution that is used to fortify 1-L samples yielding a final concentration of 1 – 3.5 ng/mL in the 1 mL extracts.

Procedure

1. Sample Preparation

- Add a 100-µL aliquot of SUR Fortification Standard to each 1-L sample for a final concentration of 25 ng/L bisphenol A-*d*16
- Fortify LFBs, LFSMs, or LFSMDs with an appropriate volume of Analyte

Fortification Standard

- c) Cap and invert each sample several times to mix

2. Cartridge Preparation

- a) Assemble a vacuum manifold. Automated extraction equipment may also be used
- b) Place **EUC18156** cartridge(s) on the manifold
- c) Add a 10 mL aliquot of methanol to the cartridge and draw through the cartridge until dry
- d) Add another 5 mL aliquot of methanol and draw through the cartridge until dry

Note: Do Not Let the Cartridge Go Dry after Starting the Following Steps

- e) Add approximately 10 mL of methanol to each cartridge
- f) Draw about 1 mL of solvent through the cartridge and turn off the vacuum temporarily
- g) Let the cartridge soak for about one minute then draw most of the remaining solvent through the cartridge leaving a thin layer of methanol on the surface of the cartridge
- h) Add 10 mL of reagent water to each cartridge and draw through leaving a thin layer of liquid on the surface of the cartridge
- i) Add another 10 mL aliquot of reagent water
- j) Draw the water through the cartridge keeping the water level above the cartridge surface
- k) Turn off the vacuum

3. Sample Extraction

- a) Add the sample to the extraction reservoir containing the conditioned cartridge and turn on the vacuum (approximately 10 to 15 in. Hg). Flow of sample through the cartridge should be a fast drip. Adjust vacuum if necessary
- b) Do not let the cartridge go dry before the entire sample volume is extracted
- c) After the entire sample has been drawn through the cartridge, add a 10 mL aliquot of 15% methanol to the sample container and wash the cartridge with the rinsate from the container
- d) Using full vacuum, draw air through the cartridge by maintaining full vacuum for 10–15 minutes

- e) After drying, turn off and release the vacuum

4. **Sorbent Elution**

- a) Insert collection tubes into the manifold to collect the cartridge extracts. The collection tube should fit around the drip tip of the base to ensure collection of all the eluent
- b) Add 5 mL of methanol to the cartridge and draw the methanol into the cartridge to soak the sorbent
- c) Allow the cartridge to soak for about one minute
- d) Using vacuum, draw the remaining methanol slowly through the cartridge into the collection tube
- e) Elute with an additional 2 x 5 mL aliquots of methanol

Note: The methanol can be eluted directly into a 15 mL tube (or larger). Otherwise 5 mL MeOH can be eluted into a 5 mL culture tube and evaporated to near dryness prior to adding the second 5 mL eluate. This process is repeated until all 15 mL MeOH has been added to the 5 mL culture tube.

5. **Extract Concentration**

- a) Concentrate the extract to dryness under a gentle stream of N₂ in a warm water bath (~45° C)
- b) Rinse the collection tube with 950 µL of 50% methanol, vortex for 2 min and transfer the rinse into a HPLC vial
- c) Add 50 µL IS Fortification Standard (ISFS) and vortex for an additional 1 min

6. **Sample Filtration**

- a) It is highly recommended that extracts be filtered with at least a 0.45 micron syringe filter prior to analysis to remove the particulates in a sample.
- b) If filtering is incorporated as part of the sample preparation, the first lot of syringe filters should be included in the procedure to document the potential interferences that are introduced or analytes are retained on the filter.
- c) Subsequent lots of syringe filters can be verified by examining CAL standards.

Sample Analysis

HPLC Conditions

HPLC instrument: Thermo Scientific Dionex UltiMate 3000 System**		
Column: Thermo Scientific Accucore C18, 100 x 2.1 mm, 2.6 µm with 10 mm guard column		
Column Temperature: 35° C		
Column Flow Rate: 0.200 mL/min		
Injection Volume: 20 µL		
Gradient		
Time (min)	Moblie phase A H ₂ O + 0.02% NH ₄ OH	Moblie phase B MeOH +0.02% NH ₄ OH
0	70	30
1	35	65
9	35	65
9.1	15	85
11	15	85
11.1	70	30
15	70	30

ESI-MS/MS Method Conditions

MS Parameters	
MS instrument	Thermo Scientific TSQ Vantage**
Polarity	HESI ⁺ & HESI ⁻
Spray Voltage V	+4500 / -3500 V
Vaporizer Temperature	350°C
Ion Transfer Capillary	300 °C
Sheath Gas Pressure	45 arbitrary units
Auxiliary Gas Pressure	40 arbitrary units
Q1 and Q3 Peak Width (FWHM)	0.4 and 0.7 Da
Collision Gas and Pressure	Ar at 1.5 mTorr
Scan Type	SRM
Cycle Time	0.75 Sec

** Alternative LC-MS/MS systems may be used

LC-ESI-MS/MS Analyte Retention Times, Precursor and Product Ions, S-lens, and Collision Energy

Analyte	Ret. Time (min.)	ESI Mode	Precursor Ion	Product Ion	S-lens V	Collision Energy eV	Internal Standard
Estriol	4.99	ESI-	286.75	144.96	100	39	Estriol-d2
Bisphenol A-d ₁₆	6.20	ESI-	240.91	223.02	60	18	¹³ C ₆ -Estradiol
Equilin	7.03	ESI-	266.80	142.97	67	37	¹³ C ₆ -Estradiol
17β-Estradiol	7.14	ESI-	270.93	144.98	67	40	¹³ C ₆ -Estradiol
Androstenedione	7.18	ESI+	286.87	96.88	76	20	Testosterone-d3
17α-Ethynylestradiol	7.19	ESI-	294.79	244.97	80	39	¹³ C ₂ -Ethynylestradiol
Estrone	7.24	ESI-	268.78	144.97	109	39	¹³ C ₂ -Ethynylestradiol
Testosterone	7.76	ESI+	288.96	96.88	76	20	Testosterone-d3

LC-ESI-MS/MS Internal Standard Retention Times, Precursor and Product Ions, S-lens, and Collision Energy

Internal Standard	Ret. Time (m)	Precursor Ion	Product Ion	S-Lens V	Collision Energy eV
Estriol-d2	5.00	288.75	146.96	100	39
¹³ C ₆ -Estradiol	7.13	276.93	146.98	67	40
¹³ C ₂ -Ethynylestradiol	7.16	296.80	144.97	80	39
Testosterone-d3	7.73	291.96	96.88	70	20

Analyte Recovery

n = 5	Estriol	Equilin	Estradiol	Androstenedione	Testosterone	Ethinylestradiol	Estrone
Fortified Conc (ng/mL)	20	20	25	10	20	10	35
Mean	97.61	82.81	96.98	92.88	90.76	96.49	93.95
SD	6.45	5.90	5.17	2.55	3.47	1.31	2.02
RSD	6.60	7.13	5.33	2.75	3.82	1.36	2.15

n = 5	Estriol	Equilin	Estradiol	Androstenedione	Testosterone	Ethinylestradiol	Estrone
Fortified Conc (ng/mL)	2	2	2.5	1	1	3.5	2
Mean	83.05	96.34	93.68	98.31	104.63	109.45	106.10
SD	4.02	8.92	7.02	3.90	4.95	6.96	6.19
RSD	4.84	9.26	7.49	3.96	4.73	6.36	5.84

*Based on EPA Method 539, Version 1.0, November 2010, Glynda A. Smith (U.S. EPA, Office of Ground Water and Drinking Water) Alan D. Zaffiro, (Shaw Environmental, Inc.) M. L. Zimmerman (Shaw Environmental, Inc.) D. J. Munch (U.S. EPA, Office Of Ground Water And Drinking Water), Technical Support Center , Standards And Risk Management Division, Office Of Ground Water And Drinking Water, U. S. Environmental Protection Agency, Cincinnati, Ohio 45268



Determination of Diesel Range Non-halogenated Organics Using GC/FID by Method 8015D*

UCT Product Number:

ECUNIPAH (2000 mg unendcapped C18, 83 mL cartridge)

February 2011

Method Summary

Method 8015D may be used to determine the concentrations of several nonhalogenated volatile organic compounds and semivolatile organic compounds using GC with flame ionization detection.

- Diesel Range Organics (**DRO**) corresponds to the range of alkanes from C₁₀ to C₂₈ and covering a boiling point range of approximately 170° C - 430° C

Sample Collection, Preservation, and Storage

- See sample collection options EPA Method 5035A

Interferences

- Solvents, reagents, glassware, and other sample processing hardware may yield artifacts and/or interferences to sample analysis
- All these materials must be demonstrated to be free from interferences under the conditions of the analysis by analyzing method blanks
- Specific selection of reagents and purification of solvents by distillation in all-glass systems may be necessary
- Contamination by carryover occurs whenever high-concentration and low-concentration samples are analyzed in sequence. To reduce the potential for carryover, the sample syringe or purging device must be rinsed out between samples with an appropriate solvent
- All glassware must be scrupulously cleaned. Clean all glassware as soon as possible after use by rinsing with the last solvent used followed by detergent washing with hot water, and rinses with tap water and organic-free reagent water
- Drain the glassware and dry it in an oven (except volumetric glassware) at 130° C for several hours or rinse it with methanol
- FID is a non-selective detector. There is a potential for many non-target compounds present in samples to interfere with this analysis and potential for analytes to be resolved poorly, especially in samples that contain many analytes.

Reagents and Standards

- Reagent grade chemicals must be used in all tests
- Organic-free reagent water
- Methanol, CH₃OH - Pesticide quality or equivalent
- Methylene Chloride – unpreserved, pesticide residue quality or equivalent
- Fuels, e.g., diesel - Purchase from a commercial source
- Alkane standard - Standard contains a homologous series of *n*-alkanes for establishing retention times (e.g., C₁₀-C₂₈ for diesel)

Stock Standards

- Prepared from pure standard materials or purchased as certified solution
- Standards must be replaced after 6 months

Secondary Dilution Standards

- Using stock standard solutions, prepare secondary dilution standards, as needed either singly or mixed together
- Secondary dilution standards should be stored with minimal headspace for volatiles
- Check frequently for signs of degradation or evaporation, especially just prior to preparing calibration standards

Recommended GC Columns

The choice of GC column will depend on the analytes of interest, expected concentrations, and intended use of the results. The capillary columns are necessary for petroleum hydrocarbon analyses and are recommended for all other analyses. Other columns may be employed if the analyst can demonstrate acceptable performance.

- Establish the GC operating conditions appropriate for the GC column being utilized and the target analytes specified in the project plan
- Optimize the instrumental conditions for resolution of the target analytes and sensitivity
- Suggested operating conditions and GC programs are given below for the recommended columns:

(DRO) 30-m x 0.53-mm ID fused-silica capillary column chemically bonded with 5% methyl silicone (DB-5, SPB-5, RTx, or equivalent), 1.5-μm film thickness

- Carrier gas (He) flow rate: 5-7 mL/min
- Makeup gas (He) flow rate: 30 mL/min
- Injector temperature: 200° C
- Detector temperature: 340° C
- Temperature program: Initial temperature: 45° C, hold for 3 minutes
- Program: 45° C to 275° C, at 12° C/min
- Final temperature: 275° C, hold for 12 minutes

Procedure**

1. Condition Cartridge

- a) Place a **ECUNIPAH** cartridge(s) on the vacuum manifold*
- b) With vacuum off add 10 mL of methylene chloride to the cartridge
- c) Let it soak for 1 minute
- d) Turn on vacuum and draw through to waste
- e) Draw vacuum through the cartridge to remove all methylene chloride
- f) Add 10 mL of methanol to the cartridge
- g) Let it soak for 1 minute
- h) Draw the methanol to the level of the frit
- i) Add 10 mL of deionized water to the cartridge
- j) Let it soak for 1 minute
- k) Draw most of the water to waste but do not allow the sorbent to dry

Note: Do not let the cartridge go completely dry after addition of methanol otherwise repeat starting at step 1.f)

2. Sample Addition

- a) Add required surrogates to the sample
- b) Shake
- c) Adjust the pH of the sample to 2 or less using 5 mL of 1:1 HCl
- d) Shake
- e) Add the sample to the cartridge under vacuum. Draw the sample through the cartridge no faster than 20 – 30 minutes per liter
- f) Allow the cartridge to dry under full vacuum for 5 minutes***

3. Extract Elution


- a) Place a collection tube or vial containing approximately 1 cm of pre-baked sodium sulfate anhydrous in the vacuum manifold
- b) Rinse sample bottle with 10 mL of methylene chloride to remove any analyte from the glass
- c) Add the methylene chloride rinse to the cartridge

- d) Allow to soak for 2 minutes then draw through
- e) Repeat this procedure three more times using 10 mL aliquots of methylene chloride
- f) Dry the extract by passing it through 5 grams of **ECSS25K** anhydrous sodium sulfate
- g) Thoroughly rinse the collection device with methylene chloride and add this solvent to the sodium sulfate


4. Concentration

- a) Concentration of final extract using a micro Kaderna-Danish flask equipped with a 3-ball Snyder column is recommended to avoid loss of low boiling alkanes. A micro-Snyder column can be used to evaporate the extract to a 1 mL final volume
- b) Extracts should not be stored with headspace for extended periods of time as low recoveries will result

Example Chromatogram of a 30 PPM Diesel Standard



**Example Chromatogram of a 30 Ppm Diesel Standard With The
Baseline Projected Between C₁₀ and C₂₈**



*Summarized from Method 8015D Nonhalogenated Organics Using GC/FID, Revision 4, June 2003, and associated appropriate methods

**(steps 1- 4 added to this method summary for use with UCT PAH/DRO cartridges)

***Faster drying results can be obtained by removing the cartridge during drying and shaking or tapping the excess moisture from the bottom of the cartridge. Drying times are approximate. Do not over dry as low recoveries may result



Polycyclic Aromatic Hydrocarbons (PAH) from a Water Matrix

UCT Product Number:

ECUNIPAH (2000 mg, unendcapped C18, 83 mL cartridge)

ECSS25K (anhydrous sodium sulfate)

January 2011

Polynuclear Aromatic Compounds Recovered in this Method

PAH	CASRN	PAH	CASRN
Acenaphthene	83-32-9	Chrysene	218-01-9
Acenaphthylene	208-96-8	Dibenzo(a,h)anthracene	53-70-3
Anthracene	120-12-7	Fluoranthene	206-44-0
Benzo(a)anthracene	56-55-3	Fluorene	86-73-7
Benzo(a)pyrene	50-32-8	Indeno(1,2,3-cd)pyrene	193-39-5
Benzo(b)fluoranthene	205-99-2	Naphthalene	91-20-3
Benzo(g,h,i)perylene	191-24-2	Phenanthrene	85-01-8
Benzo(k)fluoranthene	207-08-9	Pyrene	129-00-0

1. Condition Cartridge

- Place **ECUNIPAH** cartridge(s) on the vacuum manifold*
- With vacuum turn off add 10 mL of methylene chloride to the cartridge
- Let it soak for 1 minute
- Turn on vacuum and draw through to waste
- Add 10 mL of methanol to the cartridge
- Let it soak for 1 minute
- Draw the methanol to the level of the frit
- Add 10 mL of deionized water to the cartridge
- Let it soak for 1 minute
- Draw most of the water to waste but do not allow the sorbent to dry

Note: Do not let the cartridge go dry after addition of methanol otherwise repeat starting at step 1. e)

2. Sample Addition

- a) Add surrogates to the sample
- b) Shake
- c) Add the sample to the cartridge under vacuum. Draw the sample through the cartridge in 20 – 30 minutes (10 – 50 ml/min)
- d) Allow the cartridge to dry under full vacuum for 10 minutes**

Optional: Before proceeding to the drying step use UCT Zero-Blank™ Filter (**ECBLANK**) to reduce potential background contamination

3. Extract Elution

- a) Place a collection tube or vial in the vacuum manifold
- b) Rinse sample bottle with 10 mL of methylene chloride
- c) Add the methylene chloride rinse to the cartridge
- d) Allow to soak for 1 minute then draw through
- e) Repeat this procedure 3 more times using 10 mL aliquots of methylene chloride
- f) Dry the extract by passing it through 5 grams of **ECSS25K** anhydrous sodium sulfate
- g) Thoroughly rinse the collection device with methylene chloride
- h) Add the methylene chloride rinse to the sodium sulfate

4. Concentration and Analysis

- a) Using a standard analytical evaporator with gentle N₂ flow and low temperature (40° C) carefully concentrate the extract to a final volume for GC/MS analysis.
- b) Solvent exchange into acetonitrile and bring to a 1 ml final volume for HPLC analysis
- c) Sample is now ready for analysis

Note: Most analysis errors are caused by poor concentration technique. Do not concentrate below 0.5 mL as low recoveries may result

*The ENVIRO-CLEAN® Universal PAH cartridge can be used on standard vacuum manifolds (# VMFF016GL) and standard disk manifolds (#ECUCTVAC6) (with adapter part # ECUCTADP). The cartridge is designed to fit all automated extraction systems

**Faster drying results can be obtained by removing the cartridge during drying and shaking or tapping the excess moisture from the bottom of the cartridge. Drying times are approximate. Do not over dry as low recoveries may result



Organochlorine Pesticides and Polychlorinated Biphenyls by Solid-Phase Extraction

UCT Product Number:

EEC181M6 (1000 mg C18, 6 mL Cartridge) or

ECUNIC18 (1100mg C18, Universal Cartridge)

July 2011

Compounds Recovered Using This Method

Analyte	CASRN
1,2-Dibromo-3-chloropropane (DBCP)	96-12-8
4,4'-DDD	72-54-8
4,4'-DDE	72-55-9
4,4'-DDT	50-29-3
Alachlor	15972-60-8
Aldrin	309-00-2
Captafol	2425-06-1
Carbophenthion	786-19-6
Chlordane - not otherwise specified (n.o.s.)	57-74-9
Chlorobenzilate	510-15-6
Chloroneb	2675-77-6
Chloropropylate	5836-10-2
Chlorothalonil	1897-45-6
Dacthal (DCPA)	1861-32-1
Diallate	2303-16-4
Dichlone	117-80-6
Dichloran	99-30-9
Dicofol	115-32-2
Dieldrin	60-57-1
Endosulfan I	959-98-8
Endosulfan II	33213-65-9
Endosulfan sulfate	1031-07-8
Endrin	72-20-8
Endrin Aldehyde	7421-93-4
Endrin Ketone	53494-70-5
Etridiazole	2593-15-9
Halowax-1000	58718-66-4
Halowax-1001	58718-67-5
Halowax-1013	12616-35-2
Halowax-1014	12616-36-3
Halowax-1051	2234-13-1
Halowax-1099	39450-05-0
Heptachlor	76-44-8
Heptachlor epoxide	1024-57-3
Hexachlorobenzene	118-74-1
Hexachlorocyclopentadiene	77-47-4
Isodrin	465-73-6
Methoxychlor	72-43-5
Mirex	2385-85-5

Nitrofen	1836-75-5
Pentachloronitrobenzene (PCNB)	82-68-8
Permethrin (<i>cis</i> + <i>trans</i>)	52645-53-1
Perthane	72.56-0
Propachlor	1918-16-7
Strobane	8001-50-1
Toxaphene	8001-35-2
Trans-Nonachlor	39765-80-5
Trifluralin	1582-09-8
α -BHC	319-84-6
α -chlordane	5103-71-9
β -BHC	319-85-7
γ -BHC (Lindane)	58-89-9
γ -chlordane	5103-74-2
δ -BHC	319-86-8

PCB's Recovered Using This Method

Compound	CASRN	IUPAC #
Aroclor 1016	12674-11-2	-
Aroclor 1221	11104-28-2	-
Aroclor 1232	11141-16-5	-
Aroclor 1242	53469-21-9	-
Aroclor 1248	12672-29-6	-
Aroclor 1254	11097-69-1	-
Aroclor 1260	11096-82-5	-
2-Chlorobiphenyl	2051-60-7	1
2,3-Dichlorobiphenyl	16605-91-7	5
2,2',5-Trichlorobiphenyl	37680-65-2	18
2,4',5-Trichlorobiphenyl	16606-02-3	31
2,2',3,5'-Tetrachlorobiphenyl	41464-39-5	44
2,2',5,5'-Tetrachlorobiphenyl	35693-99-3	52
2,3',4,4'-Tetrachlorobiphenyl	32598-10-0	66
2,2',3,4,5'-Pentachlorobiphenyl	38380-02-8	87
2,2',4,5,5'-Pentachlorobiphenyl	37680-73-2	101
2,3,3',4',6-Pentachlorobiphenyl	38380-03-9	110
2,2',3,4,4',5'-Hexachlorobiphenyl	35065-28-2	138
2,2',3,4,5,5'-Hexachlorobiphenyl	52712-04-6	141
2,2',3,5,5',6-Hexachlorobiphenyl	52663-63-5	151
2,2',4,4',5,5'-Hexachlorobiphenyl	35065-27-1	153
2,2',3,3',4,4',5-Heptachlorobiphenyl	35065-30-6	170
2,2',3,4,4',5,5'-Heptachlorobiphenyl	35065-29-3	180
2,2',3,4,4',5',6-Heptachlorobiphenyl	52663-69-1	183
2,2',3,4',5,5',6-Heptachlorobiphenyl	52663-68-0	187
2,2',3,3',4,4',5,5',6-Nonachlorobiphenyl	40186-72-9	206

Procedure

1. Condition Cartridge

- a) Assemble a suitable vacuum manifold system
- b) Place cartridge(s) in the bulkhead fittings or cartridge adapters of the vacuum manifold. If nylon fittings are used, UCT CLEAN-THRU[®] Tips (CLTTP050) are recommended to prevent damage to the fittings by sulfuric acid
- c) Attach adapters and reservoirs to the cartridge(s) if necessary
- d) Rinse cartridge with 10 mL of methylene chloride (MeCl_2)
- e) Let the MeCl_2 soak for 1.5 min.
- f) Using vacuum draw the MeCl_2 to waste
- g) Add 10 mL of acetone. Let the acetone soak for 1.5 minutes
- h) Draw the acetone to waste
- i) Dry the cartridge using full vacuum for a few seconds
- j) Add 10 mL of methanol and allow the methanol to soak for 1.5 min.

Note: Do not allow the cartridge to go dry otherwise repeat starting with step 1) j)

- k) Draw some of the methanol through leaving a layer just covering the frit
- l) Add 20 mL of deionized water. Draw most of the water through to waste but do not allow the sorbent to completely dry

2. Sample Addition

- a) Adjust sample pH to ≤ 2 using 1:1 sulfuric acid

Note: Use a pH meter to accurately determine sample pH in deionized water. Do not use pH paper

- b) Mix thoroughly
- c) Start vacuum and add the sample. Draw sample through the cartridge at a rate ≤ 50 mL/minute (1 L should pass through in 20 minutes or longer)
- d) Allow the cartridge to dry under full vacuum for 10 minutes**

3. Extract Elution

- a) Place a collection tube or vial under the cartridge
- b) Add 5 mL of acetone to the sample bottle and swirl to remove any residue
- c) Add the acetone to the cartridge. Allow the solvent to soak for 1 minute then draw into collection vial
- d) Repeat this procedure 3 more times using a 10 mL portions of MeCl_2
- e) Prepare a 5 gram bed of anhydrous sodium sulfate anhydrous in a funnel using glass wool
- f) Dry the extract by passing it through the funnel of sodium sulfate anhydrous
- g) Carefully rinse the collection vial with MeCl_2 , then add to the sodium sulfate, rinsing the sodium sulfate and collect

4. Concentration and Analysis

- a) Carefully concentrate the extract. Solvent exchange if necessary

Note: Most extraction errors are caused by poor concentration technique

K-D Concentration Technique

Sample extracts may be concentrated to the final volume necessary by using the K-D technique or nitrogen evaporation.

- a) Assemble a Kuderna-Danish (K-D) concentrator by attaching a 10 mL concentrator tube to an evaporation flask
- b) Collect the dried extract in the K-D concentrator
- c) Rinse the collection tube and drying funnel then quantitatively transfer into the K-D flask with an additional 20 mL portion of solvent
- d) Add boiling chips to the flask then attach a three-ball Snyder column
- e) Attach the solvent vapor recovery glassware (condenser and collection device to the Snyder column of the K-D apparatus)
- f) Pre-wet the Snyder column by adding about 1 mL of methylene chloride or acetone
- g) Place the K-D apparatus on a hot water bath (15 - 20°C) above the boiling point of the solvent) so that the concentrator tube is partially immersed in the hot water

- h) Adjust the vertical position of the apparatus and the water temperature as necessary to complete the concentration in 10 - 20 min. At the proper rate of distillation the boiling chips of the column will actively chatter, but the chambers should not flood
- i) When the apparent volume of liquid reaches 1 mL, remove the K-D apparatus from the water bath and allow it to drain and cool for at least 10 min.
- j) If a solvent exchange is needed quickly remove the Snyder column, add 50 mL of the exchange solvent and a new boiling chip
- k) Reattach the Snyder column. Concentrate the extract increasing the temperature of the water bath to maintain a proper distillation rate
- l) Remove the Snyder column. Rinse the K-D flask and the lower joints of the Snyder column into the concentrator tube with 1 - 2 mL of solvent
- m) Adjusted to a final volume of 5.0 - 10.0 mL using an appropriate solvent

Note: If further concentration is necessary, use either the micro-Snyder column technique or a N₂ evaporation technique described below

Micro-Snyder Column Technique

- Add a fresh clean boiling chip to the concentrator tube and attach a two-ball micro-Snyder column directly to the concentrator tube
- Attach the solvent vapor recovery glassware (condenser and collection device) to the micro-Snyder column of the K-D apparatus, following the manufacturer's instructions
- Pre-wet the Snyder column by adding 0.5 mL of methylene chloride or the exchange solvent
- Place the micro-concentration apparatus in a hot water bath so that the concentrator tube is partially immersed in the hot water
- Adjust the vertical position of the apparatus and water temperature, as necessary, to complete the concentration in 5 - 10 min. At the proper rate of distillation the balls of the column will actively chatter, but the chambers will not flood
- When the apparent volume of liquid reaches 0.5 mL, remove the apparatus from the water bath and allow it to drain and cool for at least 10 min.
- Remove the Snyder column and rinse its lower joints into the concentrator tube with 0.2 mL of solvent
- Adjust the final extract volume to 1.0 - 2.0 mL

Nitrogen Evaporation Technique

- Place the concentrator tube in a warm bath (30° C) and evaporate the solvent to 0.5 mL using a gentle stream of clean, dry N₂ (filtered through a column of activated carbon)

Note: New plastic tubing must not be used between the carbon trap and the sample as phthalate interferences may be introduced

- Rinse down the internal wall of the concentrator tube several times with solvent during the concentration
- Position the concentrator tube to avoid condensing water into the extract
- Do not allow extract to become dry. If the volume of solvent is reduced below 1 mL, some analytes may be lost
- The extract may now be cleaned up or analyzed for analytes using the appropriate technique(s)
- If the sample is not analyzed immediately, cover the concentrator tube and store in a refrigerator. If the extract will be stored longer than 2 days, transfer to a vial with a PTFE-lined screw-cap and store in a refrigerator

****Faster drying results can be obtained by removing the cartridge during drying and shaking or tapping the excess moisture from the bottom of the cartridge. Drying times are approximate. Do not over dry. Low recoveries could result**

*For complete details on Method 8081B "Organochlorine pesticides by Gas Chromatography/Mass Spectrometry," December 1996, and 8082A, Polychlorinated Biphenyls by Gas Chromatography," the analyst is referred to Environmental Monitoring Systems Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268



Solid-Phase Extraction of Pesticides in Water using Graphitized Carbon Black (GCB)

UCT Part Numbers:

EUCARB1M6 (1000 mg GCB (non-porous, 120/400 mesh), 6 mL)

AD0000AS (cartridge adaptor)

RFV0075P (reservoirs, 75 mL)

May 2013

Graphitized carbon black (GCB) is a reverse phase and anion exchange sorbent. GCB retains non-polar compounds, such as organochlorine pesticides, and some very polar compounds, such as surfactants, which are difficult to retain by other reverse phase sorbents. This simple SPE method uses UCT's proprietary, treated GCB for the determination of pesticides in water providing excellent recovery.

Procedure

1. Cartridge Preparation

- a) Transfer 100 mL of aqueous sample to a glass container
- b) Adjust pH to less than 2 using 6N HCl
- c) Spike as necessary
- d) Connect **RFV0075P** reservoirs to the top of the **EUCARB1M6** cartridges using **AD0000AS** adaptor
- e) Wash cartridges with 10 mL dichloromethane (DCM)
- f) Draw full vacuum to remove all DCM
- g) Add 10 mL methanol and draw down to top of frit
- h) Add 10 mL reagent water and draw down to top of frit
- i) Do not let cartridges go dry after step g)

2. Extraction

- a) Add samples to the reservoirs adjusting vacuum to give a dropwise flow, about 10 mL/min
- b) Rinse sample containers using 10 mL reagent water and add rinsate to cartridges
- c) Dry cartridges using full vacuum for 10 min

3. Elution

- a) Insert test tubes in the manifold then elute cartridges using 5 mL ethyl acetate dropwise followed by 5 mL of DCM dropwise


- b) Dry extracts by passing through anhydrous Na₂SO₄
- c) Rinse test tubes with DCM and add to Na₂SO₄
- d) Concentrate extracts to 1 mL using a gentle stream of N₂ at 35 °C
- e) Add IS prior to GC/MS analysis

4. Analysis

Parameters

GCMS: Agilent 6890N GC coupled with 5975C MSD
MSD Injector: 1 µL splitless injection at 250 °C
Injection Vol: 1 µL
Liner: 4 mm splitless gooseneck liner with deactivated glass wool (UCT: GCLGN4MM)
Column: Restek Rxi®-5sil MS 30m x 0.25mm x 0.25µm
Guard Column: 10 m
Column Flow Rate: 1.2 mL/min
Carrier Gas: He
Full Scan: 45-500 amu
Temperature Program: Initial T 55 °C hold for 1 min; ramp at 10 °C/min to 200 °C; ramp at 7 °C/min to 300 °C; hold for 0.21 min.

Detail of Reservoir, Adaptor, and Cartridge Setup



Accuracy and Precision Data

Compound	Intra-day (n=4)		Inter-day (n=17)	
	Rec%	RSD	Rec%	RSD
alpha Lindane	93	2.1	89	9.3
beta Lindane	96	1.9	91	8.8
gamma Lindane	93	1.7	92	8.3
delta Lindane	95	3.3	89	11.7
Heptachlor	97	3.2	91	11.1
Aldrin	95	1.5	84	12.9
Heptachlor epoxide	102	2.4	97	12.0
trans-Chlordane	93	3.8	90	8.8
Endosulfan I	94	5.0	91	8.4
cis-Chlordane	96	3.3	91	9.7
p,p'-DDE	91	3.5	89	8.8
Dieldrin	98	1.9	93	9.4
Endrin	100	2.1	95	11.8
Endosulfan II	105	1.4	97	10.3
p,p'-DDD	98	2.2	92	9.8
Endrin aldehyde	95	5.4	92	9.3
Endosulfan sulfate	102	3.8	97	10.2
p,p'-DDT	99	3.0	94	9.6
Endrin ketone	106	2.1	99	10.9

Methoxychlor	105	2.7	99	10.5
Dichlofluanid	107	2.8	98	10.8
Dicofol	95	0.7	86	11.6
Tolyfluanide	106	3.1	98	11.6
Captan	119	4.2	105	13.4
Folpet	107	3.9	95	10.0
Overall average	99	2.8	93	10.4



Determination of Carbonyl Compounds In Water by Dinitrophenylhydrazine Derivatization and HPLC/UV*

EPA Method 8315A

UCT Part Number:

EUC1812M15 (Unendcapped C18 - 2000 mg/15 mL cartridge)

March 2013

Method Summary

This method provides procedures for the determination of free carbonyl compounds and aldehydes in water. A measured volume of sample is buffered to pH 3 and the analytes derivatized at 40°C for one hour using 2,4-dinitrophenylhydrazine (DNPH) then extracted through SPE cartridges containing 2000 mg of C18. The cartridge(s) are eluted with 10 mL of ethanol and the derivatives are determined by absorbance at 360 nm after separation by HPLC.

The following compounds can be determined using this method:

Analyte	CASRN
2,5-Dimethylbenzaldehyde	5779-94-2
Acetaldehyde	75-07-0
Acetone	67-64-1
Acrolein	107-02-8
Benzaldehyde	100-52-7
Butanal	123-72-8
Crotonaldehyde	123-73-9
Cyclohexanone	108-94-1
Decanal	112-31-2
Formaldehyde	50-00-0
Heptanal	111-71-7
Hexanal (hexaldehyde)	66-25-1
Hexanal	66-25-1
Isovaleraldehyde	590-86-3
m-Tolualdehyde	620-23-5
Nonanal	124-19-6
Octanal	124-13-0
o-Tolualdehyde	529-20-4
Pentanal (propionaldehyde)	110-62-3
Propanal	123-38-6
p-Tolualdehyde	104-87-0

Note: Do not rinse glassware with acetone or methanol. These solvents react with DNPH to form interferences

Stock Standard Solutions

Stock standard formaldehyde solution approximately 1 mg/mL--

- Prepare by diluting 265 μ L of formalin to 100 mL with reagent water

Standardization of formaldehyde stock solution –

- Transfer 25 mL of a 0.1M Na₂SO₃ solution to a beaker
- Record the pH
- Add a 25.0 mL aliquot of the formaldehyde stock solution and record the pH
- Titrate this mixture back to the original pH using 0.1 N HCl
- Calculate the formaldehyde concentration using the following equation:

$$\text{Concentration (mg/mL)} = 30.03 \times (\text{N HCl}) \times (\text{mL HCl}) / 25.0 \text{ where:}$$

N HCl = Normality of HCl solution used

mL HCl = mL of standardized HCl solution used

30.03 = MW of formaldehyde

Note: The pH value of the 0.1 Na₂SO₃ should be 10.5 \pm 0.2 when the stock formaldehyde solution and the 0.1 M Na₂SO₃ solution are mixed together. The pH should be 11.5 \pm 0.2. It is recommended that new solutions be prepared if the pH deviates from this value

Stock aldehyde(s) and ketone(s) –

- Prepare by adding an appropriate amount of the analyte to 90 mL of methanol
- Dilute to 100 mL to give a final concentration of 1.0 mg/mL. Replace after six weeks, or sooner, if comparison with check standards indicates a problem

Reaction Solutions

2,4-Dinitrophenylhydrazine (DNPH [2,4-(O₂N)₂ C₆ H₃]NHNH₂) (3.00 g/L) –

- Dissolve 428.7 mg of 70% (w/w) reagent in 100 mL absolute ethanol. Sonication may be needed

Note: If the DHPH does not completely dissolve, filter the solution

Citrate buffer pH 3 (1 M) –

- Prepare by adding 80 mL of 1 M citric acid solution to 20 mL 1 M sodium citrate solution. Mix thoroughly. Adjust pH with 6N NaOH or 6N HCl

Sodium Chloride Solution (saturated) –

- Prepare by mixing an excess of the reagent grade solid with reagent water

Reducing agent, ammonium chloride (100 mg/L) –

- Add to all samples containing residual chlorine. The ammonium chloride may be added as a solid with stirring until dissolved, to each volume of water
- Sodium thiosulfate is not recommended because it may produce a residue of elemental sulfur that can interfere with some analytes

Procedure

1. Derivatization Procedure

- a) Quantitatively transfer 100 mL of sample into a 250 mL Florence or Erlenmeyer flask

Note: Other volumes may be used depending on the expected concentration of the analytes. If less than 100 mL is used adjust volume to 100 mL using reagent water

- b) Add 4 mL acetate buffer and adjust pH to 5.0 ± 0.1 with 6M HCl or 6M NaOH (if only formaldehyde is being analyzed)
- c) Add 4 mL of citrate buffer and adjust the pH to 3.0 ± 0.1 with 6M HCL or 6M NaOH (if other aldehydes are being analyzed)
- d) Add 6 mL of DNPH reagent, seal the flask and place in a heated orbital shaker at 40° C for 1 hour
- e) Adjust the agitation to produce a gentle swirling of the reaction solution

2. Cartridge Conditioning

- a) Assemble a vacuum manifold
- b) Place **EUC1812M15** cartridge(s) on manifold
- c) Condition the cartridge by adding 15 mL of acetonitrile (ACN)
- d) Draw through then rinse using a 15 mL portion of reagent water
- e) While cartridge is still wet add 10 mL of dilute citrate buffer (10 mL of 1M citrate buffer dissolved in 250 mL of water)
- f) Remove the reaction vessel from the shaker after 1 hour and add 10 mL of

saturated NaCl solution to the flask

- g) Quantitatively transfer the reaction solution to the SPE cartridge and apply a vacuum to draw the solution through at 3-5 mL/min
- h) Continue to apply vacuum for about 1 minute after the liquid sample has passed through the cartridges

3. Cartridge Elution

- a) Turn off vacuum and place 10 mL volumetric flasks under cartridges
- b) Turn on vacuum and elute the cartridges with approximately 9 mL of acetonitrile directly into flasks
- c) Bring the eluate to 10 mL volume with ACN, mix thoroughly, then transfer aliquot to analysis vial and seal

Note: Because this method uses an excess of DNPH, the cartridges will retain a yellow color after this step. The color does not indicate incomplete recovery of the analyte derivatives

HPLC Analysis (suggested)

- **HPLC Column:** C18, 250 mm x 4.6 mm i.d., 5 µm particle size
- **Mobile Phase:** 70/30 methanol/water (v/v)
- **Injection Vol:** 20 µL
- **Flow Rate:** 1.2 mL/min
- **UV Detector:** 360 nm
- **Flow Program:**
 - 70/30 methanol/water (v/v) for 20 min,
 - to 100% acetonitrile in 15 minutes
 - hold at 100% acetonitrile for 10 minutes

Representative Retention Times in Reagent Water Using Specified Analysis Conditions

Analyte	Retention Time (<i>min</i>)
Formaldehyde	5.3
Acetaldehyde	7.4
Propanal	11.7
Crotonaldehyde	16.1
Butanal	18.1
Cyclohexanone	27.6
Pentanal	28.4
Hexanal	34.1
Heptanal	35.0
Octanal	40.1
Nonanal	40.4
Decanal	44.1

*Adapted from Method 8315a, Determination Of Carbonyl Compounds By High Performance Liquid Chromatography (HPLC), and EPA Method 554, Revision 1.0, James W. Eichelberger, W.J. Bashe (Technology Applications, Incorporated), Environmental Monitoring Systems Laboratory, Office Of Research And Development, U. S. Environmental Protection Agency, Cincinnati, Ohio 45268



EPA Method 1664 Hexane Extractable Material (HEM)

UCT Product Number:

ECUNIOGXF (2000 mg endcapped C18, 83 mL cartridge)

August 2013

Procedure

1. Assemble

- a) Assemble cartridge adapters on the vacuum manifold
- b) Place **ECUNIOGXF** cartridge(s) on each vacuum station needed
- c) Connect the manifold to a suitable trap and attach the trap to a vacuum system capable of attaining a minimum of 25" Hg (635 mm) of vacuum
- d) Insert waste collection vial in manifold (optional)

2. Prepare Water Sample

- a) Adjust the pH of the sample to less than 2 by adding 5 mL of 6N HCl or 2.5 mL of concentrated H₂SO₄. The pH of deionized water cannot be accurately determined using pH paper
- b) If acid was added to the sample in the field, do not add more unless the pH is greater than 2

Note: Gloves are recommended as skin oils may affect final sample weight

3. Condition the Cartridge

- a) Rinse the sides of the cartridge and bottle holder with 10 mL of hexane
- b) Allow cartridge to soak for 1 minute
- c) Draw the hexane through the cartridge using vacuum
- d) Discard the hexane from the waste vial
- e) Draw full vacuum through the cartridge for 2 minutes to dry
- f) Add 10 mL of methanol to the cartridge
- g) Slowly draw the methanol through leaving a layer on the cartridge frit
- h) Soak for one minute
- i) Draw methanol through the cartridge
- j) Remove waste vial and discard methanol
- k) Immediately, add 80 mL (fill cartridge) of DI water to the cartridge

- l) Soak for 1 minute
- m) Draw all of the water through the cartridge to waste

4. Sample Addition

- a) Add the 1 liter water sample directly to the cartridge
- b) Draw the sample through the cartridge under low vacuum. This may take several minutes depending on the solids in the sample. (Note 1) Increase vacuum pressure if necessary. Do not exceed 50 mL/minute for optimum recoveries. This is a fast drip, but not a stream
- c) Remove the cartridge and tap any excess water from the bottom of the cartridge
- d) Replace cartridge and allow to dry under full vacuum for 10 minutes
- e) Remove any water remaining in the bottom support of the cartridge with a paper towel if necessary

5. Elution

- a) Prepare an extract collection vial containing about 8 mm (0.3 inch) of anhydrous sodium sulfate
- b) Place the vial in the manifold station under the cartridge
- c) Add a thin layer of anhydrous sodium sulfate to the top of the cartridge
- d) Rinse the water sample bottle with 10 mL of hexane
- e) Add the hexane to the cartridge
- f) Soak cartridge for 3 minutes. A slow drip of hexane is permissible
- g) Turn on vacuum and slowly draw the hexane through the cartridge and into the collection vial
- h) Turn off vacuum then repeat steps 5 d) – g) 2 additional times with 10 mL of hexane
- i) Do not allow the solvent to splash into the collection vial
- j) Add another 10 mL of hexane to the cartridge, rinsing the bottle holder
- k) Soak cartridge for 3 minutes
- l) Draw the hexane through the cartridge and collect

6. Dry the Extract

- a) Remove the collection vial from the manifold and cover with a screw cap

- b) If water is still present in the extract (no free moving sodium sulfate is visible) shake the extract to form a water/hexane emulsion and immediately pour the extract through a sodium sulfate funnel or column containing approximately 30 g of anhydrous sodium sulfate held in place with a glass wool plug or frit. **Do not use filter paper**
- c) Collect the extract in a clean, pre-weighed vessel
- d) Rinse the collection vial with hexane and add it to the sodium sulfate and collect. This will rinse the vial and the sodium sulfate. Poor rinsing of the sodium sulfate will result in low recoveries

7. Gravimetric analysis

- a) Carefully evaporate the hexane using a nitrogen evaporator at 40° C until the extract just reaches dryness

Note: Do not over dry or low recoveries will result

- b) Allow to cool to room temperature in a desiccator before weighing
- c) Record this weight as the mass per unit volume of HEM and report as mg/L

Notes

- 1) If very high solids are present, add a small plug of glass wool to the cartridge prior to extraction to prevent clogging and improve flow. The glass wool must be thoroughly rinsed with hexane as part of the cartridge during the elution step.
- 2) Stearic acid must be in solution in the spiking solution or low recoveries will result. If small crystals are present in the spiking solution, sonicate or shake until dissolved.
- 3) If white crystals are present in the sample bottle after elution, the sample pH was not low enough prior to extraction. Repeat with lower pH
- 4) HCl will lose strength over time. Sulfuric acid is a good substitute
- 5) Any residue that does not rinse from the bottle or elute from the cartridge is not HEM
- 6) Placing a small amount of sodium sulfate on the top frit of the cartridge during the sorbent drying process can reduce the amount of water passed through the cartridge to the final extract



Solid-Phase Extraction and Determination of Organotin by Micro-Liquid Chromatography Electrospray Ion Trap MS (Part A Water Samples)

UCT Product Number:

ECUNIC18 (C18 endcapped, 1100 mg/83 mL cartridge)

EPA Method 8323

March 2010

Analytes Determined Using This Method

Analyte	CAS No	LOD
Tributyltin chloride	1461-22-9	780 pg
Dibutyltin dichloride	683-18-1	970 pg
Monobutyltin trichloride	1118-46-3	1 ng
Triphenyltin chloride	668-34-8	NA
Diphenyltin dichloride	1135-99-5	920 pg
Monophenyltin trichloride	1124-19-2	NA

Note: Organotins can bond to glass surfaces, glassware must be specially treated. All glassware used in the extraction and analysis of organotins must be acid washed using the following procedure.

- Wash glassware in hot soapy water then rinse with DI water
- Prepare a pH 2 acid bath using 12 N HCl and soak glassware in acid for 24 hours
- Remove glassware from bath, then rinse with DI water followed by a methanol rinse
- Place in a 60°F oven until dry

Procedure

1. Initial Preparation

- a) Fill a 2 liter volumetric flask with sample water
- b) Adjust pH to 2.5 by adding about 600 µL of 12N HCl
- c) Stopper flask and invert several times to mix acid

2. Cartridge Conditioning

- a) Add 10 mL of methanol to the cartridge to activate

- b) Briefly turn on vacuum to draw through a small amount to top of frit
- c) Wait 1-2 minutes
- d) Add 10 mL of a methanol/1% acetic acid solution
- e) Draw about 2 mL through the cartridge then turn off vacuum
- f) Let solution sit for 1-2 minutes then draw through
- g) Add 10 mL of reagent water to cartridge and partially draw through

Note: Do not let the cartridge dry out after start of activation otherwise start over at step 2. a)

3. Sample Extraction

- a) Add the 2 liter sample to the cartridge and draw through at approximately 50 mL/minute (fast drip)
- b) Rinse volumetric flask and cartridge sides with 100 mL of reagent water and draw through

4. Elution

- a) Dry cartridge by drawing full vacuum for 10 minutes
- b) Place a clean, treated collection tube in the manifold
- c) Add a **first** portion of 10 mL of methanol/1% acetic acid solution to the cartridge, rinsing the sides during addition, then slowly draw through cartridge
- d) Add a **second** 10 mL portion of methanol/1% acetic acid solution to the cartridge, then slowly draw through
- e) Add a **third** 10 mL portion of methanol/1% acetic acid solution to the cartridge, then slowly draw through

Micro-concentration by TurboVap[®] Nitrogen Evaporation

1. Place the concentrator tube in the TurboVap[®] or other analytical evaporator in a lukewarm water bath at 30° C
2. Evaporate the solvent volume to 0.5 mL using a gentle stream of clean, dry N₂
3. The internal wall of the tube must be rinsed down several times with the final solvent (methanol/1% acetic acid) during the evaporation
4. Do not allow the extract to become dry
5. Transfer the extract to a 2 mL glass vial with a PTFE-lined screw-cap or crimp-top vial and store refrigerated at 4° C
6. Sample is ready for μ -LC-ES-ITMS analysis



Nitroaromatics and Nitramines by High Performance Liquid Chromatography (HPLC) Aqueous Matrices

UCT Product Number:
ECDVB156 (500 mg DVB, 6 mL cartridge)

EPA Method 8330

September 2009

RCRA Compounds Using This Method

Table 1

Analyte	CAS	Abbreviation	% Recovery n=3
Octahydro-1,3,5,7-tetranitro-1,3,5,7-tetrazocine	2691-41-0	HMX	100
Hexahydro-1,3,5-trinitro-1,3,5-triazine	121-82-4	RDX	110
1,3,5-Trinitrobenzene	99-35-4	1,3,5-TNB	100
1,3-Dinitrobenzene	99-65-0	1,3-DNB	100
1,4-Dinitrobenzene	10025-4	1,4-DNB	97
Methy-2,4,6-trinitrophenylnitramine	47945-8	Tetryl	85
Nitrobenzene	98-95-3	NB	100
2,4,6-Trinitrotoluene	118-96-7	TNT	94
4-Amino-2,6-dinitrotoluene	19406-51-0	4-Am-DNT	120
2-Amino-2,6-dinitrotoluene	35572-78-2	2-Am-DNT	110
2,4-Dinitrotoluene	121-14-2	2,4-DNT	98
2,6-Dinitrotoluene	606-20-2	2,6-DNT	110
2-Nitrotoluene	88-72-2	2-NT	90
3-Nitrotoluene	99-08-1	3-NT	91
4-Nitrotoluene	99-99-0	4-NT	92
Nitroglycerin	55-63-0	NG	100
Pentaerythritol tetranitrate	78-11-5	PETN	100
3,5-Dinitroaniline	618-87-1	3,5-DNA	68
1-Nitronaphthalene	86-57-7	NN	97
o-Dinitrobenzene	528-29-0	o-NB	100

Safety

- Extra caution that should be taken when handling the analytical standard neat material to prevent detonation

Interferences

- 2, 4-DNT and 2, 6-DNT elute at similar retention times on C18 columns using method separation conditions. If it is not apparent that both isomers are present or are not detected an isomeric mixture should be reported
- Tetryl is thermally labile (decomposed with heat at temperature above room temperature) and decomposes in methanol/water solutions. All aqueous samples expected to contain tetryl should be diluted with acetonitrile and acidified with sodium bisulfate to pH <3 prior to filtration
- Degradation products of tetryl appear as a shoulder on the 2,4,6-TNT peak when using C18 columns

Note: All samples should be stored at 2° to 4° C prior to extraction and should be extracted within 14 days of collection

Sample Preparation--Aqueous matrices, (e.g. water)

(from Method 3535)

Important Notes:

- Any particulate matter in the original sample must be included in the sample aliquot that is extracted.
- The sample container must be rinsed with solvent as the majority of organic analytes are hydrophobic and may adhere to the sample container surfaces.
- Do not concentrate explosives residue to dryness as they may DETONATE
- For explosives and nitramines or nitroaromatics the extraction pH should be as received in the sample
- Using a graduated cylinder, measure 1 liter of sample water. A smaller sample size may be used when analytical sensitivity is not a concern
- Add methanol (ACN if tetryl is being analyzed) so that the sample is 0.5% v/v. Add surrogate standards to all samples and blanks
- Add matrix spikes standards to representative sample replicates
- Adjustment of sample pH may result in precipitation or flocculation reactions and potentially remove analytes from the aqueous portion. Transfer the precipitate with rinses to the SPE extraction cartridge
- Do not let the cartridge dry out after cartridge conditioning with acetonitrile (ACN)

Glass Apparatus and SPE Cartridge Washing

Analyte	1 st solvent wash	2 nd solvent wash	3 rd solvent wash
Explosives	5 mL acetone	15 mL isopropanol	15 mL methanol
Nitramines, Nitroaromatics	5 mL ACN	15 mL ACN	

1. Cartridge Conditioning

- a) Follow the 4 steps in the table below for solvent quantities and soak times

Analyte	Condition Step 1	Step 2	Step 3	Step 4
Explosives	20 mL ACN, 3 min*	20 mL ACN	50 mL DI water	50 mL DI water
Nitramines, Nitroaromatics	15 mL ACN, 3 min*	30 mL DI water	----	----

- b) Draw solvents through the cartridge under low vacuum

Do not let the cartridge dry out once cartridge is conditioned. This will affect analyte recovery

2. Sample Extraction

- Add the contents of the sample bottle to the cartridge
- Adjust vacuum to about 10-15 mm Hg to obtain a uniform flow rate of about 10 ml per minute.
- After all the sample is drawn through, draw air through the cartridge for 15 minutes to dry it
- Do not dry for longer than 20 minutes as lower recovery may result

3. Cartridge Elution

- a) Insert a collection tube in the base of the vacuum manifold

For Explosives

- Add 4 mL of ACN and soak for 3 minutes
- Draw through using gravity flow or very low vacuum into a collection tube
- Store extract in freezer until analysis

For Nitramines and Nitroaromatics

- Add 5 mL ACN, soak for 3 minutes
- Draw through using a gravity flow or very low vacuum into a collection tube
- Store extract in freezer until analysis

4. Extract Concentration (if necessary)

- a) Concentrate the extract to 0.7 mL under a gentle stream of nitrogen in a warm bath at 40° C
- b) Transfer the extract to a 1 mL volumetric flask
- c) Add internal standard for an extract concentration of 5 µg/mL
- d) Extract is now ready for analysis by HPLC

RP-HPLC Columns for the Analysis of Explosive Residues	
Primary Columns	C-18 reversed-phase HPLC column, 25-cm x 4.6-mm, 5 µm C8 reversed-phase HPLC column, 15-cm x 3.9-mm, 4 µm
Secondary Columns	CN reversed-phase HPLC column, 25-cm x 4.6-mm, 5 µm Luna Phenyl-Hexyl reversed-phase HPLC column, 25-cm x 3.0-mm, 5 µm

Injection volume: 100 µL

UV Detector: Dual 254 & 210 nm or Photodiode Array

Mobile phase: For C18 & CN column, 50:50 methanol:water

Flow Rate: 1.5 mL/minute

Retention Times and Capacity Factors

Analyte	LC-18 RT minutes	LC-CN RT minutes
HMX	2.44	8.35
RDX	3.78	6.15
1,3,5-TNB	5.11	4.05
1,3-DNB	6.16	4.18
3,5-DNA	6.90	NA
Tetryl	6.93	7.36
NB	7.23	3.81

NG	7.74	6.00
2,4,6-TNT	8.42	5.00
4-Am-DNT	8.88	5.10
2-Am-DNT	9.12	5.65
2,6-DNT	9.82	4.61
2,4-DNT	10.05	4.87
2-NT	12.26	4.37
4-NT	13.26	4.41
PETN	14.10	10.10
3-NT	14.23	4.45

*For complete details on Method 8330 "Nitroaromatics and Nitramines by High performance Liquid Chromatography" Revision 2 October 2006, the analyst is referred to: National Exposure Research Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268 and Method 3550 Revision 0, December 1996



Alternative Test Procedure for the Measurement of Organochlorine Pesticides and Polychlorinated Biphenyls in Wastewater

UCT Part Numbers:

ECUNIC18 (1100 mg endcapped C18, 83 mL) or

EEC181M6 (1000 mg endcapped C18, 6 mL)

Optional:

ECSS15M6 (5000 mg Na₂SO₄, 6 mL)

ECCU01K (1 kg activated copper granules)

EUFLS1M6 (1000 mg PR Grade Florisil[®], 6 mL)

EPA Method 608 ATP*

May 2013

This is a gas chromatographic (GC) method for determination of compounds listed below in municipal and industrial discharges. The EPA has approved the use of C18 cartridges for this method.

Analytes Recovered Using Method 608ATP

Analyte	CAS
Aldrin	309-00-2
α-BHC	319-84-6
β-BHC	319-85-7
γ-BHC (Lindane)	58-89-9
δ-BHC	319-86-8
α-chlordane	5103-71-9
γ-chlordane	5103-74-2
4,4'-DDD	72-54-8
4,4'-DDE	72-55-9
4,4'-DDT	50-29-3
Dieldrin	60-57-1
Endosulfan I	959-98-8
Endosulfan II	33213-65-9
Endosulfan sulfate	1031-07-8
Endrin	72-20-8
Endrin Aldehyde	7421-93-4
Endrin Ketone	53494-70-5
Heptachlor	76-44-8
Heptachlor epoxide	1024-57-3
Methoxychlor	72-43-5
Toxaphene	8001-35-2
PCB-1016	12674-11-2
PCB-1221	1104-28-2
PCB-1232	11141-16-5
PCB-1242	53469-21-9
PCB-1248	12672-29-6
PCB-1254	11097-69-1
PCB-1260	11096-82-5

Procedure

1. Condition Cartridge

- a) Insert cartridge(s) into the manifold or automated extraction system
- b) Wash with 10 mL of methylene chloride (MeCl_2)
- c) Soak for 1 minute then draw through to waste
- d) Draw air under full vacuum to dry cartridge
- e) Add 10 mL of methanol (MeOH) then slowly draw through to top of frit
- f) Soak for 1 minute

Note: Do not let the cartridge go dry after addition of methanol otherwise repeat at step1) e)

- g) Rinse the cartridge with 10 mL of reagent water
- h) Draw through leaving a thin layer on top of frit

2. Sample Extraction

- a) Adjust 1 liter of sample pH to < 2 using sulfuric acid
- b) If sample water is high in suspended solids, allow particulates to settle then slowly decant the water in the bottle. Once most of the water passes through the cartridge add the solids portion
- c) Draw the sample water through the cartridge over a 20-30 minute time period (fast drip) by adjusting the vacuum
- d) Dry the cartridge by drawing air under full vacuum through for 10 minutes

3. Extract Elution

- a) Insert a collection tube into the vacuum manifold
- b) Add 5 mL of acetone to the sample bottle then swirl
- c) Add this to the cartridge
- d) Soak for 1 minute and slowly collect eluate
- e) Add 20 of methylene chloride to the sample bottle cover and shake. Add this to the cartridge
- f) Soak for 2 minutes and slowly collect eluate
- g) Rinse the inside walls of the sample bottle using 10 mL of methylene chloride then transfer solvent to the cartridge using a disposable pipette rinsing the inside of the cartridge
- h) Soak for 2 minutes then collect eluate

4. Sample Drying

- a) Pour the combined elutes together through a drying tube **ECSS15M6**.
Alternatively, use 5 grams of sodium sulfate over a bed of glass wool in a funnel
- b) Rinse the drying tube or sodium sulfate bed with 2 x 3 mL portions of 1 methylene chloride
- c) Concentrate sample using a Kuderna-Danish (KD) concentrator while performing solvent exchange into hexane
- d) Concentrate sample under a gentle stream of N₂ while gently heating in a water bath. **Other drying techniques may be used**
- e) Rinse the inside walls of the concentrator tube two or three times with hexane during the evaporation
- f) Adjust the final volume of the extract to 10 mL

Florisil PR[®] or Copper Granule Clean-up Procedure (if needed)

Clean-up procedures may not be needed for relatively clean samples. If required, the following procedure can be used to remove polar interferences from organochlorine pesticide and PCB extracts in hexane prior to analysis.

5. Florisil PR[®] Clean-Up

- a) Place a cartridge in a vacuum manifold
- b) Pre-rinse the Florisil[®] column with 10 mL of 90:10 hexane/acetone using gravity flow (a low vacuum may be necessary to start flow)
- c) Discard solvent
- d) Add a collection tube under the column
- e) Add a 2 mL aliquot of the sample extract (in hexane) to the cartridge
- f) Collect extract by gravity
- g) Add 10 mL of 90:10 hexane/acetone to the cartridge
- h) Continue to collect by gravity or low vacuum
- i) Gently evaporate the extract to a volume of 1 mL
- j) Adjust eluate to a final volume of 2 mL with hexane
- k) Sample is now ready for analysis

6. Sulfur Clean-up

- a) Place 4 grams of **ECCU01K** copper granules in a glass vial
- b) Add 2 mL of liquid sample extract to the vial
- c) Seal the glass vial and mix sample with copper for 2 minutes
- d) Allow to stand for approximately 10 minutes
- e) If sample contains high levels of sulfur, repeat process with 4 grams of fresh copper granules

Note: For the analysis of PCB type analytes, copper may reside in the extract

7. Analysis--GC/ECD

- a) Transfer clean extract to autosampler vial
- b) Sample is now ready for analysis

*The EPA has accepted the use of C18 bonded phases in packed cartridge format expanding the method from a disk only approach. For complete details on Method 607ATP, the analyst is referred to: "An alternative test procedure for the measurement of organochlorine pesticides and polychlorinated biphenyls in waste water", Federal register/Vol.60, No.148, August 2, 1995, Environmental Monitoring Systems Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268

Florisil® is a registered trademark of U.S. Silica

DCN-216111-171



Extraction of Metals

UCT Part Number:

EUTAX15Z (Tri-acetic Acid 500 mg/10 mL)

February 2009

Metals: Tin, Nickel, Mercury, Copper, Chromium, Ruthenium

Matrix: Water, Blood, Biological Fluids, Organic Solvents & Tissue Homogenates

1. Sample Pre-treatment

- It is important when using ion exchangers to adjust the pH of the analyte of interest so that it is totally ionized. Information about the analyte pKa is important

Aqueous or Organic Solvent Samples:

- Adjust the sample to pH 7 with 100 mM dibasic sodium phosphate buffer or ammonium hydroxide then vortex

Whole Blood, Serum or Plasma:

- To 1 ml of sample add 4 ml of D.I. H₂O and vortex
- Let stand 5 minutes and centrifuge for 10 minutes at 2000 rpm and discard pellet
- Adjust to pH 9.0 with 100 mM dibasic sodium phosphate buffer or ammonium hydroxide

2. Column Conditioning

- Add 3 ml of methanol and draw to waste
- Add 3 ml of water and draw to waste
- Add 3 ml of buffer pH 7.0 and draw to waste leaving sorbent wet

3. Sample Application

- Apply the sample to the column at a rate of 1 ml per minute. A faster rate of application may exceed the rate of ion-exchange

4. Analyte Purification

- Wash the column with 2 ml of pH 7.0 buffer used in column equilibration.

5. Elution

- Elute with 3 ml of acidic methanol (2% HCl, pH 2.0).

Alternative eluants:

- Elute with 3 ml of acidic methanol (Formic acid to pH 2.0).
- Elute with 3 ml of 0.1 M Nitric Acid (pH 2.0).



Ion Exchange Sorbents for Metals Extraction-Analysis & Sorbent Use Selection Guide

UCT Product Number:

UCT ENVIRO-CLEAN® (Ion-Exchange Cartridges)*

January 2010

The determination of trace metals in aqueous environmental samples or other matrices often require sample pretreatment and clean-up procedures prior to analysis by using specific ion-exchange sorbents. The sorbents are used to eliminate matrix interferences and achieve high concentrations of metal ions for good analytical accuracy. They are important when using such techniques as AA, IES and ICP-AES.

The use of ion-exchange sorbents for the preconcentration, separation, and determination of metal ions for trace analysis is well established in the literature. Selection of an appropriate sorbent ensures both high efficiency in metal chelating while minimizing the mass of sorbent required for a particular analytical task. A high efficiency sorbent means that a smaller bed mass may be used thereby reducing the quantity of solvent required for elution yielding greater analytical sensitivity.

Recommendations in this application note include the following metal ions:

Zinc (II)	Arsenic (V)	Tin (IV)	Selenous IV)
Mercury (II)	Chromium (III)	Copper (II)	Platinum (0)

Other metal ions may be extracted by the use of these ion-exchange sorbents

Sorbent Selection for Metals Extraction

Solid-phase sorbents have differing capacity and selectivity for various metal ions due to the specific nature of the ion-exchange functional group, the metal species and the valence state of the metal of interest. Depending on the specific metal ion of interest, elution of the cartridge may be most efficient using both the acid followed by the base elution procedure. This can be determined by looking at the following Extraction Protocol Tables. For example, when eluting Hg(II) from PSA the highest recovery is obtained using acid elution (green box) followed by base elution (yellow box).

Sample Analysis

1. Sample Extraction

- a) Place a UCT ion exchange cartridge on a SPE manifold

Note: Cartridge selection will depend on the volume of sample or the concentration of metal to be extracted

- b) Condition 1 mL cartridge by adding 3 mL of methanol. (Larger cartridges will require a larger volume of solvent and water wash volume in steps c) and d))
- c) Add 3 mL of reagent water and allow to drip through the cartridge

Note: Do not allow the cartridge to dry out after addition of water, otherwise repeat step d)

- d) Add 10-50 mL of sample water to the cartridge. A larger sample volume may be used depending on metal concentration or suspended solids content
- e) Adjust vacuum setting so that the water flows at 1-3 mL/minute until sample has passed completely through the cartridge
- f) Allow the cartridge to air dry for about 1 minute under full vacuum

2. Elution--Acid

- a) Prepare a 100 mM nitric acid elution solution
- b) Place a collection vial in the vacuum manifold
- c) Add 3 mL of the nitric acid solution to the cartridge
- d) Adjust flow rate for a flow of 1-3 mL/minute
- e) Dilute eluant to an appropriate volume for detection using reagent water
- f) Sample is ready for analysis

3. Elution--Base

- a) Prepare a 100 mM triethylamine elution solution
- b) Place a collection vial in the vacuum manifold
- c) Add 3 mL of the triethylamine solution to the cartridge
- d) Adjust flow rate for a flow of 1-3 mL/minute

- e) Dilute eluant to an appropriate volume for detection using reagent water
- f) Sample is ready for analysis

4. Analysis

- a) Prepare calibration curves for use with atomic absorption (AA) or Inductively Coupled Plasma Atomic Emission Spectroscopy (ICP-AES) using appropriate metallic ion standards

Extraction Protocol Tables

How to use these tables: When choosing an ion-exchange sorbent to capture arsenic (V) for example, all ion-exchange sorbents will capture a small quantity of metal ions, however, only a base extraction would elute metal ions from these sorbents. For extraction of zinc ions, all sorbents would have moderate to high capacity but elution could only occur from the sorbent using acidic elution conditions.

Acid Extraction Protocol

Sorbent	Cu (II)	Zn (II)	As (V)	Sn (IV)	Se (IV)	Hg (II)	Cr (III)	Pt (0)
PSA	Good to High Capacity	Good to High Capacity	Little or No Capacity	Little or No Capacity	Good to High Capacity	Good to High Capacity	Moderate Capacity	Moderate Capacity
BCX-HL	Good to High Capacity	Moderate Capacity	Little or No Capacity	Little or No Capacity	Little or No Capacity	Good to High Capacity	Moderate Capacity	Little or No Capacity
CCX	Moderate Capacity	Moderate Capacity	Little or No Capacity	Moderate Capacity	Little or No Capacity	Good to High Capacity	Little or No Capacity	Good to High Capacity
TAX	Good to High Capacity	Good to High Capacity	Little or No Capacity	Moderate Capacity	Little or No Capacity	Good to High Capacity	Moderate Capacity	Moderate Capacity
THX	Good to High Capacity	Moderate Capacity	Little or No Capacity	Little or No Capacity	Moderate Capacity	Moderate Capacity	Little or No Capacity	Little or No Capacity
NAX	Good to High Capacity	Good to High Capacity	Little or No Capacity	Little or No Capacity	Good to High Capacity	Good to High Capacity	Little or No Capacity	Little or No Capacity

Base Extraction Protocol

Sorbent	Cu (II)	Zn (II)	As (V)	Sn (IV)	Se (IV)	Hg (II)	Cr (III)	Pt (0)
PSA	Little or No Capacity	Little or No Capacity	Moderate Capacity	Moderate Capacity	Good to High Capacity	Moderate Capacity	Little or No Capacity	Little or No Capacity
BCX-HL	Little or No Capacity	Little or No Capacity	Moderate Capacity	Good to High Capacity	Moderate Capacity	Good to High Capacity	Moderate Capacity	Little or No Capacity
CCX	Little or No Capacity	Little or No Capacity	Moderate Capacity	Good to High Capacity	Moderate Capacity	Moderate Capacity	Little or No Capacity	Little or No Capacity
TAX	Little or No Capacity	Little or No Capacity	Moderate Capacity	Good to High Capacity	Moderate Capacity	Good to High Capacity	Little or No Capacity	Little or No Capacity
THX	Little or No Capacity	Little or No Capacity	Moderate Capacity	Moderate Capacity	Moderate Capacity	Good to High Capacity	Little or No Capacity	Little or No Capacity
NAX	Little or No Capacity	Little or No Capacity	Moderate Capacity	Moderate Capacity	Moderate Capacity	Moderate Capacity	Little or No Capacity	Little or No Capacity

Good to High Capacity	Good to High Capacity
Moderate Capacity	Moderate Capacity
Little or No Capacity	Little or No Capacity

Ion-Exchange Sorbent Key

PSA	Primary secondary amine
BCX-HL	Benzene sulfonic acid –high load
CCX	Carboxylic acid
TAX	Triacetic acid
THX	Sulfhydryl (thiopropyl)
NAX	Aminopropyl

Primary Secondary Amine (PSA)

The PSA ion-exchange sorbent has a significant capacity for Hg(II) Se(IV) followed by a lesser capacity for Sn(IV), Cu(II), Zn(II) and Cr(III). Metal ions are readily eluted from PSA by the use of weak acid solutions such as 100 mM nitric acid solution. Additional recovery for selenium can be obtained by following the acid elution by the use of 100 mM triethylamine solution

Benzenesulfonic Acid-High Load (BCX-HL)

The BCX-HL ion-exchange sorbent is the least selective of all ion-exchange sorbents and has significant capacity for Hg(II) and Sn(IV) thus ensuring high extraction efficiency for trace analysis. It is also a strong sorbent for Cu(II), Zn(II), Cr(III) and small amounts of Pt. In most cases, metal ions are readily eluted from BCX-HL by the use of 100 mM nitric acid solution. Improvement in Hg(II) recovery yield and Sn(IV) can be achieved when eluting with 100 mM triethylamine solution.

Carboxylic Acid (CCX)

The CCX sorbents have high selectivity for Sn(IV) and Hg(II). Metal ions are readily eluted from CCX by the use of weak acid solutions such as 100 mM nitric acid solution. Sn(IV) is eluted using 100 mM triethylamine solution. Additional Hg(II) is released under basic elution.

Triacetic Acid (TAX)

TAX sorbents have the highest affinity for Sn(IV) and Hg(II) followed by lesser amounts of Cu(II) and Zn(II). Metal ions are readily eluted from TAX by the use of weak acid solutions such as 100 mM nitric acid solution. Sn(IV) is eluted using 100 mM triethylamine solution. Additional Hg(II) is released under basic elution.

Sulfhydryl THX (thiopropyl)

THX sorbents have the highest affinity for Hg(II) and Sn(IV), and approximately equal weights of Sn(IV) and Cu(II). Metal ions are readily eluted from THX by the use of weak acid solutions such as 100 mM nitric acid solution. Sn(IV), Se(IV) and Hg(II) are eluted using 100 mM triethylamine solution.

Aminopropyl (NAX)

The NAX ion-exchange sorbent has a significant capacity for Hg(II) followed by Se(IV). Metal ions are readily eluted from NAX by the use of weak acid solutions such as 100 mM nitric acid solution. Additional Hg(II) is released under basic elution.

*UCT ENVIRO-CLEAN[®] Ion-exchange cartridges are available in a variety of cartridge sizes, sorbent mass and particle size to most analytical requirements. For further information, contact UCT.

Agricultural + Food Safety





Pesticides in Fatty Matrices Extraction

UCT Part Numbers:

ECPSAC1856 (500 mg endcapped C18, 500 PSA, 6 mL cartridge)

CUMPSC18CT (150 mg MgSO₄, 50 mg PSA and 50 mg C18 in a 2 mL centrifuge tube)

ECMAG00D (500 g organic free MgSO₄ anhydrous)

ECNACL05K (5 kg NaCl)

January 2011

Procedure

1. Sample Preparation

- a) Weigh 20.0 ± 0.10 grams (g) of homogenized sample into a 250 mL plastic centrifuge bottle, tared on a balance capable of weighing to 0.01 grams
- b) Fortify each sample with process control spiking (PCS) solution
- c) Add 50 mL of ethyl acetate (EtOAc) to each tube
- d) Fortify each sample with internal standard (ISTD) spiking solution
- e) Reduce sample material particle size by using a high speed disperser for approximately 1 minute
- f) Add 2 g of anhydrous MgSO₄ (**ECMAG00D**) and 0.5 g anhydrous NaCl (**ECNACL05K**)

Note: Carefully add the reagents to the tube to avoid contaminating the threads or rims of the tubes otherwise leaks may result

- g) Seal the tube and shake vigorously for approximately 1 minute either mechanically or by hand. Make sure the solvent interacts well with the entire sample and that crystalline agglomerates are broken up
- h) Cool the sample in a -20 °C freezer for approximately 30 minutes
- i) Centrifuge at 10,000 RCF for 5 minutes
- j) Decant at least 50 mL of the EtOAc layer into a 50 mL glass graduated centrifuge tube using a funnel and filter paper. Allow the extract to come to room temperature and adjust the volume with EtOAc to 50 mL using a Pasteur pipette
- k) Concentrate the extract under a stream of nitrogen with a 70° C water bath until the volume remains constant (this will be ~ 3 mL and will take about 1 hour)

- l) Dilute to 20 mL with acetonitrile (MeCN) and cap with a glass stopper, vortex for 1 minute
- m) Freeze at -70 °C for 30 minutes
- n) Centrifuge the extract while frozen for 3 minutes (The MeCN will thaw during centrifugation)
- o) Directly after centrifugation in step n), filter > 15 mL of the MeCN layer of the extract with a 0.45 µm syringe filter into a 15 mL glass centrifuge tube
- p) Allow the extract to come to room temp, adjust the volume to 15 mL, and concentrate to 2.25 mL under a stream of nitrogen with a 70 °C water bath

2. LC-MS/MS Analysis

- a) Transfer 1 mL of extract to a 2 mL mini-centrifuge tube **CUMPSC18CT**
- b) Vortex for 1 minute and centrifuge
- c) Transfer to auto sampler vial. Sample is now ready for analysis

3. GC Analysis

- a) For GC analyses, use the dual layer cartridge **ECPSAC1856**
- b) Add approximately 0.75 – 0.80 grams (~ 0.6 cm = 0.25 inches) of anhydrous MgSO₄ added to the top of the cartridge
- c) Condition the SPE cartridge by adding one cartridge volume (4.0 mL) of MeCN using a UCT positive pressure SPE manifold
- d) Elute to waste
- e) Place a labeled 15 mL graduated disposable plastic centrifuge tube below the cartridge in the positive pressure SPE manifold
- f) Quantitatively transfer 1 mL of the sample extract from step 15 to the SPE cartridge
- g) Elute SPE cartridge in a dropwise manner (Regulated Flow Pressure = 35 psi) into a labeled 15 mL graduated glass centrifuge tube using MeCN
- h) Collect the eluate while washing the SPE cartridge **three times** with **4 mL of eluant**.
- i) After the last 4 mL portion of eluant has passed through the cartridge move the switch of the positive pressure SPE manifold from “Regulated Flow” to “Full Flow/Dry” to dry the SPE cartridge for approximately 1 minute

- j) Using an N-Evap (or equivalent) with the water bath set at 50°C and N₂ flow set at <10 liters per minute (LPM) (typical setting is 2 – 6 LPM), evaporate the sample to approximately 0.5 mL
- k) Add 3 mL of toluene to the centrifuge tube containing the sample
- l) Evaporate again to < 0.5 mL. (This is to insure all other solvents have been removed from the sample.)
- m) Bring the volume to 1.0 mL with toluene and vortex to mix solvent into sample
- n) Analyze by GCMS-EI and GCMS-NCI



Analysis of Tobacco Alkaloids

UCT Part Number:

EUBCX1HL2Z (200 mg benzenesulfonic acid high load, 10 ml)

October 2012

Tobacco alkaloids are extracted with a strong cation exchange sorbent using an acidic buffer, filtered and further acidified. Extraneous compounds are removed by washing the sorbent, yielding clean chromatography without loss of target analytes.

Procedure

1. Sample Preparation

- a) Weigh 100 mg of finely ground tobacco in a screw cap vial
- b) Add 6 mL 0.1M sodium acetate buffer (pH 4.5) and 100 μ L IS (d4- nornicotine, 1 μ g/ μ L)
- c) Mix on rotating shaker for 10 minutes
- d) Filter extract through 20 micron frit filter column
- e) Add 300 μ L glacial acetic acid and mix

2. Sample Extraction

- a) Condition column **EUBCX1H2Z** with 3 mL of MeOH:1.0M acetic acid (80:20)
- b) Pour sample onto column, draw through at 1-2 mL/min
- c) Wash column with 3 mL of MeOH:1.0M acetic acid (80:20)
- d) Dry column for 5-10 min using full vacuum

3. Elution


- a) Elute alkaloids with 3 mL CH_2Cl_2 /isopropanol/ NH_4OH (70:26:4) by gravity
- b) Evaporate eluant to dryness with N_2 and low heat ($< 40^\circ \text{C}$)
- c) Reconstitute with 200 μ L ethyl acetate
- d) Analyze on GC/FID/NPD or GC/MSD

Results

Alkaloid	pKa	Flue Cured Tobacco	Burley Lamina Tobacco
	n=15/mean	mg/gram/CV	mg/gram/CV
Mysomine	NA	48/6.2	189/7.9
Nicotine	7.94	39406/6.2	39119/8.8
Nornicotine	9.46	1381/3.5	5429/5.2
Anatasine	9.20	229/5.2	183/8.7
Anatabine	8.23	1932/2.3	1774/2.3
2,3'-dipyridyl	4.25	54/8.9	30/11.2
Cotinine	4.88	20/11.2	52/12.4
Formylornicotine	NA	31/11.9	145/12.4

The average of 15 separate runs show excellent reproducibility

GC/MSD CHROMATOGRAMS



(1) nicotine, (2) nornicotine, (3) mysomine, (4) anabasine,
(5) anatabine, (6) 2,3'-dipyridyl, (7) cotinine, (8) formylornicotine

Instrument: Agilent 5890GC/5971MSD
GC column: Rtx-5 Amine, 30 m x 0.25 mm i.d. x 1.0 µm film
Injector: 1 µL sample at 10:1 split, 250° C
Temp program: Initial 120° C, hold 1 min, ramp 2.5° C/min to 200° C,
ramp 20° C/min to 280° C, hold 1 min.
MSD conditions: SIM monitoring, EI mode, 295° C

DCN-216111-60



Analysis of Glyphosate and Glufosinate by Solid-Phase Anion Exchange Extraction with GC/MS or LC/MS Analysis

UCT Part Number:

EUQAX2M6 (1000 mg, QAX2, 6 mL cartridge)

September 2009

Analyte	CASRN	Common Name
Glyphosate ¹	1071-83-6	Roundup®,
Glufosinate ²	51276-47-2	Basta®, Challenge®)

1. Sample Preparation

- Adjust water sample pH to 6 or higher with buffer

2. Cartridge Conditioning

- Place **EUQAX2M6** cartridges(s) on manifold
- Add 5 mL of methanol to the cartridge
- Slowly draw methanol through leaving enough to cover cartridge frit
- Rinse using 10 mL of pH 6 or higher buffer leaving a layer of buffer on frit

Note: Do not let the cartridge dry out after addition of methanol otherwise repeat

3. Extraction Protocol

- Draw a known volume of sample water through the cartridge, usually 100-500 mL

Note: Sample volume is determined by the analytical quantitation limit

- Adjust vacuum so that flow is approximately 1 - 3 mL per minute
- Wash sorbent using 10 mL of pH 6 buffer
- Dry the cartridge by drawing full vacuum for 10 minutes

4. Analyte Elution

- Elute using 5 mL of 1 mol/L HCl/methanol solution (4/1)
- Add eluant to the cartridge then draw through at 1 mL/minute
- Evaporate to dryness with N₂ flow in a water bath heated to 50°C

5. GC Analysis

- a) Add 50 μL of **MTBSTFA** (N-methyl-N-(tert-butyldimethylsilyl) trifluoroacetamide) and 50 μL of dimethylformamide for derivatization
- b) Sonicate at room temperature for 2 minutes (critical)
- c) Quantitatively transfer to GC vial and cap
- d) Heat to 80°C for a minimum of 30 minutes
- e) Cool to room temperature
- f) Sample is ready for analysis

6. LC Analysis (**Alternative Analysis Procedure**—no derivatization needed)

- a) After step 4) c) dissolve the dry residue using 100 μL of methanol
- b) Quantitatively transfer to an LC vial then cap
- c) Sample is ready for LC analysis

¹N-phosphonomethyl glycine

²RS-2-amino-4-(hydroxyl-methyl-phosphoryl)butanoic acid

DCN-218030-118



Extraction of Phenoxyacetic Acid Herbicides From Soil By LC-MS/MS

UCT Part Number:

EEC181M6 (endcapped C18 - 1000 mg/6 mL)

ECUNIC18 (Universal C18 1100mg/83 mL)

April 2009

1. Sample Pretreatment

- a) Prepare an acid washed beaker*
- b) Add 10-100 grams of soil sample
- c) Add enough DI H₂O to form a loose slurry
- d) Insert a magnetic stir bar and extract for 15 minutes
- e) Adjust pH to 2 using 50% aqueous sulfuric acid (H₂SO₄)
- f) Continue extraction for 15 minutes adjusting pH as needed
- g) Filter sample through previously acidified filter media

***Note:** Acid washed glassware must be used in this procedure. Soda lime glassware must be avoided as it may interfere with the analysis

2. Condition C18 SPE Cartridge

- a) Add 5 mL CH₃OH and wait 1 minute
- b) Add 5 mL DI H₂O

Note: Aspirate at low vacuum setting. Do not let cartridge dry out otherwise repeat steps a) and b)

3. Add Sample

- a) Adjust vacuum and load cartridge at 10 mL/minute flow rate

4. Dry Cartridge

- a) Dry cartridge for 10 minutes at full vacuum

5. Elute Phenoxyacetic acid Herbicides

- a) Place a clean collection vial in manifold
- b) Add 5 mL of CH₃OH and wait 1 minute
- c) Add a second 5 mL volume of CH₃OH
- d) Adjust vacuum and collect at 1-2 mL/ minute

6. Dry Eluate

- a) Evaporate to dryness at < 40°C using N₂
- b) Reconstitute in 100 µL of mobile phase for **LC-MS/MS**
 - Inject 10-100 µL

HPLC Analysis and Instrumentation Requirements

Guard Column: C18 10mm x 2.6mm with 0.5 µm frit

Analytical Column:

- C18 100 mm x 2 mm 5 µm particle ODS-Hypersil
- C18 100 mm x 2 mm 3 µm particle MOS2-Hypersil or equivalent

HPLC/MS Interface:

- Micromixer 10-µL interface HPLC column system with HPLC post-column addition solvent

Interface:

- Thermospray ionization interface and source capable of generating both positive and negative ions and have a discharge electrode or filament

Mass Spectrometer System:

- A single quadrupole mass spectrometer capable of scanning from 1 to 1000 amu
- Scanning from 150 to 450 amu in 1.5 sec. or less using 70 volts (nominal) in positive or negative electron modes
- Capable of producing a calibrated mass spectrum for polyethylene glycol (PEG 400, 600, or 800, average mol. wts.) or other compounds used as a calibrant
- Use PEG 400 for analysis of chlorinated phenoxyacid compounds. PEG is introduced via the Thermospray interface circumventing the HPLC

Thermospray Temperatures:

Vaporizer Control: 110°C to 130°C

Vaporizer Tip: 200°C to 215°C

Jet: 210°C to 220°C

Source Block: 230°C to 265°C

Recommended HPLC Chromatographic Conditions

Chlorinated Phenoxyacid Compounds A=0.1 M ammonium acetate/methanol

Initial Mobile Phase %	Initial Time minutes	Final minutes	Final Mobile Phase %	Time minutes
75A/25	2	15	40/60	
40A/60	3	5	75/25	10

Limits of Detection in the Positive and Negative Ion Modes for HPLC Analysis of Chlorinated Phenoxyacid Herbicides and Esters

Compound	Positive Ion Mode Quantitation LOD		Negative Ion Mode Quantitation LOD	
	Ion	ng	Ion	ng
Dalapon	Not detected		141 (M ⁻ H) ⁻	11
Dicamba	238 (M ⁺ NH ₄) ⁺	13	184 (M ⁻ HCl) ⁻	3.0
2,4-D	238 (M ⁺ NH ₄) ⁺	2.9	184 (M ⁻ HCl) ⁻	50
MCPA	218 (M ⁺ NH ₄) ⁺	120	199 (M ⁻ 1) ⁻	28
Dichloroprop	252 (M ⁺ NH ₄) ⁺	2.7	235 (M ⁻ 1) ⁻	25
MCP	232 (M ⁺ NH ₄) ⁺	5.0	213 (M ⁻ 1) ⁻	12
2,4,5- T	272 (M ⁺ NH ₄) ⁺	170	218 (M ⁻ HCl) ⁻	6.5
2,4,5-TP Silvex	286 (M ⁺ NH ₄) ⁺	160	269 (M ⁻ 1) ⁻	43
Dinoseb	228 (M ⁺ NH ₄ ⁺ NO) ⁺	24	240 (M ⁻)	19
2,4-DB	266 (M ⁺ NH ₄) ⁺	3.4	247 (M ⁻ 1) ⁻	110
2,4,5-D, butoxy ethanol ester	321 (M ⁺ H) ⁺	1.4	185 (M ⁻ C ₆ H ₁₃ O ₁) ⁻	
2,4,5-T,butoxy ethanol ester	372 (M ⁺ NH ₄) ⁺	0.6	195 (M ⁻ C ₈ H ₁₅ O ₃) ⁻	
2,4,5-T, butyl ester	328 (M ⁺ NH ₄) ⁺	8.6	195 (M ⁻ C ₆ H ₁₁ O ₂) ⁻	
2,4-D, ethyl hexyl ester	350 (M ⁺ NH ₄) ⁺	1.2	161 (M ⁻ C ₁₀ H ₁₉ O ₃) ⁻	



QuEChERS Extraction and Clean-Up of Pesticides from Olive Oil

UCT Part Number:

CUMPS2CT (150 mg anhydrous MgSO_4 & 50 mg PSA)

April 2009

1. Sample Extraction

- a) In a suitable vial, add 1.5 mL of olive oil
- b) Add 1.5 mL of hexane
- c) Add 6 mL of acetonitrile
- d) Shake for 30 minutes
- e) Allow layers to phase separate for 20 minutes
- f) Collect acetonitrile layer (top layer)
- g) Repeat steps c) through f) and combine acetonitrile layers

2. Sample Clean-up

- a) Add 1 ml of combined acetonitrile to CUMPS2CT
- b) Shake for 2 minutes by hand
- c) Centrifuge at 3000 rpm for 2 minutes
- d) Remove solvent layer
- e) Analyze by HPLC using MS detection



QuEChERS Multiresidue Pesticide Method for the Determination of Multiple Pesticides in Wines*

UCT Part Number:

ECQUVIN50CT (50 mL centrifuge tube, 8.0 grams anhydrous MgSO_4 & 2 grams NaCl)

ECMPSCB15CT (900 mg anhydrous MgSO_4 , 300 mg PSA & 150 mg GCB)

February 2010

This method summary describes a multi-residue pesticide method for the determination of 72 pesticides in wines. Pesticides are extracted using acetonitrile saturated with magnesium sulfate and sodium chloride followed by a dispersive solid-phase clean-up with primary-secondary amine (PSA) and graphitized carbon black (GCB) sorbents.

Analysis is performed using ultra-performance liquid chromatography-tandem mass spectrometry (UPLC-MS/MS) in ESI (positive) mode.

1. Sample Preparation

- a) Add 20 mL acetonitrile (ACN) and *internal standard* 250 μL Fluconazole (10 $\mu\text{g/L}$) to **ECQUVIN50CT**
- b) Quantitatively add 20.0 mL of wine
- c) Shake for approximately 2 minutes
- d) Centrifuge at 4500 rpm for 5 minutes (use refrigerated centrifuge if available)
- e) Transfer 9.0 mL of top layer and add to **ECMPSCB15CT** (900 mg anhydrous MgSO_4 , 300 mg PSA & 150 mg GCB)
- f) Vortex tube for approximately 10 seconds
- g) Open tube and add 3.0 mL of toluene and shake for 1 minute
- h) Centrifuge the tube for 5 minutes @ 4500 rpm
- i) Quantitatively transfer 2.0 mL of supernatant to a glass centrifuge tube
- j) Evaporate to dryness at $< 40^\circ\text{C}$ using N_2
- k) Add 500 μL of ACN and 25 μL of *surrogate standard* (benzanilide - 20.0 $\mu\text{g/L}$) for QC and 500 μL of 20 mM ammonium acetate in 1% ACN to the dried extract
- l) Vortex for approximately 5 seconds and filter into autosampler vial using 17mm, 0.2 μm nylon membrane cartridges attached to a disposable syringe

2. UPLC/MS/MS Analysis

UPLC Conditions:

Column: Water's Acquity UPLC BEH C₁₈ column 100 x 2.1 mm, 1.7 µm particle or equivalent

Flowrate: 0.2 mL/minute

Injection volume: 3 µL

Analytical Standards: Matrix Matched

Gradient Program:

Time	% Acetonitrile	% 10mM Ammonium Acetate
0	10	90
10	90	10
14.5	90	10
14.6	10	90
20.1	10	90

Triple Quadrupole MS Conditions--electrospray ionization mode (ESI+)

Capillary Voltage: 1.5 kV

Source Temperature: 120 °C

N₂ Flow: cone 50 L/h, desolvation 800 L/h

Collision Gas: Argon

Dwell Time: 10 µS for multiple reaction monitoring (MRM) experiments

Collision Cell Pressure: 5.9 x 10⁻³ mbar

Summary of MS/MS Conditions

Pesticide	Molecular Weight	CV (V)	Quantification Transition
Acephate	183.17	20	184.0→143.0
Acetamiprid	222.67	30	223.4→126.1
Acibenzolar S-methyl	210.27	35	211.1→136.0
Aldicarb	190.27	12	208.1→116.0
Aldicarb sulfone	222.27	15	240.0→222.9
Aldicarb sulfoxide	206.26	15	224.2→206.9
Atrazine	215.69	35	215.9→173.85
Avermectin B _{1b}	873.09	20	876.6→553.4
Avermectin B _{1a}	873.09	20	890.7→567.5
Azoxystrobin	403.30	25	404.0→372.1
Benalaxyl	325.41	26	326.1→148.1
Benfuracarb	410.53	20	411.2→190.0
Benzanilide	197.24	30	198.1→105.1
Bifenazate	300.35	20	301.3→170.2
Bitertanol	337.42	20	338.2→ 99.1
Buprofezin	305.44	25	306.3→201.2
Carbaryl	201.22	22	202.1→145.1
Carbendazim	191.19	30	192.0→160.0
Carbofuran	221.26	26	222.1→123.1
Chloroxuron	290.75	35	291.0→ 72.2
Cyprodinil	225.29	45	226.1→ 93.0
Cyromazine	166.19	25	167.2→ 85.1
Diclobutrazol	328.24	30	328.1→70.2
Dimethoate	229.26	20	230.1→199.0
Dimethomorph	387.86	35	388.0→301.1
Dimoxystrobin	326.39	20	327.1→206
Dinotefuran	202.20	20	203.5→14.0
Diuron	233.10	30	233.0→72.1
Ethofumesate	286.35	30	286.9→258.9
Famoxadone	374.39	-32	373.2→ 282
Fenamidone	311.40	25	312.2→236.2
Fenbuconazole	336.82	35	337.1→125.0
Fenhexamid	302.20	65	301.9→261.9
Fenpropimorph	304.49	40	304.4→147.1
Fluconazole	306.27	30	307.2→220
Fludioxinil	248.19	-45	247.0→180.0
Furathiocarb	382.48	30	383.2→195.1
Hexaconazole	314.21	35	314.0→ 70.2
Imazalil	297.18	35	297.1→159.0
Imidacloprid	255.65	25	256.1→175.0
Ipconazole	333.86	35	334.1→70.2

Iprovalicarb	320.43	24	321.2→119.0
Kresoxim-methyl	313.35	20	314.1→116.0
Mepanipyrin	223.28	30	224.4→ 77.3
Metalaxyl	279.34	25	280.1→220.1
Methamidophos	141.13	22	142.0→ 94.0
Methomyl	162.21	20	163.0→ 88.0
Methoxyfenozide	368.47	15	369.5→149.0
Mevinphos	224.15	22	225.1→192.8
Myclobutanil	288.78	35	289.1→ 70.2
Omethoate	213.14	20	214.1→183.0
Oxadixyl	278.31	20	279.1→219.1
Piperonyl butoxide	338.45	17	356.2→177.0
Prochloraz	376.67	20	376.1→308.0
Propamocarb	188.27	30	189.1→102.1
Propargite	350.48	20	368.1→231.0
Propiconazole	342.22	35	342.0→159.0
Propoxur	209.24	20	210.0→111.0
Pyraclostrobin	387.83	23	388.0→194.0
Pyridaben	364.94	22	365.3→309.1
Pyrimethanil	199.25	40	200.1→107.0
Quinoxifen	308.14	50	307.8→196.8
Rotenone	394.42	40	395.3→213.2
Simazine	201.66	30	202.2→131.4
Spinosyn A	731.97	40	732.6→142.2
Spinosyn D	746.00	30	746.6→142.2
Spiroxamine	297.48	30	298.2→144.0
Tebuconazole	307.82	30	308.2→ 70.2
Thiabendazole	201.25	35	202.0→175.0
Triadmimefon	293.75	30	294.0→197.1
Trifloxystrobin	408.38	25	409.0→186.0
Triflumizole	345.75	20	346.0→278.1
Vamidotion	287.34	20	288.1→146.0
Zoxamide	336.54	35	336.0→187.0

Schematic Diagram of Sample Preparation Steps




Table of Average Pesticide Recoveries at 100 g/L Spike
(average values with SD, n=4)

	Pesticide Recovery	
	Red Wine @ 100 µg/L	White Wine @ 100 µg/L
Acephate	84±4	79±3
Acetamiprid	83±8	97±7
Acibenzolar S-methyl	80±15	45±5
Aldicarb	92±5	82±5
Aldicarb sulfone	91±7	83±4
Aldicarb sulfoxide	83±8	80±1
Atrazine	92±5	83±5
Avermectin B _{1b}	94±12	107±13
Avermectin B _{1a}	82±8	80±6
Azoxystrobin	93±5	86±4
Benalaxyl	92±5	84±4
Benfuracarb	ND	ND
Benzanilide	69±7	70±8
Bifenazate	86±4	86±11
Bitertanol	92±5	86±4
Buprofezin	91±4	88±6
Carbaryl	77±4	76±4
Carbendazim	126±7	106±7
Carbofuran	90±4	86±4
Chloroxuron	75±5	72±2
Cyprodinil	38±2	56±5
Cyromazine	89±6	83±5
Diclobutrazol	89±6	82±4
Dimethoate	88±7	84±4
Dimethomorph	95±5	85±4
Dimoxystrobin	85±5	74±6
Dinotefuran	88±4	78±5
Diuron	74±12	90±1
Ethofumesate	92±10	95±14
Famoxadone	87±3	86±5
Fenamidone	88±5	80±5
Fenbuconazole	133±21	90±11
Fenhexamid	91±5	83±4
Fenpropimorph	86±6	84±4
Fluconazole	112±4	101±2
Fludioxinil	91±2	87±8
Furathiocarb	81±4	77±4
Hexaconazole	90±2	77±7
Imazalil	89±5	83±5

Imidacloprid	94±6	87±4
Ipconazole	89±5	83±5
Iprovalicarb	94±6	87±4
Kresoxim-methyl	85±5	86±5
Mepanipyrin	76±6	94±12
Metalaxyl	94±5	85±5
Methamidophos	82±6	74±5
Methomyl	90±4	81±4
Methoxyfenozide	102±5	89±5
Mevinphos	84±5	71±4
Myclobutanil	96±8	90±4
Omethoate	82±4	75±4
Oxadixyl	94±3	88±4
Piperonyl butoxide	94±5	87±4
Prochloraz	84±3	84±5
Propamocarb	80±3	80±5
Propargite	93±6	86±2
Propiconazole	94±4	86±5
Propoxur	89±5	82±4
Pyraclostrobin	77±6	76±4
Pyridaben	85±4	83±4
Pyrimethanil	79±6	75±4
Quinoxifen	70±5	68±3
Rotenone	81±3	85±9
Simazine	85±9	88±7
Spinosyn A	88±7	83±4
Spinosyn D	87±4	80±3
Spiroxamine	92±5	84±4
Tebuconazole	90±4	83±5
Thiabendazole	71±3	75±5
Triadmimefon	89±8	84±7
Trifloxystrobin	90±8	84±4
Triflumizole	88±6	86±3
Vamidothion	86±4	83±6
Zoxamide	86±4	80±4

*Adapted from Kai Zhang, Jon W. Wong et al, Multiresidue Pesticide Analysis of Wines by Dispersive Solid-phase Extraction and Ultra-High Performance Liquid Chromatography-Tandem Mass Spectrometry *Journal of Agricultural and Food Chemistry*



Determination of Pesticides in Red Wine by QuEChERS Extraction, Quick QuEChERS Clean-up, and LC/MS/MS Detection

UCT Part Numbers:

RFV0050CT (50 mL polypropylene centrifuge tube)

ECQUUS2-MP (Mylar Pouch contains: 4000 mg MgSO₄, 2000 mg NaCl)

ECPURMPSCMC (Quick QuEChERS cartridge, 110 mg MgSO₄, 180 mg PSA)

March 12012

The analysis of pesticide residues in red wines is challenging due to the complexity of the matrix, which contains organic acids, sugars, phenols, and pigments, such as anthocyanins. A simple, faster, and easy to use method is developed for the determination of pesticide residues in red wines.

Eight pesticides with a wide range of polarities (LogP from -0.779 to 5.004) were selected as target analytes. Excellent accuracy and precision data were achieved using this method. Recoveries of planar pesticides, such as Carbendazim and Thiabendazole were not affected since PSA was used for clean-up instead of GCB. PSA removed organic acids, sugars and pigments from the red wine extract. Six red wine samples were extracted using this method. Cyprodinil and Carbendazim were detected in the red wine samples tested, with minimum reporting limits of 1.5 ng/mL.

Procedure

1. Extraction

- a) Add 10 mL of red wine sample to a 50 mL polypropylene centrifuge tube
(RFV0050CT)
- b) Spike with the appropriate amount of target analytes for fortified samples
- c) Vortex 30 sec, then equilibrate for 15 min
- d) Add 10 mL of acetonitrile, vortex 30 sec
- e) Add salts in Mylar pouch **(ECQUUS2-MP)**
- f) Shake vigorously for 1 min
- g) Centrifuge at 5000 rpm for 5 min at 20° C
- h) Supernatant is ready for clean-up

2. Quick QuEChERS Clean-up


- a) Draw 1 mL of supernatant into a disposable polypropylene syringe
- b) Pass the supernatant slowly through the Quick QuEChERS cartridge
(ECPURMPSCMC)

- c) Collect 0.5 mL of the cleaned extract into a 2 mL auto-sampler vial
- d) Add 10 μ L 5 ppm TPP as internal standard (IS)
- e) Samples are ready for LC/MS/MS analysis

Clean-up red wine extract with Quick QuEChERS



Quick QuEChERS before (left) and after (right) clean-up of 1 mL red wine extract



3. LC/MS/MS Detection

LC: Thermo Accela 1250 pump with PAL auto-sampler

LC Conditions

Column	Guard column: Restek C18, 2.1 x 20 mm Column: Sepax HP-C18, 2.1 x 100 mm, 3 μ m, 120 Å
Column Temperature	Ambient
Injection Volume	10 μ L at 15° C
Mobile Phase	A: 0.1% formic acid in Milli-Q-water B: 0.1% formic acid in methanol
Flow Rate	200 μ L/min

LC Gradient Program

Time	%A	%B
0	95	5
1	95	5
3	50	50
8	5	95
14.2	95	5
16	95	5

MS/MS: Thermo TSQ Vantage tandem MS

MS Conditions

Ion source:	Heated ESI
Ion polarity:	ESI +
Spray voltage:	3000 V
Sheath gas pressure:	N ₂ @ 40 psi
Auxiliary gas pressure:	N ₂ @ 10 psi
Ion transfer capillary temperature:	350 °C
Scan type:	SRM (0-16 min)
CID conditions:	Ar @ 1.5 mTorr

SRM transitions

Compound	Parent	Product ion 1	CE	Product ion 2	CE	S-Lens	Dwell time (s)
Methamidophos	142.044	94.090	14	125.050	16	59	0.15
Carbendazim	192.093	132.080	29	160.080	17	81	0.10
Thiabendazole	202.059	131.060	31	175.070	31	103	0.10
Pyrimethanil	200.116	107.060	23	183.140	22	66	0.10

Cyprodinil	226.122	77.030	40	93.050	33	88	0.10
TPP (IS)	327.093	77.020	37	152.070	33	98	0.10
Diazinon	305.135	153.090	15	169.08	14	89	0.10
Pyrazophos	374.103	194.060	20	222.130	20	104	0.10
Chlorpyrifos	349.989	96.890	32	197.940	17	69	0.10

Matrix matched calibration, LOD and LOQ


Compound	Linearity range (ng/mL)	R ²	LOD (ng/mL)	LOQ (ng/mL)
Methamidophos	2-400	0.9991	0.15	0.49
Carbendazim	2-400	0.9981	0.40	1.33
Thiabendazole	2-400	0.9940	0.09	0.31
Pyrimethanil	2-400	0.9990	0.01	0.05
Cyprodinil	2-400	0.9995	0.17	0.57
Diazinon	2-400	0.9982	0.06	0.21
Pyrazophos	2-400	0.9976	0.08	0.27
Chlorpyrifos	2-400	0.9981	0.10	0.32

Accuracy and Precision Data

Compound	Fortified at 10 ng/mL		Fortified at 50 ng/mL		Fortified at 100 ng/mL	
	Recovery%	RSD% (n=4)	Recovery%	RSD% (n=4)	Recovery%	RSD% (n=4)
Methamidophos	93.7	3.4	81.6	5.8	84.2	3.5
Carbendazim	105.7	10.8	100.1	10.6	90.5	7.6
Thiabendazole	91.2	4.9	87.9	6.8	85.0	4.0
Pyrimethanil	112.2	2.7	107.0	3.2	102.8	4.9
Cyprodinil	104.3	3.2	99.9	6.1	100.2	4.9
Diazinon	104.9	5.6	102.0	6.6	99.2	6.8
Pyrazophos	99.9	4.0	96.6	5.6	91.3	4.1

Chlorpyrifos	91.7	4.6	99.5	5.2	97.2	3.8
--------------	------	-----	------	-----	------	-----

Chromatograms of Red Wine Sample Fortified With 10 ng/mL Pesticides



Pesticides detected in red wine samples (ng/mL)

Pesticide	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	Sample 6
Methamidophos	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5
Carbendazim	< 1.5	< 1.5	< 1.5	10.2	8.7	2.3
Thiabendazole	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5
Pyrimethanil	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5
Cyprodinil	1.7	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5
Diazinon	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5
Pyrazophos	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5
Chlorpyrifos	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5	< 1.5



Analysis of Flukicides/ Anthelmintics in Animal Tissue Using QuEChERS and LC/MS/MS

UCT Part Number:

ECMSSC50CT-MP (4000 mg MgSO₄, 1000 mg NaCl)

ECMSC1850CT (1500 mg MgSO₄, 500 mg endcapped C18)

Revised March 27 2013

Procedure

1. Extraction

- a) To 10 g of homogenized/hydrated sample in a 50 mL centrifuge tube add 10 mL acetonitrile
- b) Add internal standard (100 ng/g Cyprodinil + 2,4-D)

Note: Isotopically labeled internal standards are now commercially available.

- c) Shake for 1 minute
- d) Add contents of ECMSSC50CT-MP pouch (4 g of anhydrous magnesium sulfate and 1g sodium chloride) to the centrifuge tube
- e) Immediately shake for 1 minute
- f) Centrifuge for 5 minutes at 3450 rcf

2. Sample Clean-Up

- a) Add a 3 mL aliquot of supernatant from (step 1 f) to **ECMSC1850CT**
- b) Shake for 1 minute
- c) Centrifuge for 1 minute at 3450 rcf

3. Analysis

- a) Transfer 0.5 mL of cleaned extract into a autosampler vial
- b) Add QC spike (100 ng/mL TPP).

Note: Isotopically labeled internal standards are now commercially available.

- c) Inject onto LC-MS/MS
- d) Use ESI⁺ and/or ESI⁻ mode depending upon specific analytes of interest

Note:

Abamectin, doramectin and ivermectin form sodium adducts ($[M+23]^+$) when acids are used as mobile phase additive in MS analysis. It is advisable to use ammonium formate or ammonium acetate as mobile phase buffer and monitor the ammonium adduct ($[M+18]^+$) for these three compounds. It is essential to use ammonium buffer in the organic mobile phase as the avermectins elute at 100% organic content. In addition, ammonium formate is more soluble in organic solvent than ammonium acetate.

MS amenable acids can be used for the aqueous mobile phase, which should be at a low pH (≤ 4) to get the best results. The aqueous mobile phase may also include ammonium buffer, although it is not an essential requirement.

Albendazole-sulfone and hydroxy-mebendazole are prone to isobaric interference as they have similar precursor and product ions that can't be distinguished using triple quadrupole instruments. It is therefore necessary to chromatographically separate these two compounds.

39 Flukicides/Anthelmintics

ESI+		ESI-
ISTD Triphenylphosphate	QC Spike Cyprodinil	ISTD 2,4D
Abamectin	Albendazole	Bithionol
Doramectin	Albendazole-sulfoxide	Clorsulon
Emamectin	Albendazole-sulfone	Closantel
Eprinomectin	Albendazole-amino-sulfone	Niclosamide
Moxidectin	Cambendazole	Nitroxynil
Ivermectin	Flubendazole	Oxyclozanide
Selamectin	Flubendazole, amino	Rafoxanide
Dichlorvos	Flubendazole, hydroxy	Triacclabendazole-sulfoxide
Coumaphos	Mebendazole	
Coumaphos-oxon	Mebendazole, amino	Albendazole-amino-sulfone
Haloxon	Mebendazole, hydroxy	
Morantel	Oxibendazole	
Levamisole	Thiabendazole	
Fenbendazole	Thiabendazole, 5-hydroxy	
Fenbendazole-sulfone	Triclabendazole	
Fenbendazole-sulfoxide (oxfendazole)		

Adapted from Kinsella, Lehotay et al, "New method for the Analysis of Anthelmintics in Animal Tissue"



Streamlined Method for the Determination of More Than 100 Veterinary Drugs in Animal Tissue Using Dispersive-SPE Clean-up and LC- MS/MS Detection

UCT Part Numbers:

ECMSC1850CT (500 mg C18 and 1500 mg MgSO₄)

August 2012

More than 100 veterinary drugs can be extracted and analyzed using this fast and easy multi-class, multi-residue method.

Procedure

1. Extraction

- a) Weigh 2 g of homogenized tissue sample into a 50 mL centrifuge tube
- b) Add 10 mL of MeCN/water (4/1 v/v)
- c) Shake or vortex for 5 min
- d) Centrifuge for 5 min at >3700 rcf
- e) Supernatant is ready for clean-up

2. dSPE Clean-up

- a) Transfer the supernatant into product **ECMSC1850CT**
- b) Add 10 mL hexane that has been pre-saturated with MeCN
- c) Shake or vortex for 30 sec
- d) Centrifuge for 5 min at >3700 rcf
- e) Aspirate hexane to waste
- f) Evaporate 5 mL of the extract under nitrogen at 45°C to <0.7 mL
- g) Add 0.1% formic acid to reach a final volume of 1 mL (1 g/mL sample equivalent)
- h) Transfer sample to a HPLC vial (filter with PVDF if desired)
- i) Sample is ready for LC-MS/MS analysis

3. LC-MS/MS analysis

MS: Waters TQD

HPLC: Waters Acquity UHPLC

LC Parameters

Guard	Agilent Eclipse Plus C18, 5 µm, 4.6 x 12.5 mm
Column	Waters Acquity HSS T3 (C18), 1.8 µm, 2.1 x 100 mm
Flow	0.5 ml/min
Injection Vol	20 µl
Oven	40 °C
Equilibration Time	3.3 min
Autosampler Temperature	4 °C

Mobile Phase

Aqueous A: 0.1 % formic acid in water/MeCN (95/5 v/v)

Organic B: 0.1 % formic acid in MeCN

Time	%B
0	0.2
0.1	0.2
8.0	99.8
9.5	99.8
9.6	0.2
13	0.2

For analysis of late eluting compounds, 50 µL/min of 27 mM ammonium formate in MeOH:MeCN (75:25) is infused from 5.05 to 9.45 min using the instrument's infusion syringe to enhance the signal of the late-eluting anthelmintics.

MS Instrument Settings

Capillary voltage	3000 V
Extractor voltage	3 V
Desolvation temperature	450°C
Source temperature	150°C
Dwell time	5 msec

Analyte	Drug class	RT (min)	Precursor ion	Cone V	Product 1	Collision energy (V)	Product 2	Collision energy (V)
Desacetyl Cephapirin	β-Lactam	0.69	382.1	32	152	28	124.2	48
Florfenicol Amine	Phenicol	0.68	248.1	25	230.2	10	130.1	35
Sulfanilamide	Sulfonamide	1.19	173	40	92.9	20	75.9	36
Amoxicillin	β-Lactam	1.47	366.1	20	114	22	349.3	10
Salbutamol	β -Agonist	1.46	240.2	20	148.2	20	222.3	10

Zilpaterol	β -Agonist	1.46	262.3	27	244.3	12	185.2	30
Cimaterol	β -Agonist	1.51	220	16	143	24	115.9	34
DCCD	β -Lactam	1.72	549.1	40	183	30	241.1	20
Lincomycin	Lincosamide	1.87	407.3	20	126.1	30	359.2	20
Sulfadiazine	Sulfonamide	2	251.1	30	156.1	15	108	20
Ampicillin	β -Lactam	2.01	350.1	26	106.1	24	114	30
Desethylene Ciprofloxacin	Fluoroquinolone	2.06	306.2	35	288.2	20	245.2	20
Sulfathiazole	Sulfonamide	2.1	256.1	25	156.1	15	108	25
Sulfapyridine	Sulfonamide	2.18	250.1	32	156.1	18	108.1	28
Norfloxacin	Fluoroquinolone	2.16	320.2	36	276.2	18	233.1	26
Tulathromycin	Macrolide	2.17	806.8	38	72	56	577.5	24
Oxytetracycline	Tetracycline	2.21	461.2	25	426.4	20	443.4	15
Ciprofloxacin	Fluoroquinolone	2.22	332.2	35	245.2	25	288.4	20
Ractopamine	β -Agonist	2.27	302.2	26	164	16	107	32
Sulfamerazine	Sulfonamide	2.3	265.1	28	91.9	28	155.9	16
Danofloxacin	Fluoroquinolone	2.31	358.1	28	96	26	314.2	18
Tetracycline	Tetracycline	2.35	445.2	30	154.1	30	410.2	20
Enrofloxacin	Fluoroquinolone	2.38	360.2	35	316.4	20	245.3	25
2-Quinoxalinecarboxylic Acid	Other	2.43	175	22	129	16	131	16
Sulfamethizole	Sulfonamide	2.55	271.1	28	156.1	16	92	30
Sulfamethazine	Sulfonamide	2.54	279.1	35	186.1	20	156.1	20
Sulfamethazine-13C6 (IS)		2.54	285.2	32	186.1	18	124.1	26
Cefazolin	Cephalosporin	2.56	455.1	20	156	16	323.2	12
Sulfamethoxypyridazine	Sulfonamide	2.58	281.1	30	156.1	20	126.2	20
Difloxacin	β -Lactam	2.62	400.3	35	356.4	20	299.2	30
Sarafloxacin	Fluoroquinolone	2.58	386.1	20	342.2	20	299.2	30
Clenbuterol	β -Agonist	2.56	277.2	25	259.2	10	132.1	30
Pirlimycin	Lincosamide	2.74	411.3	30	112.2	40	363.3	20
Chlortetracycline	Tetracycline	2.84	479.2	30	154.1	30	444.3	20
Clindamycin	Lincosamide	2.89	425.3	45	126.2	40	377.4	20
Gamithromycin	Macrolide	2.91	777.8	62	83	54	116	50
Sulfachloropyridazine	Sulfonamide	2.95	285	28	156.1	16	108	26
Tilmicosin	Macrolide	3.06	869.8	45	174.2	35	696.6	35
Sulfadoxine	Sulfonamide	3.1	311.2	35	156.1	20	108.1	30
Sulfamethoxazole	Sulfonamide	3.11	254	26	92.1	30	156	18
Sulfaethoxypyridazine	Sulfonamide	3.14	295.1	30	156.1	20	140.2	20
Florfenicol	Phenicol	3.15	358.1	24	241	18	206	28
Chloramphenicol	Phenicol	3.36	323.1	16	275	16	165	26
Erythromycin	Macrolide	3.49	734.8	30	158.2	36	115.9	54
Sulfadimethoxine	Sulfonamide	3.57	311.1	35	156.1	20	108	30
Sulfaquinoxaline	Sulfonamide	3.59	301.1	34	156.1	18	108	28
Prednisone	Corticosteroid	3.67	359.2	22	341.1	10	146.9	26
Tylosin	Macrolide	3.66	916.8	45	174.2	35	101.1	35
Penicillin G-d7 (IS)		3.86	342.1	46	183.1	26	160.1	24
Penicillin G	β -Lactam	3.86	335.1	18	176	16	160.1	18
Beta/Dexa-methasone	Corticosteroid	4.11	393.2	20	373.2	10	147.1	28

Sulfanitran	Sulfonamide	4.16	336.2	26	156	12	134.1	28
Sulfabromomethazine	Sulfonamide	4.21	357.1	35	92	30	156.1	25
Zeranol (±-Zearalanol)	Other	4.37	323.2	16	305.2	10	189.1	24
Oxacillin	β -Lactam	4.39	402.1	22	160	20	243.1	18
Atrazine (QC)		4.49	216.1	34	174	18	103.9	30
Cloxacillin	β -Lactam	4.66	436.2	22	160.1	12	277.1	16
Nafcillin	β -Lactam	4.79	415.2	20	199.1	14	171.1	38
Oxyphenylbutazone	NSAID	4.83	325.2	26	120.1	24	148.2	30
Flunixin	NSAID	4.86	297.1	42	279.1	22	109	50
Flunixin-d3 (IS)		4.82	300.1	40	282.1	24	112	54
Dicloxacillin	β -Lactam	5.03	470.2	22	160.1	14	311.1	16
Phenylbutazone	NSAID	5.93	309.1	28	120	20	91.8	30
Melengesterol Acetate	Other	6.3	397.4	30	279.3	20	337.5	15
2-thiouracil	Thyreostat	0.85	128.9	32	111.9	12	69.9	18
2-mercapto-1-	Thyreostat	1.14	114.9	40	87.9	16	73.9	16
6-methyl-2-thiouracil	Thyreostat	1.22	142.9	32	83.9	18	125.9	14
Metronidazole-OH	Nitroimidazole	1.42	188	22	123	14	126	18
Dipyrone	Tranquilizer	1.6	218.1	24	96.9	12	187	10
Dimetridazole-OH	Nitroimidazole	1.63	158	22	140	12	93.9	22
Metronidazole	Nitroimidazole	1.63	172	26	127.9	14	81.9	24
5-hydroxythiabendazole	Anthelmintic	1.7	218	50	190.9	26	147	32
Albendazole 2-amino-sulfone	Anthelmintic	1.85	240	36	133	28	198	20
Ronidazole	Nitroimidazole	1.85	201	18	139.9	10	54.8	20
Levamisole	Anthelmintic	1.86	205	40	178	22	90.9	34
Dimetridazole	Nitroimidazole	1.86	142	32	95.9	16	80.9	24
Thiabendazole	Anthelmintic	1.94	202	44	174.9	26	130.9	32
6-propyl-2-thiouracil	Thyreostat	2.15	171	38	154	18	112	20
2-mercaptobenzimidazole	Thyreostat	2.3	150.9	42	92.8	20	118	22
Azaperone	Tranquilizer	2.34	328.3	34	165	20	122.9	36
Orbifloxacin	Fluoroquinolone	2.39	396.2	36	352.2	18	295.1	24
Albendazole sulfoxide	Anthelmintic	2.44	282.1	28	240	14	207.3	24
Xylazine	Tranquilizer	2.48	221.1	42	164	26	147	24
Iprnidazole-OH	Nitroimidazole	2.54	186.1	22	168	14	121.8	20
Morantel	Anthelmintic	2.6	221.1	50	122.9	36	163.9	28
2-amino-Mebendazole	Anthelmintic	2.63	238.1	50	104.9	26	132.9	36
6-phenyl-2-thiouracil	Thyreostat	2.73	205	38	103	26	187.9	18
2-amino-Flubendazole	Anthelmintic	2.77	256	50	122.9	28	94.9	38
Cambendazole	Anthelmintic	2.83	303.1	34	261.1	18	217	28
Bacitracin	Other	2.87	475.3	26	85.9	24	199.1	30
Carazolol	Tranquilizer	2.9	299.3	34	116	20	97.9	22
Doxycycline	Tetracycline	2.91	445.3	28	428.2	20	97.9	46
Oxibendazole	Anthelmintic	2.95	250.1	34	218.1	18	175.9	28
Oxfendazole	Anthelmintic	3.01	316.1	40	158.9	32	191	22
Albendazole sulfone	Anthelmintic	3.02	298.1	38	266	20	159	36
Iprnidazole	Nitroimidazole	3.2	170.1	34	124	18	109	24
Clorsulon	Flukicide	3.39	377.7	24	341.8	12	241.9	20

Haloperidol	Tranquilizer	3.53	376.2	40	165	24	122.9	42
Acetopromazine	Tranquilizer	3.55	327.2	32	86	20	254	22
Promethazine	Tranquilizer	3.58	285.2	24	85.9	16	198	20
Fenbendazole sulfone	Anthelmintic	3.65	332.1	40	300	22	158.9	38
Albendazole	Anthelmintic	3.65	266.1	34	234	20	191.1	32
Mebendazole	Anthelmintic	3.7	296.1	36	264.1	20	104.9	36
Flubendazole	Anthelmintic	3.9	314.1	38	282	22	94.9	50
Propionylpromazine	Tranquilizer	3.91	341.2	32	85.9	22	268.1	24
Chlorpromazine	Tranquilizer	4.04	319.2	32	86	20	246	22
Triflupromazine	Tranquilizer	4.26	353.2	34	85.9	22	280	28
Fenbendazole	Anthelmintic	4.33	300.1	38	268	20	158.9	36
Oleandomycin triacetate	Macrolid	4.37	814.7	38	200.1	30	98	48
Nitroxylin	Flukicide	4.41	288.8	40	126.8	20	115.9	34
Virginiamycin M1	Other	4.49	526.4	26	508.3	12	108.9	44
Ketoprofen	Tranquilizer	4.71	255.1	28	104.9	24	209	14
Haloxon	Anthelmintic	5.28	415.1	44	272.9	34	210.9	36
Triclabendazole sulfoxide	Flukicide	5.37	372.8	36	357.8	18	212.9	30
Emamectin benzoate	Anthelmintic	5.49	886.8	52	158	40	126	46
Diclofenac	Tranquilizer	5.55	296	20	214.9	20	250	12
Triclabendazole	Flukicide	5.99	359	52	343.9	26	274	38
Novobiocin	Other	6.05	613.5	30	189	28	132.9	64
Oxyclozanide	Flukicide	6.08	399.6	38	363.8	14	175.9	24
Niclosamide	Flukicide	6.2	325	36	170.9	30	289	16
Tolfenamic acid	Tranquilizer	6.23	262.1	22	244	14	180	40
Bithionol	Flukicide	6.76	352.9	36	160.8	24	191.2	28
Eprinomectin	Anthelmintic	7.44	914.8	18	186.1	20	154	40
Abamectin	Anthelmintic	7.94	890.8	16	305.3	28	145	42
Closantel	Flukicide	8.07	660.9	70	126.8	54	344.8	32
Doramectin	Anthelmintic	8.3	916.9	22	331.3	26	113	56
Moxidectin	Anthelmintic	8.32	640.5	16	528.3	8	498.2	10
Rafoxanide	Flukicide	8.5	623.9	62	126.1	48	344.8	30
Selamectin	Anthelmintic	8.62	770.7	36	145	30	112.9	40
Ivermectin [M+Na] ⁺	Anthelmintic	8.77	897.8	82	183	58	329.2	56

Accuracy and Precision

A multi-day, multi-analyst validation demonstrated that the final method is suitable for screening of 113 analytes, identifying 98 and quantifying 87 out of the 127 tested drugs at or below US regulatory tolerance levels in bovine muscle. Overall, the method demonstrated reasonably good quantitative performance with recoveries ranging between 70–120% for 87 out of 127 analytes, and recovery of < 50% for only 20 analytes. 85 analytes gave RSDs ≤ 20% and 100 analytes gave RSDs ≤ 25%.

Adapted from L. Geis-Asteggianti, S.J. Lehotay, A.R. Lightfield, T. Dutko, C. Ng, L. Bluhm, Ruggedness testing and validation of a practical analytical method for >100 veterinary drug residues in bovine muscle by ultrahigh performance liquid chromatography tandem mass spectrometry, *Journal of Chromatography A* 1258 (2012) 43-54.



Acrylamide by QuEChERS Extraction with LC/MS/MS Detection

UCT Product Number:

ECMSSC50CT-MP (4000 mg MgSO₄, 1000 mg NaCl, in Mylar pouch)

CUMPS15C18CT (150 mg, MgSO₄, 150 mg PSA and 50 mg endcapped C18, 2 mL centrifuge tube)

August 2012

Acrylamide is a neurotoxin and classified as a probable human carcinogen and genotoxicant

Procedure

1. Extraction

- a) Add 5 g of homogenized sample to a 50 mL centrifuge tube
- b) Fortify with ISTD
- c) Add 10 mL of reagent water, vortex
- d) Allow >15 minutes for hydration
- e) Add 10 mL of acetonitrile, vortex
- f) Add salts from Mylar pouch **ECMSSC50CT-MP**
- g) Shake vigorously for 1 min
- h) Centrifuge at 5000 rpm for 10 min
- i) Supernatant is ready for clean-up

2. Clean-up

- a) Add 1 mL of supernatant to the 2 mL centrifuge tube **CUMPS15C18CT**
- b) Vortex for 30 sec
- c) Centrifuge at 5000 rpm for 10 min
- d) Transfer 500 µL of extract into LC vial for analysis

LC/MS/MS Instrumentation

- **LC:** Thermo Accela 1250 pump
- **Column:** Sepax C18, 150 mm x 2.1 mm, 3µm
- **Guard:** Restek C18 2.1 x 20mm
- **Column Temperature:** Ambient
- **Injection:** 20 µL at 15° C
- **Mobile Phase:** A: water; B: methanol
- **Flow Rate:** 200 µL/min

Mobile Phase Program

Time	Mobile Phase
0 min	100% A
3 min	100% A
5 min	100% B
6 min	100% B
7 min	100% A
12 min	100% A*

* divert to waste for 0-1 min & 3-12 min

MS/MS: Thermo TSQ Vantage

- **Ion source:** APCI
- **Ion polarity:** positive mode
- **Discharge Current:** 22V
- **Declustering voltage:** 11 V
- **Sheath gas:** N₂ at 10 arbitrary units
- **Auxiliary gas:** N₂ at 15 arbitrary units
- **Vaporizer Temp:** 380 °C
- **Ion Transfer Cap. Temp:** 250 °C
- **Scan Type:** SRM
- **Dwell Time:** 150 ms
- **CID Pressure:** 0.5 mTorr

SRM Transitions				
Analyte	Parent ion	Product ion 1	Collision Energy	S-lens
Acrylamide	72.0	55.0	9	43
	72.0	27.0	20	43
¹³ C ₃ -Acrylamide	75.1	58.0	10	50

Acrylamide Recovery

Matrix n = 5	Analyte Conc. ng/mL	Recovery %	Mean Conc. ng/mL	SD (%)	RSD (%)
French Fries	50	106	53	6.1	11.56
	250	106	265	9.8	3.71
Potato Chips	50	111	56	9.4	16.84
	250	103	257	32.2	12.54
Multigrain Cereal	50	98	49	5.9	11.93
	250	93	232	9.2	3.97

DCN-211380-242



Trichothecene Type A & B Analysis in Wheat and Corn Using the QuEChERS Approach*

UCT Part Number:

ECMSSC50CT-MP (50 mL centrifuge tube, 4 g anhydrous magnesium sulfate, 1 g NaCl)

CUMPS2CT (150 mg anhydrous magnesium sulfate and 50 mg PSA)

January 2010

An extraction and purification method for the simultaneous LC-MS determination of five mycotoxins is described including three type A, diacetoxyscirpenol (DAS), T-2 toxin and HT-2 toxin, and two type B trichothecenes, deoxynivalenol (DON) and nivalenol (NIV). These mycotoxins are responsible for a wide range of disorders in animals. They have been found to inhibit proteins synthesis and to have immunosuppressive and cytotoxic effects. Health risks associated with human exposure to *Fusarium* toxins are recognized worldwide and depend on concentration in a particular diet. The major dietary sources of trichothecenes are cereal products wheat and corn. The analysis has been optimized using a modified QuEChERS approach.

Procedure

1. Sample Preparation

- a) Thoroughly homogenize a sample of grain products using a laboratory mill
- b) Weigh 5 g of sample into the 50 ml centrifuge tube
- c) Add 10 mL of methanol:acetonitrile (85:15) into 50 mL centrifuge tube
- d) Shake to disperse solvent
- e) Add the contents of the **ECMSSC50CT-MP** pouch containing 4 g anhydrous magnesium sulfate, 1 g sodium chloride to the centrifuge tube
- f) Vortex for 1 minute then centrifuge @ 4,000 rpm for 10 minutes

2. Sample Clean-up


- a) Transfer a 1 mL aliquot to a 2 mL **CUMPS2CT** tube (150 mg anhydrous magnesium sulfate and 50 mg PSA)
- b) Shake for 1 minute
- c) Centrifuge for 10 minutes @ 4,000 rpm
- d) Filter extract through a 0.45 µm filter into an LC injection vial if supernatant is not clear
- e) Sample is now ready for analysis

3. Analysis

- MSD detection with atmospheric pressure ionization (API) configured for electrospray positive ion mode
- Analytical column: Luna C18 (250mm x 4.6 mm x 5 µm) or equivalent may be used but may change elution times
- Mobile phase A: 1% formic acid in water, B: 1% formic acid in methanol
- Gradient, Flow 0.5 mL/minute, Initial 40%B, 10 minutes 90% B until 25 minutes

Mass Ions for Mycotoxins [Na+M]	
Ion	M/Z
NIV	355
DON	319
DAS	389
HT2	447
T2	489

Chromatogram Showing Elution of Mycotoxins
Peaks in order of elution: NIV, DON, DAS, HT-2, T-2



*Modified from Sospedra et al, "Use of the Modified Quick, Easy, Cheap, Effective, Rugged and Safe Sample Preparation Approach for the Simultaneous Analysis of Type A and B Trichothecenes in Wheat Flour," J of Chromatography A

DCN-102201-182



Multiresidue Analysis in Cereal Grains Using Modified QuEChERS Method with UPLC-MS/MS and GC-TOFMS*

UCT Part Number:

ECMSSC50CT-MP (50 mL centrifuge tube, 4 g anhydrous. magnesium sulfate, 1 g NaCl)

CUMPS15C18CT (150 mg anhydrous magnesium sulfate, 150 mg PSA and 50 mg C18)

February 2010

This QuEChERS procedure is specifically developed for cereal grains (corn, oats, rice and wheat) using ultra pressure liquid chromatography-tandem mass spectrometry UPLC MS/MS and automated direct sample introduction GC-TOFMS to achieve good recoveries of over 150 analytes

Pesticide Reference Standards (Chemservice (West Chester, PA)

- Prepare individual pesticide stock solutions (2000 - 5000 µg/mL) in ethyl acetate or acetonitrile (MeCN) and store at -18° C
- Prepare two composite pesticide stock solutions, MIX-1 and MIX-2 at 10 µg/mL in MeCN
- Add 0.1% acetic acid to prevents degradation of base-sensitive analytes in MeCN

Isotopically Labeled Internal Standards (Cambridge Isotope Laboratories, Inc. (Andover, MA))

Prepare 5 µg/mL in acetone

- atrazine (ethylamine-d5)
- carbofuran (ring-¹³C6)
- dimethoate (o,o-dimethyl-d6)
- 2,4-DDT (ring-¹³C6)
- α-HCH (¹³C6)
- parathion (diethyl-d10)

QC Working Solution

- trans-permethrin (phenoxy-¹³C6) (1 and 5 µg/mL in acetone)

Procedure

1. Sample Preparation

- a) Thoroughly homogenize a sample of grain products using a laboratory mill to a flour-like consistency
- b) Place appropriate weight** of sample into the 50 ml centrifuge tube
- c) Add 10 mL of deionized water (15 mL for rice) and 10 mL of acetonitrile
- d) Add 200 µL of ISTD standard solution
- e) Vortex tube to disperse sample and standard for 1 hour using a wrist action shaker
- f) Add the contents of the **ECMSSC50CT-MP** pouch into the centrifuge tube
- g) Immediately seal tube and vortex for 1 minute
- h) Centrifuge @ rcf >3,000 for 10 minutes

2. Sample Clean-up

- a) Transfer a 1 mL aliquot to a 2 mL **CUMPS15C18CT** tube
- b) Vortex for 30 seconds
- c) Centrifuge for 5 minutes
- d) Transfer 300 µL of the supernatant into the chamber of a Mini-UniPrep syringeless filter vial (Whatman) and add 30 µL 1 µg/mL QC solution*
- e) Mix thoroughly
- f) Transfer 125 µL of the extract in the Mini-UniPrep vial into a deactivated glass insert placed in a GC autosampler vial and cap the vial with a heat treated septum (overnight at 250° C)
- g) Press the 0.2 µm polyvinylidene fluoride (PVDF) filter of the Mini-UniPrep to filter the extract for the UPLC-MS/MS analysis
- h) Add 30 µL of QC standard solution
- i) Sample is now ready for analysis

3. Analysis UPLC-MS/MS

- Acquity UPLC interfaced to a Quattro Premier triple-quad mass spectrometer (Water's Corp.) MassLynx software v 4.1 or equivalent
- **Column:** Acquity UPLC BEH C18 (50 x 2.1 mm, 1.7 µm particle size, 130 Å pore size) or equivalent

- **Temperature:** 40°C
- **Injection Volume:** 2 µL

Binary Mobile Phase:

- **A** 10 mM ammonium formate in water (pH 3, adjusted with formic acid)
- **B** 10 mM ammonium formate in methanol

Gradient:

Flowrate: 450 µL/minute

Time minutes	% B
0	30
4	30
7.5	60
8.5	60
10.5	100
12.5	100
12.6	30
15.0	30

MS Determination

- Electrospray (ESI) positive mode combined with monitoring of the two most abundant MS/MS (precursor f product) ion transitions.

The MS source conditions:

- capillary voltage of 1.7 kV
- extractor voltage of 4.0 V
- RF lens at 0.9 V
- source temperature of 130° C
- desolvation temperature of 350° C
- collision gas (argon) pressure of 4.31 x 10⁻³ mbar
- desolvation gas (N₂) flow of 600 L/h
- cone gas (N₂) flow of 100 L/h

4. For GC amenable pesticides use automated DSI-GC-TOF Mass Analyzer

GC Column: Use a combination of a 20 m x 0.25 mm id x 0.25 µm film thickness RTX-5 ms column and a 1m x 0.1 mm id x 0.1 µm film thickness RTX-pesticide 2 column (Restek). This translates into a 1.68 m x 0.1 mm id “virtual” column setting in the ATAS Evolution software or equivalent

Oven Temperature Program (start after a 4.5 minutes solvent vent period):

- 60° C, hold for 4 minutes then ramped to 180° at 20° C/minutes, then ramp 5° C/minutes to 230° C, then 20° C/minutes to 280° C, and finally ramp to 300° C at 40° C/minutes, and hold for 12 minutes. The total run time is 35 minutes.

Automated DSI-GC-TOFMS Analysis.

- Agilent 6890 GC equipped with a secondary oven and nonmoving quad-jet dual stage modulator for two-dimensional comprehensive GC/GC chromatography or equivalent
- Pegasus 4D (Leco Corp., St. Joseph, MI) TOF mass spectrometer or equivalent
- Inject using CombiPAL autosampler (Leap Technologies, Carrboro, NC) or equivalent
- Automated DSI accessory (LINEX) with an Optic 3 programmable temperature vaporizer (PTV) inlet (ATAS-GL International, Veldhoven, The Netherlands) or equivalent
- Leco Chroma TOF (version 3.22) software for GC TOFMS control and data acquisition/processing or equivalent
- CombiPAL Cycle Composer with macro editor (version 1.5.2) and ATAS Evolution software (version 1.2a) to control the automated DSI process and PTV (including column flow) or equivalent


Automated DSI Injection:

- Inject 10 μL into a disposable microvial (1.9 mm i.d., 2.5 mm o.d., 15 mm, (Scientific Instrument Services, Ringoes, NJ), Siltek deactivated (Restek Bellefonte, PA) or equivalent
- Wash with acetone heated at 250° C
- Place in a LINEX DMI tapered liner
- The liner is then transferred into the Optic inlet

Optic 3 PTV Conditions:


- Solvent vent at an injector temperature of 100° C for 4.5 minutes
- Initial column flow of 0.8 mL/minutes and a split flow of 50 mL/minutes,
- Follow by a splitless transfer of analytes for 4 minutes. The injector temperature was ramped to 280° C (at 16° C/s) Column flow changed to 1.5 mL/minutes (kept constant for the entire GC run). After the splitless period, the split flow adjusted 50 mL/minutes for 6 minutes. After 6 minutes reduce split flow to 25 mL/minutes and decrease injector temperature to 250° C

Shown Below are the UPLC-MS/MS Extracted Ion Chromatograms of Selected Pesticides Spiked at 25 ng/g in Wheat Extract



Total Ion Chromatogram

DSI-LVI-GC-TOFMS analysis of a corn extract prepared using 5 g of sample, original QuEChERS (with 10 mL of water addition for swelling), and 50 mg of PSA in the dispersive SPE step. The highlighted region of the chromatogram is saturated with fatty acids. The dotted trace represents optimized analysis using 2.5 g of corn sample using dispersive SPE with 150mg of PSA and 50 mg of C₁₈



*Summarized from Mastovska et al, "Pesticide Multiresidue Analysis in Cereal Grains Using Modified QuEChERS Method Combined with Automated Direct Sample Introduction GC-TOFMS and UPLC-MS/MS Techniques" , " J of Agricultural and Food Chemistry, Full article may be found at <http://forums.unitedchem.com/>

** Corn 2.5 g, oat 3.5 g, rice 5.0 g, wheat 5.0 g

Listing of chemical suppliers and instrument manufacturers does not constitute endorsement by UCT



Extraction of Pesticides from Tomato Using the QuEChERS Approach

(This method is applicable to all pigmented fruit and vegetables)

UCT Product Number:

ECQUEU750CT-MP (4000 mg magnesium sulfate anhydrous, 1000 mg sodium chloride, 500 mg sodium citrate dibasic sesquihydrate, 1000 mg sodium citrate tribasic dihydrate)

ECQUEU32CT (2 mL micro-centrifuge tube with 150 mg magnesium sulfate anhydrous, 25 mg primary secondary amine bonded phase (PSA) and 2.5 mg graphitized carbon black)

ECQUEU515CT (15 mL centrifuge tube with 900 mg magnesium sulfate anhydrous, 150 mg primary secondary amine (PSA) bonded phase and 15 mg graphitized carbon black)

March 2010

Procedure

1. Sample Preparation

- a) Add 15g of homogenized and hydrated tomato product (> 80% moisture) to a centrifuge tube
- b) Add 15 mL acetonitrile including internal standard
- c) Shake or vortex for 30 seconds
- d) Add contents of a package of **ECQUEU750CT-MP** to centrifuge tube
- e) Immediately, shake vigorously for 2 minutes
- f) Centrifuge for 2 minutes at 3450 rcf
- g) Draw 1 or 6 mL of supernatant for clean-up

2. Clean-Up

- a) For 1 mL of supernatant, use product **ECQUEU32CT**
- b) For 6 mL of supernatant, use product **ECQUEU515CT**
- c) Add supernatant to centrifuge tube and shake vigorously for 1 minute
- d) Centrifuge for 2 minutes at 3450 rcf

3. **Analysis by GC** (suggested)

- a) Transfer an aliquot of supernatant from step 2 to a centrifuge tube
- b) Add TPP solution and 1 mL of toluene
- c) Evaporate using nitrogen at 50°C to approximately 0.3 to 0.6 mL.
- d) Bring to 1 mL final volume with toluene
- e) Inject 8 µL on LVI/GC/MS

4. **Analysis by LC** (suggested)

- a) Transfer 0.25 mL of supernatant from step 2 to a LC vial.
- b) Add TPP solution and 0.86 mL of 6.7 mM formic acid
- c) Analyze by LC/MS/MS

References:

QuEChERS Method EN 15662

Anastassiades, et al (2003) "Fast and Easy Multiresidue method employing acetonitrile extraction partitioning and dispersive solid-phase extraction for the determination of pesticide residues in product" Journal of AOAC International Vol 86 no. 2



Extraction of Polycyclic Aromatic Hydrocarbons from Fish Using the QuEChERS Approach

UCT Product Number:

ECMSSC-MP (4000 mg Magnesium Sulfate, 1000 mg sodium chloride)

ECMPSC1815CT (900 mg MgSO₄, 300 mg PSA and 150 mg endcapped C18)

ECPAHFR50CT (ENVIRO-CLEAN® PAH 50 mL centrifuge tubes)

June 2010

The QuEChERS approach is used for the extraction of 16 PAH compounds from fish with analysis by LC with fluorescence detection

1. Extraction: Protect samples from exposure to light to avoid degradation

- Transfer 5 grams of homogenized fish to a 50 mL centrifuge tube
- Add 10 mL of acetonitrile
- Mix by shaking
- Add the contents of pouch containing 4 grams MgSO₄ and 1 gram of NaCl to the centrifuge tube
- Immediately vortex the mixture for 3 minutes
- Centrifuge samples for 3 minutes at 3400 rpm
- Recover the clear supernatant for clean-up

2. Clean-up, Dispersive Solid-phase (dSPE)

- Use 3 mL of the supernatant for clean-up
- Add supernatant to centrifuge tube **ECMPSC1815CT**
- Shake for 1 minute
- Centrifuge for 1 minute at 3400 rpm
- Filter supernatant through 0.20 µm PTFE membrane filter
- Samples are ready for analysis

3. Analysis Conditions

Separation of the compounds is performed in a C18 column (CC 150/4Nucleosil 100-5 C18 PAH, 150 x 4.0 mm; 5 µm particle size; Macherey-Nagel, Duren, Germany) or equivalent maintained at room temperature

- a) Inject 15 μ L
- b) The initial composition of the mobile phase is 50% of ACN and 50% water
- c) Program A linear gradient to 100% in 15 minutes
- d) Hold 13 minutes
- e) Initial conditions are achieved within 1 minute and maintained for 6 min before next run
- f) Total run time 40 minutes
- g) Flow rate 0.8 mL/minute

Fluorescence wavelength program:

Each compound is detected at its optimum excitation/emission wavelength pair:

315/260 nm naphthalene, acenaphthene and fluorene

366/260 nm Phenanthrene

430/260 nm anthracene, fluoranthene, pyrene, benz(a)anthracene, chrysene, benzo(b)fluoranthene, benzo(k)fluoranthene, benzo(a)pyrene, dibenz(a,h)anthracene, benzo(g,h,i)perylene and dibenzo(a,l)pyrene

505/290 nm (indeno(1,2,3-cd)pyrene)

***Adapted from**, Ramalhosa, Maria Joao et al, "Analysis of polycyclic aromatic hydrocarbons in fish: evaluation of a quick, easy, cheap, effective, rugged, and safe extraction method", J. Sep. Sci. 2009, 32, 3529 – 3538



Determination of Polycyclic Aromatic Hydrocarbons (PAHs) in Seafood Using GC/MS

UCT Part Numbers:

ECQUUS2-MP (4 g of muffled anh. MgSO_4 and 2 g of NaCl)

ECPAHFR50CT (50 mL centrifuge tubes, PAHs removed)

EUSILMSSM26 (6 mL, 1g silica gel cartridge with 200 mg of muffled anhydrous sodium sulfate on top)

January 2012

This method is used for the determination of polycyclic aromatic hydrocarbons (PAHs), polychlorinated biphenyls (PCBs) and polybrominated diphenyl ethers (PBDEs) in fish and seafood--oyster, shrimp, and mussel. Benzo[a]pyrene is the main analyte of interest. GC/MS instrumentation is used for analysis.

PAH Analytes Covered in this Method

PAH	Abbreviation	PAH	Abbreviation
Anthracene	Ant	Indeno[1,2,3- <i>cd</i>]pyrene	IcdP
Benz[a]anthracene	BaA	Naphthalene	Naph
Benzo[a] pyrene	BaP	Phenanthrene	Phe
Benzo[b]fluoranthene	BpF	Pyrene	Pyr
Benzo[k]fluoranthene	BkF	3-Methylchrysene	3-MC
Benzo[g,h,i]perylene	BghiP	1-Methylnaphthalene	1-MN
Chrysene	Chr	1-Methylphenanthrene	1-MP
Dibenz[a,h]anthracene	DBahA	2,6 Dimethylnaphthalene	2,6-DMN
Fluoranthene	Flt	1,7-Dimethylphenanthrene	1,7-DMP
Fluorene	Fln		

Procedure

1. Extraction

- a) To the 50 mL polypropylene centrifuge tube add 10 ± 0.1 g of homogenized seafood sample
- b) Add 50 μ L of 1 μ g/mL ^{13}C -PAHs solution
- c) Vortex sample for 15 sec and then equilibrate for 15 min
- d) Add 5 mL of reagent water and 10 mL of ethyl acetate (EtOAc)
- e) Shake for 1 min

2. Partition

- a) Add the contents of pouch **ECQUUS2-MP**. Tightly seal the tube to ensure that salts do not get into the screw threads
- b) Shake for 1 min
- c) Centrifuge at $> 1,500$ rcf for 10 min
- d) Remove 5 mL aliquot of the upper ethyl acetate layer, add 50 μ L of isooctane as a keeper
- e) Evaporate all ethyl acetate until only isooctane and co-extracted sample fat remain
- f) Reconstitute with 1 mL of hexane

3. Clean-Up

- a) Condition a silica SPE column **EUSILMSSM26 (Note 1)** (1 g of silica gel with approx. 0.2 g of muffled anh. sodium sulfate on the top) with 6 mL of hexane:dichloromethane (3:1 v/v) and 4 mL of hexane
- b) Apply the reconstituted extract to the silica SPE cartridge (**Note 2**)
- c) Elute with 10 mL of hexane:dichloromethane (3:1 v/v) and collect the eluent
- d) Add 0.5 mL isooctane to the eluent as a keeper and evaporate to 0.5 mL to remove hexane and dichloromethane from the final extract
- e) Transfer the final extract into an autosampler vial for GC/MS analysis

Notes:

1. The fat capacity of the 1-g silica gel SPE column is approx. 0.1 g. If the ethyl acetate extract aliquot contains more than 0.1 g of fat, use a smaller aliquot to avoid sample breakthrough
2. Ethyl acetate should not be present in the extract applied to the silica cartridge as it affects extract polarity and potential retention of fat and analytes on the silica gel.

GC Conditions for the Analysis of PAHs

Column	BPX-50 (30 m x 0.25 mm i.d. x 0.25 µm film thickness)
Oven Temperature Program	80°C (hold for 4.3 min), 30°C/min to 220°C, 2°C/min to 240°C, and 10°C/min to 360°C (hold for 17 min)
He Flow Rate	1.3 mL/min (hold for 19 min), then 50 mL/min to 2 mL/min (hold for 16 min)
Injection Technique	PTV solvent vent
Injection Volume	1 x 8 µL
Vent Time	2.3 min
Vent Flow	50 mL/min
Vent Pressure	50 psi
Inlet Temperature Program	50°C (hold for 2.3 min), then 400°C/min. to 300°C

MS Conditions

Any GC-MS instrument (single quadrupole, triple quadrupole, time-of-flight or ion trap) with electron ionization (EI) may be used

MS Ions (*m/z*) for Quantification and Identification of Target PAHs for Single-stage MS Instruments

Analyte PAH's	Abbreviation	Confirmation Ions (<i>m/z</i>)	Quantification Ions (<i>m/z</i>)
Anthracene	Ant	177	178
Benz[<i>a</i>]anthracene	BaA	226	228
Benzo[<i>a</i>] pyrene	BaP	253	252
Benzo[<i>b</i>]fluoranthene	BpF	253	252
Benzo[<i>k</i>]fluoranthene	BkF	253	252
Benzo[<i>g,h,i</i>]perylene	BghiP	277	276
Chrysene	Chr	226	228
Dibenz[<i>a,h</i>]anthracene	DBahA	276	278
Fluoranthene	Flt	200	202


Fluorene	Fln	165	166
Indeno[1,2,3-cd]pyrene	IcdP	277	276
Naphthalene	Naph	127	128
Phenanthrene	Phe	177	178
Pyrene	Pyr	200	202
3-Methylchrysene	3-MC	241	242
1-Methylnaphthalene	1-MN	115	142
1-Methylphenanthrene	1-MP	189	192
2,6 Dimethylnaphthalene	2,6-DMN	141	156
1,7-Dimethylphenanthrene	1,7-DMP	191	206

**MS Ions (*m/z*) for Quantification and Identification of Target ¹³C-PAHs
for Single-stage MS Instruments**

Analyte PAH's	Abbreviation	Confirmation Ions (<i>m/z</i>)	Quantification Ions (<i>m/z</i>)
Naphthalene (¹³ C ₆)	Naph- ¹³ C ₆	133	134
Fluorene (¹³ C ₆)	Fln- ¹³ C ₆	171	172
Phenanthrene (¹³ C ₆)	Phe- ¹³ C ₆	183	184
Anthracene (¹³ C ₆)	Ant- ¹³ C ₆	183	184
Fluoranthene (¹³ C ₆)	Flt- ¹³ C ₆	205	208
Pyrene (¹³ C ₆)	Pyr- ¹³ C ₆	208	205
Benz[a]anthracene (¹³ C ₆)	BaA- ¹³ C ₆	232	234
Chrysene (¹³ C ₆)	Chr- ¹³ C ₆	232	234
Benzo[b]fluoranthene (¹³ C ₆)	BbF- ¹³ C ₆	259	258
Benzo[k]fluoranthene (¹³ C ₆)	BkF- ¹³ C ₆	259	258
Benzo[a]pyrene (¹³ C ₄)	BaP- ¹³ C ₄	257	256
Indeno[1,2,3-cd]pyrene (¹³ C ₆)	IcdP- ¹³ C ₆	283	282

Dibenz[<i>a,h</i>]anthracene ($^{13}\text{C}_6$)	DBahA- $^{13}\text{C}_6$	282	284
Benzo[<i>g,h,i</i>]perylene ($^{13}\text{C}_{12}$)	BghiP- $^{13}\text{C}_{12}$	289	288

An Example Chromatogram of A GC Separation of PAH's and Their Alkyl Homologues In A Standard Solution Mixture At 25 ng/mL In Isooctane



1 – naphthalene, 2 – 2-methylnaphthalene, 3 – 1-methylnaphthalene, 4 – acenaphthylene, 5 – acenaphthene, 6 – fluorene, 7 – dibenzothiophene, 8 – phenanthrene, 9 – anthracene, 10 – fluoranthene, 11 – pyrene, 12 – 1-methylpyrene, 13 – benzo[*c*]fluorene, 14 – benz[*a*]anthracene, 15 – cyclopenta[*c,d*]pyrene, 16 – chrysene, 17 – 1-methylchrysene, 18 – 5-methylchrysene, 19 – 3-methylchrysene, 20 – benzo[*b*]fluoranthene, 21 – benzo[*k*]fluoranthene, 22 – benzo[*j*]fluoranthene, 23 – benzo[*a*]pyrene, 24 – dibenz[*a,h*]anthracene, 25 – indeno[1,2,3-*cd*]pyrene, 26 – benzo[*g,h,i*]perylene, 27 – dibenzo[*a,l*]pyrene, 28 – dibenzo[*a,e*]pyrene, 29 – dibenzo[*a,h*]pyrene, 30 – dibenzo[*a,i*]pyrene

*The analyst should refer to Katerina Mastovska, Wendy R. Sorenson, Covance Laboratories Inc
Jana Hajslova, Institute of Chemical Technology, Prague "Determination of Polycyclic Aromatic Hydrocarbons (PAHs) in Seafood using Gas Chromatography-Mass Spectrometry: A Collaborative Study"

References

Lucie Drabova, Kamila Kalachova, Jana Pulkrabova, Tomas Cajka, Vladimir Kocourek and Jana Hajslova. "Rapid Method for Simultaneous Determination of Polycyclic Aromatic Hydrocarbons (PAHs), Polychlorinated Biphenyls (PCBs) and Polybrominated Diphenyl Ethers (PBDEs) in Fish and Sea Food Using GC-TOFM," ICT document, Prague, Czech Republic, 2010.



QuEChERS Analysis of Miticides and Other Agrochemicals in Honey Bees, Wax or Pollen*

UCT Part Number:

ECMSSA50CT-MP (6000 mg MgSO₄ and 1500 mg sodium acetate)

CUMPSC18CT (150 mg MgSO₄, 50 mg PSA, 50 mg C18)

ECPSACB256 (dual layer cartridge 250 mg GCB, 500 mg PSA)

ECMAG00D (organic free magnesium sulfate anhydrous)

March 2012

An analytical method using QuEChERS type procedures for 121 different pesticide residues is described. Extracts of wax, beebread, and adult bees or brood can also be analyzed for metabolites of primary miticides and insecticides using this method. This includes the oxon and phenolic metabolites of coumaphos, chlorferone, the sulfoxide and sulfone metabolites of aldicarb, and the toxic olefin and 5-hydroxy metabolites of imidacloprid.

Sample Collection and Preservation

- Wrap in aluminum foil and store on dry ice until placement in a -80° C freezer as soon as possible
- Beebread and brood can be removed from the combs at room temperature and then stored along with the remaining beeswax at -20° C until processed

1. Sample Preparation

- a) Weigh 3 grams beebread (or comb wax) into a 50 mL centrifuge tube
- b) Add 100 µL of a process control spiking solution
- c) Add 27 mL of extraction solution*
- d) *44% DI water, 55% acetonitrile & 1% glacial acetic acid

Brood and adults are extracted without using DI water

- e) Add 100 µL of an internal standard
- f) For beebread, reduce particle size by using a high speed disperser for 1 minute
- g) For comb wax melt the sample at 80° C in a water bath followed by cooling to RT
- h) Add the contents of **ECMSSA50CT-MP** pouch to the mixture

- i) Seal tube and shake vigorously for 1 minute
- j) Centrifuge for 1 minute

2. **Clean-Up for LC/MS-MS**

- a) Transfer 1 mL of supernatant to **CUMPSC18CT** micro centrifuge tube
- b) Vortex for 1 minute and centrifuge
- c) Transfer supernatant to an autosampler vial for LC analysis

3. **Clean-Up for GC/MS**

- a) Prepare a dual layer solid-phase extraction cartridge **ECPSACB256** by adding about 80 mg of anhydrous magnesium sulfate (ECMAG00D) to the top of the cartridge
- b) Condition cartridge by adding 4.0 mL of acetone/toluene (7:3 v:v)
- c) Using a positive pressure or vacuum manifold, elute solvent to waste
- d) Add 2 mL of supernatant to the cartridge
- e) Elute cartridge using 3 x 4 mL of acetone/toluene & 7:3, (v:v) into a 15 mL graduated glass centrifuge tube
- f) Using an analytical evaporator @ 50°C, dry eluate to a final volume of 0.4 mL
- g) Sample is ready for analysis

4. **Analysis—by LC or GC**

LC analysis is necessary for neonicotinoids, other polar pesticides and their metabolites

For LC analysis

- a) Analysis by LC/MS-MS use a 3.5 µm 2.1 X 150 mm Agilent Zorbax SB-C18 (or equivalent)
- b) Agilent 1100 LC with a binary pump interfaced to a Thermo-Fisher TSQ Quantum Discovery triple quadrupole MS 9 (or equivalent)

For GC analysis

- a) For analysis use Agilent 6890 (or equivalent) GC equipped with a 0.25 mm, 30 m J&W DB-5MS (2 µm film) capillary column
- b) Interface to an Agilent 5975 triple quadrupole MS (or equivalent)
- c) Use GC/MS in the electron impact and negative chemical ionization modes

Pesticides Representatives Found in Wax Samples	
Aldicarb sulfoxide	Flutolanil
Aldicarb sulfone	Fluvalinate
Allethrin	Heptachlor
Atrazine	Heptachlor epoxide
Azinphos methyl	Hexachlorobenzene
Azoxystrobin	Imidacloprid
Bendiocarb	Iprodione
Bifenthrin	Malathion
Boscalid	Metalaxyl
Captan	Methidathion
Carbaryl	Methoxyfenozide
Carbendazim	Metribuzin
Carbofuran	Norflurazon
Carbofuran, 3-hydroxy	Oxyfluorfen
Carfentrazone ethyl	Parathion methyl
Chlorfenapyr	p-Dichlorobenzene
Chlorferone (coumaphos)	Pendimethalin
Chlorothalonil	Permethrin
Chlorpyrifos	Phosmet
contrast	Piperonyl butoxide
Coumaphos	Prallethrin
Coumaphos	Pronamide
Cyfluthrin	Propiconazole
Cyhalothrin	Pyraclostrobin
Cypermethrin	Pyrethrins
Cyprodinil	Pyridaben
DDE p,p'	Pyrimethanil
Deltamethrin	Pyriproxyfen
Diazinon	Quintozone
Dicofol	Spirodiclofen
Dieldrin	Tebufenozide
Dimethomorph	Tebuthiuron
DMA (amitraz)	Tefluthrin
DMPF (amitraz)	Tetradifon
Endosulfan I	Thiabendazole
Endosulfan II	Thiacloprid
Endosulfan sulfate	Triadimefon
Esfenvalerate	Tribufos
Ethion	Trifloxystrobin
Ethofumesate	Trifluralin
Fenamidone	Vinclozolin
Fenbuconazole	
Fenhexamid	
Fenpropathrin	
Fipronil	

*Summarized and adapted from: Mullin CA, Frazier M, Frazier JL, Ashcraft S, Simonds R, et al. (2010) **High Levels of Miticides and Agrochemicals in North American Apiaries: Implications for Honey Bee Health**. PLOS ONE 5(3): e9754. doi:10.1371/journal.pone.0009754



QuEChERS Pesticide Analysis for Fresh Produce by GC/MS/MS*

UCT Product Number:

ECMSSC50CTFS-MP (6 g MgSO₄, 1.5 g NaCl)

ECQUEU1115CT (1.2 g MgSO₄, 0.4 g PSA, 0.4 g GCB)

ECMSC1850CT (1500 mg MgSO₄, 500 mg endcapped C18)

ECMAG00D (organic free MgSO₄ anhydrous)

October 12, 2010

This modified QuEChERS procedure uses GC-MS/MS for analysis of organohalogen, organophosphorus, and pyrethroid pesticides in produce. It is an improvement over the traditional QuEChERS procedure since the sample extracts are in toluene instead of acetonitrile and cleaner due to additional clean-up procedures. In addition, the method uses smaller sample sizes and less solvent than standard multiresidue procedures, and the solid-phase dispersive steps involving GCB/PSA/C18 provide sufficient clean-up for GC-MS/MS analysis.

1. Sample Extraction

- a) Combine 15 g of cryo-ground sample with 15 mL acetonitrile
- b) Add contents of **ECMSSC50CTFS-MP**
- c) Shake by hand for 2 minutes
- d) Add IS (500 µL of 3.4 µg/mL solution of tris(1,3-dichloroisopropyl) phosphate)
- e) Centrifuge 4500 rpm for 5 minutes

2. Clean-Up

- a) Transfer upper layer (12 mL) to a clean centrifuge tube **ECMSC1850CT** containing 0.5 grams C₁₈ and 1.2 g MgSO₄
- b) Shake for 1 minute and centrifuge @ 4500 rpm for 5 minutes
- c) Transfer 9 mL of supernatant to extraction tube containing **ECQUEU1115CT**
- d) Vortex 15 seconds
- e) Add 3 mL toluene
- f) Shake the centrifuge tube for 2 minutes
- g) Centrifuge @ 4500 rpm for 5 minutes
- h) Transfer extract to clean tube
- i) Reduce 6 mL volume to < 100 µL using N₂ in an evaporator (35°C)
- j) Add 1.0 mL toluene and QC standard (20 µg/mL deuterated polycyclic hydrocarbons) along with 50 mg anhydrous MgSO₄
- k) Centrifuge @ 1500 rpm for 5 minutes

- I) Transfer 1.0 mL of extract to ALS vials for analysis

Note:

- Use matrix-matched calibration standards in toluene rather than standards prepared in solvent. This will compensate for matrix enhancement effects
- Coextractives in the sample matrix have been shown to cause an enhancement of the pesticide peak response in the matrix compared to that of the same amount of the pesticide in the matrix-free solvent

GC-MS/MS Tandem Mass Spectrometry

Varian CP-3800 series gas chromatograph coupled with a Varian 1200 L triple-quadrupole mass spectrometer with a CTCOMBI PAL autosampler (Varian Inc., Palo Alto, CA)

- Column: Deactivated guard column (5 m x 0.25 mm i.d., Restek Corp.) Varian 30 m x 0.25 mm x 0.25 μ m, VF-5 fused silica capillary analytical column
- Head pressure 13.2 psi with 1.2 mL/min flow rate
- He carrier gas
- Column temperature programmed as follows:
 - initial temperature 105° C for 6 min
 - increased to 130° C at 10° C/min
 - ramp to 230° C at 4° C/min and to 290° C at 1° C/min
 - Hold for 5.5 min.
 - Total run time 45 min.
- Injector temperature: 280° C
- injection volume: 1.0 μ L in splitless mode
- Ion source and transfer line temperatures are 240° and 300° C, respectively
- Set Electron multiplier voltage to 1400V by automatic tuning
- Use argon collision gas for all MS/MS
- Pressure in the collision cell 1.8 mTorr

Table of Analytes Covered in this Method


Analytes		
acenaphthene-d ₁₀	Diamidafos (nellite)	<i>p,p'</i> -methoxychlor
acrinathrin	Diazinon	metolachlor
akton	Dibutyl chlorenate	mevinphos
alachlor	Dicaphon	mirex
aldrin	Dichlobenil	naphthalene-d ₈
allethrin	Dichlofenthion	<i>cis</i> -nonachlor
atrazine	Dichlofluanid	<i>trans</i> -nonachlor
azamethiophos	3,4'-dichloroaniline	parathion
aziphos-ethyl	4,4'-dichlorobenzophenone	parathion-methyl
aziphos-methyl	Dichlorvos	pentachloroaniline
α-BHC	Dicloran	pentachlorobenzene
β-BHC	Dieldrin	pentachlorobenzonitrile
δ-BHC	Dimethachlor	Pentachlorophenyl methyl ester
benfluralin	dioxabenzofos	pentachlorothioanisole
bifenthrin	Dioxathion	<i>cis</i> -permethrin
bromophos	Disulfoton	<i>trans</i> -permethrin
bromophos-ethyl	Ditalimfos	phenanthrene-d ₁₀
bromopropylate	Edifenphos	phenothrin
captafol	α-endosulfan	phorate
captan	β-endosulfan	phosalone
carbophenothion	Endosulfan ether	phosmet
<i>cis</i> -chlordane	Endosulfan sulfate	phenthoate
<i>trans</i> -chlordane	Endrin	pirimiphos-ethyl
α-chlordene	Endrin aldehyde	pirimiphos-methyl
β-chlordene	Endrin ketone	procymidone
γ-chlordene	EPN	profenofos
β-chlorfenvinphos	Ethalfuralin	propachlor
chlorobenzilate	Ethion	propazine
chloroneb	Ethoprop	propetamphos
chlorothalonil	Etridazole	propyzamide
chlorpyrifos	Famphur	prothiophos
chlorpyrifos-methyl	Fenamiphos (ronnel)	pyraclofos
chlorthiophos	Fenarimol	pyrazophos
chrysene-d ₁₂	Fenchlorphos	pyridaphenthion
coumaphos	Fenitrothion	quinalphos
cyanazine	Fensulfothion	quintozone
cyanophos	Fenthion	resmethrin
Cyfluthrin 1	Fenvalerate 1	simazine
Cyfluthrin 2	Fenvalerate 2	sulfotep-ethyl
Cyfluthrin 3	Fluchloralin	sulprofos
Cyfluthrin 4	Flucythrinate 1	tebupirimfos
λ-cyhalothrin	Flucythrinate 2	propachlor
Cypermethrin 1	Fluridone	propazine
Cypermethrin 2	Fluvalinate 1	Tecnazene (TCNB)
Cypermethrin 3	Fluvalinate 2	tefluthrin
Cypermethrin 4	Folpet	temephos
Dacthal (DCPA)	Fonophos	terbufos
<i>o,p'</i> -DDD	Heptachlor	terbuthylazine
<i>p,p'</i> -DDD	Heptachlor epoxide	2,3,5,6-tetrachloroaniline
<i>o,p'</i> -DDE	hexachlorobenzene	tetrachlorvinphos
<i>p,p'</i> -DDE	Iprobenfos (IBP)	tetramethrin
<i>o,p'</i> -DDT	Iprodione	thiometon
<i>p,p'</i> -DDT	Isazophos	tolclofos-methyl
DEF (tribufos)	Isofenfos	tolyfluanid
deltamethrin	Jodfenphos (iodofenphos)	triallate
demeton-S	Leptophos	triazophos
demeton-S-methyl	Lindane (BHC)	trifluralin
dialifor	Malathion	triphenyl
Diallate 1	methidathion	tris(1,3-dichloroisopropyl) phosphate
Diallate 2	<i>o,p'</i> -methoxychlor	vinclozolin

Problems with pesticides with low (<70%) recoveries or large variances (SD > 20%) may be attributed to the following issues:

- early eluting analytes
- sensitivity to pH changes
- prone to volatility loss (i.e., 3,4'-dichloroaniline, dichlorvos, diclobenil, and etridazole),
- strongly adsorbed to the PSA or GCB sorbents (i.e., chlorothalonil, endrin aldehyde, hexachlorobenzene, pentachlorobenzene, pentachlorobenzonitrile, and tachlorothioanisole)
- difficult to ionize by mass spectrometric detection (i.e., captafol, captan, dichlofluanid, folpet, and tolylfluanid).
- Highly nonpolar or late-eluting pesticides such as temephos and fluridone may also be problematic

For recovery data, target, qualifier and transition ions please reference original paper*

Reconstructed GC-MS/MS chromatograms of various commodities containing various pesticides including chlorpyrifos (A), o,p'-DDE (B), p,p'-DDT (C), and p,p'-DDE (D) present in carrot; bifenthrin (E) and endosulfan sulfate (F) present in bell pepper; β - (G) and R- (H) endosulfan, endosulfan sulfate (I) and chlorothalonil (J) present in tomato; and phosmet (K) in peach. Included are the transitions from precursor to product ions and the relative ion ratios between the two transitions, primary (top) and secondary (bottom), which are used for pesticide identification



Reagents and Materials

Pesticide standards may be obtained from:

- U.S. Environmental Protection Agency National Pesticide Standard Repository (U.S. EPA, Ft. Meade, MD)
- ChemServices (West Chester, PA), Sigma/Aldrich/Fluka Chemicals (St. Louis, MO),
- Crescent Chemicals (Islandia, NY)
- tris(1,3-dichloroisopropyl) phosphate from TCI America (Portland, OR)
- Quality control standards, naphthalene-d₈, acenaphthalened₁₀, phenanthrene-d₁₀, and chrysene-d₁₂ (Sigma/Aldrich/Fluka Chemicals (Milwaukee, WI).

*Adapted and used by permission from Jon W. Wong, Kai Zhang, "Multiresidue Pesticide Analysis In Fresh Produce By Capillary Gas Chromatography-Mass Spectrometry/Selective Ion Monitoring (GC-MS/SIM) and -Tandem Mass Spectrometry", (GC-MS/MS), J Agric. Food Sci., DOI: 10.1021/Jf903854n

Listing of instrument manufacturers and standards suppliers does not constitute endorsement by UCT. Equivalent systems may be used



Analysis of Cyromazine in Poultry Feed Using a QuEChERS Approach

UCT Product Number:

ECMSSA50CT-MP (6 g anhydrous MgSO₄ and 1.5 g Na Acetate)

EEC18156 (500 mg endcapped C18, 6 mL cartridge)

October 2010

Introduction

This summary outlines a QuEChERS procedure for the analysis of the insecticide cyromazine (Trigard or Larvadex) in poultry feed by LC-MS/MS. Processing time is significantly faster than EPA method AG-555 and uses less solvent. Modifications include adding glacial acetic acid to the acetonitrile to increase extraction efficiency.

Procedure

1. Sample Preparation

- a) Homogenize 2 grams of poultry feed and add to a 50 mL centrifuge tube
- b) Add 10 mL of acetonitrile/acetic acid (75:25)
- c) Sonicate at 50/60 Hz for 15 minutes
- d) Add the contents of **ECMSSA50CT-MP** pouch and shake for 1 minute
- e) Centrifuge at 3400 rpm for 10 minutes
- f) Transfer 1 mL of supernatant to a calibrated test tube and add 9 mL of water: acetonitrile (95:5) with 0.1% acetic acid

2. Sample Clean-up

- a) Add the 10 mL from 1) f) above to a **EEC18156** cartridge and elute dropwise
- b) Filter eluant using a 0.45 µm Teflon filter (Millipore, Billerica, MA) or equivalent
- c) Transfer 2 mL of eluant to an HPLC vial for analysis by LC-MS/MS

3. Analysis LC-MS/MS

Waters Alliance 2695 HPLC (Waters) coupled with a micromass Quattro Micro triple-quadrupole mass spectrometer (Micromass, Manchester, U.K.) or equivalent

HPLC conditions:

Guard column (Alltima, C18, 5 µm, 2.1 x 7.5 mm, Deerfield, IL) or equivalent

Analytical column (Alltima, C18, 5 µm, 2.1 x 250 mm, Waters) or equivalent

Mobile phase: (A) acetonitrile with 0.1% formic acid and (B) water with 0.1% formic acid

Gradient:

- 0-2 min, 5%A
- 2-5 min from 5 to 10% A
- 5-5.5 min from 10 to 90% A
- 5.5-8 min 90 to 5% A
- 8-10 min, from 90 to 5% A
- 10-12 min, 5% A

Flow rate 0.2 mL/minute

Injection volume: 25 µL

- **Mass Spectrometer**

- Positive ion mode electrospray ionization
- Monitor the ion transition of the parent ion (m/z 167) to the product ion (m/z 85) in multiple reaction monitoring (MRM)

Mass Spectrometry Conditions for Cyromazine Quantitation

Capillary Voltage	3.1 kV
Cone Voltage	65 V
Collision Energy	21-24 V
Source Temperature	120° C
Desolvation Temperature	350° C
Cone Gas Flow	135 L/h
Desolvation Gas Flow Rate	750 L/h
Collision Gas	Argon
Parent Ion	(m/z) 167
Product Ion	(m/z) 85

*Summarized with permission from Xia, Kang, Atkins, Jack et al, "Analysis of Cyromazine in Poultry Feed Using the QuEChERS Method Coupled with LC-MS/MS" J. Agric. Food Chem, DOI:10.1021/jf9034282

Listing of instrument manufacturers does not constitute endorsement by UCT



Determination of Anthelmintic Drug Residues in Milk Using Ultra High Performance Liquid Chromatography- Tandem Mass Spectrometry*

UCT Products:

ECMSSC50CT-MP (4000 mg anhydrous MgSO₄, 1000 mg NaCl)

ECMSC1850CT (1500 mg anhydrous MgSO₄ and 500mg C18)

Revised March 27 2013

Introduction

A modified QuEChERS-based method is used with an additional concentration step to detect 38 anthelmintic residues (nematicides, flukicides, endectocides) in milk at $\leq 1\mu\text{g/kg}$ using UHPLC-MS/MS detection. The drugs covered by this method include benzimidazoles, avermectins and flukicides.

Procedure

1. Sample Preparation

- a) Weigh 10.0 g milk into a 50 mL centrifuge tube
- b) Add IS and allow to sit for 15 minutes
- c) Add 10 mL acetonitrile (MeCN) and the contents of **ECMSSC50CT-MP** pouch
- d) Shake vigorously, then centrifuge for 12 minutes @ $\geq 3,500$ rcf

2. Dispersive Sample Clean-up

- a) Add the supernatant to **ECMSC1850CT**
- b) Vortex sample for 30 seconds
- c) Centrifuge for 10 minutes @ ≥ 3000 rcf
- d) Transfer 5 mL of supernatant to an evaporation tube
- e) Add 0.25 mL DMSO (keeper solvent) and vortex briefly
- f) Evaporate the MeCN @ 50° C using nitrogen evaporation to 0.25 mL
- g) Filter extract using 0.2 μm PTFE syringe filter
- h) Sample is ready for UHPLC-MS/MS analysis

3. Analysis UHPLC-MS/MS

- Waters Acquity UPLC system (Milford MA; USA) or equivalent
- Analytical column HSS T3 C18 (100 × 2.1 mm, particle size 1.8 µm) (or equivalent) with appropriate guard column
- Column temperature: 60° C
- Pump flow rate of 0.6 mL/min
- Binary gradient:
mobile phase A 0.01% formic acid in water:MeCN (90:10, v/v)
mobile phase B 5mM ammonium formate in MeOH:MeCN (75:25 v/v)
Gradient profile:
0 – 0.5 min, 100% A
5 min, 50% A
7 min, 10% A
8.5 min, 10% A
8.51 min, 0% A
9.5 min, 0% A
9.51 min, 100% A
13 min, 100% A
- Injection volume 5 µL
- Waters Quattro Premier XE triple quadrupole mass spectrometer
- Electrospray ionization (ESI) interface using fast polarity switching
- System controlled by MassLynx™ software and data was processed using TargetLynx™ Software (Waters)

Note:

Ammonium formate is used in the organic mobile phase because abamectin, doramectin and ivermectin form sodium adducts ([M+23]⁺) when acids are used. In this case, the ammonium adducts ([M+18]⁺) should be monitored for these three compounds and not the protonated precursor ions.

MS amenable acids can be used for the aqueous mobile phase, which should be at a low pH (≤4) to get the best results. It is essential to use ammonium buffer in the organic mobile phase as the avermectins elute at 100% organic content. The aqueous mobile phase may also include ammonium buffer, although it is not an essential requirement. Additionally, ammonium formate is more soluble in organic solvent than ammonium acetate.

Albendazole-sulfone and hydroxy-mebendazole are prone to isobaric interference as they have similar precursor and product ions that can't be distinguished using triple quadrupole instruments. It is therefore necessary to chromatographically separate these two compounds.

Standards, Internal Standards, Stock Solutions & Suppliers

Sigma-Aldrich

Analyte*	Abbreviation	Analyte	Abbreviation
Abamectin	ABA	Ivermectin	IVER
Albendazole	ABZ	Levamisole	LEVA
Bithionol	BITH	Morantel	MOR
Clorsulon	CLOR	Niclosamide	NICL
Closantel	CLOS	Nitroxynil	NITR
Coumaphos	COUM	Oxfendazole	OFZ
Doramectin	DORA	Oxyclozanide	OXY
Enamectin	EMA	Rafoxanide	RAF
Fenbendazole	FBZ	Thiabendazole	TBZ
Haloxon	HAL		

Witega Laboratories Berlin-Aldershof GmbH (Berlin, Germany)

Analyte**	Abbreviation
Albendazole-2-amino-sulfone	ABZ-NH ₂ -SO ₂
Albendazole sulfone	ABZ-SO ₂
Albendazole-sulfoxide	ABZ-SO
Amino-oxibendazole	OXI-NH ₂
5-hydroxy-thiabendazole	5-OH-TBZ
Fenbendazole-sulfone	FBZ-SO ₂
Triclabendazole	TCB
Triclabendazole-sulfone	TCB-SO ₂
Triclabendazole sulfoxide	TCB-SO

Deuterated forms of these standards are available from Witega & QUCHEM (Belfast, UK)

Janssen Animal Health (Beerse, Belgium)

Analyte**	Abbreviation
Amino-flubendazole	FLU-NH ₂
Amino-mebendazole	MBZ-NH ₂
Hydroxy-flubendazole	FLU-OH
Hydroxy-mebendazole	MBZ-OH
Flubendazole	FLU
Mebendazole	MBZ

Greyhound Chromatography and Allied Chemicals, (Merseyside, UK)

Analyte**	Abbreviation
Coumaphos-oxon	COUM-O

QMX Laboratories (Essex, UK)

Analyte**	Abbreviation
Cambendazole	CAM
Oxibendazole	OXI

Meril Animal Health (Lyon, France)

Analyte**	Abbreviation
Eprinomectin	EPR

Fort Dodge Animal Health (Princeton, NJ, USA)

Analyte**	Abbreviation
Moxidectin	MOXI

Non-Isotopically Labeled Internal Standards Used

Internal Standard	Abbreviation & Source
Selamectin	SELA (Pfizer Animal Health)
Salicylanide	SALI (Sigma-Aldrich)
4-nitro-3-(trifluoromethyl)phenol	TFM (Sigma-Aldrich)
Ioxynil	IOX (Sigma-Aldrich)

Primary Stock Standard Solutions:

- 4,000 µg/mL from the certified standard materials-- ABZ, ABZ-SO, ABZ-SO₂, ABZ-NH₂-SO₂, FBZ, OFZ, FBZ-SO₂, EPR, CLOS, OXY, NITR, CLOR, BITH and MOR
- The remaining standards are prepared at concentrations of 2,000 µg/mL
- All internal standards are prepared at concentration of 1,000 µg/mL
- Avermectins were prepared in MeCN
- Flukicides, CAM, LEVA and TCB metabolites are prepared in MeOH
- Benzimidazoles are prepared in DMSO

Intermediate working standard mix solutions:

- 100 µg/mL for OXY, CLOR, BITH and MOR
- 50 µg/mL in MeOH for the remaining analytes

Prepare working IS as follows:

- 20 µg/mL for SELA and TCB-NH₂, 4 µg/mL for LEVA-D₅, TBZ164 D₃ and IOX
- 2 µg/mL for the remaining analytes in MeOH- D

Primary, intermediate and working standard solutions are stable for at least six months when stored at -20°C

*Adapted and used with permission from Whelen, M., Kinsella, B., "Determination Of Anthelmintic Drug Residues In Milk Using Ultra High Performance Liquid Chromatography-Tandem Mass Spectrometry With Rapid Polarity Switching", doi:10.1016/j.chroma.2010.05.007, CHROMA 351049, J. of Chromatography A

**Listing of instrument manufacturers and standards suppliers does not constitute endorsement by UCT. Equivalent systems may be used



EURL-FV Multiresidue Method Using QuEChERS by GC-QqQ/MS/MS & LC-QqQ/MS/MS for Fruits & Vegetables

UCT Product Number:

ECQUEU750CT 50 mL centrifuge tube contains: (4 g MgSO₄, 1 g NaCl, 0.5 grams Na Citrate Dibasic Sesquihydrate, 1 g Na Citrate Tribasic Dihydrate)

ECMPS15CT 15 mL centrifuge tube contains: (900 mg MgSO₄ & 150 mg PSA) *(other configurations are available)*

January 2011

This summary of the European Union Reference Laboratory Residue method describes a QuEChERS approach for the analysis of 138 pesticides included in the Coordinated Multiannual Community Control Programme for 2010, 2011 and 2012 (Commission Regulation (EC) No 901/2009). Analysis is developed for **avocado, carrot, orange and pepper**.

Samples are prepared according to the Quality Control procedure established in the "Method Validation and Quality Control Procedures for Pesticide Residues Analysis in Food and Feed" (Document No. SANCO/10684/2009)

Procedure

1. Sample Preparation (for pesticides analyzed by **HPLC-MS/MS**)

- a) Homogenize the sample using a food processor according the typical QuEChERS procedures
- b) Weigh 10 g \pm 0.1 g of sample into a 50 mL centrifuge tube
- c) Add 10 mL of acetonitrile
- d) Shake vigorously or vortex for 1 minute to disperse contents
- e) Centrifuge for 5 minutes @ 4000 rpm

2. Clean-Up

- a) Transfer 6 mL of the supernatant to product **ECMPS15CT**
- b) Shake vigorously or vortex for 1 minute
- c) Centrifuge @ 6000 rpm for 2 minutes
- d) Transfer 1 mL of extract to a test tube and add 220 μ L of acetonitrile
- e) Using a 0.45 μ Teflon syringe filter, transfer extract to an LC injection vial

- f) Sample is ready for HPLC analysis

1A. Sample Preparation (For pesticides analyzed by GC-MS/MS)

- a) Procedure is the same as for HPLC analysis through steps 2) c)
- b) Transfer 1 mL of extract to a test tube
- c) Evaporate to dryness
- d) Add 1 mL of cyclohexane:acetone (9:1) to the dried extract
- e) Shake or vortex until completely dissolved
- f) Filter extract using a 0.45µ Teflon syringe filter into a GC vial
- g) Sample is ready for analysis by GC-MS/MS

**1) Instrumentation and Analytical Conditions for the LC/QqQ (MS/MS)
System**

- LC-MS/MS System 3200 Q TRAP, Applied Biosystems
- Column: Atlantis T3 2.1x100 mm, 3 µm
- Column temperature: 40 °C
 - Mobile phase A: H₂O, 2 mM ammonia formate, 0.1 % formic acid
 - Mobile phase B: methanol
 - Injection volume: 10µL
 - Autosampler temperature: 10° C
 - Analysis time: 18 min.

HPLC Flow Rate and Elution Gradient Table

Time (min)	A (%)	B (%)	Flow (µL/min)
0.0	95	5	300
1.0	95	5	300
1.1	70	30	300
10.0	0	100	300
13.0	0	100	300
13.1	95	5	300

3. Instrumentation and Analytical Conditions for the GC/QqQ (MS/MS)

- GC: Agilent 7890 Series or equivalent
- Autosampler: Agilent 7683 Injector and sample tray
 - Inlet: Splitless
 - Carrier gas: He
 - Inlet pressure: 22.73 psi
 - Inlet temperature: 250°C

- Injection volume: 1 µL
- Analytical column: Agilent J&W HP-5ms 30 m x 250 µm x 0.25 µm or equivalent
- Retention time locking: Chlorpyrifos methyl locked to 16.596 min
- Spectrometer: Agilent 7000B Series
- Source temperature: 280°C
- Quadrupole temperature: Q1 and Q2 = 150°C
- Collision gas flows: N₂ at 1.5 mL/min, He at 25 mL/min

GC Oven Temperature Program

	Rate (°C/min)	Value (°C)	Hold Time (min)	Run Time (min)
Initial		70	2	2
Ramp 1	25	150	0	5.2
Ramp 2	3	200	0	21.9
Ramp 3	8	280	10	41.9

Spike Level with Method Validation Results

Pesticide	Mean 0.01 mg/Kg	RSD 0.01 mg/Kg	Mean 0.1 mg/Kg	RSD 0.1 mg/Kg	Technique w MS/MS
Acephate	88	7	89	3	HPLC
Acetamiprid	84	10	102	4	HPLC
Acrinathrin	103	12	91	13	GC
Aldicarb(RD)	105	13	98	5	HPLC
Amitraz(RD)	89	8	80	7	HPLC
DMPF	84	7	92	5	HPLC
DMF	75	5	103	7	HPLC
Azinphos-methyl	86	17	100	5	HPLC
Azoxystrobin	87	6	104	4	HPLC
Bifenthrin	95	9	92	10	GC
Bitertanol	118	11	103	8	HPLC
Boscalid	98	8	100	4	HPLC
Bromopropylate	95	11	94	12	GC
Bupirimate	96	6	94	10	GC
Buprofezin	87	8	89	8	GC
Cadusafos	88	16	96	6	HPLC
Captan	76	11	85	9	GC
Carbaryl	110	9	97	5	HPLC
Carbendazim(rd)	96	6	95	5	HPLC
Carbofuran(rd)	91	14	103	4	HPLC
Chlorfenvinphos	99	8	86	27	GC

Chlorothalonil	56	4	67	3	GC
Chlorpropham(RD)	99	7	95	7	GC
Chlorpyrifos	98	6	92	7	GC
Chlorpyrifos-methyl	96	5	93	7	GC
Clofentezin(RD)	93	27	88	9	HPLC
Clothianidin	109	8	98	7	HPLC
Cyfluthrin(RD)	104	13	97	16	GC
Cypermethrin(RD)	109	11	109	16	GC
Cyproconazole	97	13	95	5	HPLC
Cyprodinil	97	7	91	9	GC
Deltamethrin	106	14	95	16	GC
Diazinon	99	7	91	7	GC
Dichlofluanid	68	5	71	2	GC
Dichlorvos	85	9	98	6	HPLC
Dicloran	111	17	97	8	GC
Difenoconazole	94	9	99	5	HPLC
Dimethoate(RD)	86	11	99	5	HPLC
Dimethomorph	107	10	97	5	HPLC
Diphenylamine	101	8	86	7	GC
Endosulfan(rd)	92	10	89	10	GC
Epoxiconazole	86	8	101	3	HPLC
Ethion	114	9	102	10	GC
Etofenprox	95	13	94	14	GC
Ethoprophos	98	4	92	6	GC
Fenarimol	98	12	90	13	GC
Fenazaquin	89	9	88	11	GC
Fenbutatin oxide		-	78	4	HPLC
Fenbuconazole	106	12	100	5	HPLC
Fenhexamid	79	9	94	5	HPLC
Fenitrothion	116	10	103	7	GC
Fenoxycarb	105	24	98	8	HPLC
Fenpropathrin	99	10	96	11	GC
Fenpropimorph	81	7	99	3	HPLC
Fenthion(RD)	92	7	90	8	GC
Fenthion sulfoxide	86	15	92	10	HPLC
Fenvalerate	96	15	93	15	GC
Fipronil(RD)	98	5	91	11	GC
Fludioxinil	132	3	96	2	HPLC
Flufenoxuron	83	18	110	10	HPLC
Fluquinconazole	94	21	99	7	HPLC
Flusilazole	98	8	93	10	GC
Flutriafol	99	8	102	5	HPLC
Folpet	59	21	65	13	GC
Formetanate(RD)	87	10	95	4	HPLC
Fosthiazate	79	20	114	8	HPLC

Hexaconazole	95	12	96	4	HPLC
Hexythiazox	104	15	97	12	HPLC
Imazalil	88	6	97	6	HPLC
Imidacloprid	96	17	100	6	HPLC
Indoxacarb(RD)	90	40	113	9	HPLC
Iprodione	106	13	94	13	GC
Iprovalicarb	105	8	99	4	HPLC
Kresoxim-methyl	103	14	107	5	HPLC
Lambda-cyhalothrin(RD)	108	12	96	13	GC
Linuron	97	29	99	6	HPLC
Lufenuron	132	16	110	12	HPLC
Malathion(RD)	107	6	100	9	GC
Mepanipyrim(RD)	103	12	95	12	GC
Metalaxyl(rd)	111	11	94	8	GC
Metconazole	101	9	95	4	HPLC
Methamidophos	91	13	89	6	HPLC
Methidathion	115	6	101	9	GC
Methiocarb(RD)	109	33	104	15	HPLC
Methomyl(RD)	109	7	106	5	HPLC
Methoxyfenozide	125	12	99	13	HPLC
Monocrotophos	79	10	98	6	HPLC
Myclobutanil	89	11	91	11	GC
Oxadixyl	94	12	88	9	GC
Oxamyl	96	8	96	5	HPLC
Oxydemeton-methyl(RD)	94	7	97	4	HPLC
Paclobutrazole	91	13	100	5	HPLC
Parathion	119	6	101	9	GC
Parathion-methyl(RD)	109	6	100	7	GC
Pencycuron	100	18	101	6	HPLC
Penconazole	97	11	97	5	HPLC
Pendimethalin	116	8	96	8	GC
Permethrin(rd)	98	12	93	13	GC
Phenthoate	109	7	99	8	GC
Phosalone	111	14	99	14	GC
Phosmet(RD)	100	10	99	7	HPLC
Pyraclostrobin	95	11	110	4	HPLC
Pirimicarb(RD)	92	6	94	7	GC
Pirimiphos-methyl	111	8	94	9	GC
Prochloraz(RD)	88	8	95	5	HPLC
Procymidone	94	8	93	8	GC
Profenofos	100	10	97	10	GC
Propamocarb(RD)	68	9	69	6	HPLC
Propargite	112	10	104	7	HPLC
Propiconazole	94	9	91	10	GC
Propyzamide	99	5	96	7	GC

Prothioconazole	78	24	33	12	HPLC
Pyridaben	95	14	92	13	GC
Pyrimethanil	119	13	95	6	GC
Pyriproxyfen	103	16	95	13	GC
Quinoxifen	114	8	87	10	HPLC
Spinosad(RD)	97	11	98	4	HPLC
Spiroxamine	124	15	80	15	HPLC
Taufluvallinate	102	13	96	17	GC
Tebuconazole	113	11	95	11	GC
Tebufenozide	124	45	97	13	HPLC
Tebufenpyrad	92	11	94	12	GC
Teflubenzuron	96	31	103	14	HPLC
Tefluthrin	89	5	89	7	GC
Tetraconazole	107	10	91	8	GC
Tetradifon	87	12	90	13	GC
Thiabendazole	92	9	93	8	HPLC
Thiamethoxam(RD)	84	16	101	5	HPLC
Thiacloprid	86	8	105	4	HPLC
Thiophanate-methyl	69	13	104	6	HPLC
Tolclofos-methyl	91	4	97	5	GC
Tolyfluanid(RD)	69	9	72	10	GC
Triadimenol(RD)	-	-	105	29	HPLC
Triazophos	117	7	102	11	GC
Trichlorfon	75	19	106	7	HPLC
Trifloxystrobin	93	13	103	7	HPLC
Triflumuron	-	-	121	6	HPLC
Trifluralin	92	3	88	6	GC
Triticonazole	104	14	97	5	HPLC
Vinclozolin(RD)	97	7	95	7	GC
Zoxamide	79	17	112	5	HPLC

Summarized and adapted from "EURL-FV Multiresidue Method using QuEChERS followed by GC-QqQ/MS/MS and LC-QqQ/MS/MS for Fruits and Vegetables," European Reference Laboratory in Pesticide Residue, 2009



Extraction of Pyrethrin and Pyrethroid Pesticides from Fish Using the QuEChERS Approach

UCT Product Number:

EC4MSSA50CT-MP (4000 mg MgSO₄ and 1000 mg sodium acetate)

CUMPSC18CT (150 mg MgSO₄, 50 mg PSA and 50 mg endcapped C18)

January 2011

The QuEChERS approach is used for the determination of trace levels of natural pyrethrins and synthetic pyrethroids (cypermethrin & deltamethrin) in fish.

1. Extraction

- a) Weigh 10 grams of homogenized fish into a 50 mL centrifuge tube
- b) Add 500 ng *cis*-permethrin (phenoxy-¹³C₆) surrogate standard
- c) Add 10 mL 1% acetic acid in acetonitrile
- d) Add the contents of pouch **EC4MSSA50CT-MP**
- e) Shake vigorously for 1 minute then centrifuge

2. Clean-up, Dispersive Solid-phase (dSPE)

- a) Transfer 1 mL of supernatant to a 2 mL micro-centrifuge tube
CUMPSC18CT
- b) Shake for 1 minute then centrifuge
- c) Transfer 0.5 mL of extract to a graduated tube then evaporate to near dryness
- d) Add 50 ng *trans*-permethrin (phenoxy-¹³C₆) and bring to exactly 0.5 mL with trimethyl phosphate (TMP)
- e) Add MgSO₄ to the 0.2 mL mark then vortex
- f) Transfer supernatant to injection vial for analysis

3. Analysis

- a) Use GC/MS in CI mode
- b) Column: HP-5, 30m X 0.32 mm with 0.25 µm film (or equivalent)
- c) Splitless mode @ 240°

GC Oven program:

- Initial 80°C, hold 1 minute
- 50°C/min to 200°C
- 5°C/min to 285°C
- 50°C/min to 325°C, hold 5 minutes
- Transfer line 250°C

MS Conditions:

- Source 150°C
- Methane reagent gas
- Selected Ion Monitoring Mode

Calibration using matrix matching may be required

***Adapted from** Roscoe, Veronica, Judge, Judy, Rawn, Dorothea F.K., "Application of the QuEChERS Extraction Method for the Analysis of Pyrethrin and Pyrethroid Pesticides in Fin and non-Fin Fish", Health Products and Food Program, Winnipeg, Manitoba and Bureau of Chemical Safety, Food Research Division, Ottawa, Ontario, Canada, Florida Pesticide Residue Workshop, July 2009



QuEChERS-Based LC/MS/MS Method for Multiresidue Pesticide Analysis in Fruits and Vegetables*

UCT Product Number:

EC4MSSA50CT-MP (4 g anhydrous MgSO₄, 1.0 g Sodium Acetate)

ECMS12CPSA415CT (1.2 g anhydrous MgSO₄, 400 mg PSA)

July 2011

A high-throughput, QuEChERS analytical method (LC-MS/MS) is described for the part per trillion (ppt) determination of 191 pesticides in orange, peach, spinach and ginseng. Pesticide classes include carbamates, polar organophosphates, phenylureas, anilides, benzoyl phenylureas, conazoles, macrocyclic lactone, neonicotinoids, strobilurines, and triazines. This method was validated by the U.S. Food and Drug Administration (FDA).

Analytes Covered in this Method

Table 1

Analyte	CASRN	Analyte	CASRN
Acephate	30560-19-1	Imazalil	35554-44-0
Acetamiprid	135410-20-7	Imidacloprid	138261-41-3
Acibenzolar-S-	135158-54-2	Indoxacarb	173584-44-6
Alanycarb	83130-01-2	Ipconazole	125225-28-7
Aldicarb	116-06-3	Iprovalicarb	140923-17-7
Aldicarb sulfone	1646-88-4	Isoprocarb	2631-40-5
Aldicarb sulfoxide	1646-87-3	Isoproturon	34123-59-6
Ametryn	834-12-8	Isoxaflutole	141112-29-0
Aminocarb	2032-59-9	Ivermectin	70288-86-7
Amitraz	33089-61-1	Kresoxim-methyl	143390-89-0
Avermectin B _{1a}	65195-55-3	Linuron	330-55-2
Avermectin B _{1b}	65195-56-4	Lufenuron	103055-07-8
Azoxystrobin	131860-33-8	Mefenacet	73250-68-7
Benalaxyl	71626-11-4	Mepanipyrim	110235-47-7
Bendiocarb	22781-23-3	Mepronil	55814-41-0
Benfuracarb	82560-54-1	Mesotrione	104206-82-8
Benzoximate	29104-30-1	Metalaxyl	57837-19-1
Bifenazate	149877-41-8	Metconazole.1	125116-23-6
Bitertanol	55179-31-2	Methabenzthiazuron	18691-97-9
Boscalid	188425-85-6	Methamidophos	10265-92-6
Bromuconazole 46	116255-48-2	Methiocarb	2032-65-7
Bromuconazole 47	116255-48-2	Methomyl	16752-77-5
Bupirimate	41483-43-6	Methoprotrotryne	841-06-5
Buprofezin	953030-84-7	Methoxyfenozide	161050-58-4
Butafenacil	134605-64-4	Metobromuron	3060-89-7

Butocarboxin	34681-10-2	Metribuzin	21087-64-9
Butoxycarboxin	34681-23-7	Mevinphos-E	813-78-5
Carbaryl	63-25-2	Mevinphos-Z	7786-34-7
Carbendazim	10605-21-7	Mexacarbate	315-18-4
Carbetamide	16118-49-3	Monocrotophos	6923-22-4
Carbofuran	1563-66-2	Monolinuron	1746-81-2
Carbofuran, 3OH-	16655-82-6	Moxidectin	113507-06-5
Carboxin	5234-68-4	Myclobutanil	88671-89-0
Carfentrazone-ethyl	128639-02-1	Neburon	555-37-3
Chlorfluazuron	71422-67-8	Nitenpyram	150824-47-8
Chlorotoluron	15545-48-9	Novaluron	116714-46-6
Chloroxuron	1982-47-4	Nuarimol	63284-71-9
Clethodim	99129-21-2	Omethoate	1113-02-6
Clofentezine	74115-24-5	Oxadixyl	77732-09-3
Clothianidin	210880-92-5	Oxamyl	23135-22-0
Cyazofamid	120116-88-3	Paclobutrazol	76738-62-0
Cycluron	2163-69-1	Penconazole	66246-88-6
Cymoxanil	57966-95-7	Phenmedipham	13684-63-4
Cyproconazole A	94361-06-5	Picoxystrobin	117428-22-5
Cyproconazole B	94361-07-6	Piperonyl butoxide	51-03-6
Cyprodinil	121552-61-2	Pirimicarb	23103-98-2
Desmedipham	13684-56-5	Prochloraz	67747-09-5
Diclobutrazol	75736-33-3	Promecarb	2631-37-0
Dicrotophos	141-66-2	Prometon	1610-18-0
Diethofencarb	87130-20-9	Prometryn	7287-19-6
Difenoconazole	119446-68-3	Propamocarb	24579-73-5
Diflubenzuron	35367-38-5	Propargite	2312-35-8
Dimethoate	60-51-5	Propham	122-42-9
Dimethomorph A	110488-70-5	Propiconazole	60207-90-1
Dimethomorph B	2274-67-1	Propoxur	114-26-1
Dimoxystrobin	149961-52-4	Pymetrozine	123312-89-0
Diniconazole	83657-24-3	Pyracarbolid	24691-76-7
Dioxacarb	6988-21-2	Pyraclostrobin	175013-18-0
Diuron	330-54-1	Pyridaben	96489-71-3
Doramectin	117704-25-3	Pyrimethanil	53112-28-0
Emamectin B_{1a}	155569-91-8	Pyriproxyfen	95737-68-1
Epoxiconazole	133855-98-8	Quinoxifen	124495-18-7
Eprinomectin B_{1a}	123997-26-2	Rotenone	83-79-4
Etaconazole	60207-93-4	Secbumeton	372137-35-4
Ethiofencarb	29973-13-5	Siduron	26259-45-0
Ethiprole	181587-01-9	Simetryne	1014-70-6
Ethofumesate	26225-79-6	Spinosyn A	168316-95-8
Etoxazole	153233-91-1	Spirodiclofen	148477-71-8
Famoxadone	131807-57-3	Spiromefesin	283594-90-1
Fenamidone	161326-34-7	Spiroxamine	118134-30-8

Fenarimol	60168-88-9	Sulfentrazone	122836-35-5
Fenazaquin	120928-09-8	Tebuconazole	107534-96-3
Fenbuconazole	114369-43-6	Tebufenozide	112410-23-8
Fenhexamid	126833-17-8	Tebufenpyrad	119168-77-3
Fenoxycarb	79127-80-3	Tebuthiuron	34014-18-1
Fenpropimorph	67564-91-4	Teflubenzuron	83121-18-0
Fenpyroximate	134098-61-6	Terbumeton	33693-04-8
Fenuron	134098-61-6	Terbutryn	886-50-0
Fludioxinil	131341-86-1	Tetraconazole	112281-77-3
Flufenacet	142459-58-3	Thiabendazole	148-79-8
Flufenoxuron	101463-69-8	Thiacloprid	111988-49-9
Fluometuron	2164-17-2	Thiamethoxam	153719-23-4
Fluoxastrobin	361377-29-9	Thidiazuron	51707-55-2
Fluquinconazole	136426-54-5	Thiobencarb	28249-77-6
Flusilazole	85509-19-9	Thiofanox	39196-18-4
Flutolanil	66332-96-5	Thiophanate-methyl	23564-05-8
Flutriafol	76674-21-0	Triadimefon	43121-43-3
Forchlorfenuron	68157-60-8	Triadimenol	55219-65-3
Formetanate HCl	22259-30-9	Tricyclazole	41814-78-2
Fuberidazole	3878-19-1	Trifloxystrobin	141517-21-7
Furalaxyl	57646-30-7	Triflumizole	99387-89-0
Furathiocarb	65907-30-4	Triflumuron	64628-44-0
Hexaconazole	79983-71-4	Triticonazole	131983-72-7
Hexythiazox	78587-05-0	Vamidotion	2275-23-2
Hydramethylnon	67485-29-4	Zoxamide	156052-68-5

Deuterium Isotope Internal Standards	
D10-Diazinon	D6-diuron
D6-Dichlorvos	D6-Linuron
D6-Dimethoate	D6-Malathion

CDN-Isotopes (Montreal, QC, Canada)

Analytical Stock Solutions

Prepare separate stock solutions of analytical standards, including the isotope labeled internal standards (ILIS) for individual compounds.

- Weigh 10-75 mg each and dissolve in 10 or 25 mL of acetonitrile, methanol, or methanol/water (50:50 v/v) in volumetric flasks
- Prepare intermediate solutions in 100mL volumetric flasks by mixing stock solutions

- Prepare five levels of matrix-matched calibration standards from intermediate solutions by using sample matrix extract and matrix buffer (20 mM ammonium formate) in concentrations of 1, 5, 10, 50, and 100 ppb
- Add the ILIS solution prior to sample preparation and use as an internal standard in the quantitative analysis

Procedure

1. Sample Preparation--orange, peach, spinach

- Weigh 10 ± 0.1 g of cryoground sample into 50 mL centrifuge tube
- Add 10 mL of 1% acetic acid in acetonitrile and contents of **EC4MSSA50CT-MP** pouch
- Shake by hand then add 200 μ L of surrogate solution and a steel ball
- Place on a Geno/Grinder shaker (or equivalent) for 1 min @ 1000 strokes/minute
- When shaking is complete centrifuge @ 4500 rpm for 5 min
- Transfer 9 mL of supernatant to a 15 mL centrifuge tube containing **ECMS12CPSA415CT**
- Shake on Geno/Grinder for 1 min @ 500 strokes/min
- Centrifuge @ 4500 rpm for 5 min
- Transfer 2.0 mL of supernatant to injection vials for analysis. Filter cloudy extracts using 0.2 nylon or PTFE membrane filter directly into the LC autosampler vials

2. Calibration Standards-- orange, peach, spinach

- Prepare matrix-matched calibration standards by mixing 300 μ L of 0.0167, 0.033, 0.067, 0.167, and 0.333 ppm standard solutions. Use 200 μ L of matrix blank extracts and 500 μ L of 20 mM ammonium formate sample buffer
- Add 500 μ L of sample buffer just prior to sample analysis
- Filter cloudy extracts using 0.2 nylon or PTFE membrane filter directly into the LC autosampler vials
- Filtered samples should be clear and can be stored in a freezer until analysis

1a. **Sample Preparation--ginseng**

- a) Prepare ginseng samples by using 1.0 ± 0.05 g of ginseng
- a) Add 10 mL of HPLC-grade water and a steel ball bearing
- b) Shake on a GenoGrinder at 1000 strokes/min for 1 minute
- c) Add 10 mL of 1% acetic acid in acetonitrile, 200 μ L of surrogate solution and contents of **EC4MSSA50CT-MP** pouch
- d) Shake by hand
- e) Place on a Geno/Grinder shaker (or equivalent) for 1 min @ 1000 strokes/minute
- f) When shaking is complete centrifuge @ 4500 rpm for 5 min
- g) Transfer 9 mL of supernatant to a 15 mL centrifuge tube containing **ECMS12CPSA415CT**
- h) Shake on Geno/Grinder for 1 min @ 500 strokes/min
- i) Centrifuge @ 4500 rpm for 5 min
- j) Transfer 2.0 mL of supernatant to injection vials for analysis. Filter cloudy extracts using 0.2 nylon or PTFE membrane filter directly into the LC autosampler vials

2a. **Calibration Standards--ginseng**

- a) Prepare matrix-matched calibration standards by adding 100 μ L of 0.033, 0.067, 0.167, 0.333, 0.8, 1.6 ppm standard solutions to 400 μ L of ginseng blank extracts
- b) Add 500 μ L of sample matrix buffer just prior to analysis to achieve matrix-matched calibration standards of 1.67, 3.33, 6.67, 16.7, 33.3, 80, and 160 ppb, respectively
- c) Filter using 0.2 m Nylon or PTFE membrane filters
- d) Filtered samples should be clear and can be stored in a freezer until analysis

3. **Sample Analysis**

- a) HPLC analysis with Shimadzu Prominence/20 series (Columbia, MD) or equivalent interfaced to an ABSciex (Forest City, CA) 4000QTrap mass spectrometer through an ESI interface (IonSpray)
- b) Acquire MRM data in positive ion mode

- c) Identify target pesticides using two specific MRM transitions for each pesticide to achieve an identification point (IP) of 4
- d) Quantify using either external standard calibration (NRCG) or internal standard calibration (FDA and MOE) with $^2\text{H}_{10}$ -diazinon as IS
- e) Use N_2 of 99% purity from a nitrogen generator (Parker Balston, Haverhill, MA) in the ESI source and the collision cell
- f) Restek LC column (Bellefonte, PA; Ultra Aqueous, C-18, 100 x 2.1 mm, 3 μm) and guard column (Ultra Aqueous, C-18 cartridges, 10 x 2.1 mm in guard cartridge holder) or equivalent
- g) Curtain, collision, nebulizer, auxiliary gases, and source temperature of the ESI source were set at 15, 6, 35, and 45 psi and 450° C, respectively
- h) Ion spray voltage: 5200
- i) Declustering potential (DP), collision energy (CE), and collision cell exit potential (CXP) are optimized by direct infusion. The two most intense ion pairs of each analyte are chosen for the analysis. Values of DP, CE, and CXP and the two specific, most intense MRM pairs are listed in Table 3. Principal component analysis (PCA) is carried out using Infometrix Pirouette 4 (Bothell, WA)
- j) Table 2 lists mobile phases, column temperatures, injection volume, flow rate, and LC gradient parameters

Table 2

HPLC Gradient Elution Parameters	
Mobil Phase	A: 5 mM ammonium formate, 0.1% formic acid in water
	B: 5 mM ammonium formate, 0.1% formic acid in MeOH
Column Temperature	35° C
Flow rate	0.3 mL/min
Total run time	14.0 min
Gradient program	10% B at 0 min, hold for 1 min 5% B at 0 min 20% B at 0 min to 98%
Injection volume	20 μL

Table 3

DP: declustering potential, V; CE: collision energy, V; CXP: collision cell exit potential

Pesticide	Formula	Mol Wt	MRM	DP	CE	CXP
			Transitions #1 & #2			
Carbofuran, 3OH-	C ₁₂ H ₁₅ NO ₄	237	238→163 / 181	66	21	16
Acephate	C ₄ H ₁₀ NO ₃ PS	183	184→143 / 49	61	13	4
Acetamiprid	C ₁₀ H ₁₁ N ₃ CIN ₄	223	223→126 / 99	61	29	12
Acibenzolar-S-methyl	C ₈ H ₆ N ₂ OS ₂	210	211→136 / 140	46	39	9
Alanycarb	C ₁₇ H ₂₅ N ₃ O ₄ S ₂	400	400→238 / 91	35	14	4
Aldicarb sulfoxide	C ₇ H ₁₄ N ₂ O ₃ S	206	207→132 / 89	30	10	8
Aldicarb	C ₇ H ₁₄ N ₂ O ₂ S	190	208→116 / 89	36	11	10
Aldicarb sulfone	C ₇ H ₁₄ N ₂ O ₄ S	222	223→86 / 148	52	21	5
Ametryn	C ₉ H ₁₇ N ₅ S	227	209→152 / 137	71	21	8
Aminocarb	C ₁₁ H ₁₆ N ₂ O ₂	208	209→152 / 137	71	21	8
Amitraz	C ₁₉ H ₂₃ N ₃	293	294→163 / 107	46	21	4
Avermectin B1a	C ₄₈ H ₇₂ O ₁₄	873	895→751 / 449	176	61	20
Avermectin B1b	C ₄₈ H ₇₀ O ₁₄	859	890→567 / 305	76	23	18
Azoxystrobin	C ₂₂ H ₁₇ N ₃ O ₅	403	404→372 / 344	51	19	4
Benalaxyl	C ₂₀ H ₂₃ NO ₃	325	326→148 / 294	71	31	8
Bendiocarb	C ₁₁ H ₁₃ NO ₄	223	224→109 / 167	61	27	20
Benfuracarb	C ₂₀ H ₃₀ N ₂ O ₅ S	411	411→195 / 252	50	30	4
Benzoximate	C ₁₈ H ₁₈ CINO ₅	364	364→199 / 105	51	13	14
Bifenazate	C ₁₇ H ₂₀ N ₂ O ₃	300	301→170 / 198	61	29	10
Bitertanol	C ₂₀ H ₂₃ N ₃ O ₂	337	338→70 / 269	51	31	12
Boscalid	C ₁₈ H ₁₂ Cl ₂ N ₂ O	343	343→307 / 140	91	27	4
Bromuconazole 46	C ₁₃ H ₁₂ BrCl ₂ N ₃ O	377	378→159 / 70	61	37	14
Bromuconazole 47	C ₁₃ H ₁₂ BrCl ₂ N ₃ O	377	378→159 / 70	61	37	14
Bupirimate	C ₁₃ H ₂₄ N ₄ O ₃ S	316	317→166 / 108	86	33	12
Buprofezin	C ₁₆ H ₂₃ N ₃ OS	305	306→201 / 116	46	17	4
Butafenacil	C ₂₀ H ₁₈ ClF ₃ N ₂ O ₆	475	492→331 / 349	61	35	20
Butocarboxin	C ₇ H ₁₄ N ₂ O ₂ S	190	213→75 / 116	50	20	5
Butoxycarboxin	C ₇ H ₁₄ N ₂ O ₄ S	222	223→106 / 166	45	15	8
Carbaryl	C ₁₂ H ₁₁ NO ₂	201	202→145 / 127	56	15	10
Carbendazim	C ₉ H ₉ N ₃ O ₂	191	192→160 / 132	80	24	10
Carbetamide	C ₁₂ H ₁₆ N ₂ O ₃	236	237→192 / 118	56	13	12
Carbofuran	C ₁₂ H ₁₅ NO ₃	221	222→123 / 165	66	31	22
Carboxin	C ₁₂ H ₁₃ NO ₂ S	235	484→452 / 285	66	23	14
Carfentrazone-ethyl	C ₁₃ H ₁₀ Cl ₂ F ₃ N ₃ O ₃	412	412→346 / 366	81	31	4
Chlorfluazuron	C ₂₀ H ₉ Cl ₃ F ₅ N ₃ O ₃	541	540→158 / 383	91	27	4
Chlorotoluron	C ₁₀ H ₁₃ CIN ₂ O	213	213→72 / 46	61	31	4
Chloroxuron	C ₁₅ H ₁₅ CINO ₂	291	291→72 / 218	65	30	4
Clethodim	C ₁₇ H ₂₆ CINO ₃ S	360	360→164 / 268	61	29	10
Clofentezine	C ₁₄ H ₈ Cl ₂ N ₄	303	303→138 / 102	61	23	8
Clothianidin	C ₆ H ₈ CIN ₅ O ₂ S	250	250→169 / 132	51	17	4
Cyazofamid	C ₁₃ H ₁₃ CIN ₄ O ₂ S	325	325→108 / 261	61	21	10
Cycluron	C ₁₁ H ₂₂ N ₂ O	198	199→89 / 72	50	21	4
Cymoxanil	C ₇ H ₁₀ N ₄ O ₃	198	199→128 / 111	60	13	4
Cyproconazole A	C ₁₅ H ₁₈ CIN ₃ O	292	292→70 / 125	66	39	12
Cyproconazole B	C ₁₅ H ₁₈ CIN ₃ O	292	292→70 / 125	66	39	12
Cyprodinil	C ₁₄ H ₁₅ N ₃	225	226→93 / 77	101	51	16
Desmedipham	C ₁₆ H ₁₆ N ₂ O ₄	300	318→182 / 136	41	19	12
Diclobutrazol	C ₁₅ H ₁₉ Cl ₂ N ₃ O	328	328→70 / 158	81	49	12
Dicrotophos	C ₈ H ₁₆ NO ₅ P	237	238→112 / 193	66	19	8
Diethofencarb	C ₁₄ H ₂₁ NO ₄	267	268→226 / 124	61	15	14

Difenoconazole	C ₁₉ H ₁₇ Cl ₂ N ₃ O ₃	406	406→251 / 253	81	37	16
Diffubenzuron	C ₁₄ H ₉ Cl ₂ FN ₂ O ₂	311	311→158 / 141	71	23	10
Dimethoate	C ₅ H ₁₂ NO ₃ PS ₂	229	230→199 / 125	50	14	15
Dimethomorph A	C ₂₁ H ₂₂ CINO ₄	388	388→301 / 165	66	25	4
Dimethomorph B	C ₂₁ H ₂₂ CINO ₄	388	388→301 / 165	66	25	4
Dimoxystrobin	C ₁₉ H ₂₂ N ₂ O ₃	326	327→205 / 116	40	15	4
Diniconazole	C ₁₅ H ₁₇ Cl ₂ N ₃ O	326	326→70 / 158	86	51	12
Dioxacarb	C ₁₁ H ₁₃ NO ₄	223	224→167 / 123	51	13	10
Diuron	C ₉ H ₁₀ Cl ₂ N ₂ O	233	233→72 / 72	56	33	4
Doramectin	C ₅₀ H ₇₄ O ₁₄	899	921→777 / 449	71	65	15
Fenpyroximate	C ₂₄ H ₂₇ N ₃ O ₄	422	422→366 / 135	56	23	4
Emamectin B _{1a}	C ₄₉ H ₇₅ NO ₁₃	886	886→158 / 82	111	51	10
Epoxiconazole	C ₁₇ H ₁₃ ClFN ₃ O	330	330→121 / 101	66	29	10
Eprinomectin B _{1a}	C ₅₀ H ₇₅ NO ₁₄	914	914→186 / 154	76	27	12
Etaconazole	C ₁₄ H ₁₅ Cl ₂ N ₃ O ₂	328	328→159 / 205	46	37	10
Ethiofencarb	C ₁₁ H ₁₅ NO ₂ S	225	226→106 / 164	41	21	4
Ethiprole	C ₁₃ H ₉ Cl ₂ F ₃ N ₄ OS	397	397→350 / 255	81	29	24
Ethofumesate	C ₁₃ H ₁₈ O ₅ S	286	287→121 / 259	81	23	8
Etoxazole	C ₂₁ H ₂₃ F ₂ NO ₂	359	360→141 / 57	76	45	4
Famoxadone	C ₂₂ H ₁₈ N ₂ O ₄	374	392→331 / 238	31	15	4
Fenamidone	C ₁₇ H ₁₇ N ₃ OS	311	312→92 / 236	66	39	16
Fenarimol	C ₁₇ N ₁₂ Cl ₂ N ₂ O	331	331→268 / 81	61	31	4
Fenazaquin	C ₂₀ H ₂₂ N ₂ O	306	307→161 / 147	71	25	12
Fenbuconazole	C ₁₉ H ₁₇ CIN ₄	337	337→124 / 70	81	41	8
Fenhexamid	C ₁₄ H ₁₇ Cl ₂ NO ₂	302	302→97 / 55	66	35	18
Fenoxycarb	C ₁₇ H ₁₉ NO ₄	301	302→88 / 116	66	31	6
Fenpropimorph	C ₂₀ H ₃₃ NO	303	304→147 / 117	66	39	4
Fenuron	C ₉ H ₁₂ N ₂ O	164	165→72 / 46	56	25	4
Fludioxinil	C ₁₂ H ₆ F ₂ N ₂ O ₂	248	266→229 / 227	41	23	14
Flufenacet	C ₁₄ H ₁₃ F ₄ N ₃ O ₂ S	363	364→152 / 194	51	29	10
Flufenoxuron	C ₂₁ H ₁₁ ClF ₆ N ₂ O ₃	489	489→158 / 141	86	29	10
Fluometuron	C ₁₀ H ₁₁ F ₃ N ₂ O	232	233→72 / 46	71	37	12
Fluoxastrobin	C ₂₁ H ₁₆ ClFN ₄ O ₅	459	459→427 / 188	55	28	4
Fluquinconazole	C ₁₆ H ₈ Cl ₂ FN ₅ O	376	376→307 / 349	71	33	4
Flusilazole	C ₁₆ H ₁₅ F ₂ N ₃ Si	315	316→247 / 165	81	27	16
Flutolanil	C ₁₇ H ₁₆ F ₃ NO ₂	323	324→262 / 242	76	27	16
Flutriafol	C ₁₆ H ₁₃ F ₂ N ₃ O	301	302→70 / 123	66	37	12
Forchlorfenuron	C ₁₂ H ₁₀ CIN ₃ O	248	248→129 / 93	52	25	4
Formetanate HCl	C ₁₁ H ₁₅ N ₃ O ₂	221	222→165 / 120	60	21	12
Fuberidazole	C ₁₁ H ₈ N ₂ O	184	185→157 / 65	81	33	14
Furalaxyl	C ₁₇ H ₁₉ NO ₄	301	302→95 / 242	56	41	18
Furathiocarb	C ₁₈ H ₂₆ N ₂ O ₅ S	382	383→195 / 252	76	27	12
Hexaconazole	C ₁₄ H ₁₇ Cl ₂ N ₃ O	314	314→70 / 159	56	41	12
Hexaflumuron	C ₁₆ H ₈ Cl ₂ F ₆ N ₂ O ₃	461	461→158 / 141	56	25	4
Hexythiazox	C ₁₇ H ₂₁ CIN ₂ O ₂ S	353	353→228 / 168	61	23	14
Hydramethylnon	C ₂₅ H ₂₄ F ₆ N ₄	494	495→323 / 151	146	45	20
Imazalil	C ₁₄ H ₁₄ Cl ₂ N ₂ O	297	297→159 / 201	66	33	14
Imidacloprid	C ₉ H ₁₀ CIN ₅ O ₂	256	256→209 / 175	61	23	12
Indoxacarb	C ₂₂ H ₁₇ ClF ₃ N ₃ O ₇	528	528→203 / 218	86	55	12
Ipconazole	C ₁₈ H ₂₄ CIN ₃ O	334	334→70 / 125	76	55	12
Iprovalicarb	C ₁₈ H ₂₈ N ₂ O ₃	320	321→119 / 203	66	29	8
Isoprocarb	C ₁₁ H ₁₅ NO ₂	193	194→95 / 137	61	23	16
Isoproturon	C ₁₂ H ₁₈ N ₂ O	206	207→72 / 46	66	29	4
Isoxaflutole	C ₁₅ H ₁₂ F ₂ NO ₄ S	359	377→251 / 360	36	41	16
Ivermectin	C ₄₈ H ₇₄ O ₁₄	875	897→754 / 610	65	65	8
Kresoxim-methyl	C ₁₈ H ₁₉ NO ₄	313	314→116 / 206	51	21	4
Linuron	C ₉ H ₁₀ Cl ₂ N ₂ O ₂	249	249→160 / 182	61	23	4

Lufenuron	C ₉ H ₁₀ Cl ₂ N ₂ O ₂	511	511→158 / 141	61	27	4
Mefenacet	C ₁₆ H ₁₄ N ₂ O ₂ S	298	299→148 / 120	56	21	10
Mepanipirim	C ₁₄ H ₁₃ N ₃	223	224→106 / 77	86	37	8
Mepronil	C ₁₇ H ₁₉ NO ₂	269	270→119 / 228	76	33	8
Mesotrione	C ₁₄ H ₁₃ NO ₇ S	339	357→228 / 288	60	31	9
Metalaxyl	C ₁₅ H ₂₁ NO ₄	279	280→220 / 192	61	21	14
Metconazole.1	C ₁₇ H ₂₂ ClN ₃ O	319	320→70 / 125	81	51	12
Methabenzthiazuron	C ₁₀ H ₁₁ N ₃ OS	221	222→165 / 150	51	21	4
Methamidophos	C ₂ H ₆ NO ₂ PS	141	142→94 / 125	55	20	4
Methiocarb	C ₁₁ H ₁₅ NO ₂ S	225	226→169 / 121	61	13	12
Methomyl	C ₅ H ₁₀ N ₂ O ₂ S	162	163→88 / 106	35	12	5
Methoprotrotyne	C ₁₁ H ₂₁ N ₅ OS	271	272→240 / 198	50	27	4
Methoxyfenozide	C ₂₂ H ₂₈ N ₂ O ₃	368	369→149 / 313	56	25	10
Metobromuron	C ₉ H ₁₁ BrN ₂ O ₂	259	259→170 / 148	56	23	4
Metribuzin	C ₈ H ₁₄ N ₄ OS	214	215→84 / 187	71	29	4
Mevinphos-Z	C ₇ H ₁₃ O ₆ P	224	225→127 / 193	55	20	8
Mevinphos-E	C ₇ H ₁₃ O ₆ P	224	225→127 / 193	55	20	8
Mexacarbate	C ₁₂ H ₁₈ N ₂ O ₂	222	223→166 / 151	66	23	12
Monocrotophos	C ₇ H ₁₄ NO ₅ P	223	224→127 / 98	51	23	12
Monolinuron	C ₉ H ₁₁ ClN ₂ O ₂	215	215→126 / 99	51	23	4
Moxidectin	C ₃₇ H ₅₃ NO ₈	640	662→549 / 467	90	45	16
Myclobutanil	C ₁₅ H ₁₇ ClN ₄	289	289→70 / 125	71	37	12
Neburon	C ₁₂ H ₁₆ Cl ₂ N ₂ O	275	275→88 / 114	56	23	4
Nitenpyram	C ₁₁ H ₁₅ ClN ₄ O ₂	271	271→225 / 126	51	17	14
Novaluron	C ₁₇ H ₉ ClF ₈ N ₂ O ₄	493	493→158 / 141	71	27	4
Nuarimol	C ₁₇ H ₁₂ ClFN ₂ O	315	315→252 / 81	81	31	16
Omethoate	C ₅ H ₁₂ NO ₄ PS	213	214→124 / 182	46	29	4
Oxadixyl	C ₁₄ H ₁₈ N ₂ O ₄	278	279→219 / 132	61	17	14
Oxamyl	C ₇ H ₁₃ N ₃ O ₃ S	219	237→72 / 90	36	25	4
Paclobutrazol	C ₁₅ H ₂₀ ClN ₃ O	294	294→70 / 125	66	49	12
Penconazole	C ₁₃ H ₁₅ Cl ₂ N ₃	284	284→159 / 70	71	39	10
Phenmedipham	C ₁₆ H ₁₆ N ₂ O ₄	300	301→136 / 168	50	26	4
Picoxystrobin	C ₁₈ H ₁₆ F ₃ NO ₄	367	368→145 / 205	56	27	4
Piperonyl butoxide	C ₁₉ H ₃₀ O ₅	338	356→177 / 119	51	19	10
Pirimicarb	C ₁₁ H ₁₈ N ₄ O ₂	238	239→72 / 182	66	35	12
Prochloraz	C ₁₅ H ₁₆ Cl ₃ N ₃ O ₂	377	376→308 / 70	46	17	10
Promecarb	C ₁₂ H ₁₇ NO ₂	207	208→109 / 151	36	23	8
Prometon	C ₁₀ H ₁₉ N ₅ O	225	226→142 / 86	76	33	10
Prometryn	C ₁₀ H ₁₉ N ₅ S	241	242→200 / 158	71	19	4
Propamocarb	C ₉ H ₂₀ N ₂ O ₂	188	189→102 / 144	61	25	8
Propargite	C ₁₉ H ₂₆ O ₄ S	350	368→231 / 175	46	15	14
Propham	C ₁₀ H ₁₃ NO ₂	179	180→138 / 120	36	13	10
Propiconazole	C ₁₅ H ₁₇ Cl ₂ N ₃ O ₂	342	342→159 / 69	61	39	10
Propoxur	C ₁₁ H ₁₅ NO ₃	209	210→111 / 168	39	19	6
Pymetrozine	C ₁₀ H ₁₁ H ₅ O	217	218→105 / 78	71	27	4
Pyracarbolid	C ₁₃ H ₁₅ NO ₂	217	218→125 / 97	61	27	8
Pyraclostrobin	C ₁₉ H ₁₈ ClN ₃ O ₄	388	388→194 / 163	31	19	4
Pyridaben	C ₁₉ H ₂₅ ClN ₂ OS	365	365→147 / 309	46	31	4
Pyrimethanil	C ₁₂ H ₁₃ N ₃	199	200→107 / 82	71	33	4
Pyriproxyfen	C ₂₀ H ₁₉ NO ₃	321	322→96 / 185	46	21	4
Quinoxifen	C ₁₅ H ₈ Cl ₂ FNO	308	308→162 / 197	81	65	10
Rotenone	C ₂₃ H ₂₂ O ₆	394	395→213 / 192	91	33	14
Secbumeton	C ₁₀ H ₁₅ N ₅ O	225	226→170 / 100	50	35	4
Siduron	C ₁₄ H ₂₀ N ₂ O	232	233→137 / 94	66	21	4
Simetryne	C ₈ H ₁₅ N ₅ S	213	214→124 / 144	51	27	4
Spinosyn A	C ₄₁ H ₆₅ NO ₁₀	732	748→142 / 98	86	45	8
Spirodiclofen	C ₂₁ H ₂₄ Cl ₂ O ₄	411	411→313 / 71	71	17	8

Spiromefesin	C ₂₃ H ₃₀ O ₄	370	371→273 / 255	71	19	8
Spiroxamine	C ₁₈ H ₃₅ NO ₂	297	298→144 / 100	76	29	12
Sulfentrazone	C ₁₁ H ₁₀ Cl ₂ F ₂ N ₄ O ₃ S	387	387→307 / 146	81	27	4
Tebuconazole	C ₁₆ H ₂₂ CIN ₃ O	308	308→70 / 125	81	49	12
Tebufenozide	C ₂₂ H ₂₈ N ₂ O ₂	352	353→133 / 297	51	25	10
Tebufenpyrad	C ₁₈ H ₂₄ CIN ₃ O	334	334→117 / 145	71	47	4
Tebuthiuron	C ₉ H ₁₆ N ₄ OS	228	229→172 / 116	46	21	4
Teflubenzuron	C ₁₄ H ₆ Cl ₂ F ₄ N ₂ O ₂	381	381→141 / 158	66	53	4
Terbumeton	C ₁₀ H ₁₉ N ₅ O	225	226→170 / 100	76	27	12
Terbutryn	C ₁₀ H ₁₉ N ₅ S	241	242→186 / 68	71	27	12
Tetraconazole	C ₁₃ H ₁₁ Cl ₂ F ₄ N ₃ O	372	372→159 / 70	76	45	10
Thiabendazole	C ₁₀ H ₇ N ₃ S	201	202→175 / 131	85	35	12
Thiacloprid	C ₁₀ H ₉ CIN ₄ S	253	253→126 / 99	71	31	10
Thiamethoxam	C ₈ H ₁₀ CIN ₅ O ₃ S	292	292→211 / 181	61	19	12
Thidiazuron	C ₉ H ₈ N ₄ OS	220	221→102 / 127	66	21	4
Thiobencarb	C ₁₂ H ₁₆ CINOS	258	258→125 / 89	56	27	8
Thiofanox	C ₉ H ₁₈ N ₂ O ₂ S	218	219→76 / 57	36	20	8
Thiophanate-methyl	C ₁₂ H ₁₄ N ₄ O ₄ S ₂	342	343→151 / 311	61	29	14
Triadimefon	C ₁₄ H ₁₆ CIN ₃ O ₂	294	294→197 / 225	66	23	14
Triadimenol	C ₁₄ H ₁₈ CIN ₃ O ₂	296	296→70 / 227	46	31	12
Tricyclazole	C ₉ H ₇ N ₃ S	189	190→163 / 136	81	33	10
Trifloxystrobin	C ₂₀ H ₁₉ F ₃ N ₂ O ₄	408	409→186 / 206	31	23	4
Triflumizole	C ₁₅ H ₁₅ ClF ₃ N ₃ O	346	346→278 / 73	51	15	8
Triflumuron	C ₁₅ H ₁₀ ClF ₃ N ₂ O ₃	359	359→156 / 139	51	23	4
Triticonazole	C ₁₇ H ₂₀ CIN ₃ O	318	318→70 / 125	66	45	12
Vamidothion	C ₈ H ₁₈ NO ₄ PS ₂	287	288→146 / 118	61	19	10
Zoxamide	C ₁₄ H ₁₆ Cl ₃ NO ₂	337	336→187 / 159	45	35	15
D10-Diazinon	C ₁₂ D ₁₀ H ₁₁ N ₂ O ₃ PS	314	315→170	50	29	4
D6-Dimethoate	C ₅ D ₆ H ₆ NO ₃ PS ₂	235	236→131	50	30	4
D6-diuron	C ₉ D ₆ H ₄ Cl ₂ N ₂ O	239	239→78	90	30	4
D6-Linuron	C ₉ D ₆ H ₄ Cl ₂ N ₂ O ₂	255	255→166	90	30	4
D6-Dichlorvos	C ₄ D ₆ H ₁ Cl ₂ O ₄ P	227	227→115	70	27	4
D6-Malathion	C ₁₀ D ₆ H ₁₃ O ₆ PS ₂	330	337→291	55	12	4

*Summarized with permission from Wong, Jon, Hao, Chunyan, Zhang, Kai, et al., J. Agric. Food Chem. 2010, 58, 5897–5903 5897, DOI:10.1021/jf903849n

Listing of instrument manufacturers does not constitute endorsement by UCT



A Summary of US FDA LIB 4465: Collaboration of the QuEChERS Procedure for the Multiresidue Determination of Pesticides in Raw Agricultural Commodities by LC/MS/MS

UCT Product Numbers:

ECMSSC50CTFS-MP (6000 mg anhydrous magnesium sulfate, 1500 mg sodium chloride)

CUMPS2CT (150 mg anhydrous magnesium sulfate, 50 mg PSA)

ECMS12CPSA415CT (1200 mg anhydrous magnesium sulfate, 400 mg PSA)

March 2013

Method Summary

The analysis of fruits and vegetables for 173 pesticides using a single level calibration standard has been demonstrated to be an effective screening tool and can be completed in less than 20 minutes with overall accuracy of 105% and precision of 3% RSD. Pesticides are selected from a broad range of classes representing *carbamates*, *mectins*, *azoles*, *neonicotinoids*, *benzimidazoles*, *phenylureas*, *strobilurins*, *organophosphorus*, *anilides*, *tetrazines*, *anilides*, *benzoylphenylureas*, and *others*.

Procedure

Sample Preparation

Samples are composited by grinding in a vertical cutter mixed with dry ice

1. Sample Extraction

- a) Weigh 15 g of hydrated sample into the 50 ml centrifuge tube
- b) Add 15 mL acetonitrile (ACN)

Note: Adjust ACN volume of spike samples to account for spike solution volume to maintain ratio of 1g sample/mL of ACN, e.g. for 5 ml spike volume add 10 mL ACN to 15 g sample

- c) Shake for 1 min
- d) Add internal standard
- e) Add spike standard if needed
- f) Add the contents of pouch **ECMSSC50CTFS**

- g) Shake 1 min
- h) Centrifuge @ ~4500 rpm for 5 min

2. **PSA Clean-up**

- a) Transfer 1.0 mL of extract to **CUMPS2CT** (or alternative, step b)
- b) Transfer all extract to **ECMS12CPSA415CT**
- c) Vortex and centrifuge
- d) Dilute 0.5 mL extract to 5.0 mL with LC-MS aqueous buffer
- e) Filter through 0.2 or 0.45 µm Nylon filter
- f) Sample is ready for analysis

LC-MS/MS---Instrumentation

- AB Sciex 4000 QTrap: scheduled MRM in the positive ionization mode
- Shimadzu High Pressure HPLC System
- LC-20AD Pump
- Sil-20AC Autosampler
- CTO-20AC Column oven

HPLC Columns

- Ultra Aqueous C18, 3 µm, 100 x 2.1 mm with 10 x 2.1 mm guard column (Restek)

HPLC Instrument Parameters

Equilibration time (min)	1.5
Injection volume (µL)	20
Total Flow (mL/min)	0.5
Rinsing volume (µL)	200
Rinsing speed (µL/sec)	35
Sampling speed (µL/sec)	15
Cooler temperature (°C)	15
Column oven temp (°C)	40

Standards

Pesticide standard mixes may be purchased from AccuStandards and consist of 9 mixes of 20-25 analytes (total of 196 compounds)

The following injection and spiking standards were prepared in acetonitrile from the 3.0 µg/mL mixture of all standards:

Injection Standard: 200 ng/mL

Internal Standard: 200 ng/mL BDMC

Spike standards: 3000, 1200, 300, and 60 ng/mL

HPLC Mobile Phase Composition

Pump A: Water with 4 mM ammonium formate and 0.1 % formic acid

Pump B: Methanol with 4 mM ammonium formate and 0.1 % formic acid

Time	Parameter
Min	% B
0.0	5
1.0	5
9.0	95
11.3	95
12.0	5
13.4	5
13.5	stop

Mass Spectrometer Parameters

Typical MS Settings

MRM Detection Window (sec)	60
Target Scan Time (sec)	.5
Resolution Q1	unit
Resolution Q2	unit
MR Pause (msec)	5
Collision gas	med
Curtin gas (mL/min)	30
Exit Potential (volts)	10
Ion Source gas 1 (mL/min)	50
Ion Source gas 2 (mL/min)	50
Interface heater	on
Ion Spray (Volts)	5000
Turbo Spray T (°C)	400

MS/MS Transition Parameters

Compound	Transition 1					Transition 2				
	Q1	Q2	DP	CE	EXP	Q1	Q2	DP	CE	EXP
3-Hydroxycarbofuran	238.1	163	66	21	15	238.1	181	66	16	11
Acephate	184.1	143	61	13	5	184.1	49	61	33	6
Acetamiprid	223	126	60	29	10	223	99	60	51	14
Acibenzolar-S-methyl	211	136	46	39	8	211	140	46	31	8
Alanycarb	400.1	238.2	35	14	5	400.1	91.1	35	40	5
Aldicarb+NH4	208.1	116	35	11	10	208.1	89	35	23	16
Aldicarb Sulfoxide	207.1	132.1	30	10	8	207.1	89.1	30	19	6
Aldoxycarb	223.1	86.1	52	21	5	223.1	148	52	13	9
Aminocarb	209.1	152	71	21	8	209.1	137.1	71	35	10
Amitraz	294.2	163.2	46	21	4	294.2	107.1	46	57	4
AvermectinB1a+NH4	890.9	567.7	75	23	18	890.9	305.4	72	35	22
AvermectinB1b+Na	876.5	291	41	35	4	876.5	145	41	43	4
Azoxystrobin	404.1	372.1	51	19	5	404.1	344.1	51	27	5
BDMC	260	122	52	34	5	260	107	52	54	5
Benalaxyl	326.2	148.1	71	31	8	326.2	294.1	71	17	10
Bendiocarb	224.1	109	61	27	20	224.1	167.1	61	15	12
Benfuracarb	411.2	195.1	50	30	5	411.2	252.1	50	19	5
Bentazon	241	199	76	19	8	241	107	76	39	8
Benzoximate	364	199	51	13	13	364	105	51	35	4
Bifenazate	301.1	170.1	59	30	9	301.1	198.1	59	21	10
Bitertanol	338.2	70	51	31	12	338.2	269.2	48	13	14
Boscalid	343	307	90	27	7	343	140	90	27	6
BromuconazoleA	378	159	61	39	12	378	70	61	43	12
BromuconazoleB	378.1	159.1	61	39	12	378.1	70.1	61	43	12
Bupirimate	317	166.1	86	33	12	317	108	86	37	10
Buprofezin	306.2	201.1	46	17	5	306.2	116.2	46	21	5
Butafenacil+NH4	492.1	331	58	33	16	492.1	349	61	21	12
Butocarboxim+Na	213.1	75	50	21	6	213.1	116	50	13	6
Butoxycarboxin	223.1	106	45	15	8	223.1	166	45	11	5
Carbaryl	202.1	145	57	15	9	202.1	127	54	41	8
Carbendazim	192.2	160.2	80	24	10	192.2	132.1	80	41	7
Carbetamide	237.1	192	55	13	10	237.1	118.1	56	19	10
Carbofuran	222.1	123	66	31	19	222.1	165.1	66	19	11
Chlorantraniliprole	484	452.9	66	23	14	484	285.9	66	19	16
Chlorfluazuron	540	158	91	27	4	540	383	91	28	4
Chlorotoluron	213.1	72.2	61	31	5	213.1	46.2	61	27	5
Chloroxuron	291.1	72.4	65	34	5	291.1	218.1	65	30	5
Clethodim	360.1	164	61	28	9	360.1	268.1	61	17	8
Clofentezine	303	138	65	22	8	303	102	65	51	14
Clothianidin	250	169	51	17	4	250	132	51	21	10
Cyazofamid	325	108	60	20	9	325	261.1	60	15	13
Cycluron	199.1	89.1	50	21	5	199.1	72.2	50	21	4
Cyflufenamid	413.1	295.1	56	23	8	413.1	223.1	56	33	14
Cymoxanil	199	128	60	13	5	199	111	60	25	5
CyproconazoleA	292	70	63	37	10	292	125	63	43	8

CyproconazoleB	292.1	70.1	63	37	10	292.1	125.1	63	43	8
Cyprodinil	226	93	95	49	13	226	77	95	64	12
Cyromazine	167.1	85.1	62	27	15	167.1	125.1	62	27	8
Desmedipham+NH4	318.1	182	42	19	10	318.1	136	39	34	9
Diclobutrazol	328.1	70	81	49	12	328.1	158.9	81	49	10
Dicrotophos	238.1	112.1	66	19	8	238.1	193	66	15	13
Diethofencarb	268.1	226.1	60	15	12	268.1	124	61	45	8
Difenoconazole	406.1	251.1	80	37	13	408.2	253.1	76	33	5
Diflubenzuron	311	158.2	71	23	10	311	141.1	71	45	10
Dimethoate	230	199	49	16	12	230	125	50	27	8
DimethomorphA	388.1	301	66	25	5	388.1	165.1	66	45	5
DimethomorphB	388.2	301.1	66	25	5	388.2	165.2	66	45	5
Dimoxystrobin	327.1	205	40	15	5	327.1	116	40	35	5
Dinotefuran	203.1	129.2	51	19	8	203.1	157.2	51	13	14
Dioxacarb	224.1	167	51	13	10	224.1	123	51	23	21
Diuron	233.1	72	56	33	5	235.1	72.1	56	38	10
Doramectin+NH4	916.9	593.6	68	20	16	916.9	331.5	65	33	22
Eamectin	886.5	158.1	111	51	10	886.5	82.1	111	127	13
Eprinomectin	914.5	186.2	77	27	12	914.5	154.2	77	58	10
Ethaboxam	321	183.1	86	33	12	321	200.1	86	39	12
Ethiofencarb	226.1	106.9	41	21	5	226.1	164.1	41	11	5
Ethiprole	397.3	350.9	81	29	24	397.3	255.2	81	49	16
Ethirimol	210.2	140.1	81	31	8	210.2	98.1	81	39	18
Etoxazole	360.1	141	76	45	5	360.1	57.2	76	45	5
Famoxadone+NH4	392	331	32	15	6	392	238	37	23	6
Fenamidone	312.1	92	66	39	16	312.1	236.1	66	21	14
Fenazaquin	307.1	161.1	68	27	10	307.1	147	68	28	9
Fenbuconazole	337	124.9	81	41	8	337	70	81	39	12
Fenhexamid	302	97	75	34	14	302	55	75	67	9
Fenobucarb	208.1	95.1	61	21	18	208.1	152.1	61	13	10
Fenoxycarb	302.1	88	65	30	6	302.1	116.1	65	17	7
Fenpyroximate	422	366.1	56	23	5	422	135.1	56	43	5
Fenuron	165.1	72.1	56	25	5	165.1	46	56	29	5
Flonicamid	230.1	203.1	55	35	4	230.1	174	55	35	4
Flubendiamide	683	408	56	17	12	683	274	56	43	16
Fludioxinil+NH4	266	229	41	23	14	266	227.1	41	13	14
Flufenoxuron	489	158	86	29	10	489	141.1	86	63	8
Fluometuron	233.1	72.1	71	37	12	233.1	46	71	35	4
Fluoxastrobin	459.2	427.2	55	28	5	459.2	188	55	35	5
Flusilazole	316.1	247.1	78	27	14	316.1	165.1	78	38	9
Flutolanil	324.1	262.1	74	26	14	324.1	242.1	74	34	12
Flutolanil+NH4	341.1	242.1	61	35	4	341.1	262.1	61	35	4
Flutriafol	302.1	70.1	66	37	12	302.1	123	66	41	8
Forchlorfenuron	248	129.1	52	25	5	248	93.1	52	48	5
Formetanate	222.1	165	71	22	9	222.1	93	76	53	14
Fuberidazole	185	157	81	33	13	185	65	81	67	11
Furathiocarb	383.1	195.1	74	26	10	383.1	252.1	74	19	14
Halofenozide	331.1	275	41	11	16	331.1	105.1	41	25	8

Hexaflumuron	461.1	158.2	56	25	5	461.1	141.1	56	65	5
Hexythiazox	353.1	228	63	23	12	353.1	168	63	36	9
Hydramethylnon	495.2	323.2	146	45	18	495.2	151.1	146	95	8
Imazalil	297	159	65	34	12	297	201	65	29	10
Imidacloprid	256	209.1	61	23	10	256	175.1	61	28	10
Indoxacarb	528	203	89	54	10	528	218	86	33	14
Ipconazole	334.2	70	74	52	10	334.2	125	74	50	17
Iprovalicarb	321.2	119	66	29	8	321.2	203.1	66	13	13
Isoprocarb	194.1	95	60	23	13	194.1	137	60	13	10
Isoproturon	207.2	72.1	66	29	5	207.2	46.1	66	31	5
Isoxaflutole	360.1	251.1	62	24	9	360.1	220.1	62	50	9
Isoxaflutole+NH4	377	251.1	56	29	14	377	69	56	35	12
Ivermectin+NH4	892.8	569.7	70	21	20	892.8	713.8	71	15	24
Kresoxim-methyl	314	116	51	21	4	314	206	51	13	4
Linuron	249.1	160	60	23	5	249.1	182.1	60	21	5
Lufenuron	511.1	158.1	61	27	5	511.1	141.2	61	67	5
Malathion	331	127	71	19	8	331	285	71	11	16
Mandipropamide	412.1	328.1	81	21	10	412.1	356.1	81	17	10
Mepanipyrim	224	106	86	37	8	224	77	86	59	14
Metaflumizone	507.1	178.1	101	39	12	507.1	287.1	101	37	16
Metalaxyl	280.1	220.2	60	20	12	280.1	192.2	60	26	10
Metconazole	320.1	70	81	51	12	320.1	125	81	59	10
Methamidophos	142	94	54	20	5	142	125	54	19	7
Methiocarb	226.1	169.1	61	13	11	226.1	121.1	61	27	8
Methomyl	163.1	88.1	35	12	6	163.1	106	35	13	6
Methoxyfenozide	369.1	149.1	56	24	9	369.1	313.2	56	13	10
Metobromuron	259	170.2	56	23	4	259	148.2	56	21	4
Mevinphos-E	225.1	127.1	51	20	7	225.1	193.2	51	10	10
Mevinphos-Z	225	127	51	20	7	225	193.1	51	10	10
Mexacarbate	223.2	166.1	64	23	10	223.2	151	64	35	9
Monocrotophos	224.1	127.1	53	23	10	224.1	98	53	17	5
Monolinuron	215.1	126.1	51	23	5	215.1	99	51	41	5
Moxidectin	640.5	528.5	61	12	16	640.5	498.5	61	17	16
Myclobutanil	289	70	71	37	12	289	125	71	47	8
Novaluron	493	158.1	71	27	5	493	141.1	71	69	5
Nuarimol	315	252.1	75	31	13	315	81	75	44	12
Omethoate	214	124.9	46	29	5	214	182.8	46	17	5
Oxadixyl	279.1	219.1	61	17	13	279.1	132.1	61	43	21
Oxamyl+NH ₄	237.1	72.1	36	25	5	237.1	90.1	36	12	6
Paclobutrazol	294	70	62	46	10	294	125	58	49	8
Pencycuron	329.1	125	76	37	22	329.1	218.1	76	25	14
Phenmedipham	301.1	136	50	26	5	301.1	168.1	50	14	4
Phorate Sulfone	293.1	97.1	36	41	5	293.1	171.1	36	17	5
Picoxystrobin	368	145	56	27	4	368	205	56	15	4
PiperonylButox+NH ₄	356.2	177.2	49	22	9	356.2	119.1	49	46	8
Pirimicarb	239.2	72.1	64	35	10	239.2	182.1	64	23	10

Prochloraz	376	308	45	17	10	376	70	45	44	12
Promecarb	208.1	109	37	23	8	208.1	151	37	13	10
Propamocarb	189.2	102	60	25	8	189.2	144	61	19	13
Propargite+NH ₄	368.2	231.1	46	15	13	368.2	175.1	46	23	12
Propiconazole	342.1	159	62	40	9	342.1	69	62	36	10
Propoxur	210.1	111	39	19	6	210.1	168.1	39	11	10
Pymetrozine	218	105	71	27	5	218	78	71	47	5
Pyracarbolid	218.1	125	59	27	8	218.1	97	59	40	14
Pyraclostrobin	388	194	31	19	5	388	163	31	29	5
Pyridaben	365	147	46	31	5	365	309	46	19	5
Pyrimethanil	200	107	71	33	5	200	82	71	35	5
Pyriproxyfen	322	96	45	21	5	322	185	45	29	5
Rotenone	395.1	213.1	90	32	12	395.1	192.1	90	34	10
Siduron	233.3	137.2	66	21	5	233.3	94	66	31	5
Spinetoram A	748.5	142.2	86	45	8	748.5	98.1	86	109	18
Spinetoram B	760.5	142.2	96	41	10	760.5	98.1	96	101	18
SpinosynA	732.5	142.2	111	43	10	732.5	98.1	111	103	16
Spirodiclofen	411.3	313.3	72	23	8	411.3	71.3	71	33	10
Spiromesifen	371.2	273.2	73	16	6	371.2	255.2	74	33	4
Spiromesifen+NH ₄	388.2	273.2	41	19	12	388.2	255.2	41	39	16
Spirotetramat	374.2	330.2	66	23	8	374.2	302.2	66	25	20
Spiroxamine	298.2	144.2	72	28	10	298.2	100.1	72	46	14
Sulfentrazone	387	307.1	81	27	5	387	146	81	57	5
Tebuconazole	308.2	70	81	49	11	308.2	125	81	51	8
Tebufenozide	353.2	133	54	24	9	353.2	297.2	54	14	9
Tebuthiuron	229.1	172.4	46	21	5	229.1	116.1	46	35	5
Teflubenzuron	381.1	141.2	66	52	5	381.1	158.2	66	23	5
Temephos	467	419.1	101	29	12	467	405	101	23	12
Thiabendazole	202.1	175.1	84	35	10	202.1	131.2	84	45	8
Thiacloprid	253	126	68	30	9	253	99	68	60	14
Thiamethoxam	292	211	64	18	10	292	181	64	32	10
Thidiazuron	221.1	102.1	57	28	6	221.1	128.2	57	22	7
Thiophanate-methyl	343	151.1	61	29	14	343	311	61	17	10
Triadimefon	294	197.1	63	22	12	294	225	63	19	8
Triadimenol	296.1	70	46	31	12	296.1	227.1	46	19	14
Trichlorfon	256.9	109.1	66	25	20	256.9	127	66	25	8
Tricyclazole	190	163	81	33	10	190	136	81	41	11
Trifloxystrobin	409	186	31	23	5	409	206	31	21	5
Triflumizole	346.1	278.1	51	15	8	346.1	73	51	27	6

Triflumuron	359.1	156.2	52	23	6	359.1	139	52	44	6
Triticonazole	318.1	70	63	42	10	318.1	125	63	41	8
Vamidothion	288	146	61	19	10	288	118	61	33	10
Zoxamide	336.1	187	55	33	11	336.1	159	53	39	12

Adapted from: Sack, Chris*, Smoker, Michael, KAN, Lenexa, Chamkasem, Narong, SRL, Thompson, Richard, Satterfield, Greg, ARL, MacMahon, Shaun, Masse, Claude NERL, Mercer, Greg, Neuhaus, Barbara, PRL-NW, Cassias, Irene, Chang, Eugene, Lin, Yi, PRL-SW, Wong, Jon, Zhang, Kai, CFSAN, *Development and Validation of a Multiresidue Determination for Pesticides by LC-MS/MS* DFS/ORA/FDA No. 4464 Pesticides and *Collaboration of the QuEChERS Procedure for the Multiresidue Determination of Pesticides by LC-MS/MS In Raw Agricultural Commodities*, DFS/ORA/FDA, No. 4465 Pesticides



Modified QuEChERS Procedure for Analysis of Bisphenol A in Canned Food Products

UCT Product Number:

ECQUEU750CT-MP (pouch contains 4000 mg MgSO₄, 1000 mg NaCl, 500 mg Na citrate dibasic sesquihydrate, 1000 mg Na citrate tribasic)

ECPSAC1856 (500 mg of PSA and 500 mg endcapped C18, 6 mL cartridge)

ECPSACB6 (400 mg PSA and 200 mg GCB, 6 mL cartridge)

SMTBSTFA-1-1 (MTBSTFA w/1% TBDMCS)

ECSS10K (Sodium sulfate, anhydrous, ACS Grade, Granular 60 Mesh)

July 2011

Procedure

(Developed for tuna, baby food, pineapple and tea)

1. Sample Preparation

- Weigh 10 g homogenized sample into a 50 mL QuEChERS tube
- Spike with 25 ng/g internal standard: Bisphenol A d16 (BPA d16)
- Add 10 mL MeCN
- Add the contents of **ECQUEU750CT-MP**
- Shake vigorously for 2 min
- Centrifuge at 3500 rpm for 3 min

2. Clean-up (products without pigments)

- Attach **ECPSAC1856** to the vacuum manifold
- Add 3 g muffled Na₂SO₄ to the cartridge
- Rinse with 2 x 2 mL MeCN
- Insert test tubes into the manifold
- Load 5 mL of the supernatant to the cartridge
- Turn on vacuum to collect extract dropwise
- Transfer 1 mL of cleaned extract into a 5 mL test tube
- Concentrate to dryness by N₂ at 35 °C

3. Clean-up (products with pigments)

- Attach **ECPSACB6** to the vacuum manifold
- Add 3 g muffled Na₂SO₄ to the cartridge

- c) Rinse with 2 x 2 mL MeCN
- d) Insert test tubes into the manifold
- e) Load 5 mL of the supernatant to the cartridge
- f) Turn on vacuum to collect extract dropwise
- g) Transfer 1 mL of cleaned extract into a 5 mL test tube
- h) Concentrate to dryness by N₂ at 35 °C

4. **Derivatization**

- a) Add 50 µL pyridine to the dried test tube and vortex
- b) Add 50 µL MTBSTFA/1%TBDMCS and cap the test tube
- c) Heat at 75° C for 30 min
- d) Cool then concentrate to dryness by N₂ at 35 °C
- e) Reconstitute with 75 µL toluene and 25 µL surrogate: 2 ppm triphenyl phosphate (TPP)
- f) Vortex
- g) Transfer to auto-sampler vial with 100 µL insert
- h) Inject 1 µL to GC/MS

5. **GC/MS Analysis** **GC/MS**

- Agilent 6890N GC coupled with 5975C MSD, equipped with 7683 auto sampler. Chemstation software for data acquisition and analysis. Equivalent instrumentation may be used

Injector

- 1 µL splitless injection at 250 °C, split vent of 30 mL/min at 1 min

Liner

- 4 mm splitless gooseneck, 4mmID x 6.5mmOD x 78.5mm (UCT #: **GCLGN4MM**)

Glass wool for liner

- Restek® Deactivated Wool

GC capillary column

- Restek® Rxi-5sil MS 30m x 0.25mm x 0.25µm

Oven temperature program

- Initial oven temperature of 100 °C, hold for 1 min
 - Ramp at 20 °C/min. to 300 °C, hold for 1 min
 - Ramp at 40 °C/min. to 320 °C, hold for 2.5 min
- Total run time 15 min. Data acquisition begins at 9 min

Carrier Gas

- He constant flow 1.2 mL/min

MSD Conditions

- Aux temperature: 280 °C, MS Source: 230 °C, MS Quad: 150 °C

Simultaneous Scan/SIM:

Scan range: 50-500

SIM:

Group 1: 9.0 min.: 326.1, 325.1 (Triphenyl phosphate)

Group 2: 10.5 min.: 441.3, 456.3, 442.3 (derivatized: Bisphenol A-2TBDMS)
452.4, 470.4, 453.4 (derivatized: Bisphenol A d16-2TBDMS)

Dwell time: 100 ms for all ions

Matrix matched curves are generated with the adjustment of the concentrations in the blank and baby food sample

Experimental Data


Calculations

C_{blank} and C_{sample} were calculated by equations 1 and 2:


Equation 1: $A_{\text{sample}}/A_{\text{blank}} = (C_{\text{sample}} + C_{\text{blank}}) * R_1\% / C_{\text{blank}} * R_1\%$

Equation 2: $A_{\text{spiked sample}}/A_{\text{sample}} = (C_{\text{sample}} + C_{\text{blank}} + C_{\text{spike}}) * R_2\% / (C_{\text{sample}} + C_{\text{blank}}) * R_2\%$

Tuna: Tuna in pure olive oil (3oz can)



BPA: Linear dynamic range: 1-500 ng/mL (conc. per 1 mL extract) $R^2=0.9933$




BPA d16: Linear dynamic range: 1-500 ng/mL (conc. per 1 mL extract) $R^2=0.9983$


Bisphenol A in Canned Tuna

	BPA in blank	BPA in tuna	Tuna fortified with 25 ng/g BPA	Tuna fortified with 50 ng/g BPA
Analyte	Conc. (ng/g) (ng/mL)	Conc. (ng/g)	Recovery% \pm RSD% (n=3)	Recovery% \pm RSD% (n=3)
BPA d16	0	0	74.0 \pm 6.4	76.5 \pm 5.9
BPA	0.56	6.64	86.6 \pm 7.5	100 \pm 8.0


Chromatogram of Canned Tuna




Chromatogram of canned tuna sample fortified with 25 ng/g BPA d16 and BPA



Baby Food: Peaches, Stage 2



BPA: Linear dynamic range: 1-500 ng/mL (conc. per 1 mL extract) $R^2=0.9998$




BPA d16: Linear dynamic range: 1-500 ng/mL (conc. per 1 mL extract) $R^2=0.9999$


Bisphenol A in Baby Food

	BPA in blank (ng/g)	BPA in baby food (ng/g)	Baby food fortified with 3 ng/g BPA	Baby food fortified with 10 ng/g BPA
Analyte	Conc. (ng/mL)	Conc. (ng/g)	Recovery% \pm RSD% (n=3)	Recovery% \pm RSD% (n=3)
BPA d16	0	0	98.0 \pm 4.7	98.8 \pm 7.7
BPA	0.33	< 1 ng/g (0.33 ng/g)	99.2 \pm 2.0	95.3 \pm 11


Chromatogram of Baby Food




Chromatogram of baby food sample fortified with 3 ng/g BPA d16 and BPA



Pineapple: crushed pineapple in 100% pineapple juice (canned), pH=3 .



BPA: Linear dynamic range: 1-500 ng/mL (conc. per 1 mL extract) $R^2=0.9989$




BPA d16: Linear dynamic range: 1-500 ng/mL (conc. per 1 mL extract) $R^2=0.9999$


Bisphenol A in Canned Pineapple

	BPA in blank	BPA in pineapple	Pineapple fortified with 5 ng/g BPA	Pineapple fortified with 25 ng/g BPA
Analyte	Conc. (ng/g) (ng/mL)	Conc. (ng/g)	Recovery% \pm RSD% (n=3)	Recovery% \pm RSD% (n=3)
BPA d16	0	0	112 \pm 2.3	93.4 \pm 6.1
BPA	0.33	1.65	112 \pm 5.7	96.1 \pm 5.7


Chromatogram of Canned Pineapple




Chromatogram of canned pineapple sample fortified with 5 ng/g BPA d16 and BPA



Tea Sample: Green Tea with Ginseng and Honey (canned), pH=3.5



BPA: Linear dynamic range: 1-500 ng/mL (conc. per 1 mL extract) $R^2=0.9998$




BPA d16: Linear dynamic range: 1-500 ng/mL (conc. per 1 mL extract) $R^2=0.9993$


Bisphenol A in Canned Tea

	BPA in blank (ng/g)	BPA in tea (ng/g) (ng/g)sample	Tea fortified with 3 ng/mL BPA	Tea fortified with 10 ng/mL BPA
Analyte	Conc. (ng/mL)	Conc. (ng/mL)	Recovery% \pm RSD% (n=3)	Recovery% \pm RSD% (n=3)
BPA d16	0	0	120 \pm 5.9	107 \pm 2.7
BPA	0.46	2.28	104 \pm 8.2	90.0 \pm 5.8

Chromatogram of Canned Tea Sample



Chromatogram of canned tea sample fortified with 3 ng/mL BPA d16 and BPA



DCN-119270-225



QuEChERS Sample Preparation For The Analysis Of Pesticide Residues In Olives*

UCT Product Number:

ECMSSC50CT-MP (4 g MgSO₄, 1.0 g NaCl)

ECQUEU122CT (2 mL centrifuge tube, 150 mg MgSO₄, 50 mg PSA, 50 mg C18 and 50 mg GCB)

CUMPSC1875CB2CT – For better recovery of planar pesticides
(2 mL centrifuge tube, 150 mg MgSO₄, 50 mg PSA, 50 mg C18, 7.5 mg GCB)

April 1, 2011

Summary

This application is a summary of the original paper “Evaluation of the QuEChERS sample preparation approach for the analysis for pesticide residues in olives”*. It describes the use of QuEChERS for the extraction and clean-up of 16 pesticide residues contained in olives. LC-MS/MS with positive ESI was used for pesticides that are difficult to detect by GC-MS. Matrix matched calibration standards were used to compensate for matrix effects. The method achieves acceptable quantitative recoveries of 70–109% with RSDs <20% for DSI-GC-MS and 88–130% with RSDs <10% for LC-MS/MS, and LOQ at or below the regulatory maximum residue limits. Analyte protectants were used with DSI to improve analyte peak shapes and intensities.

Analytes Covered in this Method

Analyte	CASRN
Ometholate	1113-02-6
Dimethoate	60-51-5
Simazine	122-34-9
Diazinon	65863-03-8
p,p'-DDE	82413-20-5
Diuron	56449-18-4
Carbaryl	63-25-2
Malathion	121-75-5
Fenthion	55-38-9
Methidathion	950-37-8
Napropamide	15299-99-7
Oxyfluorfen	42874-03-3
Carfentrazone-ethyl	128639-02-1
Phosmet	732-11-6

Pyriproxyfen	95737-68-1
Deltamethrin	64121-95-5

1. Sample Extraction

- Weigh 10 g of homogenized sample into a 50 mL centrifuge tube
- Add 10 mL of acetonitrile (MeCN)
- Add contents of **ECMSSC50CT-MP**
- Shake vigorously by hand for 1 minute
- Centrifuge @ 3450 rcf for 1 minute

2. Dispersive Clean-Up

- Transfer 1 mL the supernatant to a micro-centrifuge tube **ECQUEU122CT**
or CUMPSC1875CB2CT
- Mix for 20 s
- Centrifuge @ 3450 rcf for 1 minute
- Transfer 400 µL of supernatant to an autosampler vial
- Add 25 µL of TPP solution (10 g/mL triphenylphosphate on MeCN with 1.6% formic acid)
- Shake for 5 s
- Extract is ready for analysis

3. Automated DSI-GC-MS Analysis

GC-MS was performed using an Agilent (Little Falls, DE, USA) 5890 Series II GC and 5972 MS instrument. Injection was performed using a Combi-PAL autosampler (CTC Analytics, Zwingen, Switzerland) using second generation automated DSI accessory (Linex) in combination with an Optic 3 PTV (Atas-GL International BV, Veldhoven, NL)

Note: Equivalent instrumentation and analytical columns can be used

Analyte Protectant Solution

(95% or better purity, prepare at 10:1:1 mg/mL in 7:3 water/MeCN, Sigma or Fluka)

- 3-ethoxy-1,2-propanediol
- D-sorbitol
- L-gulonic acid
- γ -lactone

A quality check standard solution of 16 µg/mL triphenylphosphate (TPP) is prepared in MeCN containing 1.6% formic acid (FA)

For analysis by DSI-GC-MS, 20 µL of the analyte protectant solution was added to all the final extracts and matrix-matched calibration standards by transfer of 400 µL of extract into an autosampler vial and adding 25 µL of TPP solution

Conditions:

- Injection volume 10 µL
- 100° C (held 3.5 min with 50:1 split ratio)
- Ramp at 5° C/s to 280°C (use splitless for 3.5 min, then 50:1 split until 9 min, then change split flow to 20:1 and cool injector temperature to 150° C)

GC Separation:

- Varian VF-5 EZ-guard column (30 m x 0.25 mm id x 0.25 µm film thickness) with an integrated retention gap (5 m x 0.25 mm) at the inlet and an additional 1 m of uncoated capillary at the MS entrance
- He carrier gas @ 1 mL/ min

Oven temperature program:

- Start at 3.5 min (after sample introduction)
- 80° C hold for 3.5 min
- Ramp to 230° C at 108° C/min
- Then ramp to 300° C at 45° C/min, hold for 10 min.
- MS transfer line temperature at 290° C
- Electron ionization (EI) at –70 eV in SIM and full-scan (50–600 m/z) modes in different experiments

Agilent Chemstation for data acquisition/processing and GC-MS control, and Cycle Composer and Atas Evolution software are used to control the automated DSI process and PTV. The pesticide analytes in GC-MS and SIM ions are shown in the table below.

GC-MS SIM Conditions for the Monitored Pesticides

Pesticide	Start time (min)	t _R (min)	m/z (% relative abundance)	
			Quantitation ion	Qualifier ions
Dimethoate	4.5	15.89	87 (100)	125 (45), 93 (54), 58 (19)
Simazine		16.00	201 (78)	173 (41), 186 (51), 158 (25)
Diazinon	16.09	16.18	179 (100)	137 (98), 304 (47), 152 (70)
Diuron		16.52	72 (100)	232 (38), 234 (26), 187 (11)
Carbaryl	17.49	17.70	144 (100)	115 (33), 116 (26), 145 (15)
Malathion		18.03	173 (94)	125 (100), 93 (93), 127 (75)
Fenthion	18.1	18.27	278 (100)	125 (37), 109 (33), 79 (19)
Methidathion	19.05	19.29	145 (88)	93 (40), 125 (27), 302 (19)
Napropamide	13.39	19.58	271 (26)	72 (100), 128 (63)
p,p'-DDE		19.367	318 (64)	246 (100), 248 (64), 316 (56)
Oxyfluorfen		19.71	361 (38)	252 (100), 300 (35), 280 (14)
Carfentrazone-ethyl	20	20.28	312 (100)	330 (65), 340 (63), 376 (31)
TPP	20.38	20.96	326 (100)	325 (87), 77 (88), 215 (20)
Phosmet		21.17	160 (100)	133 (15), 104 (15), 193 (4)
Pyriproxyfen	21.30	21.50	136 (100)	226 (12), 185 (6)
Deltamethrin	22.8	23.59	253 (85)	181 (100), 251 (44), 152 (20)

4. LC-MS/MS Analysis

Suggested Instrumentation: Agilent 1100 HPLC (consisting of vacuum degasser, autosampler Model WPALS, and a binary pump) equipped with a Prodigy ODS-3 (150 mm x 3 mm) and 5 µ particle size analytical column coupled to a ODS-C18 (4 mm x 2 mm and 5 µ particle size) guard column from Phenomenex (Torrance, CA, USA).

- Column temperature: 30° C
- Injection volume 5 µl.
- Mobile phase A water, B MeCN, both with 0.1% FA
- Gradient program:
 - Flow rate 0.3 mL/min
 - 25% solvent B linear gradient to 100% over the first 5 min
 - Hold for 7 min until 12 min
 - 11-min post run column wash

The LC system is connected to an API 3000 triple-quadrupole mass spectrometer (Applied Biosystems, Toronto, Canada) operated in ESI positive mode. Optimizations of the mass analyzer parameters were done by infusion of 1 µg/mL analyte solutions at 10 µL/min with a

syringe pump (Harvard Apparatus, Holliston, MA, USA) using the autotune function.

Note: Equivalent instrumentation and analytical columns can be used

Final MS/MS conditions include:

- N₂ pressure 55 psi
- nebulizer gas setting 14
- curtain gas setting 11
- collision gas setting 12
- 4200 V ionspray voltage
- ESI temperature 525° C
- focusing potential 100 V
- entrance potential 10 V
- 0.15 s dwell time

The pesticide analytes by LC-MS/MS are shown in the table below with respective analytical ions

LC-MS/MS Conditions for the Monitored Pesticides

(Quantitation ion is shown as first mass)

Pesticide	Start time (min)	t _R (min)	Precursor ion (m/z)	Product ions (m/z)
Omethoate	2.5	2.68	214.0	183.2, 125.2
Dimethoate	5	6.83	230.0	199.1, 125.1
Simazine	7.6	7.98	202.0	124.2, 132.2
Carbaryl		8.48	202.2	145.1, 127.1
Diuron		8.67	233.1	72.2, 160.1
Phosmet	9	9.27	318.0	160.2, 133.2
Methidathion		9.28	303.0	145.1, 85.1
Malathion		9.64	331.0	127.2, 285.2
TPP	9.8	10.18	327.0	77.2, 152.0

*Adapted and used by permission from Cunha, Sara C., Lehotay, Steven J., Mastovska, Katerina, Fernandes, Jos O., Beatriz, Maria, Oliveira, P. P., Sep. Sci. 2007, 30, 620 – 632, DOI 10.1002/jssc.200600410

Listing of instrument manufacturers and standards suppliers does not constitute endorsement by UCT.
Equivalent systems may be used



Determination of Carbendazim in Orange Juice Using QuEChERS with LC/MS/MS Detection

UCT Part Numbers:

ECQUEU750CT-MP: Pouch contains: 4000 mg MgSO_4 , 1000 mg NaCl, 500 mg Na citrate dibasic sesquihydrate and 1000 mg Na citrate tribasic dihydrate

CUMPSC18CT: 2 mL centrifuge tube contains: 150 mg MgSO_4 , 50 mg PSA, 50 mg endcapped C18

February 1, 2012

Introduction

The planar fungicide carbendazim (CASRN 10605-21-7) can be used to control mold on citrus crops but is not approved for use in the US or on imported products. Concentrations in citrus products can be rapidly and accurately determined using a QuEChERS extraction with dSPE clean-up. LOD and LOQ for this method are 0.4 and 1.4 ng/mL, respectively.

Procedure

1. Extraction

- a) Add 10 mL of orange juice to a 50 mL centrifuge tube
- b) Add 10 mL acetonitrile then vortex
- c) Add the contents of **ECQUEU750CT-MP**
- d) Shake vigorously for 1 min
- e) Centrifuge at 5000 rpm for 5 min at 20° C
- f) Supernatant is ready for clean-up

2. Dispersive Clean-up

- a) Add 1 mL of supernatant to **CUMPSC18CT** tube
- b) Shake sample(s) for 1 min
- c) Centrifuge at 10,000 rpm for 5 min
- d) Transfer 0.5 mL to 2 mL autosampler vial
- e) Add 25 μL 1 ppm TPP, vortex
- f) Samples are ready for analysis

3. LC/MS//MS Analysis

Instrumentation: Thermo Accela HPLC with autosampler

LC Conditions

Column	Guard column: Restek C18, 2.1 x 20 mm Column: Sepax HP-C18, 2.1 x 100 mm, 3 μ m, 120 Å
Column Temperature	Ambient
Injection Volume	10 μ L at 15° C
Mobile Phase	A: 0.1% formic acid in water B: 0.1% formic acid in methanol
Flow Rate	200 μ L/min

LC Gradient Program

Time	%A	%B
0	50	50
3	0	100
8	0	100
9	50	50
14	50	50

MS Conditions

Thermo TSQ Vantage MS

Ion source:	Heated ESI
Ion polarity:	ESI +
Spray voltage:	3000 V


Sheath gas pressure:	N ₂ @ 40 psi
Auxiliary gas pressure:	N ₂ @ 10 psi
Ion transfer capillary temperature:	350 °C
Scan type:	SRM (0-10 min)
CID conditions:	Ar @ 1.5 mTorr

SRM transitions

Compound	Parent	Product ion 1	CE	Product ion 2	CE	S-Lens	Dwell time (s)
Carbendazim	192.093	132.080	29	160.080	17	81	0.20
TPP (IS)	327.093	77.020	37	152.070	33	98	0.10

Matrix Matched Calibration Curve

Dynamic linearity range is from 2 to 500 ng/mL with $R^2=0.9992$




Accuracy and Precision of Carbendazim Data from Spiked Orange Juice Sample

Fortified (ng/mL)	Recovery%	RSD % n=4
10	96.6	4.5
50	100.2	3.4
250	103.7	2.1


Chromatogram of Orange Juice Sample and Fortified at 50 ng/mL

RT: 1.96 - 3.96 SM: 7G



50 ng/mL Fortified

RT: 1.96 - 3.96 SM: 7G



DCN-212020-230



Determination of Pesticide Residues in Marijuana and Tea by QuEChERS and LC/MS/MS

UCT Part Numbers:

ECPAHR50CT (50 mL centrifuge tubes)

CUMPSGG2CT (2 mL dSPE tube with 150 mg MgSO₄, 50 mg PSA and 50 mg ChloroFiltr[®])

ECQUUS2-MP (pouch containing 4 g MgSO₄ and 2 g NaCl)

April 2013

Summary

Various pesticides can be rapidly determined in dried tea leaves or marijuana using this simple method. Samples are hydrated in water and then extracted by QuEChERS and dSPE clean-up using ChloroFiltr[®]. Results are determined by LC/MS/MS.

Internal standard

Prepare a 10 ppm TPP solution by mixing 20 µL of the 5000 ppm TPP solution with 10 mL of MeCN

Transfer all standards to amber glass vials and store at -20°C until needed

Procedure

1. QuEChERS Extraction

- a) Weigh 2 g of the homogenized tea or marijuana into a 50 mL centrifuge tube
- b) Add 10 mL of reagent water to each tube and hydrate the samples for 1 hour using a horizontal shaker
- c) Add 100 µL of the 10 ppm TPP solution to all samples
- d) Add 10 mL of acetonitrile (MeCN). Vortex for 1 min
- e) Add contents of **ECQUUS2-MP** pouch then shake vigorously for 1 min
- f) Centrifuge at 5000 rpm for 5 min

2. dSPE clean-up

- a) Transfer 1 mL of the supernatant to 2 mL dSPE tube **CUMPSGG2CT**
- b) Shake for 30 seconds
- c) Centrifuge at 10,000 rpm for 5 min
- d) Transfer 0.3 mL of the cleaned extract into a 2 mL auto-sampler vial
- e) Add 0.3 mL of reagent water
- f) Vortex, then filter using a 0.45 µm syringe filter
- g) The samples are ready for LC/MS/MS analysis

Instrumentation

Thermo Scientific Dionex Ultimate 3000 LC system coupled to a TSQ Vantage[®] triple quadrupole mass spectrometer

HPLC Parameters

HPLC: Thermo Scientific Dionex UltiMate 3000 [®] LC System		
Column: Thermo Scientific, Accucore aQ [®] , 100 x 2.1 mm, 2.6 µm		
Guard Column: Thermo Scientific, Accucore aQ [®] , 10 x 2.1 mm, 2.6 µm		
Column Temperature: 40 °C		
Column Flow Rate: 0.200 mL/min		
Auto-sampler Temperature: 10 °C		
Injection Volume: 10 µL		
Gradient Program:		
Mobile Phase A: 0.3 % formic acid and 0.1 % ammonia formate in water		
Mobile Phase B: 0.1 % formic acid in MeOH		
Time (min)	Mobile Phase A (%)	Mobile Phase B (%)
0	99	1
1.5	99	1
3.5	20	80
10	10	90
12	0	100
15	0	100
15.2	99	1
20	99	1
Divert mobile phase to waste from 0 - 0.5 and 15 - 20 min to prevent ion source contamination		

MS Parameters

Polarity	ESI +
Spray voltage V	4000 V
Vaporizer Temperature	300 °C
Ion transfer capillary	200 °C
Sheath gas pressure	50 arbitrary units
Auxiliary gas pressure	25 arbitrary units
Q1 and Q3 peak width	0.2 and 0.7 Da
Collision gas and pressure	Ar at 1.5 mTorr
Scan type	SRM
Cycle time	1 sec
Acquisition method	EZ Method

SRM Transitions

Name	Rt (min)	Precursor ion	Product ion 1	CE 1	Product ion 2	CE 2	S-lens (V)
Methamidophos	1.24	142.007	124.57	14	111.6	5	60
Carbendazim	6.37	192.093	132.08	29	160.08	17	81
Dicrotophos	6.46	238.009	126.58	17	108.60	33	73
Acetachlor	6.48	269.417	111.86	15	71.69	33	72
Thiabendazole	6.61	202.059	131.06	31	175.07	24	103
DIMP	7.30	181.283	96.60	13	78.62	32	44
Tebuthiuron	7.32	228.946	115.59	26	171.63	17	72
Simazine	7.34	201.400	67.68	33	103.60	24	85
Carbaryl	7.41	201.956	126.63	30	144.63	7	40
Atrazine	7.68	215.957	67.65	35	173.60	16	79
DEET	7.72	191.947	118.63	15	90.66	28	92
Pyrimethanil	8.10	200.116	107.06	23	183.14	22	66
Malathion	8.12	331.011	98.57	23	126.86	12	60
Bifenazate	8.22	300.925	169.82	15	197.62	5	51
Tebuconazole	8.73	308.008	69.66	29	124.56	35	97
Cyprodinil	8.81	226.122	77.03	40	93.05	33	88
TPP (IS)	8.81	327.093	77.02	37	152.07	33	98
Diazinon	8.87	305.135	153.09	15	169.08	14	89

Zoxamide	8.90	335.807	186.50	20	158.51	38	102
Pyrazophos	8.99	374.103	194.06	20	222.13	20	104
Profenofos	9.59	372.300	302.37	19	143.48	35	104
Chlorpyrifos	10.23	349.989	96.89	32	197.94	17	69
Abamectin	11.20	890.486	304.40	18	306.68	15	102
Bifenthrin	12.77	440.039	165.21	39	180.42	11	66

**Accuracy and Precision Data Obtained from the Fortified
Tea Samples**

Analyte	Spiked at 2 ng/mL		Spiked at 10 ng/mL		Spiked at 50 ng/mL	
	Rec%	RSD (n=6)	Rec%	RSD (n=6)	Rec%	RSD (n=6)
Methamidophos	112.5	7.3	100.7	1.5	85.9	10.1
Carbendazim	nd	nd	87.8	14.6	79.5	13.2
Dicrotophos	114.2	3.3	102.3	3.4	93.9	11.0
Acetachlor	108.3	10.4	111.3	4.6	105.5	6.5
Thiabendazole	86.7	6.0	84.0	2.9	73.5	10.8
DIMP	111.7	6.7	109.3	2.7	102.4	7.8
Tebuthiuron	113.3	3.6	108.3	2.4	100.4	8.1
Simazine	110.8	6.0	118.8	2.1	105.6	10.2
Carbaryl	115.8	5.0	122.0	2.1	111.2	10.4
Atrazine	124.2	7.8	117.5	2.6	105.2	9.8
DEET	149.2	12.8	125.5	3.1	106.3	10.1
Pyrimethanil	99.2	5.9	98.0	5.0	90.3	2.6
Malathion	143.3	9.5	125.0	6.7	110.2	4.4
Bifenazate	114.2	12.2	106.8	6.7	98.2	3.8
Tebuconazole	72.5	7.2	80.2	5.4	79.1	4.7
Cyprodinil	90.8	5.4	77.5	2.7	74.1	3.2
Diazinon	108.3	7.5	99.5	1.9	97.7	4.3
Zoxamide	95.0	7.4	92.2	1.3	90.1	3.4
Pyrazophos	90.0	5.0	91.7	2.4	88.3	5.1
Profenofos	96.7	7.8	80.8	4.5	75.0	5.3
Chlorpyrifos	80.0	7.9	82.5	2.7	81.6	5.4
Abamectin	99.2	8.1	89.8	4.0	82.2	9.4
Bifenthrin	89.2	7.5	119.8	7.6	126.5	19.0

nd: < LOQ, not determined

**Crude Tea and Marijuana Extracts Before, and After dSPE Clean-up Using
150 mg MgSO_4 , 50 mg PSA, and 50 mg ChloroFiltr®**



Figure 1a Tea extracts before (left) and after dSPE clean-up (right)




Figure 1b Marijuana (Sample #4) extracts before (left) and after dSPE clean-up (right)

DCN-316140-262



Determination of Pesticides in Strawberries by QuEChERS Extraction, Quick QuEChERS Clean-up, and GC/MS Detection

UCT Part Numbers:

ECQUEU750CT-MP: (4000 mg magnesium sulfate, 1000 mg sodium chloride, 500 mg sodium citrate dibasic sesquihydrate, 1000 mg sodium citrate tribasic dihydrate)

ECPURMPSMC: (Quick QuEChERS push thru cartridge contains: 110 mg MgSO_4 , 180 mg PSA)

July 2012

Procedure

1. Extraction

- a) Add homogenized and hydrated strawberry sample (10 g) to a 50 mL centrifuge tube
- b) Add 10 mL acetonitrile, vortex 30 sec
- c) Add the contents of pouch (**ECQUEU750CT-MP**)
- d) Shake vigorously for 1 min
- e) Centrifuge at >1500 rcf for 1 min at 20° C
- f) Supernatant is ready for clean-up

2. Quick QuEChERS Clean-up

- a) Load 1 mL of supernatant into a disposable syringe
- b) Pass the supernatant slowly through Quick QuEChERS cartridge (**ECPURMPSMC**)
- c) Collect 0.5 mL cleaned extract into a GC vial
- d) Add triphenyl phosphate as internal standard (200 ng/mL)
- e) Samples are ready for GC/MS analysis

Clean-up of Strawberry Extract with Quick QuEChERS



3. GC/MS Detection

Thermo TRACE GC Ultra gas chromatograph coupled with a Thermo ISQ single quadrupole mass spectrometer and TriPlus autosampler

GC/MS Conditions (Using a matrix matched calibration)

Column	Rtx-5MS, 30 m x 0.25 mm x 0.25 μ m
Carrier Gas	Helium
Flow Rate	1.2 mL/min
Ramp	55°C for 1 min, 20°C/min to 300°C, hold for 4 min
Injector Temperature	220°C
Injection Volume	1 μ L in splitless mode
Ion Source Temperature	200°C
Transfer Line Temperature	250°C
MS Operation	SIM and Full Scan

Accuracy and Precision Data

Compound	Fortified at 10 ng/mL		Fortified at 50 ng/mL		Fortified at 100 ng/mL	
	Recovery%	RSD% (n=4)	Recovery%	RSD% (n=4)	Recovery%	RSD% (n=4)
Methamidophos	93.7	3.4	81.6	5.8	84.2	3.5
Carbendazim	105.7	10.8	100.1	10.6	90.5	7.6
Thiabendazole	91.2	4.9	87.9	6.8	85.0	4.0
Pyrimethanil	112.2	2.7	107.0	3.2	102.8	4.9
Cyprodinil	104.3	3.2	99.9	6.1	100.2	4.9
Diazinon	104.9	5.6	102.0	6.6	99.2	6.8
Pyrazophos	99.9	4.0	96.6	5.6	91.3	4.1
Chlorpyrifos	91.7	4.6	99.5	5.2	97.2	3.8



Multi-residue Pesticide Analysis of Botanical Dietary Supplements using SPE Clean-up and GC-Triple Quadrupole MS/MS*

UCT Part Numbers:

ECPSACB256 (500 mg PSA, 250 mg GCB, 6 mL cartridge)

ECMSSC50CT-MP (4000mg MgSO₄, 1000mg NaCl)

May 2013

Summary

A screening method for the analysis of 310 pesticides, isomers of organohalogen, organophosphorus, organonitrogen and pyrethroid pesticide metabolites in a variety of dried botanical dietary supplements, spices, medicinal plants, herbals, teas, and phyto-medicines is described. Acetonitrile/water is added to the dried botanical along with anhydrous MgSO₄ and NaCl for extraction. This is followed by clean-up using a tandem SPE cartridge consisting of graphitized carbon black (GCB) and primary- secondary amine sorbent (PSA). Pesticides in the study were spiked at 10, 25, 100 and 500 µg/kg. Mean pesticide recoveries were 97%, 91%, 90% and 90%. Percent RSDs were 15%, 10%, 8%, and 6% respectively.

Some Pesticides Screened by this Method

Azoxystrobin	Chlorpyrifos	DDT
Diazinon	Dimethomorph	Hexachlorobenzene
Hexachlorocyclohexanes	methamidophos	Pentachloroaniline
Pentachloroanisole	Pentachlorobenzene	Pentachlorothioanisole
Quinoxifen	Quintozone	Tecnazene
Tetraconazole	Tetramethrin	

Prepare stock solutions of individual standards by dissolving 25–100 mg of pesticide in 25 mL of toluene.

Procedure

1. Botanical Preparation

- a) Add dry botanical powder (1.00 ± 0.02 g) to the 50 mL centrifuge tube
- b) Add 10 mL water and 10 mL extraction solvent (60 µg/L of the internal standard, tris-(1,3-dichloroisopropyl)phosphate in acetonitrile)
- c) Shake vigorously to insure the botanical is completely wetted
- d) Allowed to stand for 15 minutes
- e) Add the contents of **ECMSSC50CT-MP** pouch to each centrifuge tube
- f) Shake vigorously after addition to disperse the salts
- g) Shake samples vigorously for 1 minute
 - a. Centrifuge at 4500 rpm (4200g) x 5 min

2. Solid-phase Clean-up

- a) Condition **ECPSACB256** cartridge(s) on a manifold using 3 x 6 mL acetone
- b) Do not let cartridge go to dryness after last acetone wash
- c) Insert 15 mL disposable centrifuge tubes in the vacuum manifold
- d) Add a layer of anhydrous sodium sulfate to the top of each cartridge
- e) Add a 1.25 mL aliquot of the extract to the cartridge
- f) Allow to percolate through the cartridge. Apply low vacuum if needed
- g) Rinse cartridge with 1 mL of acetone and continue to collect

3. Cartridge Elution

- a) Elute cartridge with 12 mL of 3:1 acetone:toluene
- b) Reduce extract to approximately 100 µL with a gentle N₂ stream in a water bath at 50-55 °C
- c) Add 0.5 mL toluene, QC standards (50 µL of deuterated polycyclic aromatic hydrocarbons mixture, 500 µg/L), and 25 mg of magnesium sulfate
- d) Centrifuge at 3500 rpm x 5 min
- e) Divide the toluene extracts between two GC vials with 250 µL vial inserts keeping one vial as a reserve spare

4. GC-MS/MS Analysis

GC-MS/MS Parameters

(Equivalent equipment may be used)

GC: TRACE Ultra Gas Chromatograph
MS: TSQ Quantum triple quadrupole
Autosampler: TriPlus (Thermo Fisher Scientific)
Column: 30 m x 0.25 188 mm id HP-5MS fused silica capillary column (Agilent Technologies, Santa Clara, CA, USA)
Guard Column: deactivated 5 m x 0.25 mm I.D, Restek Corp., Bellefonte, PA
Oven Temperature: Program, initial 105° C for 3 min, 130° C/ @ 10° C/min, 200° C @ 4°C/min, 290° C @ 8° C/min. Hold 6 min.
Column Flow Rate: 1.4 mL/min He
Injector: PTV 100° C for 0.05 min, ramp 12° C/sec to 280° C
Autosampler: TriPlus Thermo Fisher Scientific
Auto-sampler Temperature: 10 °C
Injection Volume: 2.0 µL splitless mode
Injection Liner: 2 mm id x 120 mm open baffled fused silica deactivated
Ion Source & Transfer T: 250°C and 280°C, respectively
Electron Multiplier V: auto-tune approx. 1400 V
Ar Collision gas: 1.5 mTorr
Cycle Time: 0.5 sec
Q1 entrance mass width (FWHM): 0.7 amu.
Stock pesticide standards: Full scan 50-550 m/z

There is not complete agreement over which transitions for a given pesticide are optimal for foods or dietary supplements. Reference information on SRM transitions for these analytes is provided in references.¹⁻⁴

Representative Recoveries (RSD) and Percent LOQ's in Each Botanical Matrix

Representative Recoveries (mean, n = 4) ± percent relative standard deviation (RSD) for pesticides by botanical, at 10 and 500 µg/kg and the number not detected (ND) at each fortification concentration

Botanical		10 µg/kg	ND	500 µg/kg	ND
Astragalus	Astragalus membranaceus	94 ±13	68	92 ±3	15
Bitter Orange Peel	Citrus aurantium	112 ±15	63	90 ±5	23
Black Cohosh Root	Cimicifuga racemosa	84 ±11	39	82 ±4	14
Chamomile	Matricaria chamomilla	87 ±11	68	91 ±4	29
Cinnamon	Cinnamomum verum	63 ±26	149	101 ±7	9
Comfrey Root	Symphytum officinale	89 ±18	69	83 ±10	15
Dong Quai	Angelica sinensis	107 ±19	156	97 ±8	16
Echinacea	Echinacea purpurea	97 ±16	61	101 ±8	11
Fenugreek	Trigonella foenum	99 ±14	82	81 ±7	11
Garlic	Allium sativum	98 ±18	78	87 ±6	15
Ginger	Zingiber	103 ±14	211	104 ±6	59
Ginkgo Biloba	Ginkgo biloba	99 ±16	89	80 ±7	14
Ginseng	Panax quinquefolius	88 ±11	64	86 ±6	8
Green Tea		91 ±13	43	79 ±6	11
Hoodia	Hoodia gordonii	104 ±19	94	93 ±5	20
Hops	Humulus lupulus	111 ±10	233	102 ±6	53
Jasmine	Jasminum odoratissimum	100 ±14	65	84 ±4	10
Kava Kava	Piper methysticum	111 ±10	164	100 ±4	59
Licorice Root	Glycyrrhiza glabra	93 ±14	43	87 ±6	15
Milk Thistle	Silybum marianum	90 ±13	73	77 ±10	17
Psyllium	Plantago psyllium	99 ±11	39	95 ±4	16
Saw Palmetto	Serenoa serrulata	103 ±13	111	98 ±7	13
St. John's Wort	Hypericum perforatum	93 ±10	100	83 ±6	16
Valerian Root	Valeriana wallichii	101 ±19	68	94 ±10	13

* Adapted from, Douglas G. Hayward, Jon W. Wong, Feng Shi, Kai Zhang, Nathaniel S. Lee, Alex L. DiBenedetto, & Mathew J. Hengel. "Multi-residue Pesticide Analysis of Botanical Dietary Supplements using Salt-out Acetonitrile Extraction, Solid-phase extraction clean-up column and Gas Chromatography-Triple Quadrupole Mass Spectrometry" DOI: 0.1021/ac400481w

References

- 1) Pang, G. F.; Fan, C. L.; Liu, Y. M.; Cao, Y. Z.; Zhang, J. J.; Li, X. M.; Li, Z. Y.; Wu, Y. P.; Guo, T. T. J. AOAC Int. 2006, 89, 740–771
- 2) Walorczyk, S.; Gnusowski, B. J. Chromatogr. A. 2006, 1128(1-2), 236-243
- 3) Okihashi, M.; Takatori, S.; Kitagawa, Y.; Tanaka, Y. J. AOAC Int. 2007, 90(4), 1165-1179
- 4) Frenich, A. G.; González-Rodríguez, M. J.; Arrebola, F. J.; Vidal, J. L. M. Anal. Chem. 2005, 77, 4640-4648



ChloroFiltr[®]: A Novel Sorbent for Chlorophyll Removal using QuEChERS

UCT Part Numbers:

ECPAHFR50CT (50 mL polypropylene centrifuge tube)

ECQUUS2-MP (Mylar pouch with 4000 mg MgSO₄ and 2000 mg NaCl)

CUMPSGGC182CT (2 mL centrifuge tube with 150 mg MgSO₄, 50 mg PSA, 50 mg C18 and 50 mg ChloroFiltr[®])

January 2013

Spinach and other highly pigmented vegetables contain chlorophylls, carotenoids, xanthophylls, and anthocyanins. Chlorophylls have the greatest adverse effect on GC systems due to their non-volatile characteristics. This QuEChERS procedure uses ChloroFiltr[®] to significantly reduce chlorophylls without sacrificing the recoveries of planar pesticides.

Procedure

1. QuEChERS Extraction

- a) Homogenize 500 g of spinach in a food processor for 1-2 minutes
- b) Weigh 10 grams of homogenized spinach sample into 50 mL centrifuge tube
- c) Spike with 100 µL of 50 ppm triphenyl phosphate* as internal standard (IS)
- d) Add 10 mL of acetonitrile then shake for 1 min
- e) Add contents of Mylar pouch **ECQUUS2-MP** then shake vigorously for 1 min
- f) Centrifuge at 5,000 rpm for 5 min
- g) Supernatant is ready for clean-up

*50 ppm TPP solution: mix 50 µL of 5000 ppm TPP solution with 4.95 mL of MeCN

2. dSPE Clean-up


- a) Transfer 1 mL of the extract to the 2 mL **CUMPSGGC182CT** (ChloroFiltr) dSPE micro centrifuge tube
- b) Shake for 30 sec
- c) Centrifuge at 3,000 rpm for 5 min
- d) Transfer 0.4 mL of the supernatant to a 2 mL autosampler vial
- e) Sample is ready for LC/MS/MS analysis

MS Parameters	
MS instrument	Thermo TSQ Vantage triple
Polarity	ESI +
Spray Voltage	3000 V
Vaporizer Temperature	350 °C
Ion Transfer Capillary	300 °C
Sheath Gas Pressure	40 arbitrary units
Auxiliary Gas Pressure	10 arbitrary units
Q1 and Q3 Peak Width	0.7 Da
Collision Gas and Pressure	Argon at 1.5 mTorr

HPLC Conditions	
HPLC system	Thermo Accela 1250 LC equipped with PAL auto-sampler
LC Column	Sepax HP-C18, 2.1*100 mm, 3 µm
Guard column	Restek C18, 2.1*20 mm, 3 µm
Column temperature	ambient
Auto-sampler Temp.	15 °C
Injection volume	10 µL
Mobile phase A	0.1% formic acid in Milli-Q water
Mobile phase B	0.1% formic acid in methanol
Flow rate	200 µL/min

Gradient Program		
Time (min)	Mobile phase A	Mobile phase B
0	95	5
1	95	5
3	50	50
8	5	95
14	5	95
14.2	95	5
16	95	5

SRM Transitions							
Compound	Precursor ion	Product ion 1	CE 1	Product ion 2	CE 2	S-Lens	Dwell time
Carbendazim	192.093	132.080	29	160.080	17	81	0.10
Thiabendazole	202.059	131.060	31	175.070	31	103	0.10
Pyrimethanil	200.116	107.060	23	183.140	22	66	0.10
Cyprodinil	226.122	77.030	40	93.050	33	88	0.10
TPP (IS)	327.093	77.020	37	152.070	33	98	0.10
Diazinon	305.135	153.090	15	169.08	14	89	0.10
Pyrazophos	374.103	194.060	20	222.130	20	104	0.10
Chlorpyrifos	349.989	96.890	32	197.940	17	69	0.10



Extract cleaned with ChloroFiltr[®] (Middle) is less green than that cleaned with GCB (Right), indicating that ChloroFiltr[®] is slightly more efficient in Chlorophyll removal.

Comparison of Pesticide Recoveries and RSDs Obtained by dSPE Clean-up of Spinach Sample using ChloroFiltr [®] and GCB (n=4)				
Pesticide	ChloroFiltr [®]		GCB (7.5 mg)	
	Recovery%	RSD%	Recovery%	RSD%
Carbendazim	87.1	1.0	71.2	4.0
Thiabendazole	93.2	1.9	55.9	2.6
Pyrimethanil	97.3	1.2	85.0	1.2
Cyprodinil	91.2	0.5	79.3	3.1
Diazinon	104.5	2.3	100.0	0.6
Pyrazophos	92.0	0.9	92.7	1.6
Chlorpyrifos	95.6	2.5	96.3	2.1



Analysis of 136 pesticides in avocado using a modified QuEChERS method with LC-MS/MS and GC-MS/MS*

UCT Part Numbers:

ECMSSA50CT-MP (6 g of MgSO₄ and 1.5 g anhydrous sodium acetate)

CUMPSC18CT (2 mL dispersive cleanup tubes containing 150 mg of anhydrous MgSO₄, 50 mg of PSA, and 50 mg endcapped C-18)

July 2013

Summary

A simple, high-throughput modified QuEChERS screening method for the analysis of 136 pesticides in highly fat rich avocado is described. The average recoveries for 79 pesticides by LC-MS/MS at 10, 50, and 200 ng/g fortifying levels were 86% or better (with maximum RSD at 9.2%). GC-MS/MS analysis demonstrated 70% recovery or better (RSD < 18%) from 57 pesticides at the same spike levels.

Table of Pesticides Evaluated for this Method

	Name	Class
Fungicides	Pyrachlostrobin	Strobilurin
	Chlorothalonil	OC
	Pyrimethanil	Anilnopyrimidine
	Imazalil	Imidazole
	o-Phenylphenol	Phenol
	Procymidone	Dicarboximide
	Tebuconazole	Triazole
	Thiabendazole	Benzimidazole
	Tolyfluand	N-Trihalomethylthio
	Hexachlorobenzene	OC
Insecticides	Bifenthrin	Pyrethroid
	Aminocarb	Carbamate
	Chlorpyrifos	Pyridine OP
	Chlorpyrifos-methyl	Pyridine OP
	Dichlorvos	OP
	DDT	OC
	DDE	OC
	Endosulfan	OC
	Ethion	OP
	Methamidophos	OP
	Acephate	OP
	Permethrin	Pyrethroid
	Acetamiprid	Neonicotinoid
	Prometryn	Triazine
Herbicides	Linuron	Phenylurea
	Trifluralin	Dinitroaniline

OC=organochlorine OP=organophosphate

Procedure

1. Sample Preparation

- a) Add 3 g of homogenized sample to a 50 mL centrifuge tube
- b) Add fortification and/or internal standards
- c) Add 5 mL of reagent water and 25 mL of 1% acetic acid in acetonitrile (MeCN) to each sample tube
- d) Cap tube and shake for 10 minutes with an SPEX 2000 Geno grinder (or equivalent) @ 1000 stroke/min
- e) Add one **ECMSSA50CT** packet to each sample tube and shake for additional 10 min @ 1000 strokes/min
- f) Centrifuge @ 3000 rpm for 10 min

2. Sample Clean-up for LC

- a) Transfer 1 mL of supernatant to an autosampler vial
- b) Sample is ready for LC-MS/MS analysis (*if sample clean-up is desired, see Sample Clean-up for GC below*)

3. Sample Clean-up for GC

- a) Pipette 1 mL of supernatant into **CUMPSC18CT** tube
- b) Vortex for 1 min
- c) Centrifuge @ 2000 rpm for 10 min
- d) Sample is ready for GC analysis

Note: Extract a clean matrix and clean-up with the steps above. This extract must be used to prepare matrix-matched calibration standards. Matrix-matching is necessary for this procedure.

LC-MS/MS Parameters
(Equivalent instrumentation may be used)

HPLC Conditions	
LC:	Shimadzu with two LC 20AD pumps
MS:	4000 Q-TRAP mass spectrometer AB Sciex
Autosampler:	Sil-20AC autosampler
Column:	Ultra Aqueous C18 column (3 µm, 100 x 2.1 mm) Restek
Guard Column:	(10 x 2.1 mm) Restek
Column Oven:	CTO-20AC column oven (Shimadzu)
Separation Temp:	50 °C
Software:	Analyst software version 1.4
Mobile Phase:	A 4 mM ammonium formate and 0.1% formic acid in water, B 4 mM ammonium formate and 0.1% formic acid in methanol
Mobile Phase Program:	Gradient start at 5% B (0.0 - 0.4 min); flow rate of 0.5 mL/min. 60% B at 5 min, then 95% B at 12.5 min, hold until 14.5 min, and concluded by column equilibration at initial condition for 3 min. Total run time 18 min.
Injection Volume:	1.0 µL
MS/MS Conditions	
Electrospray:	positive ion
Ion Transition:	60 sec each analyte
Curtain gas (CUR):	30 psi
Ion Spray V:	4500 volts
Nebulizer Gas (GSI):	60 psi
Heater Gas (GS2)	60 psi
Source Temp (TEM):	350 °C

GC/MS Parameters
(Equivalent instrumentation may be used)

GC Conditions	
GC:	Agilent 7890A GC,
MS:	7000 triple-quadrupole MS, MassHunter software (version B.05.00412)
Autosampler:	7693 autosampler
Column:	two HP-5ms Ultra Inert capillary columns from Agilent (0.25 mm ID x 15m, 0.25 µm film thickness) connected at backflush union
Column Head Pressure:	12.772 psi
Oven Temperature:	initial 60° C for 1 min, 40°/min to 170° C , then 10°/min to 310° C. Hold 1.2 min. Total run time 19 min.
Column Flow Rate:	1.335 mL/min He
Injector:	60° C for 0.2 min, ramp to 280° C @ 600° C/min
Autosampler:	TriPlus Thermo Fisher Scientific
Back Flush:	column 1 for 2 minutes at 310°
Injection Volume:	1.0 µL splitless mode
MS Parameters	
Ion Source & Transfer Temp:	300 °C
Electron Multiplier V:	1400 V by auto tune
Collision gas:	He & N ₂ @ 1.5 and 2.25 mTorr, respectively

Retention Time (RT) and MRM Conditions for LC-MS/MS Analysis

Compound dependent parameters:

DP = declustering potential, CE = collision energy, EP = entrance potential, CXP = collision cell exit potential

Q1	Q3	RT (min)	Analyte	DP	EP	CE	CXP
184.1	143	2.4	Acephate. 1	61	10	13	4
184.1	49	2.4	Acephate. 2	61	10	33	4
223	126	5.2	Acetamiprid. 1	61	10	29	12
223	99	5.2	Acetamiprid. 2	61	10	53	18
228.1	186.1	7	Ametryn. 1	71	10	21	4
228.1	96	7	Ametryn 2	71	10	35	4
209.1	152	3.1	Aminocarb.1	71	10	21	8
209.1	137.1	3.1	Aminocarb.2	71	10	35	10
318	160.1	7.1	Azinphos-methyl	41	10	13	10
318	132	7.1	Azinphos-methyl	41	10	21	10
224.1	109	5.8	Bendiocarb 1	61	10	27	20
224.1	167.1	5.8	Bendiocarb 2	61	10	15	12
440.1	181.2	13.6	Bifenthrin 1	51	10	39	14
440.1	166.1	13.6	Bifenthrin 2	51	10	65	10
343	307	7.8	Boscalid.1	91	10	27	4
343	140	7.8	Boscalid.2	91	10	27	4
197	117.2	4.4	Chlordimeform	81	10	41	18
197	89	4.4	Chlordimeform	81	10	71	14
350	198	12.3	Chlorpyrifos	56	10	25	10
350	97	12.3	Chlorpyrifos	56	10	47	10
362.8	227	10.2	Coumaphos	71	10	37	12
362.8	306.9	10.2	Coumaphos	71	10	25	18
241.1	214.2	5.7	Cyanazine	66	10	27	18
241.1	104.1	5.7	Cyanazine	66	10	47	4
199.1	89.1	7.3	Cycluron	50	10	21	4
199.1	89	7.3	Cycluron	50	10	21	4
292	70	8	Cyproconazole A1	66	10	39	12
292	125	8	Cyproconazole A2	66	10	45	8
292.1	70.1	8.4	Cyproconazole B1	66	10	39	12
292.1	125.1	8.4	Cyproconazole B2	66	10	45	8
318.1	182	6.7	Desmedipham.1	41	10	19	12
318.1	136	6.7	Desmedipham.2	41	10	33	10
305	169.1	9.9	Diazinon	86	10	31	10
305	153.1	9.9	Diazinon	86	10	29	8
350	123	8.3	Dichlorfluanid 1	21	10	41	10
350	224	8.3	Dichlorfluanid 2	21	10	21	10
220.8	127.1	5.9	Dichlorvos	71	10	27	22
220.8	109.1	5.9	Dichlorvos	71	10	25	18
238.1	112.1	4.6	Dicrotophos.1	66	10	19	8
238.1	193	4.6	Dicrotophos.2	66	10	15	14
406.1	251.1	11.6	Difenoconazole 1	81	10	37	16

408.2	253.1	11.6	Difenoconazole 2	76	10	33	4
230	199	4.6	Dimethoate.1	50	10	14	15
230	125	4.6	Dimethoate.2	50	10	27	8
388.1	301	8.1	Dimethomorph A1	66	10	25	4
388.1	165.1	8.1	Dimethomorph A2	66	10	45	4
388.2	301.1	8.4	Dimethomorph B1	66	10	25	4
388.2	165.2	8.4	Dimethomorph B2	66	10	45	4
224.1	167	4.7	Dioxacarb.1	51	10	13	10
224.1	123	4.7	Dioxacarb.2	51	10	23	24
330	121.1	9.5	Epoxiconazole. 1	66	10	29	10
330	101.1	9.5	Epoxiconazole. 2	66	10	69	18
162	119	8.4	Ethiolate. 1	106	10	23	20
162	120.1	8.4	Ethiolate. 2	106	10	19	20
384.8	199.2	12	Ethion. 1	51	10	15	18
384.8	142.9	12	Ethion. 2	51	10	39	24
287.1	121.1	7.1	Ethofumesate. 1	81	10	23	8
287.1	259.1	7.1	Ethofumesate. 2	81	10	15	16
394.2	177.3	13.6	Etofenprox NH₄ +1	46	10	21	12
394.2	107.2	13.6	Etofenprox NH₄ +2	46	10	61	18
337	124.9	9.4	Fenbuconazole.1	81	10	41	8
337	70	9.4	Fenbuconazole.2	81	10	39	12
302.1	88	9.2	Fenoxycarb.1	66	10	31	6
302.1	116.1	9.2	Fenoxycarb.2	66	10	17	8
304	147	7.2	Fenpropimorph.1	66	10	39	4
304	117	7.2	Fenpropimorph.2	66	10	71	4
266	229	7.6	Fludioxinil.1	41	10	23	14
266	227.1	7.6	Fludioxinil.2	41	10	13	14
376	307	8.5	Fluquinconazole.1	71	10	33	4
376	349	8.5	Fluquinconazole.2	71	10	25	4
324.1	262.1	7.5	Flutolanil.1	76	10	27	16
324.1	242.1	7.5	Flutolanil.2	76	10	37	14
314.1	70	10.3	Hexaconazole.1	56	10	41	12
314.1	159	10.3	Hexaconazole.2	56	10	41	14
297	159	6.5	Imazalil.1	66	10	33	14
297	201	6.5	Imazalil.2	66	10	27	12
249.1	160	7.7	Linuron.1	61	10	23	4
249.1	182.1	7.7	Linuron.2	61	10	21	4
331	127.1	7.5	Malathion. 1	46	10	17	10
331	99.1	7.5	Malathion. 2	46	10	31	10
142	94	1.7	Methamidophos.1	55	10	20	4
142	125	1.7	Methamidophos.2	55	10	19	8
284.2	252.2	8.7	Metolachlor. 1	56	10	21	10
284.2	176.2	8.7	Metolachlor. 2	56	10	33	10
166.2	109.1	5.6	Metolcarb. 1	36	10	15	10
166.2	94.2	5.6	Metolcarb. 2	36	10	37	10
225.1	127.1	4.7	Mevinphos-E.1	55	10	20	8

225.1	193.2	4.7	Mevinphos-E.2	55	10	10	13
225	127	5.2	Mevinphos-Z.1	55	10	20	8
225	193.1	5.2	Mevinphos-Z.2	55	10	10	13
224.1	127.1	4.1	Monocrotophos.1	51	10	23	12
224.1	98	4.1	Monocrotophos.2	51	10	17	4
215.1	126.1	6.4	Monolinuron.1	51	10	23	4
215.1	99	6.4	Monolinuron.2	51	10	41	4
289	70	8.3	Myclobutanil.1	71	10	37	12
289	125	8.3	Myclobutanil.2	71	10	47	8
315	252.1	7.4	Nuarimol.1	81	10	31	16
315	81	7.4	Nuarimol.2	81	10	45	14
214	124.9	3	Omethoate.1	46	10	29	4
214	182.8	3	Omethoate.2	46	10	17	4
284.1	159	10.4	Penconazole.1	71	10	39	10
284.1	70	10.4	Penconazole.2	71	10	37	12
318	160	7.1	Phosmet.1	51	10	19	10
318	133	7.1	Phosmet.2	51	10	49	10
356.2	177.2	12.1	Piperonyl butoxide 1	51	10	19	10
356.2	119.1	12.1	Piperonyl butoxide 2	51	10	51	8
239.2	72.1	5.9	Pirimicarb.1	66	10	35	12
239.2	182.1	5.9	Pirimicarb.2	66	10	23	12
376	308	10.9	Prochloraz.1	46	10	17	10
376	70	10.9	Prochloraz.2	46	10	45	12
242.2	158.1	7.8	Prometryn.1	71	10	35	4
242.2	200.1	7.8	Prometryn.2	71	10	19	4
212.2	169.9	6.6	Propachlor. 1	66	10	23	30
212.2	93.9	6.6	Propachlor. 2	66	10	39	16
368.2	231.1	12.6	Propargite.1	46	10	15	14
368.2	175.1	12.6	Propargite.2	46	10	23	12
342.1	159	10.6	Propiconazole.1	61	10	39	10
342.1	69	10.6	Propiconazole.2	61	10	37	12
210.1	111	5.8	Propoxur.1	39	10	19	6
210.1	168.1	5.8	Propoxur.2	39	10	11	11
218.1	125	6	Pyracarbolid.1	61	10	27	8
218.1	97	6	Pyracarbolid.2	61	10	41	18
388	194	10.5	Pyraclostrobin.1	31	10	19	4
388	163	10.5	Pyraclostrobin.2	31	10	29	4
365	147	13.3	Pyridaben.1	46	10	31	4
365	309	13.3	Pyridaben.2	46	10	19	4
200	107	7.7	Pyrimethanil.1	71	10	33	4
200	82	7.7	Pyrimethanil.2	71	10	35	4
308.1	162.1	12.9	Quinoxifen.1	81	10	65	10
308.1	197.1	12.9	Quinoxifen.2	81	10	45	12
226.2	170.1	6.5	Secbumeton.1	50	10	35	4
226.2	100	6.5	Secbumeton.2	50	10	35	4
298.2	144.2	7.9	Spiroxamine.1	76	10	29	12

298.2	100.1	7.9	Spiroxamine.2	76	10	47	18
323	115	8.9	Sulfotep. 1	46	10	39	10
323	97.1	8.9	Sulfotep. 2	46	10	45	10
308.2	70	9.9	Tebuconazole.1	81	10	49	12
308.2	125	9.9	Tebuconazole.2	81	10	51	8
334	117	12.1	Tebufenpyrad.1	71	10	47	4
334	145	12.1	Tebufenpyrad.2	71	10	37	4
230.3	174.2	7.7	Terbutylazine 1	41	10	27	10
230.3	68	7.7	Terbutylazine 2	41	10	59	10
372.1	159	8.8	Tetraconazole.1	76	10	45	10
372.1	70	8.8	Tetraconazole.2	76	10	47	12
202.1	175.1	4.9	Thiabendazole.1	85	10	35	12
202.1	131.2	4.9	Thiabendazole.2	85	10	45	8
364	237.9	9.5	Tolyfluand.1	6	10	19	10
364	137.1	9.5	Tolufluand.2	6	10	37	10
294	197.1	7.8	Triadimefon.1	66	10	23	14
294	225	7.8	Triadimefon.2	66	10	19	8
296.1	70	8	Triadimenol.1	46	10	31	12
296.1	227.1	8	Triadimenol.2	46	10	19	14
314	162	8.3	Triazophos 1	56	10	25	10
314	119	8.3	Triazophos 2	56	10	49	10
190	163	5.8	Tricyclazole 1	81	10	33	10
190	136	5.8	Tricyclazole 2	81	10	41	12
409	186	11.2	Trifloxystrobin. 1	31	10	23	4
409	206	11.2	Trifloxystrobin. 2	31	10	21	4
346.1	278.1	11.7	Triflumizole. 1	51	10	15	8
346.1	73	11.7	Triflumizole. 2	51	10	27	6
346.1	278.1	11.8	Triflumizole. 1	51	10	15	8
346.1	73	11.8	Triflumizole. 2	51	10	27	6

GC-MS/MS Conditions for GC-amenable Pesticides

Analyte	Precursor 1	Product 1	Collision Energy	Precursor 2	Product 2	Collision Energy	RT (min)
Amitraz	293.1	162	6	293.1	132	25	14.77
Benfluralin	292	160	22	292	206	12	7.29
BHC-alpha	219	183	7	181	145	15	7.64
BHC-beta	219	183	8	217	181	7	8.03
BHC-delta	219	183	8	217	181	7	8.51
BHC-gamma	219	183	8	217	181	7	8.04
Bromopropylate	338.9	182.9	18	342.9	184.9	18	13.89
Cadusafos	159	97	24	158	81	15	7.44
Chlorothalonil	265.9	133	53	265.9	169.9	28	8.59
Chlorpyrifos-methyl	285.9	93	24	285.9	208	15	9.13
Cypermethrin	181	152	30	163	127	4	16.56
Dacthal	298.9	164.9	54	300.9	222.9	30	10.04
DEF	202	147	2	202	113	18	11.57

Dieldrin	262.9	192.9	40	262.9	190.9	38	11.7
Dinitramine	261	195	23	261	241	10	8.4
Endosulfan Sulfate	271.9	236.9	15	271.9	116.9	48	13
Endosulfan-I	240.9	205.9	15	195	159	8	11.25
Endosulfan-II	195	159	8	240.9	205.9	15	12.25
Endrin	262.9	192.9	40	262.9	190.9	38	12.1
EPN	157	110	14	185	110.1	25	13.92
Etridiazole	210.9	182.9	9	210.9	139.9	26	5.87
Fenarimol	219	107	12	251	139	15	15.06
Fenvalerate 1	167	125	12	125	89	23	17.38
Fenvalerate 2	167	125	12	125	89	23	17.58
Fluvalinate 1	250	55	18	250	200	24	17.55
Fluvalinate 2	250	55	18	250	200	24	17.6
Heptachlor	352.8	262.8	15	352.8	281.9	18	10.6
Hexachlorobenzene	283.9	213.9	40	283.8	248.9	22	7.78
L-Cyhalothrin	197	141	13	181	152	29	14.85
Iprodione	314	56	24	314	245	10	13.68
Methyl Parathion	263	109	12	263	79	32	9.13
MGK-264	164	80	32	164	98	12	10.42
Napropamide	271.1	72	15	271.1	128	2	11.39
o,p'-DDT	235	165	30	235	199	18	12.42
o,p'-Methoxychlor	227	121	15	121	78	26	13.19
o-phenylphenol	170	115.1	45	170	141	30	6.27
Oxadixyl	163	132	10	163	117	30	12.42
p,p'-DDE	246	176	35	318	246	25	11.6
p,p'-DDT	235	165	30	235	199	18	13.01
Parathion	291	109	10	291	81	35	9.96
Pentachloroaniline	262.9	191.9	25	264.9	193.9	28	8.91
Pentachlorobenzene	249.9	214.9	21	249.9	141.9	50	6.38
Permethrin-cis	183	153	18	183	115	30	15.62
Permethrin-trans	183	153	18	183	115	30	15.74
Phosalone	182	75	36	182	111	17	14.56
Pirimiphos-methyl	290	125	24	290	233	10	9.58
Procymidone	283	96	10	283	67	37	10.83
Profenofos	336.9	266.9	14	336.9	188	32	11.53
Pronamide	173	74	50	173	109	30	8.18
Propanil	161	99	30	217	161	7	8.93
Pyriproxifen	136	41.1	18	136	78.1	32	14.6
Quinalphos	157	102	28	146	118	10	10.72
Tetradifon	353.9	159	12	353.9	227	9	14.39
Tolclofos-methyl	265	93	26	265	109	52	9.22
Triallate	268	183.9	20	268	226	12	8.56
Trifluralin	306	264	7	306	160	25	7.25
Vinclozolin	212	172	16	187	124	22	9.1

**Average Recovery and RSD of 79 Pesticides Spiked in Avocado at Three Concentrations
via LC-MS/MS Analysis**

Analyte	10 ng/g spike level N=5		50 ng/g spike level N=5		200 ng/g spike level N=5	
	Recovery %	RSD %	Recovery %	RSD %	Recovery %	RSD %
Acephate	104.9	5.0	82.6	11.8	92.6	6.3
Acetamiprid	102.7	6.7	84.6	8.9	96.4	3.9
Ametryn	99.8	3.9	84.3	11.4	91.4	6.1
Aminocarb	104.4	2.4	83.9	10.3	93.4	5.3
Azinphos-methyl	115.0	7.3	87.7	11.1	98.3	5.6
Bifenthrin	104.7	6.1	85.1	10.7	93.9	8.8
Boscalid	121.6	7.1	105.6	14.3	85.7	6.0
Chlordimeform	120.2	6.7	88.5	15.9	95.2	3.9
Chlorpyrifos	102.3	9.3	86.7	12.7	91.9	4.4
Coumaphos	99.0	5.7	81.9	11.2	91.8	4.4
Cyanazine	115.0	2.9	87.1	13.2	87.1	13.2
Cycluron	121.0	4.6	91.9	10.9	103.0	3.8
Cyproconazole A	140.0	9.2	86.7	14.3	93.9	6.0
Cyproconazole B	116.6	8.3	115.2	34.6	102.5	5.1
Desmedipham	110.8	5.9	108.3	29.9	103.4	6.8
Diazinon	112.0	3.7	87.1	11.7	95.3	4.7
Dichlorfluanid	99.8	9.1	84.1	11.0	92.3	4.6
Dichlorvos	83.2	18.8	77.2	9.3	86.8	4.0
Dicrotophos	80.8	14.8	74.7	5.7	93.8	9.5
Difenoconazole	103.6	3.2	84.1	11.9	92.6	5.2
Dimethoate	111.6	5.1	87.3	12.6	100.3	6.9
Dimethomorph A	103.3	4.6	83.9	12.3	92.8	4.2
Dimethomorph B	97.1	5.3	90.3	9.3	98.1	4.9
Dioxacarb	116.0	8.6	86.6	9.5	100.6	5.1
Epoxiconazole	97.2	4.0	83.5	12.1	92.9	5.3
Ethiolate	107.5	4.9	86.5	11.6	98.7	6.8
Ethion	102.0	8.3	88.8	15.8	94.0	8.7
Ethofumesate	98.3	6.5	83.3	11.4	92.4	5.5
Fenbuconazole	107.4	16.9	84.3	14.2	96.3	6.6
Fenoxycarb	104.1	14.4	92.0	12.3	102.3	7.2
Fenpropimorph	105.0	7.1	82.1	11.2	94.1	4.9
Fludioxinil	110.6	8.4	82.0	11.7	92.0	6.2
Fluquinconazole	118.0	13.6	83.9	16.9	102.5	8.9
Flutolanil	146.2	7.5	90.4	18.0	97.6	5.1
Hexaconazole	109.0	4.9	85.8	13.7	93.3	3.5
Imazalil	117.0	4.2	88.4	14.6	100.9	9.5
Linuron	123.4	8.6	94.5	13.9	97.7	6.3
Malathion	103.2	12.4	87.4	14.2	97.1	5.3
Methamidophos	113.0	2.3	83.1	15.9	93.3	7.7
Metolachlor	102.5	2.5	81.7	11.3	94.4	6.3
Metolcarb	100.1	5.9	83.3	13.3	93.5	4.6
Mevinphos-E	108.1	8.2	84.1	11.0	90.4	3.1
Mevinphos-Z	99.6	14.7	83.9	9.3	91.1	4.8
Monocrotophos	97.0	3.3	82.5	8.7	90.4	4.5

Monolinuron	105.0	4.8	85.1	11.8	93.1	5.4
Myclobutanil	110.4	3.1	87.0	11.6	93.0	4.5
Nuarimol	111.2	12.6	91.8	7.8	96.5	4.5
Omethoate	137.0	15.4	83.8	10.9	98.6	7.3
Penconazole	113.4	7.9	88.4	13.6	96.4	5.8
Phosmet Piperonyl	104.8	3.1	85.4	8.3	96.0	6.4
butoxide	106.0	4.1	83.0	10.4	91.5	6.6
Pirimicarb	104.3	2.9	84.5	11.1	93.0	5.6
Prochloraz	124.2	29.9	83.9	10.7	92.6	5.8
Prometryn	101.0	8.5	85.6	10.4	95.7	5.4
Propachlor	101.0	4.5	81.2	12.6	92.2	5.5
Propargite	109.2	6.7	84.2	7.0	91.8	5.5
Propiconazole	106.2	7.1	85.0	13.2	97.3	9.7
Propoxur	97.0	5.1	83.8	10.3	92.4	4.2
Pyracarbolid	101.3	3.2	82.9	13.4	93.0	5.7
Pyraclostrobin	109.6	7.6	83.8	10.9	93.0	5.3
Pyridaben	95.2	7.1	78.6	10.2	85.8	5.5
Pyrimethanil	107.0	15.4	91.2	12.0	93.3	6.5
Quinoxifen	105.6	6.5	84.6	9.3	92.0	3.1
Secbumeton	103.8	5.8	82.2	8.7	92.7	5.1
Spiroxamine	104.6	4.3	83.4	12.6	94.5	6.1
Sulfotep	108.2	7.7	84.7	11.8	91.7	5.6
Tebuconazole	110.6	5.9	88.2	9.8	102.7	9.6
Tebufenpyrad	106.8	10.9	81.9	11.6	95.3	5.9
Terbutylazine	101.4	5.9	84.0	8.8	93.4	4.3
Tetraconazole	112.4	10.7	89.4	5.6	104.0	5.9
Thiabendazole	110.6	4.2	84.9	10.2	94.4	7.3
Tolyfluanid	129.6	4.2	86.7	9.6	89.8	6.0
Triadimefon	95.9	16.4	86.8	7.0	99.9	6.1
Triazophos	102.9	25.3	89.1	7.7	102.9	6.6
Tricyclazole	104.3	4.5	84.0	9.1	93.3	4.2
Trifloxystrobin	96.8	4.2	82.7	10.6	90.9	5.9
Triflumizole	101.5	5.7	84.0	11.1	92.1	4.7
Average	107.1		86.1		94.8	
Std. Dev	9.9		5.8		4.0	
RSD %	9.2		6.7		4.2	


Average Recovery and RSD of 57 Pesticides Spiked in Avocado at Three Concentrations with GC-MS/MS Analysis

Analyte	10 ng/g spike level N=5		50 ng/g spike level N=5		200 ng/g spike level N=5	
	Recovery %	RSD %	Recovery %	RSD %	Recovery %	RSD %
Amitraz	31.8	12.7	38.3	18.0	58.0	7.2
Benfluralin	81.3	9.4	68.5	12.5	91.3	4.8
BHC-alpha	74.9	5.2	76.1	11.9	95.7	3.5
BHC-beta	93.4	12.2	73.2	20.4	103.5	2.7

BHC-delta	70.5	4.8	76.5	12.0	95.4	4.1
BHC-gamma	84.2	12.2	73.1	20.5	101.7	3.5
Bromopropylate	60.2	15.7	69.2	13.7	97.1	5.1
Cadusafos	69.8	3.4	68.8	11.4	92.0	3.1
Chlorothalonil	70.4	28.2	52.2	14.2	81.9	19.4
Chlorpyrifos-methyl	79.0	9.0	73.7	12.4	92.3	7.5
Cypermethrin	130.7	11.0	104.2	10.3	92.3	5.9
Dacthal	70.1	7.5	71.1	14.5	90.2	3.4
DEF	57.1	18.9	61.6	11.0	94.2	6.6
Dieldrin	83.0	26.3	73.8	11.3	94.4	3.6
Dinitramine	92.2	6.5	77.7	12.0	95.2	4.6
Endosulfan Sulfate	106.9	14.2	69.2	22.4	106.2	5.8
Endosulfan-I	91.4	31.7	72.6	16.2	92.2	11.3
Endosulfan-II	78.2	7.3	70.6	9.2	100.0	5.9
Endrin	99.7	12.6	73.4	11.9	100.0	5.7
EPN	66.7	26.7	68.5	13.9	107.5	4.8
Etofenprox	82.8	8.9	78.8	11.6	89.0	4.8
Etridiazole	104.7	7.0	68.7	15.1	110.4	11.2
Fenarimol	63.2	7.7	65.8	15.3	96.9	6.6
Fenvalerate 1	72.2	27.7	76.9	14.3	102.9	7.7
Fenvalerate 2	75.4	20.2	63.9	22.5	92.3	3.9
Fluvalinate 1	58.4	31.4	65.0	17.9	99.6	5.2
Fluvalinate 2	51.5	37.4	57.5	27.5	81.7	11.9
Heptachlor	65.4	17.7	69.7	13.3	95.1	6.1
Hexachlorobenzene	60.6	9.1	61.6	11.9	81.0	6.1
L-Cyhalothrin	66.3	13.9	75.2	9.3	98.0	6.2
Iprodione	37.0	82.8	68.7	14.1	92.7	16.9
Methyl Parathion	75.0	14.1	77.0	13.8	95.6	5.2
MGK-264	74.1	10.1	70.8	11.7	97.7	2.0
Napropamide	74.4	10.2	74.7	15.4	103.7	4.9
o,p'-DDT	94.2	20.3	62.1	29.8	119.2	23.1
o,p'-Methoxychlor	80.5	12.3	84.9	18.5	112.0	15.3
o-phenylphenol	105.0	17.9	76.7	11.3	83.6	5.1
Oxadixyl	64.6	8.6	73.9	13.4	76.6	6.6
p,p'-DDE	61.5	7.4	67.2	14.3	89.0	4.7
p,p'-DDT	NA	NA	NA	NA	NA	NA
Parathion	58.5	14.6	66.4	13.3	94.2	4.6
Pentachloroaniline	71.3	5.0	70.0	11.7	89.9	3.8
Pentachlorobenzene	70.5	4.6	68.2	13.0	85.4	3.8
Permethrin-cis	89.9	12.5	62.1	13.8	93.6	4.8
Permethrin-trans	98.5	14.1	74.7	34.7	111.6	9.1
Phosalone	74.4	15.0	75.6	11.0	108.0	8.5
Pirimiphos-methyl	77.7	11.5	72.2	12.7	92.5	2.1
Procymidone	76.8	5.0	75.6	11.6	98.5	13.5
Profenofos	52.2	37.2	95.1	6.5	89.6	3.7
Pronamide	71.3	8.6	71.7	15.7	93.2	5.2

Propanil	72.4	9.0	72.2	13.8	96.1	6.4
Pyriproxifen	64.8	7.4	67.9	13.4	96.1	6.4
Quinalphos	79.5	15.8	67.5	13.4	91.1	5.0
Tetradifon	66.3	5.9	72.1	11.3	88.4	8.5
Tolclofos-methyl	81.6	3.7	75.4	10.9	94.5	3.7
Triallate	70.3	4.4	67.4	17.1	92.3	4.7
Trifluralin	63.9	9.2	70.8	10.4	95.5	5.7
Vinclozolin	71.5	11.0	70.6	9.8	101.3	6.5
Average	73.9		70.2		94.3	
Std. Dev	15.0		7.9		17.0	
RSD %	20.3		11.3		18.0	

Reconstructed LC-MS/MS Chromatogram of Avocado Blank, Avocado Blank Fortified at 10 ng/g, and Avocado Blank Spiked with 50 ng/g Standard Mix Sample Concentration is 0.12 G
Sample/MI Solvent With 1 μ L Injection Volume



* Adapted from: 'Analysis of 136 Pesticides in Avocado Using a Modified QuEChERS Method with LC-MS/MS and GC-MS/MS'
Narong Chamkasem ^a, Lisa W. Ollis ^a, Tiffany Harmon ^a, Sookwang Lee ^a and Greg Mercer ^b
^a Southeast Regional Laboratory, U.S. Food and Drug Administration, 60 Eighth Street, N.E., Atlanta, GA, 30309
^b Pacific Regional Laboratory – Northwest, U.S. Food and Drug Administration, 22201 23rd Drive, S.E., Bothell, WA 98021



QuEChERS Determination of 2, 4-Dichlorophenoxyacetic Acid (2, 4-D) and Other Acidic Herbicides by LC-MS/MS

UCT Part Numbers:

ECQUEU750CT-MP (4000 mg MgSO₄, 1000 mg NaCl, 500 mg Na Citrate dibasic sesquihydrate, 1000 mg Na Citrate tribasic dihydrate)

ECMSC1850CT contains (500 mg C18 and 1500 mg MgSO₄)

August 2013

The pesticide 2,4-dichlorophenoxyacetic acid (2,4-D; Agent Orange; CASRN 94-75-7; pKa 2.73) along with other acidic pesticides can be extracted and analyzed using this QuEChERS approach.

Procedure

1. Sample Preparation

- a) Weigh 10 g of homogenized, hydrated sample into a 50 mL centrifuge tube
- b) Add ISTD
- c) Add NH₄OH dropwise with shaking until pH 12
- d) Cap and let sit for 30 min at room temperature
- e) Neutralize sample by evaporation of NH₄. Dry N₂ may be used
- f) Add 10 mL of acetonitrile
- g) Shake
- h) Add contents of **ECQUEU750CT-MP** pouch to centrifuge tube
- i) Immediately shake or vortex for 2 min
- j) Centrifuge for 5 min at >1500 rcf
- k) Remove 5 - 10 mL aliquot and add to **ECMSC1850CT** tube
- l) Shake vigorously or vortex for 2 min
- m) Centrifuge for 5 min at > 1500 rcf
- n) Pass extract through a 0.45 µm syringe filter into sample vial
- o) Analyze by LC-MS/MS

Instruments

MS: AB Sciex 5500 QTrap

HPLC: Shimadzu LC-20AD pump, Sil-20AC autosampler & CTO-20AC column oven

LC Parameters

Guard Column	Agilent Eclipse Plus C18, 5 µm, 4.6 x 12.5 mm
Column	Agilent Eclipse Plus C18, 1.8 µm, 4.6 x 75 mm
Flow	0.5 ml/min
Injection Vol	10 µl
Oven	40 °C
Equilibration Time	0 min
Rinsing Vol	200 µl
Rinsing Speed	35 µl/sec
Sampling Speed	15 µl/sec
Cooler Temperature	15 °C

Mobile Phase

Aqueous A: 0.1 % formic acid in water

Organic B: 0.1 % formic acid in methanol

Time	%B
0	10
3.0	70
8.0	95
12.0	95
12.5	10
14.5	10


MS Instrument Settings

MRM Window	40 sec
Target Scan	0.25 sec
Cycle	0.5 sec
Resolution	Unit
PosNeg Switch	50 mS
Pause	5 mS
Curtin Gas	50 ml/min
Collision	High
Ion Spray Voltage	4500 volts
Turbospray Temperature	450 °C
Ion Source Gas 1	50 ml/min
Ion Source Gas 2	80 ml/min
Interface Heater	On

Analyte	ESI	RT	Q1	Q2	DP	EP	CE	CXP
2,3,6-TBA	-	5.9	224.6	180.8	-30	-10	-12	-18
2,3,6-TBA	-	5.9	222.6	178.9	-30	-10	-12	-18
2,4,5-T	-	7.7	252.6	194.9	-35	-10	-20	-18
2,4,5-T	-	7.7	254.8	196.9	-60	-10	-18	-18
2,4,5-TB	-	8.7	282.9	196.8	-35	-10	-20	-18
2,4,5-TB	-	8.7	280.8	194.8	-35	-10	-16	-18
2,4-D	-	6.8	218.7	160.9	-20	-10	-20	-18
2,4-D	-	6.8	220.6	162.9	-35	-10	-20	-18
2,4-DB	-	7.7	246.6	160.9	-20	-10	-18	-18
2,4-DB	-	7.7	248.6	162.9	-35	-10	-18	-18
4-CPA	-	6	184.7	126.9	-15	-10	-20	-18
4-CPA	-	6	186.7	128.9	-40	-10	-20	-18
Acifluorfen	-	8.2	359.8	316	-25	-10	-14	-18
Acifluorfen	-	8.2	359.8	194.9	-25	-10	-38	-18
Bromoxynil	-	6.4	273.4	78.8	-25	-10	-62	-18
Bromoxynil	-	6.4	275.4	80.8	-20	-10	-60	-18
Chloramben	-	5.2	203.6	159.9	-35	-10	-10	-18
Chloramben	-	5.2	205.6	161.9	-35	-10	-12	-18
Dalapon	-	4.4	140.6	96.9	-50	-10	-12	-18
Dalapon	-	4.4	142.6	98.9	-40	-10	-12	-18
Dicamba	-	6	218.7	175	-45	-10	-10	-18
Dicamba	-	6	220.4	176.9	-30	-10	-10	-18
Dicamba	-	6	218.7	35	-45	-10	-40	-18
Dichloroprop	-	7.4	232.6	160.9	-20	-10	-18	-18

Dichloroprop	-	7.4	234.5	162.9	-20	-10	-18	-18
Diclofop	-	9	324.9	253	-25	-10	-22	-18
Diclofop	-	9	326.9	255	-25	-10	-22	-18
Fenac	-	7.2	192.7	35	-55	-10	-35	-18
Fenac	-	7.2	194.7	35	-55	-10	-35	-18
Haloxifop	-	8.4	359.8	288.1	-40	-10	-20	-18
Haloxifop	-	8.4	361.9	290	-35	-10	-20	-18
MCPA	-	6.8	198.6	140.9	-20	-10	-20	-18
MCPA	-	6.8	200.7	142.9	-30	-10	-20	-18
MCPB	-	7.8	226.8	141	-30	-10	-18	-18
MCPB	-	7.8	228.5	143	-20	-10	-18	-18
Mecoprop	-	7.4	212.7	141	-20	-10	-22	-18
Mecoprop	-	7.4	214.6	143	-30	-10	-24	-18
PCP	-	10	264.4	35	-50	-10	-66	-18
PCP	-	10	262.4	35	-60	-10	-68	-18
Quizalofop	-	8.4	343	271	-20	-10	-22	-18
Silvex	-	8.2	266.6	194.8	-30	-10	-18	-18
Silvex	-	8.2	266.6	158.9	-30	-10	-40	-18
Triclopyr	-	7.2	255.4	197.9	-35	-10	-18	-18
Triclopyr	-	7.2	253.3	195.8	-35	-10	-18	-18
Diflufenzopyr	-	6.5	333	128	-50	-10	-34	-18

Clean-Up and Fractionation





Clean-Up of Organochlorine Pesticides and PCB Extracts Using FLORISIL®

UCT Part Numbers:

EUFLSA1M6 (1000 mg small particle Grade A Florisil® for slower gravity flow)

EUFLS1M6 (1000 mg regular particle PR Grade Florisil® for more viscous samples)

December 2012

This Florisil® is a magnesium silicate with basic properties. It is used to remove polar interferences from organochlorine pesticide and PCB extracts in hexane prior to chromatographic analysis. Other analytes may potentially be cleaned up using this method provided that adequate performance is demonstrated.

Procedure

1. Pre-rinse a column with 9 mL of 90:10 hexane/acetone by gravity flow
(A low vacuum may be required to start flow)
2. Discard solvent
3. Add a collection tube under the column
4. Add a 2 mL aliquot of the sample extract (in hexane) to the column
5. Collect extract by gravity
6. Add 9 mL of 90:10 hexane/acetone to the column
7. Continue to collect by gravity
8. Gently evaporate the extract to a volume of 1 mL
9. Bring to a final volume of 2 mL with hexane

Florisil® is a registered trademark of U.S. Silica Co.



Fractionation of Aliphatic And Aromatic Hydrocarbons Using ENVIRO-CLEAN® EPH SILICA

(Developed in cooperation with Lancaster Laboratories, Inc.)

UCT Part Numbers:

XRSIHT13M15 (Enviro-Clean® EPH Fractionation 3000 mg/15 mL)

February 2009

The composition of petroleum is a complex mixture of hundreds of different hydrocarbon compounds. The resultant makeup of hydrocarbons released into the environment is variable and dependent on the conditions to which it is subsequently exposed. While in the environment, petroleum composition is influenced by a number of factors including volatilization, leaching, and/or biological degradation. These environmental effects yield a mixture whose toxicological properties can be vastly different than the parent product. Based on the known toxicological properties of petroleum products we can assume that:

- aromatic compounds are more toxic than aliphatic compounds
- the toxicity of aliphatic compounds is dependent upon their molecular weight with low molecular weight compounds showing relatively higher toxicity

The fractionation of the total petroleum hydrocarbon extract is necessary to determine the concentration of the aliphatic versus aromatic compounds. The Massachusetts Department of Environmental Protection (MADEP) has taken the approach of fractionating the C9-C18 aliphatics (n-nonane to n-octadecane), C19-C36 aliphatics (n-nonadecane to hexatriacontane), and the C11-C22 aromatics (naphthalene to benzo (g,h,i)perylene). These compounds are associated with the release of hydrocarbons in the environment. The aromatics are considered the most toxic form of hydrocarbon.

Procedure

1. Prepare Extract

- a) Solvent exchange the hydrocarbon extract from methylene chloride to hexane using a K-D apparatus

2. Prepare Cartridge

- a) Thoroughly rinse cartridge with two, 10 mL aliquots of pentane
- b) Add 1 mL of the extract to the cartridge
- c) Elute aliphatic fraction with pentane by gravity and collect everything in an ampoule.
A total of 10 mL should be collected
- d) Place a fresh ampoule under the cartridge and elute the aromatic fraction with methylene chloride by gravity. A total of 10 mL should be collected

- e) Concentrate each fraction separately to a final volume on a steam bath using an ampoule and micro-Snyder column combination. Other techniques may be used but the loss of C9-C18 hydrocarbons may result

It is very important to keep the silica cartridges dry and away from room air prior to use. Moisture and contaminants in the air will reduce the effectiveness of the silica and may cause contamination of the extract. Pre-rinsing the cartridges with acetone may reduce this problem.

Results

Classification	Range	Percent Recovery
Aromatics	C11-C22	88
surrogates	2-fluorobiphenyl	123
surrogates	o-terphenyl	100
Aliphatics	C9-C18	85
	C19-C36	89
surrogates	1-chlorooctadecane	58

MA EPH DATA from Lancaster Labs



Removal of Sulfur from Environmental Samples Using Copper Beads

UCT Part Numbers:

ECCU01K (Copper granules 99.5% 30 mesh - 1kg)

February 2009

1. Post Sample Extraction

- a) Place 4 g of copper beads in a glass vial
- b) Add 2 mL of liquid sample extract to the vial

2. Sulfur Removal

- a) Seal the glass vial and mix sample with copper beads for 2 minutes
- b) Allow to stand for approximately 10 minutes
- c) If sample contains high levels of sulfur, repeat process with 4 g of fresh copper beads

Note: For the analysis of PCB type analytes, copper may reside in the extract

3. Analysis by GC/MS or GC/ECD

- a) Transfer clean extract to autosampler vial
- b) Inject 1-2 μL for GC
- c) Inject 5-10 μL for LC



Determination of Phthalates Leached from Toys into Artificial Saliva

UCT Product #: ECPAHFR50CT (50 mL polypropylene centrifuge tubes)

ECQUUS2-MP (Mylar pouch with 4000 mg MgSO₄ and 2000 mg NaCl)

November 12, 2012

Phthalates are plasticizers added to polymers to make them softer and more flexible. They are commonly used in the manufacturing of toys. Since phthalates are not permanently bonded in polymers, they may be released through touching, licking, or chewing. This new method artificially reproduces the licking and chewing actions of children by shaking cut toys in artificial saliva. Following extraction, phthalate levels are determined by GC/MS.

Preparation of Artificial Saliva (AS)¹

AS was prepared by mixing 0.18 g xanthan gum, 1.2 g potassium chloride, 0.85 g sodium chloride, 0.05 g magnesium chloride, 0.13 g calcium chloride, and 0.13 g di-potassium hydrogen orthophosphate with 1 L of reagent water (1) and stirred for 4 hours.

Procedure

1. Cut plastic toys into small pieces
2. Weigh 1 to 2 g of toy samples to a 50 mL centrifuge tube **ECPAHFR50CT**
3. Add 10 mL of **AS** and two stir bars to the 50 mL tube
4. Shake for 1 hour using a horizontal shaker
5. Transfer the **AS** to a new 50 mL tube, add 10 mL of ethyl acetate and shake for 1 min
6. Add salts from pouch **ECQUUS2-MP** and shake vigorously for 1 min
7. Centrifuge at 5000 rpm for 5 min
8. Transfer 1 mL of the supernatant into a 2 mL auto-sampler vial
9. Add 10 µL of 50 ppm triphenyl phosphate (TPP) as internal standard
10. The sample is ready for GC/MS analysis

INSTRUMENT CONDITIONS:

GC/MS: Agilent 6890N GC coupled with 5975C MSD, equipped with 7683 auto sampler

Injector: 1µL splitless injection at 250 °C, 30 mL/min split vent at 1 min

Liner: 4 mm splitless gooseneck, 4mmID*6.5mmOD*78.5mm (UCT cat#: GCLGN4MM)

GC capillary column: Restek Rxi-5sil MS 30m*0.25mm*0.25µm integrated with 10m guard column

Oven temperature program: Initial temperature of 70 °C, hold for 1 min; ramp at 20 °C/min to 315 °C, hold for 4.75 min. Acquire data from 6 to 14 min.

Carrier gas: He at a constant flow of 1.2 mL/min.

MSD: Transfer line: 280 °C; MS Source (EI): 250 °C; MS Quad: 150 °C

Simultaneous Scan/SIM: Scan range: 55-350 amu

SIM Table

Compound	Abbreviation	Rt (min)	Group #	Start (min)	Dwell time (ms)	Quantify ion	Qualifier ion 1	Qualifier ion 2
Dimethyl phthalate	DMP	6.989	1	6	25	163	194	133
Diethyl phthalate	DEP	7.858	2	7.5	25	149	177	105
Dibutyl phthalate	DBP	9.865	3	9	25	149	223	150
Benzyl butyl phthalate	BBP	11.716	4	11	25	149	91	206
Triphenyl phosphate	TPP	11.964			25	326	325	
Bis(2-ethylhexyl) phthalate	DEHP	12.432	5	12.3	25	149	167	279
Di-n-octyl phthalate	DOP	13.138			25	149	279	150


Calibration Curves

Matrix matched calibration curves are constructed by analyzing matrix matched standards (0-1000 ng/mL). The responses were linear with correlation coefficient higher than 0.9973. The limit of quantification of this method is 25 ng/mL.

Accuracy and Precision Data (n=4)*		
Analyte	Fortified 200 ng/mL	Fortified 500 ng/mL
	Recovery ± RSD %	Recovery ± RSD %
Dimethyl phthalate	90.2 ± 2.6	94.6 ± 2.1
Diethyl Phthalate	91.1 ± 1.4	95.6 ± 2.0
Dibutyl phthalate	90.6 ± 3.8	97.0 ± 2.2
Benzyl butyl phthalate	85.5 ± 1.2	92.1 ± 2.5
Bis(2-ethylhexyl) phthalate	93.2 ± 2.2	92.7 ± 2.6
Di-n-octyl phthalate	88.8 ± 3.1	92.8 ± 1.3

*The control sample did not contain phthalates

Chromatogram of AS fortified with 200 ng/mL phthalates



Results

Diethyl phthalate was found leaching from toy samples into artificial saliva at a concentration of 285 ng/g (RSD=5.9%, n=3).

References

- (1) A. Preetha and R. Banerjee, Comparison of Artificial Saliva Substitutes, *Trends Biomater. Artif. Organs* **18**(2),178-186 (2005).