

# Pulsed Discharge Detector for the Thermoquest TRACE<sup>TM</sup> GC


- Non-radioactive
- Sub-picogram sensitivity
- Capillary or packed columns
- Plug-in installation on the Thermoquest TRACE GC

## **Description and Operating Principle**

The VICI pulsed discharge detector (PDD), winner of R&D Magazine's prestigious R&D 100 Award, utilizes a stable, low-powered, pulsed DC discharge in helium as an ionization source. Elutants from the column, flowing



counter to the flow of helium from the discharge zone, are ionized by photons from the helium discharge. (*See the illustration below.*) The bias electrode focuses the resulting electrons toward the collector electrode, where they cause changes in the standing current which are quantified as the detector output. The PDD's performance, verified by hundreds of users, is equal to or better than conventional detectors with radioactive sources.



## Mode Selection

### **Helium Photoionization Mode**

In the helium photoionization mode, the PDD is a universal, non-destructive, high sensitivity detector. The closeto-true-mass response to both inorganic and organic compounds is linear over a wide range. Response to fixed gases is positive (standing current increases), with an MDQ in the low ppb range. The PDD in this mode is an excellent replacement for flame ionization detectors in petrochemical or refinery environments, where the flame and the use of hydrogen can be problematic.

#### Selective Photoionization Mode

When the helium discharge gas is doped with a suitable noble gas, such as argon, krypton, or xenon (depending on the desired ionization potential cutoff point), the PDD can function as a specific photoionization detector for selective determination of aliphatics, aromatics, amines, and other species. Any problems associated with the presence of a window between the photon source and the ionization chamber are eliminated. In most applications involving current commercial PIDs, analyte condensation and decomposition on the window attenuate the lamp energy, necessitating frequent cleaning and recalibration.


The following patent numbers apply to this product: 5.153.519 5.317.271 5.394.090

5,153,519	5,317,271	5,394,090
5,394,091	5,541,519	5,532,599
5,528,150	5,594,346	5,394,092



CE

## Sample Chromatograms



Impurities in Valco Instruments Co. Inc. food grade  $CO_2$ J. Madabushi, Ph.D. cos bpm 9 qdd Acetaldehyde, m-Xylene, 18 Foluene, 11 ppb Benzene, 14 ppb 450 ppb ĉ cos. 2 ppm <sub>2</sub>S, 130 ppb Food and beverage grade CO<sub>2</sub> is routinely analyzed for trace level impurities which could impart an undesirable flavor to the product or pose a health hazard. Used in the argon photoionization mode, the PDD has little response to the CO<sub>2</sub> matrix peak but excellent sensitivity to all peaks of interest. This analysis previously required sample preconcentration and multiple detectors to achieve similar results.

## **Ordering Information**

#### **Complete Detector Kit**

Includes the detector and all the connection hardware, filters, and electronics to plug it directly into TRACE GC's

Description

Product number

Available from

Designed and manufactured by



**Valco Instruments Co. Inc.** tel: 800 367-8424 fax: 713 688-8106 valco@vici.com Printed in USA © 2000 9/00