

Dual Cell Microvolume Thermal Conductivity Detector

- Stand-alone unit
- Optimized for capillary chromatography
- Thermal stability to ±0.02°C
- Dual filaments capable of independent or referenced (differential) operation



The Valco Microvolume Thermal Conductivity Detector (TCD) is useful in a wide variety of capillary and packed column applications. Constant filament temperature control provides a linear dynamic range permitting measurement of a wide range of concentrations without the need for multiple standards or sample dilution.


Since the detector is non-destructive of the sample and contributes virtually no band spreading, it can be used in series with other detectors without affecting the performance characteristics of either.


Description

The detector consists of the cell housing and the electronics controller. The cell design permits mounting in virtually any orientation with no effect on performance. It can be installed easily on virtually any gas chromatograph, comprising a stand-alone unit requiring nothing else for operation but carrier gas flow.

Each of the two cell chambers is independent of the other, except for block temperature. Filaments can be replaced individually. Front panel controls set the temperature (up to 400°C) for the cell and for each filament. Since each detector cell can be operated separately or simultaneously, two analyses can be run using a single Valco TCD.


To insure compatibility with any system, two outputs are provided: 0-1 mV full scale attenuated output for recorders, and 0-10 V full scale unattenuated output for integrators and data acquisition systems.





Gas Standard 30 m x 530 micron PLOT column, 100 µl valve injection

Independent vs. referenced operation



Column temp: 65°C Detector temp: 100°C Filament temp setting: 5.0 Flow rate Channel A: 5.5 ml/min

Channel B: 5.42 ml/min

Specifications

Overal

Linear range 1 nanogram to 3 micrograms $nC_{\!4}$	
Minimum detectable approx. 50 picograms n-butane quantity	
Time constant < 150 milliseconds	
Cell temperature Automatic proportional control with ±0.02°C stability	
Maximum cell 400°C temperature	

Detector assembly

Dimensions	. 3.12" x 6" x 3.75" high (8 cm x 15 cm x 9 cm)
Gas connections	Valco 1/16" zero dead volume fittings
Single multi-pin	. 5 foot cable supplied

Control unit

Dimensions	12" x 8" x 5" high
	(30 cm x 20 cm x 13 cm)

	· · · · · · · · · · · · · · · · · · ·
Electrical connections	Single multi-pin connector
Operator controls	 Cell temperature control (40-400°C) 10-turn filament temperature potentiometers (A & B) 10-turn coarse and fine baseline adjustment potentiometers (A & B) 12 position recorder attenuator output switch (A, B, or A-B) Filament power on/off switch
Indicator LEDs	Detector heater "on" Filament power "on"
Power requirements	Universal 100-250 VAC 50/60 Hz, 100W maximum

Product numbers

Dual cell microvolume TCD with:	110 VAC	230 VAC
nickel/iron filaments tungsten/rhenium filaments		TCD2-NIFE-220 TCD2-WRE-220



Valco Instruments Co. Inc. tel: 800 367-8424 fax: 713 688-8106 valco@vici.com



Valco International tel: Int + 41 41 925-6200 fax: Int + 41 41 925-6201 VICI vici@vici.com Valco Instruments