

CHROMSPEC™ AUTOSAMPLER VIALS

Our new line of CHROMSPEC™ vials and caps offer the same high quality and affordability that you have come to expect. The new line also introduces products with important new features:

Bonded Cap/Septa (Available in 9-425 screw cap format only)

Septa are bonded to the cap to prevent "push through" of the septa. This type of cap works very well with Waters Alliance autosamplers and other systems that have blunt point sampling needles.

Ultra-Low Bleed (ULB) Septa (Available in Snap Cap and Screw Cap)

The septa have been conditioned at high temperatures to remove extractables. Use these for highly sensitive analyses where low background levels are a necessity.

- **40% larger opening reduces autosampler needle damage**
- **Manufactured from Type 1 borosilicate glass**
- **Compatible with robotic autosamplers**

12 x 32mm WIDE OPENING CRIMP TOP GLASS VIALS

Product Number	Description	Qty
C223682CLM	Clear Wide Opening Crimp Top Glass Vial	100
C223683MLM	Clear Wide Opening Crimp Top Glass Vial	1000
C7711800W	Clear Wide Opening Crimp Top Glass Vial with White Marking Spot	100
C7711800WM	Clear Wide Opening Crimp Top Glass Vial with White Marking Spot	1000
C7711800A	Amber Wide Opening Crimp Top Glass Vial	100
C7711800AM	Amber Wide Opening Crimp Top Glass Vial	1000
C7711800AW	Amber Wide Opening Crimp Top Glass Vial with White Marking Spot	100
C7711800AWM	Amber Wide Opening Crimp Top Glass Vial with White Marking Spot	1000



11mm CRIMP TOP SEAL WITH SEPTA

C221150C	Silver Aluminum Seal with Teflon/Orange Rubber Septa	100
C221150	Silver Aluminum Seal with Teflon/Orange Rubber Septa	1000



2 mL 9-425 WIDE OPENING SCREW THREAD VIALS

C779100	Clear 2 mL 9-425 Wide Opening Screw Thread Vial	100
C779100M	Clear 2 mL 9-425 Wide Opening Screw Thread Vial	1000
C779100W	Clear 2 mL 9-425 Wide Opening Screw Thread Vial with White Marking Spot	100
C779100WM	Clear 2 mL 9-425 Wide Opening Screw Thread Vial with White Marking Spot	1000
C779100A	Amber 2 mL 9-425 Wide Opening Screw Thread Vial	100
C779100AM	Amber 2 mL 9-425 Wide Opening Screw Thread Vial	1000
C779100AW	Amber 2 mL 9-425 Wide Opening Screw Thread Vial with White Marking Spot	100
C779100AWM	Amber 2 mL 9-425 Wide Opening Screw Thread Vial with White Marking Spot	1000



PREASSEMBLED CAPS AND SEPTA FOR 2 mL 9-425 WIDE OPENING SCREW THREAD VIALS

C779200B	9-425 Screw Cap, Blue*; Red PTFE/White Silicone Septa Non-bonded	100
C779200BM	9-425 Screw Cap, Blue*; Red PTFE/White Silicone Septa Non-bonded	1000
C779200BB	9-425 Screw Cap, Blue*; Red PTFE/White Silicone Septa Bonded	100
C779200BBM	9-425 Screw Cap, Blue*; Red PTFE/White Silicone Septa Bonded	1000
C779200XBB	9-425 Screw Cap, Blue*; Red PTFE/White Silicone Septa Bonded Pre-Slit	100
C779200XBBM	9-425 Screw Cap, Blue*; Red PTFE/White Silicone Septa Bonded Pre-Slit	1000
C779400B	9-425 Screw Cap, Blue*; Clear PTFE/White Silicone Septa Ultra Low Bleed	100
C779400BM	9-425 Screw Cap, Blue*; Clear PTFE/White Silicone Septa Ultra Low Bleed	1000

* Caps are available in a variety of colours. Please inquire.