

Product Datasheet

Goat Anti-Mouse IgG H&L, BF647 conjugated (orb1809470)

Catalog Number	orb1809470
Category	Antibodies
Description	The Goat Anti-Mouse IgG H&L, BF647 conjugated is a BF647 conjugated polyclonal antibody suitable for FC, ICC, IF. This product reacts with Mouse samples.
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	BF647
Reactivity	Mouse
Form/Appearance	Liquid
Concentration	2.0 mg/ml
Buffer/Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Purification	Affinity purified by Protein G, nonspecific adsorbed
Immunogen	Native Mouse IgG
Tested applications	FC, ICC, IF
Dilution range	IF=1:200-1000, Flow-Cyt=1:50-200, ICC/IF=1:100-1000
Antibody Type	Secondary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)