



CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

Compressed Gas Purifiers

CHROMSPEC™ He Specific Multi-Contaminant Gas Trap



CHROMSPEC's layered bed, multi-purpose trap is designed to remove oxygen, moisture, hydrocarbons and carbon dioxide from gas streams. This 750cc unit will take the place of three separate traps, simplifying installation and reducing the number of potential leak sites that contribute to the aspiration of contaminants into the system. End fittings are stainless steel Swagelok™ fitted with 40 micron stainless steel frits.

Specifications

Oxygen Reduction.....	<2ppb
Oxygen Capacity	1200mg
Moisture Reduction	<10ppb
Moisture Capacity.....	46gm
Hydrocarbon Reduction.....	<20ppb
Hydrocarbon Capacity	30gm
Carbon Dioxide Reduction	<1ppm
Carbon Dioxide Capacity	100gm

Ordering Information

Description	End Fitting	Product Number
CHROMSPEC™ 750cc Trap	1/8"	C36190
CHROMSPEC™ 750cc Trap	1/4"	C36192
Mounting Clip (2 recommended)		C36194

For more information, please contact our Technical Support Team at:

www.chromspec.com • tech@chromspec.com • 1-800-267-8103

