



CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

CHROMSPEC™ - Compressed Gas Purifiers Oxygen Moisture Trap



CHROMSPEC™ oxygen/moisture traps are ideal for use in gas chromatography and GC/MS. These traps utilize a one-piece spun aluminum body and Swagelok™ end fittings. Each fitting incorporates a 40 micron stainless steel frit for particulate control. A regeneration service is offered for the high capacity model.

Specifications

	100cc unit	750cc high capacity
Oxygen reduction	<2ppb	<2ppb
Moisture reduction	<10ppb	<10ppb
Oxygen capacity	450mg	3200mg
Moisture capacity	2gm	15gm
Maximum pressure	250psi	250psi

Ordering information

Description	End Fitting	Item Number
Trigon Oxygen/Moisture Trap	1/8"	C36060
Trigon Oxygen/Moisture Trap	1/4"	C36065
High Capacity Oxygen/Moisture Trap	1/8"	C36082
High Capacity Oxygen/Moisture Trap	1/4"	C36084
Mounting Clip for Oxygen/Moisture Trap (2 recommended)		C36078
Mounting Clip for High Capacity Trap (2 recommended)		C36194

1-800-267-8103 • sales@chromspec.com • tech@chromspec.com

