

# CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

# **ZPure™ Purifiers**

- Available with Brass or Stainless Steel Compression Fittings, 1/8" or 1/4"
- Stainless steel bodies with 10 micron stainless steel frits at both ends to protect against dust
- High-quality activated adsorbents for long filter life and efficient contaminant removal
- · Designed to reduce contaminants to low ppb levels
- All ZPure filters are helium leak checked and shipped filled with helium
- · Dual seal design protects against leaks for added safety
- Recommended flow rate is 2 L/min and maximum operating pressure is 13.8 bar (200 psi)





# **Stainless Steel Inline Purifiers**

#### ZPure™ Dual Purifier

This is a high capacity purifier for oxygen and water removal. It reduces oxygen and water to low ppb-levels from inert gases, He, Ar,  $N_2$ , and  $H_2$ . This purifier will improve baseline noise and sensitivity. The advanced adsorbent technology is exclusive to CRS.

# ZPure™ O, Purifier

This is a high capacity purifier for oxygen removal. It reduces oxygen from GC/MS carrier gas to low ppb-levels. The high performance design and advanced adsorbent technology, exclusive to CRS, provides outstanding efficiency for oxygen removal.



#### ZPure™ HC Purifier

This is a high capacity purifier for hydrocarbon removal C5 and heavier. It reduces hydrocarbons from inert gases, air, or hydrogen to low-ppb levels. It is also recommended for use with FID gas lines.

# ZPure™ H<sub>o</sub>O Purifier

This is a high capacity purifier for moisture removal. It reduces moisture from inert gases, air, or hydrogen to less than 20 ppb-levels.



# **Glass Indicating Purifiers**

## ZPure™ Glass Oxygen Purifier

This high-efficiency gas purifier is recommended as a downstream indicator for high capacity, non-indicating purifiers, such as the ZPure Dual or ZPure Oxygen. It reduces oxygen to low-ppb levels from chromatography carrier gas streams and provides a visual breakthrough indicator. The adsorbent materials are held in a sturdy glass tube for added safety.



#### **ZPure Glass Moisture Purifier**

This is an advanced-design gas purifier for applications requiring moisture removal with a visual breakthrough indicator. Up to 6 grams of water are reduced to low-ppb levels. The adsorbing materials, Molecular Sieve 5a and cobalt-free CoFree indicator, are held in a sturdy glass tube for added safety.



# ZPure Glass H,O/HC Purifier

This indicating purifier removes water and hydrocarbons C5 and heavier from inert gases or air. It contains a visual breakthrough indicator. The adsorbing materials, activated carbon, cobalt-free CoFree indicator, and Molecular Sieve 5a, are held in a sturdy glass tube for added safety.



## **ZPure Glass High Sensitivity Moisture Purifier**

This is an advanced-design gas purifier for applications requiring high-efficiency moisture removal with an extremely sensitive moisture indicator. It reduces water to low-ppb levels from gas streams and provides a visual indication of breakthrough when exposed to moisture levels of 1 ppm or more. The adsorbent materials are held in a sturdy glass tube for added safety.



## **ZPure Stainless Steel Inline Purifiers**

### **ZPure Glass Indicating Purifiers**

Description Brass Fittings	Item No.	Description Brass Fittings	Item No.
ZPure Dual, 1/8" fittings	202206-B	ZPure Indicating O <sub>2</sub> , 1/8" fittings	202236-B
ZPure Dual, 1/4" fittings	202207-B	ZPure Indicating O <sub>2</sub> , 1/4" fittings	202237-B
ZPure O <sub>2</sub> , 1/8" fittings	202216-B	ZPure Moisture, 1/8" fittings	202234-B
ZPure O <sub>2</sub> , 1/4" fittings	202217-B	ZPure Moisture, 1/4" fittings	202235-B
ZPure HC, 1/8" fittings	202232-B	ZPure H₂O/HC, 1/8" fittings	202256-B
ZPure HC, 1/4" fittings	202233-B	ZPure H <sub>2</sub> O/HC, 1/4" fittings	202257-B
ZPure H₂O, 1/8" fittings	202230-B	ZPure High Sensitivity Moisture, 1/8" fittings	202283-B
ZPure H₂O, 1/4" fittings	202231-B	ZPure High Sensitivity Moisture, 1/4" fittings	202284-B
Description Stainless Steel Fittings	Item No.	Description Stainless Steel Fittings	Item No.
Description Stainless Steel Fittings ZPure Dual, 1/8" fittings	Item No. 202206-SS	Description Stainless Steel Fittings  ZPure Indicating O <sub>2</sub> , 1/8" fittings	Item No. 202236-SS
approximate Continue Action Ac	I MANAGETT TOWNSON		TP. MARKET ME PERCONAL
ZPure Dual, 1/8" fittings	202206-SS	ZPure Indicating O <sub>2</sub> , 1/8" fittings	202236-SS
ZPure Dual, 1/8" fittings ZPure Dual, 1/4" fittings	202206-SS 202207-SS	ZPure Indicating O <sub>2</sub> , 1/8" fittings ZPure Indicating O <sub>2</sub> , 1/4" fittings	202236-SS 202237-SS
ZPure Dual, 1/8" fittings ZPure Dual, 1/4" fittings ZPure O <sub>2</sub> , 1/8" fittings	202206-SS 202207-SS 202216-SS	ZPure Indicating $O_2$ , 1/8" fittings  ZPure Indicating $O_2$ , 1/4" fittings  ZPure Moisture, 1/8" fittings	202236-SS 202237-SS 202234-SS
ZPure Dual, 1/8" fittings ZPure Dual, 1/4" fittings ZPure O <sub>2</sub> , 1/8" fittings ZPure O <sub>2</sub> , 1/4" fittings	202206-SS 202207-SS 202216-SS 202217-SS	ZPure Indicating O <sub>2</sub> , 1/8" fittings ZPure Indicating O <sub>2</sub> , 1/4" fittings ZPure Moisture, 1/8" fittings ZPure Moisture, 1/4" fittings	202236-SS 202237-SS 202234-SS 202235-SS
ZPure Dual, 1/8" fittings  ZPure Dual, 1/4" fittings  ZPure O <sub>2</sub> , 1/8" fittings  ZPure O <sub>2</sub> , 1/4" fittings  ZPure HC, 1/8" fittings	202206-SS 202207-SS 202216-SS 202217-SS 202232-SS	ZPure Indicating O <sub>2</sub> , 1/8" fittings  ZPure Indicating O <sub>2</sub> , 1/4" fittings  ZPure Moisture, 1/8" fittings  ZPure Moisture, 1/4" fittings  ZPure H <sub>2</sub> O/HC, 1/8" fittings	202236-SS 202237-SS 202234-SS 202235-SS 202256-SS

