

CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

CRS Premium Injection Port Septa



- · Center point guides the needle for easy penetration
- · Reduces needle bending
- · Precision molding assures accurate fit
- · Packaged in glass jars for high purity

Premium septa are precision-molded with a recess on the injection side, to guide the syringe needle to the same point with every injection. The CenterGuide™ design requires less

force for initial penetration for a smoother feel. Septum life may be enhanced

in some applications by minimizing tearing and coring.

The CenterGuide design are available in 5, 9, 11, 11.5, 12.7 and 17 mm sizes for Septa BTO[®], Advanced Green 3™ and Marathon™.



Septa BTO is bleed and temperature optimized for the most demanding GC and GC-MS applications: it is the finest GC septum we have tested. Septa BTO is formulated to extend low bleed and outstanding mechanical properties of premium GC septa to the highest-temperature applications. It retains remarkable softness and pierceability at high temperatures, with extremely low bleed, and has been optimized to reduce injection port adhesion.

Marathon Septa™

This is an advanced GC septum for autosampler use, with significantly longer life. The Marathon septum typically achieves 400 injections without failure when used with a rounded-tip (HP-style) needle and autosampler or needle guide. Now you can make extended autosampler runs without fear of sample loss caused by blow-back of leaking carrier gas. The Marathon septum is also ideal for standard manual injection GC and GC-MS.

Advanced Green 3 Septa™

The CRS Advanced Green 3 septum was created to combine significantly longer injection life, low bleed and low injection port adhesion. The result is a green septum you can use for all your critical analyses.

Description	Qty/Pk	Item No.	Description	Qty/Pk	Item No
BTO® CenterGuide Septa			Marathon™ CenterGuide Septa		
5 mm	50pk	298787	5 mm	25pk	239587
9 mm	50pk	298713	5 mm	50pk	239597
7/16" (11 mm)	50pk	298717	9 mm	25pk	239778
7/16" (11 mm)	100pk	298718	9 mm	50pk	239779
1/2" (12.7 mm)	50pk	298727	7/16" (11 mm)	25pk	239287
17 mm	50pk	298735	7/16" (11 mm)	50pk	239297
Shimadzu	50pk	298735	7/16" (11 mm)	100pk	239888
*For 17 mm injectors the maximum is 300°			11.5 mm	25pk	239787
AC2TM ContarCylide South			11.5 mm	50pk	239798
AG3™ CenterGuide Septa 5 mm 50pk 246525		1/2" (12.7 mm)	25pk	239387	
5 mm	CONTRACTOR AND	246525	1/2" (12.7 mm)	50pk	239397
9 mm	50pk		Shimadzu	25pk	239488
7/16" (11 mm)	50pk	246225	Shimadzu	50pk	239498
7/16" (11 mm)	100pk	246228		- 1997 C. 18 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	100001000000000000000000000000000000000
1/2" (12.7 mm)	50pk	246325			
17 mm	50pk	246424			
Shimadzu	50pk	246424	Green Control	S Chromatogra	