



CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

High Power Electronic Crimping Tool



11 mm crimper and decapper jaws



20 mm crimper and decapper jaws

Features:

- Higher-power version of electronic crimping tool with a fixed power supply and cord (No battery)
- Uses interchangeable jaw sets
- Stores separate programs for different caps and seals
- Strong enough for all steel and magnetic caps
- Optional base and mounting kit



Jaw Sets are easy to change. Simply insert the jaw set into the bushing at the bottom of the tool. Push up and twist until the set locks into position

To remove, push the outside button and rotate.



Display



Base with Mounting Kit

Description	Item No.
High Power Electronic Crimping Tool	C696AHPS0
Base with Mounting Kit	C695ABAS0
11 mm Crimper Jaw Set	C69311955
20 mm Crimper Jaw Set	C69320955
13 mm Crimper Jaw Set	C69313955
8 mm Crimper Jaw Set	C69308955

Description	Item No.
13 mm Jaw Set, for Flip Off Caps	C69313956
20 mm Jaw Set, for Flip Off Caps	C69320956
11 mm Decapper Jaw Set	C69311965
20 mm Decapper Jaw Set	C69320965
13 mm Decapper Jaw Set	C69313965



crs | Chromatography
Research Supplies

*For more information, contact our Technical Support Team at:
1-800-267-8103 • tech@chromspec.com*