



CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

High Power Electronic Crimping Tool



11 mm crimper and decapper jaws



20 mm crimper and decapper jaws

Features:

- Higher-power version of electronic crimping tool with a fixed power supply and cord (No battery)
- Uses interchangeable jaw sets
- Stores separate programs for different caps and seals
- Strong enough for all steel and magnetic caps
- Optional base and mounting kit



Jaw Sets are easy to change. Simply insert the jaw set into the bushing at the bottom of the tool. Push up and twist until the set locks into position

To remove, push the outside button and rotate.



Display



Base with Mounting Kit

Description	Item No.	Description	Item No.
High Power Electronic Crimping Tool	5AHPS0	13 mm Jaw Set, for Flip Off Caps	313956
Base with Mounting Kit	5ABAS0	20 mm Jaw Set, for Flip Off Caps	320956
11 mm Crimper Jaw Set	311955	11 mm Decapper Jaw Set	311965
20 mm Crimper Jaw Set	320955	20 mm Decapper Jaw Set	320965
13 mm Crimper Jaw Set	313955	13 mm Decapper Jaw Set	313965
8 mm Crimper Jaw Set	308955		



crs | Chromatography Research Supplies

*For more information, contact our Technical Support Team at:
1-800-267-8103 • tech@chromspec.com*

High Power Electronic Crimping Tool with External Control



11 mm crimper and decapper jaws



20 mm crimper and decapper jaws

Features:

- Higher-power programmable version of electronic crimping tool with a fixed power supply and cord (No battery)
- Serial (RS485) connection allows for external control of all crimper's functions with a programmable controller or PC
- Demonstration software for interface included
- Uses interchangeable jaw sets
- Stores separate programs for different caps and seals
- Strong enough for all steel and magnetic caps
- Optional base and mounting kit



Jaw Sets are easy to change. Simply insert the jaw set into the bushing at the bottom of the tool. Push up and twist until the set locks into position

To remove, push the outside button and rotate.



Display



Base with Mounting Kit

Description	Item No.
High Power Electronic Crimping Tool with External Control	5AHPR0
Base with Mounting Kit	5ABAS0
11 mm Crimper Jaw Set	311955
20 mm Crimper Jaw Set	320955
13 mm Crimper Jaw Set	313955
8 mm Crimper Jaw Set	308955

Description	Item No.
13 mm Jaw Set, for Flip Off Caps	313956
20 mm Jaw Set, for Flip Off Caps	320956
11 mm Decapper Jaw Set	311965
20 mm Decapper Jaw Set	320965
13 mm Decapper Jaw Set	313965



crs | Chromatography Research Supplies

*For more information, contact our Technical Support Team at:
1-800-267-8103 • tech@chromspec.com*