



CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

CHROMSPEC™ GC FERRULES

- Large inventory for quick delivery
- Low prices without sacrificing quality
- Guaranteed product performance
- Packaged in convenient dial packs
 - Easy dispensing
 - Minimize handling and contamination

CHROMSPEC™ ... The Name for Affordable Quality

Try CHROMSPEC™ Ferrules in your Lab!
Contact us at tech@chromspec.com
or call 1-800-267-8103 to request a sample.

CHROMSPEC™ Compact Ferrules for Agilent 5890 / 6890 / 6580 GCs

- High-purity, high density graphite and graphite/Polyimide®
- Pre-conditioned to eliminate out-gassing
- Smoother surface and cleaner edges than conventional ferrules



Product Number	Description	Equivalent to Agilent Product Number	Pkg Size	Qty 1-4 Price	Qty 5+ Price
C75510	Reducing Single Ferrules, 1/16" – 0.4 mm, Graphite	500-2114	10	\$30.60	\$25.60
C75502	Reducing Single Ferrules, 1/16" – 0.5 mm, Graphite	5080-8853	10	30.60	25.60
C75514	Reducing Single Ferrules, 1/16" – 0.8 mm, Graphite	500-2118	10	30.60	25.60
C75506	Reducing Single Ferrules, 1/16" – 1.0 mm, Graphite	5080-8773	10	30.60	25.60
C75872	Reducing Single Ferrules, 1/16" – 0.4 mm, Polyimide	5181-3322	10	50.40	45.40
C75874	Reducing Single Ferrules, 1/16" – 0.5 mm, Polyimide	5062-3513	10	50.40	45.40
C75876	Reducing Single Ferrules, 1/16" – 0.8 mm, Polyimide	5062-3511	10	50.40	45.40
C75522	Reducing Single Ferrules, 1/16" – 0.4 mm, Graphite/Polyimide	5181-3323	10	50.40	45.40
C75524	Reducing Single Ferrules, 1/16" – 0.5 mm, Graphite/Polyimide	5062-3514	10	50.40	45.40
C75526	Reducing Single Ferrules, 1/16" – 0.8 mm, Graphite/Polyimide	5062-3512	10	50.40	45.40

...2

CHROMSPEC™ Universal Graphite / Polyimide® Ferrules

- 15% graphite / 85% Polyimide®
- Will not fragment, can be reused many times and maintain their seal under vacuum
- Ideal for GC/MS applications as they will not contaminate the ion source with particles
- Temperature limit of 400°C

Product Number	Description	Equivalent to Agilent Product Number	Pkg Size	Qty 1-4 Price	Qty 5+ Price
C75529	Reducing Single Ferrules, 1/16" – 0.4 mm	5062-3508 8004-0212	10	\$55.00	\$50.00
C75242	Reducing Single Ferrules, 1/16" – 0.5 mm	5062-3506	10	55.00	50.00
C75244	Reducing Single Ferrules, 1/16" – 0.8 mm	5062-3538	10	55.00	50.00
C75530	Reducing Single Ferrules, 1/16" – 1.0 mm		10	55.00	50.00
C75106	Two Hole Ferrules, 1/16" – 0.4 mm	5062-3580	10	66.80	61.80
C75107	Two Hole Ferrules, 1/16" – 0.5mm	5062-3581	10	66.80	61.80
C75542	No Hole Single Ferrules, 1/16"	5181-3308	10	53.50	48.50

CHROMSPEC™ Universal Graphite Ferrules

- Produced from specially formulated 100% graphite
- Ideal for glass to metal connections
- Temperature range: subambient to 450°C

C75204	Ferrules, 1/4"	0100-1324	10	\$36.70	\$31.70
C75202	Ferrules, 1/8"	0100-1325	10	33.60	28.60
C75200	Ferrules, 1/16"	0100-1326	10	30.60	25.60
C75272	Ferrules, Reducing, 1/4" – 1/8"		10	36.70	31.70
C75270	Ferrules, Reducing, 1/4" – 1/16"		10	36.70	31.70
C75250	Ferrules, Reducing, 1/8" – 1/16"		10	33.60	28.60
C75230	Ferrules, Reducing, 1/16" – 0.4 mm		10	30.60	25.60
C75232	Ferrules, Reducing, 1/16" – 0.5 mm		10	30.60	25.60
C75234	Ferrules, Reducing, 1/16" – 0.8 mm		10	30.60	25.60
C75236	Ferrules, Reducing, 1/16" – 1.0 mm		10	30.60	25.60

CHROMSPEC™ Universal Polyimide® Ferrules

- One piece design
- Machined from 100% polyimide
- Temperature limit of 350°C

C76040	Ferrules, Straight, 1/4" – 1/4"	5080-8774	10	\$77.50	\$72.50
C75830	Ferrules, Straight, 1/8" – 1/8"	G6600-80023	10	64.20	59.20
C75850	Ferrules, Straight, 1/16" – 1/16"	0100-1333	10	55.00	50.00
C75886	Ferrules, Reducing, 1/8" – 1/16"	0100-01342	10	64.20	59.20
C76048	Ferrules, Reducing, 1/4" – 1/8"		10	77.50	72.50
C76046	Ferrules, Reducing, 1/4" – 1/16"		10	77.50	72.50
C75854	Ferrules, Reducing, 1/16" – 0.4mm		10	55.00	50.00
C75856	Ferrules, Reducing, 1/16" – 0.5mm		10	55.00	50.00
C75858	Ferrules, Reducing, 1/16" – 0.8mm		10	55.00	50.00
C75860	Ferrules, Reducing, 1/16" – 1.0mm		10	55.00	50.00

More ferrule sizes are available – please inquire.

Note: Vespel® is the trademark name for a range of polyimide products manufactured by Dupont.
Vespel® is a registered trademark of E.I. du Pont de Nemours and Company.

Pricing subject to change without notice.

