

CHROMSPEC™ ULB Septa

CHROMSPEC™ ULB (Ultra Low Bleed) Septa for Open Top Caps

Storage vial septa should provide easy access to the contained material, while maintaining sample integrity. For over 20 years, our customers have used Tegrabond™ septa for these applications. They have provided excellent value and performance. Today, we are proud to introduce an even better product: ULB septa.



| | Good: Original Tegrabond™ Septa | Better: ULB Septa |
|---------------------|--|--|
| Clean Silicone | Clean silicone for smooth baselines in standard detectors | Ultra-clean silicone for smooth baselines in ultra-sensitive MS detectors |
| Pierceable Silicone | Soft 45-50 Durometer Silicone: industry standard for easy piercing | Silicone tested with autosampler needles for easier piercing than 45-50 Durometer |
| PTFE | 0.010" skived PTFE directly bonded to silicone | 0.005" cast PTFE: same chemical resistance, lower raw material costs, easier to pierce |
| Manufacturing | North American made for over 20 years by the same manufacturer | North American made by the manufacturer of Tegrabond with 20 more years of R&D |
| Low Prices | Good | Even better: improved production methods have lowered our costs! |

Regardless of which storage septa you currently use, we invite you to upgrade the performance AND save money!


Contact our Technical Team to request a free sample.

CHROMSPEC™ ULB (Ultra Low Bleed) Septa

| Product Number | Description | Qty* |
|---------------------------|--|------|
| Original Tegrabond | ULB Equivalent | |
| C665107C | C66U0713C Clear PTFE / White ULB Silicone, 7mm x 1.3mm thick | 100 |
| C665108M | C66U0813M Clear PTFE / White ULB Silicone, 8mm x 1.3mm thick | 1000 |
| C669110C | C66U1030C Clear PTFE / White ULB Silicone, 10mm x 3mm thick | 100 |
| C665112C | C66U1213C Clear PTFE / White ULB Silicone, 12mm x 1.3mm thick | 100 |
| C669112C | C66U1230C Clear PTFE / White ULB Silicone, 12mm x 3mm thick | 100 |
| C665113M | C66U1313M Clear PTFE / White ULB Silicone, 13mm x 1.3mm thick | 1000 |
| C669113C | C66U1330C Clear PTFE / White ULB Silicone, 13mm x 3mm thick | 100 |
| C665116C | C66U1613C Clear PTFE / White ULB Silicone, 16mm x 1.3mm thick | 100 |
| C669116C | C66U1630C Clear PTFE / White ULB Silicone, 16mm x 3mm thick | 100 |
| C665118M | C66U1813M Clear PTFE / White ULB Silicone, 18mm x 1.3mm thick | 1000 |
| C669118M | C66U1830M Clear PTFE / White ULB Silicone, 18mm x 3mm thick | 1000 |
| C669119C | C66U1930C Clear PTFE / White ULB Silicone, 19mm x 3mm thick | 100 |
| C665122M | C66U2013M Clear PTFE / White ULB Silicone, 20mm x 1.3mm thick | 1000 |
| C669120M or C8812120M | C66U2030M Clear PTFE / White ULB Silicone, 20mm x 3mm thick | 1000 |
| C6613520M | C66U2A30M Clear PTFE / White ULB Silicone, 20.32mm x 3mm thick | 1000 |
| C669122M or C8812122M | C66U2230M Clear PTFE / White ULB Silicone, 22mm x 3mm thick | 1000 |
| C669126C | C66U2630C Clear PTFE / White ULB Silicone, 26mm x 3mm thick | 100 |
| C8812130C | C66U3030C Clear PTFE / White ULB Silicone, 30mm x 3mm thick | 100 |

*Packs of 100 and 1000 are available for all products.
Also available in diameters ranging from 8 to 30mm.

CHROMSPEC™ ... The Name For Affordable Quality



**CHROMATOGRAPHIC
SPECIALTIES INC.**

1-800-267-8103 • www.chromspec.com • tech@chromspec.com

