



# CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

Prices up to  
**40%**  
lower

## CHROMSPEC™ Replacement Deuterium Lamps

All CHROMSPEC™ replacement lamps are fully guaranteed and manufactured to the same standards as OEM lamps, but at a greatly reduced cost.

- **Save up to 40%** compared to the OEM pricing
- Fully guaranteed with a 1,000 or 2,000 hour pro-rated warranty
- Each lamp is carefully designed to match the OEM product specifications



### Bright Ideas for Lamps...

#### Handle with care

Deuterium lamps produce a plasma ball approximately 0.5mm in diameter. Disturbing the filaments by even tenths of a millimeter will wreck the alignment of your bulb.

#### Shelf-life is not warranty life

Lamps stored away from jostling and bumping last a long time. If you keep a spare bulb for a year before using it, we will still be there to support you.

#### Clean lamps are happy lamps

Wear cotton or thin rubber gloves to protect your lamp from fingerprints. Wipe all lamp bulbs down with IPA before installing.

#### Track your operating time

CHROMSPEC™ brand lamps are produced in accordance with RoHS rules. On-off times in a lab book or calculating bulb time based on working hours per/day and number of days are accurate for tracking bulb life.

#### Warranties

A 1,000 or a 2,000 hour lamp doesn't need to be replaced when its warranty runs out. A 1,000 hour warranty says we expect a lamp to last *at least* 1,000 hours.

#### Lamps have turn-offs too

If your instrument will be left unused for more than a few hours, turn the lamp off to preserve its useful life. If it will be less than two hours until your next run, leave it on. The thermal cycling and warm-up time will offset the advantages of turning the lamp off.

Product Number	Model	OEM Part Number	Type
<b>For Agilent</b>			
C55D401AGLL*	1100 VWD (G1314A)	G1314-60100	D2
C55D401AGLL3*	1100 VWD (G1314A), Pkg. 3	G1314-60100	D2
C55D501AG	1100 DAD (G1315), 1100 MWD (G1365)	2140-0590 / 2140-0605	D2
C55D501AGLL*	1100 DAD, 1100 MWD - Long Life	2140-0813 / 5181-1530	D2
C55D501AGLL3*	1100 DAD, 1100 MWD - Long Life, Pkg. 3	2140-0813 / 5181-1530	D2
C55D601AG	8452A	08452-60104	D2
<b>For Waters</b>			
C55D201WA	486LC Module 1	WAT080678	D2
C55D401WALL*	966 2996 DAD	WAT057760	D2
C55D401WALL3*	966 2996 DAD, Pkg. 3	WAT057760	D2
C55D501WALL*	Alliance 2487, 2488	WAS081142	D2
C55X051SWA	470, 474, 2475	WAT052665	Xe

\* Double life lamp with 2,000 hour warranty.

Select Tungsten and Xenon lamps are also available. Many lamps are also available in packages of 6.

.../cont'd

**For more information, or for tips on how to install your lamp,  
visit [www.chromspec.com/lamps](http://www.chromspec.com/lamps),  
or contact our Technical Support Team at  
1-800-267-8103 • [tech@chromspec.com](mailto:tech@chromspec.com)**



## CHROMSPEC™ Replacement Deuterium Lamps (cont'd)

Product Number	Model	OEM Part Number	Type
<b>For ABI / Perkin-Elmer</b>			
C55D001AB	757, 759, 783A, 785A, 1000S, FS980, 120A, 130A	2900-0489	D2
<b>For Beckman</b>			
C55D001BE	163	22 947029	D2
C55D301BE	166 (Pre-Aligned)	239372	D2
C55D501BE	168 (Pre-Aligned)	538711	D2
<b>For Gilson</b>			
C55D101GI	115, 116, 117, 118, 119, 151, 152, 155, 156	100326	D2
<b>For Hitachi</b>			
C55D001ME	L & U Series	AN0-0424AS	D2
<b>For Shimadzu</b>			
C55D101SHLL*	SPD10A, 10AS, 10AV, 10AVP, M10AVP	228-34016-02	D2
C55D301SHLL*	LC2010	228-37401-00	D2
<b>For TSP</b>			
C55D201TS	SP100/200, UV100/150/200/1000/2000/3000, Focus, CE	9551-0023	D2
<b>For Varian</b>			
C55D301VALL*	Vista 5000, 5500, UV50/100/200, 9050, Prostar 310	0391615691	D2
C55D401VA	Prostar 330 DAD	0393570502	D2
C55D501VALL*	Prostar 325, 335	110715400	D2

\* Double life lamp with 2,000 hour warranty.  
 Select Tungsten and Xenon lamps are also available.  
 Many lamps are also available in packages of 6.

We offer lamps for the following instruments:

Agilent	Ciba Corning	GBC	LDC	SSI
Alltech	Coleman	Gilford	Linear	Shimadzu
Beckman	Cole-Parmer	Gilson	LKB	Sigma
Biorad	Desaga	Hitachi (Merck Hitachi)	Michrom BioResources	Spectra Physics
Biotronic	Dionex (Gynkotek)	IBM	Micrometrics	TSP
Carlo-Erba	DuPont	Isco	Perkin Elmer	Unicam
Cary	EM Science	Jasco	Rainin	Varian
Cecil Instruments	Fisher	Knauer	Saitron	Waters

