
CHROMATOGRAPHIC SPECIALTIES INC.

PARTS FOR AGILENT 1100 AND 1200 HPLCs

Check Valves (with Ceramic Ball and Seat)

EXPERT ADVICE

Product Number	Description	Qty	Equivalent to Agilent Number
M275601330 I	Housing & Cartridge, OPTI-MAX [®] Outlet	1	G1311-60012 G1312-60012G G1312-60067
	Replacement Cartridges, OPTI-MAX®,1/8" Ball Replacement Cartridge for Active Inlet Valve, 600 bar	2 1	5062-8562
	DPTI-MAX [®] Outlet Cartridge for UHPLC Pressures (20,000 psi+), I/16" Ceramic Ball & Seat, SS Cartridge, Agilent 1200, 1220, &1260	1	G1312-60067
	DPTI-MAX [®] Outlet Cartridge for UHPLC Pressures (20,000 psi+), I/16" Double Ceramic Ball & Seat, SS Replacement Cartridge, Agilent 1200, 1220, & 1260	1	G1312-60067
M274705033 I	Passive Inlet Check Valve, Agilent 1220/1260	1	G1312-60066
Faster changes Free turn housir istons and Sea	•		
M273601570	Piston Seal, UHMW-PE, OPTI-SEAL®	1	0905-1420
M273601571	Piston Seal, UHMW-PE, OPTI-SEAL®	10	0905-1420
M273601572	Piston Seal, ITB [™] PTFE	1	5063-6589
M273601573	Piston Seal, ITB [™] PTFE	10	5063-6589
M273400658	Piston, Sapphire	1	5063-6586
M273602364	Wash Seal, UHMW-PE, OPTI-SEAL®	1	0905-1175
M273602851	Gasket, Wash Seal	6	5062-2484
urge and Mixi	ng Valve		
C56018AG	Frits, PTFE, Purge valve	5	01018-22707
M273102853	Seal Cap, 2-in-1 (Replaces two part Agilent cap with a single part)	1	5067-4728
UAG131167701T	Valve, Multi-Channel Gradient	1	G1311-67701
utosampler Pa	arts		
M273103212	Rotor Seal Assembly, Vespel	1	0100-1853
RK26424	Stator Face Assembly	1	0100-1851
C56028AG	Needle	1	G1313-87201
M273103329	Needle Seat Assembly with convenient 1 Piece Ferrule	1	G1329-87017
M273103146	Needle Seat Assembly with 2 Piece Ferrule	1	G1329-87017
etector Lamps	5		
C55D401AGLL	CHROMSPEC [™] Lamp, Longlife D ₂ for VWD*	1	G1314-60100
C55D501AGLL	CHROMSPEC [™] Lamp, Longlife D ₂ for DAD & MWD*	1	2140-0813
C55D701AGLL	Deuterium Lamp, CHROMSPEC, Longlife, For Agilent 1100/1200 DAD (To Fit G1315/G1365 C&D, G7115A, G7165	1 A Series	2140-0820
C55D801AGLL	Deuterium Lamp, CHROMSPEC, Longlife,	1	5190-0917
	For Agilent 1290 DAD		
Save	Lamps are also available in packs of 6		
More	*CHROMSPEC [™] brand Deuterium lamps are compatible with Agilent 1200 s an RFID chip. Installation date and serial numbers need to be entered man		
	Want to keep your 1100 running after Agi	lent's	s end of servi











M273102853 on Agilent Purge Valve Base



C55D401AGLL



1-800-267-8103 • www.chromspec.com • sales@chromspec.com

UNIVERSAL HPLC AND UHPLC FITTINGS

HPLC Tubing and Fittings For systems running at <6,000 psi. Also ideal for UHPLC systems after the column!

1/16" PEEK® Tubing — Dual Layer Colour Coded, Straight

- Undyed core & dyed shell = permanent, easy to see size information
- · No leaching concerns



· Same pressure rating, same great price as our dyed, undyed, striped tubing

Product Number	ID mm (in)	Colour Code	Pressure Limit (psi)	Qty
C446007D	0.13 (0.005)	Red	6700	10 ft
C446008D	0.18 (0.007)	Yellow	6400	10 ft
C446009D	0.25 (0.010)	Blue	6100	10 ft

More diameters and lengths are available.

One-Piece PEEK® Fingertight Fittings

· Easy to use, no separate ferrule to fumble with



 Inexpensive, seal to 5,000 psi · Contact us for colour coded part numbers

Product Number	Description	Qty
C44550205	Fingertight Fitting, Natural	5
PEEK® Nuts And Double Ferrules		

Nuts And Double Ferrules

Separate nut and ferrule for better seal

• 6,000 psi h	and or wrench tight	
C4455805	PEEK [®] Handtight Nut	5
C4455105	PEEK [®] Hex-Head Nut	5
C4450045	PEEK [®] Double-Ferrule	5

Optimize Opti-Solv® In-Line Filters

· In-line filter gives easy, inexpensive protection for standard analytical columns



- · Use Gorilla-Grip for easy hand-tight operation; wrench-tight filter is narrower for confined spaces
- · Mini-filter uses depth filtration for reliable low-volume / fast HPLC work

	OPTI-SOLV [®] In-Line Filter, Gorilla-grip* OPTI-SOLV [®] In-Line Filter, wrench-tight*	1
	OPTI-SOLV [®] In-Line Filter,	10
	Replacement Frits, 2µm	
M100400097	OPTI-SOLV [®] Mini-Filter, 2µm, SS frit	5

* Requires Frits

Mobile Phase Filters



· Your first line of defense against particles or dust entering through your mobile phase bottles

C4446772	Filter, Last Drop Economy, 2µm,	1
	22mm OD, SS	
C4490000530	Filter, Last Drop with Tripod Tubing	1
	Connector, 2µm, SS	
RK26431	Bluestem Glass Solvent Filter, 15µm	1
RK26392	PTFE Adapter for Bluestem Filter	4
M100400071	Solvent Reservoir Filter, Hastelloy C,	1
	for 1/16"ID Tubing, 2μ m	

UHPLC Tubing and Fittings

Recommended for systems operating at >6,000 psi (400 bar). Also suitable for demanding HPLC applications.

1/16" 316 Stainless Steel Capillary Tubing

- Bulk Tubing
 - High quality, low cost per foot
 - Cuts depend on user skill



- Pre-cut Tubing
 - All cuts smooth, square, free of metal fragments
 - Convenient sizes, more economical than custom cut tubing

Length	0.005" ID (17,250psi*)	0.010" ID (15,750psi*)	0.020" ID (12,750psi*)
5 cm	VT5C5D	VT5C10D	VT5C20D
10 cm	VT10C5D	VT10C10D	VT10C20D
20 cm	VT20C5D	VT20C10D	VT20C20D
30 cm	VT30C5D	VT30C10D	VT30C20D
10ft	VTSS105S10F	VTSS110S10F	VTSS120S10F

* Calculated ASME maximum working pressure for 316 grade tubing

Optimize EXP[™] Reusable Fittings



- · Handtight seal to 8,700 psi, wrench-tight to 20,000 psi
- · No permanent swage for easy repeat connections

Product Number	Description	Qty
M152003888	EXP [™] Hand-Tight Nut and Titanium Hybrid Ferrule	1
M152003833	EXP [™] 10-32 Hex Head Male Nut, Stainless Steel, Standard Length	1
M152403831 M152403832	EXP [™] Titanium Hybrid Ferrule EXP [™] Titanium Hybrid Ferrule	1 10

VICI SS Nuts & Ferrules



· Ideal for making repeated or long term connections to one location

C34200 C34202	Standard Zero Volume Nut, 1/16", SS Ferrule, 1/16", 303SS	10 10
00.202		

We also offer:

- Tubing Cutters
- Membrane Filters
- Degassers
- Vials and Septa



1-800-267-8103 • www.chromspec.com • sales@chromspec.com