

CHROMATOGRAPHIC SPECIALTIES INC.

COST SAVING PARTS FOR WATERS ALLIANCE HPLCs

Lower cost than the OEM parts, with no compromise in quality

Check Valves

OPTI-MAX® cartridges have the same external dimensions as Waters Performance PLUS™ cartridges, and can be used in the same holders.

Product Number	Description	Qty	Equivalent to Waters Number
M394602187	Alliance Conversion Kit OPTI-MAX®	1	WAT270941, 700002332
M104602004	Replacement Cartridges, OPTI-MAX®, 1/16"	2	700002332





UHMW-PE OPTI-SEAL® Piston Seals Save Money for Reversed Phase Methods

- · UHMW-PE lasts longer than PTFE in water and polar organic solvents
- · Less filter maintenance: UHMW-PE sheds strands that sit on frits rather than plug them

Pistons and Seals

M393602596	Plunger Seal, UHMW-PE, OPTI-SEAL®	1	WAT270938
M393602603	Plunger Seal, UHMW-PE, OPTI-SEAL®	2	WAT270938
M393602597	Plunger Seal, UHMW-PE, OPTI-SEAL®	10	WAT270938
M393602858	Tube Seals, Seal Wash	4	WAT270940
M393602860	Piston Seal, Seal Wash	2	WAT271018
M393402180	Sapphire Plunger	1	WAT270959
Filters M390403524	Filter Insert Assembly, 2μ m	1	WAT088084

Autosampler Parts

M393001106	Injector Wash Frit	3	WAT005106
M393001114	PTFE Washer, HPLC	6	WAT005135
H80024	Syringe, 250µL, 1725 for WISP	1	WAT073109
H80020	Syringe, 25µL, 1702 for WISP	1	-

Detector Lamps

•				
C55D401WALL	CHROMSPEC™ D₂ Lamp, LongLife, For 996, 2996	1	WAT052586	
	CHROMSPEC™ D, Lamp, LongLife, For 2487, 2488	1	WAS081142	
C55X051WA	CHROMSPEC™ Xe Lamp, For 470, 474 & 2475	1	WAT047018	
C55D801WALL	CHROMSPEC™ Deuterium Lamp, LongLife, For Acquity	1	201000281	







Lamps are also available in packs of 6

Parts for Acquity® UHPLC Systems

Individual parts also available, please contact our Technical Support Team tech@chromspec.com.

Product Number	Description	Comparable to OEM#
C56022WA C56023WA C56078WA	ACQUITY® H-Class QSM Performance Maintenance Kit ACQUITY® I2V BSM Performance Maintenance Kit ACQUITY® BSM Performance Maintenance Kit	201000233 201000197 201000173



UNIVERSAL HPLC AND UHPLC FITTINGS

HPLC Tubing and Fittings
For systems running at <6,000 psi. Also ideal for UHPLC systems after the column!

1/16" PEEK® Tubing — Dual Layer Colour Coded, Straight

- Undyed core & dyed shell = permanent, easy to see size information
- · No leaching concerns
- · Same pressure rating, same great price as our dyed, undyed, striped tubing



Product Number	ID mm (in)	Colour Code	Pressure Limit (psi)	Qty
C446007D	0.13 (0.005)	Red	6700	10 ft
C446008D	0.18 (0.007)	Yellow	6400	10 ft
C446009D	0.25 (0.010)	Blue	6100	10 ft

More diameters and lengths are available.

One-Piece PEEK® Fingertight Fittings

- · Easy to use, no separate ferrule to fumble with
- Inexpensive, seal to 5,000 psi

•	is for colour coded part number	s
Product Number	Description	Qty

C44550205 Fingertight Fitting, Natural PEEK® Nuts And Double Ferrules

- · Separate nut and ferrule for better seal
- 6,000 psi hand or wrench tight

C4455805	PEEK® Handtight Nut	5
C4455105	PEEK® Hex-Head Nut	5
C4450045	PEEK® Double-Ferrule	5

Optimize Opti-Solv® In-Line Filters

- · In-line filter gives easy, inexpensive protection for standard analytical columns
- · Use Gorilla-Grip for easy hand-tight operation; wrench-tight filter is narrower for confined spaces
- · Mini-filter uses depth filtration for reliable low-volume / fast HPLC work

	OPTI-SOLV® In-Line Filter, Gorilla-grip* OPTI-SOLV® In-Line Filter, wrench-tight*	1 1
M100400103	OPTI-SOLV® In-Line Filter,	10
	Replacement Frits, 2µm	
M100400097	OPTI-SOLV® Mini-Filter, 2µm, SS frit	5

^{*} Requires Frits

Mobile Phase Filters

· Your first line of defense against particles or dust entering through your mobile phase bottles



_				
	C4446772	Filter, Last Drop Economy, 2µm,	1	
		22mm OD, SS		
	C4490000530	Filter, Last Drop with Tripod Tubing	1	
		Connector, 2µm, SS		
	RK26431	Bluestem Glass Solvent Filter, 15µm	1	
	RK26392	PTFE Adapter for Bluestem Filter	4	
	M100400071	Solvent Reservoir Filter, Hastelloy C,	1	
		for 1/16"ID Tubing 2um		

UHPLC Tubing and Fittings

Recommended for systems operating at >6,000 psi (400 bar). Also suitable for demanding HPLC applications.

1/16" 316 Stainless Steel Capillary Tubing

- · Bulk Tubing
 - High quality, low cost per foot
 - Cuts depend on user skill
- · Pre-cut Tubing
 - All cuts smooth, square, free of metal fragments
 - Convenient sizes, more economical than custom cut tubing

Length	0.005" ID (17,250psi*)	0.010" ID (15,750psi*)	0.020" ID (12,750psi*)	(Pkg. 2)
5 cm	VT5C5D	VT5C10D	VT5C20D	
10 cm	VT10C5D	VT10C10D	VT10C20D	
20 cm	VT20C5D	VT20C10D	VT20C20D	
30 cm	VT30C5D	VT30C10D	VT30C20D	
				Each
10ft	VTSS105S10F	VTSS110S10F	VTSS120S10F	Luon
* Calculator	ASME maximum w	orking proceure for	316 grade tubing	

* Calculated ASME maximum working pressure for 316 grade tubing

Optimize EXP™ Reusable Fittings

- · Handtight seal to 8,700 psi, wrench-tight to 20,000 psi
- · No permanent swage for easy repeat connections

Product Number Description	Qty
M152003888 EXP™ Hand-Tight Nut and Titanium	1
Hybrid Ferrule M152003833 EXP™ 10-32 Hex Head Male Nut,	1
Stainless Steel, Standard Length	
M152403831 EXP™ Titanium Hybrid Ferrule	1
M152403832 EXP™ Titanium Hybrid Ferrule	10

VICI SS Nuts & Ferrules

· Ideal for making repeated or long term connections to one location



C34200	Standard Zero Volume Nut, 1/16", SS	10
C34202	Ferrule, 1/16", 303SS	10

