



Vials =

Solid Phase Microextraction (SPME) Specialty Vials and Seals

Thicker Lip Vials

Thin Seals



The new Finneran SPME Vials with a thicker lip are ideal for most SPME applications. The thicker lip vial creates a tighter seal on the vial. The Finneran SPME Vials are manufactured from high quality Type I borosilicate glass. This new thicker lip vial is available in 10mL and 20mL capacities. To complete the package – J.G. Finneran offers the 20mm gold magnetic seal, 8mm hole with a thin 0.60" PTFE/silicone seal. The thin seal is a perfect choice for SPME applications**. The thin seal extends the life of the needle and is available in blue or clear.

Shelf Pack – 100 pieces; Case Pack 1000 pieces

Cat. No.	Description

31020SP-2346 10mL SPME Vial, clear, radius bottom, 23x46mm, 20mm flat top crimp, thicker lip finish
32020SP-2375 20mL SPME Vial, clear, radius bottom, 23x75mm, 20mm flat


top crimp, thicker lip finish

Seals

Case Pack – 100 pieces

Cat. No.

Description



5145TP8-20SP 20mm Gold Magnetic Seal, 8mm hole, thin .060" blue PTFE/white silicone septa

5156TP8-20SP 20mm Gold Magnetic Seal, 8mm hole, thin .060" clear PTFE/white silicone septa

**Suitable for Gerstel MPS, CTC and Leap autosamplers with SPME options