

11mm Aluminum Seals & Snap Top CapsTM with Metal O-Rings







Our 11mm aluminum seals and snap caps are now available assembled with a septa and a round metal o-ring for use in any application where CTC Analytics PAL* instruments must transport the vial to another location using the magnets found on these instruments.

Aluminum Seals with Metal O-Rings

Cat. No.	Description	Case
C735140MR11M	11mm Silver Aluminum Seal, PTFE/Butyl Lined with Metal O-Ring	1000
C735150MR11M	11mm Silver Aluminum Seal, PTFE/Silicone Lined with Metal O-Ring	1000
C735160MR11M	11mm Silver Aluminum Seal, PTFE/Silicone/PTFE Lined with Metal O-Ring	1000

Snap Top Caps™ with Metal O-Rings

Cat. No.	Description	Case
C735810MR11M	11mm Clear Snap Cap, PTFE Lined with Metal O-Ring	1000
C735840MR11M	11mm Clear Snap Cap, PTFE/Butyl Rubber Lined with Metal O-Ring	1000
C735850MR11M	11mm Clear Snap Cap, PTFE/Silicone Lined with Metal O-Ring	1000
C735860MR11M	11mm Clear Snap Cap, PTFE/Silicone/PTFE Lined with Metal O-Ring	1000
C735870MR11M	11mm Clear Snap Cap, PTFE/Silicone Lined, with Slit and Metal O-Ring	1000









Engineered and Manufactured in compliance with ISO 9001:2008 certified guidelines



1-800-267-8103 • www.chromspec.com • tech@chromspec.com