

Solutions Roadmap

Shimadzu® Nexera Series

Produced and Illustrated by Sciencix Engineers






Shimadzu[®] Nexera Series

Autosampler Models: SIL-30A/C, SIL30A/C/M/P				
Product Description	Sciencix Part #	Comparable OEM Part #		
Plunger Seal	CTS-10694	228-35145-00		
HPV Rotor, SIL-30AC, SIL-30ACMP	11-4179	228-52139-00		
Pt Coated Needle, 30 Series	CTS-21261	228-41024-95		
Uncoate Needle Kit	CTS-21264	228-41024-96		
PEEK™ Needle Seal, SIL-30AC, SIL-30ACMP	11-4048	228-53178-91		
Needle Seal, SIL-30AC, SIL-30ACMP	11-4052	228-52253-00		
LPV Rotor, SIL-30AC, SIL-30ACMP	11-4180	225-51922-00		
Sapphire Plunger	CTS-11434	228-35010-91		
Rinse Port Cap (w/hole), 10/pk	CTS-21690	228-48331-92		
Needle Seal Kit, SIL- 30AC	CTS-21422	228-52401-92		
Needle Seal, (For pH 1- 14, for operation up to 130 Mpa)	CTS-21089	228-57729-41		



SIL-30A/C/M/P, LC-30A/D/S/F



Pump Models: LC-30A/D, LC-30A/D/S/F		
Product Description	Sciencix Part #	Comparable OEM Part #
Suction Filter Assembly, SS Improved	CTS-21352	228-45708-91
Suction Filter, Stainless Improved	CTS-A11316	228-45707-91
Inlet Check Valve	11-4060	228-52964-95
Plunger Seal and Back Up Ring	11-4126	228-52711-93
LC-30AD Plunger Holder Assembly, UHP2 (w/Diaphragm)	11-4073	228-52069-94
PTFE Diaphragm, LC- 30AD/i-Series, 2/pk	11-4122	228-55272-41
SS Line Filter	11-4129	228-35871-99
Outlet Check Valve	11-4063	228-53334-96
PTFE Diaphragm, 2/PK	11-2891	228-32784-91
Gold Plunger Seal	11-4133	228-32628-91
Inlet Check Valve	CTS-21145	228-48249-96
Outlet Check Valve	CTS-A11187	228-45705-91
Line Filter Assembly	CTS-A11421	228-35871-96
Sapphire Plunger	CTS-11297	228-35601-93
Sapphire Plunger w/Holder	11-3975	228-35281-95



Shimadzu® Nexera Series


Comparable

OEM Part #

228-34016-02

228-35097-95

228-34016-00



Disclaimer: Sciencix, Inc. is not affiliated with or an authorized distributor/official partner of Original Equipment Manufacturer (OEM) parts. All OEM part names, numbers, brands, registered trademarks are the sole property of their respective owners in the United States or any other country. All OEM part names, numbers, and company names used in this catalog are for cross reference purposes only. Use of these names and brands does not imply endorsement. Statements made by Sciencix regarding the quality and performance of our products are our own and should not be construed as an indication of OEM acceptance or approval. © 2021 Sciencix, Inc.



Shimadzu[®] Nexera Series PM Kit

DESCRIPTION Maintenance Kit, LC-30AD	Sciencix Part # CTS-21494	Comparable OEM Part # 228-53265-44
Needle Seat Kit, SIL-30AC	CTS-21422	228-52401-92, Old # 228-52228-00