

Solutions Roadmap Waters® Alliance® Series

Produced and Illustrated by Sciencix Engineers



Sildin 0 • CE

01)

2

Display Contrast

90

Step Fletz

9

0

4

Enit


1-800-267-8103 • www.chromspec.com • tech@chromspec.com

© 2021 Sciencix, Inc.



Waters[®] Alliance[®] Series

Solvent Manager: 2690/D, 2695/D, 2790			
Product Description	Sciencix Part #	Comparable OEM Part #	
Alliance® Head Plunger Seal Kit, 2/pk	CTS-10566	WAT270938, (2 of WAT270789)	
Seal Wash Plunger Seal, Alliance® 2/pk	CTS-10560	WAT271018, (2 of WAT270160)	
Sapphire Plunger, Alliance®	CTS-10476	WAT271067, WAT270959, Old# WAT270488	
Seal Wash Face Seal	CTS-10532	WAT271017	
Was Tube Seal Kit, 4/pk	CTS-10597	WAT270940 (4 of WAT270668)	
Face Seals, 4/pk	CTS-10989	WAT270939, (4 of WAT270163)	
Check Valve Cartridge, Ceramic Ball & Seat, 2/pk	11-2401	700002399	
Sparge Diffuser	CTS-10984	WAT007272	
Cartridge Check Valve Housing, 2/pk	CTS-A11431-02	700002332	
Check Valve Cartrdige, 2/pk	CTS-A094	WAT270941	
Alliance® Head Plunger Seal Kit, Black, 2/pk	CTS-10711	WAT271066	
Solvent Reservoir filter, 10µm	CTS-5531	WAT025531, 700003615	
Check Valve Cartridge, 2/pk	CTS-10600	700000254, (2 of 289000144)	
O-ring, TFE	CTS-1010T	WAT097387	
HPMV Gold Seals, 2/pk	CTS-10692	WAT045454	
Cartridge Check Vavle System Kit	CTS-21011	700000253	



Sample Manager: 2690/D, 2695/D, 2790				
Product Description	Sciencix Part #	Comparable OEM Part #		
2µm Precolumn Filter Insert	CTS-10919	WAT088084		
250μL Syringe Tip	CTS-3195	WAT073195		
250μL Syringe	CTS-3109	WAT073109		
Syringe Seal, 250µL/2500µL	CTS-3209	WAT077347		
Precolumn Filter Assembly	CTS-10527	WAT035190, Old # WAT084561		
2µm Precolumn Filter Insert, 5/pk	CTS-10919-5	WAT005139		
Needle Wash Frit	CTS-2324	WAT072324		
Performance Autosampler Needle	CTS-7550	700001247		
HPMV 1 &2 Rebuild Kit w/Spool	CTS-10761	WAT045424		
Seal Pack Rebuild Kit	CTS-10839	WAT271019		
Needle Wash Frit, 5/pk	CTS-10788-5	700001318, (5 of 415000230)		
Seal Pack Rebuild Kit	CTS-21707	700011783		
Waste Valve, V3 (2 Way Solenoid Valve)	CTS-10770	700001229		
TFE 3-Way Valve 4 (V4)	CTS-10285	WAT045407, WAT078800		
Seal Pack Assembly Kit	CTS-11222	700002791		

Disclaimer: Sciencix, Inc. is not affiliated with or an authorized distributor/official partner of Original Equipment Manufacturer (OEM) parts. All OEM part names, numbers, brands, registered trademarks or trademarks are the sole property of their respective owners in the United States or any other country. All OEM part names, numbers, and company names used in this catalog are for cross reference purposes only. Use of these names and brands does not imply endorsement. Statements made by Sciencix regarding the quality and performance of our products are our own and should not be construed as an indication of OEM acceptance or approval. © **2021 Sciencix, Inc.**

2690/D, 2695/D

2790



Waters[®] Alliance[®] Series

Detector Models: 2	487/2488, 996/29	96
Product Description	Sciencix Part #	Comparable OEM Part #
Deuterium Lamp (2000 hr)	CTS-10465	WAS081142
Lamp Housing Window Assembly	CTS-A11184	WAS081341
Deuterium Lamp (2000 hr)	CTS-10268	WAT052586
 2487/2488 996/2996 		



Waters[®] Alliance[®] Series PM Kits

DESCRIPTION	Sciencix Part #	Comparable OEM Part #
2690/2695 Performance Maintenance Kit	CTS-10988	WAT270944
Alliance® Sovlent Manager Performance Maintenance Kit	CTS-11080	700004073
2690/2695 Performance Maintenance Kit w/Seal Pack	CTS-11380	

Disclaimer: Sciencix, Inc. is not affiliated with or an authorized distributor/official partner of Original Equipment Manufacturer (OEM) parts. All OEM part names, numbers, brands, registered trademarks or trademarks are the sole property of their respective owners in the United States or any other country. All OEM part names, numbers, and company names used in this catalog are for cross reference purposes only. Use of these names and brands does not imply endorsement. Statements made by Sciencix regarding the quality and performance of our products are our own and should not be construed as an indication of OEM acceptance or approval. © **2021 Sciencix, Inc.**