



A 3-IN-ONE  
SOLUTION







## 24 Well Plate Manifold System

**Part Number: VMF24WP**


- System includes collection plate, extraction plate, manifold and vacuum gauge waste base.
- Extraction plate can accommodate (24) 6 mL SPE columns.
- UCT SPE columns can be packed in 500 mg, 1 and 2 gram formats.
- Collection plate will accommodate 16x100mm and 16x125mm test tubes.



## 48 Well Plate Manifold System

**Part Number: VMF48WP**

- System includes collection plate, extraction plate, manifold and vacuum gauge waste base.
- Extraction plate can accommodate (48) 4.5 mL SPE columns in a staggered array.
- Elute samples directly into 12 x 75 mm test tube or 2 mL analytical vials.



## 96 Deep Well Plate Manifold System

**Part Number: VMF96WP**

- System includes manifold and vacuum gauge waste base.
- Well plates are separate and not included with system.
- Accommodate most standard 96 well plates.
- Collect directly into microplate.



## Universal Vacuum Manifold

**Part Number: VMFUVWP**

- A manifold system that accommodate a 24 well, 48 and 96 well plate configuration.
- System includes collection plate, extraction plate, manifold and vacuum gauge waste base.
- Extraction and collection plates for 24 well and 48 well configurations are included.
- Erlenmeyer flask and tubing are included.




24 Well Collection Plate  
VMF24CP



24 Well Extraction Plate  
VMF24EP




48 Well Collection Plate  
VMF48CP



48 Well Extraction Plate  
VMF48EP



Base Gasket  
VMFUEG




Lid Gasket  
VMFUVNG




Manifold Top  
VMFUVMT



Manifold Base  
VMFUVMB



Half Inch Spacer  
VMFUV05SP



One Inch Spacer  
VMFUV10SP



Collection Plate Riser  
VMFUVR1



Barbed Hose Adaptor  
VMFUVST

## PRICES AND TERMS

Our prices are subject to change without notice. The price in effect when we receive your order will apply. All prices are in US Dollars and are F.O.B. Lewistown, PA 17044.

Terms of payment are net 30 days.

## MINIMUM ORDERS

We welcome all orders, therefore, we do not have a minimum order requirement. When ordering, please include your purchase order number, complete "Ship To" and "Bill To" address, catalog number, quantity, and description of product(s). Also include your name and a phone number where you can be reached should we have any questions concerning your order.

## SHIPMENTS

Normal processing is within 24 hours after receipt of an order. Unless special shipping requests have been made, our trained staff will send all orders by UPS Ground service. The appropriate shipping charges (freight & insurance costs) will be added to the invoice, unless otherwise instructed by the customer.

## SPECIAL PRICING

We offer special pricing for volume purchases and standing orders. These discounts apply to bonded phase extraction column purchases only. Please call a sales representative for more information on special pricing qualifications.

## RETURN POLICY

Our Quality Manager will handle all returns. Before returning merchandise, please call to obtain a return authorization number from the quality manager. We will need to know the reason for the return, date of purchase, purchase order number and invoice number in order to issue a return authorization number. Return merchandise must be received before a credit can be issued. Returns will not be accepted after 90 days. A restocking fee of 25% of the price paid, or a minimum of \$25.00 (whichever is greater) will be charged on all returns.

## WARRANTY

All products manufactured by UCT are guaranteed against defects in materials and workmanship for a period of 90 days after shipment. UCT will replace any items that prove to be defective during this time period.

The exclusive remedy requires the end user to first advise UCT of the defective product by phone or in writing. Secondly, the defective product must be returned within 30 days after proper approval from our Quality Manager. All returns must indicate the purchase order number, the lot number and the shipping date. UCT's total liability is limited to the replacement cost of UCT products.

This warranty does not apply to damage resulting from misuse.

## Contact Us

Phone: 215.781.9255  
800.385.3153  
Fax: 215.785.1226

UCT, Inc.  
2731 Bartram Rd.  
Bristol, PA 19007

Email: [info@unitedchem.com](mailto:info@unitedchem.com)  
Web: [www.unitedchem.com](http://www.unitedchem.com)



4101-26-04

