



Micro Valves for GC and LC

- Convenient panel mount
- Compact 1" design
- 200 psi, .060" bore
- Variety of configurations

Simplify your liquid or gas handling application with a VICI Precision Sampling Micro valve. The unique design of the fitting detail allows a leak-free seal with no potential for rotor damage from overtightening. Internal parts are PEEK and PTFE.



(Fittings not included)

"T" flow path

3 ports 660100 660110 4 ports





180° flow path

660200 2 ports 4 ports 660210





90° flow path

2 ports 660300 3 ports 660310 4 ports 660320







North America, South America, and Australia/Oceania contact:

